

Surrogate Summary

SDG No.: **P4103**

Client: **ENTACT**

Analytical Method: **8082A**

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Rec	Qual	Limits	
								Low	High
I.BLK-PO106313.D	PIBLK-PO106313.D	Tetrachloro-m-xylene	1	20	17.3	86		70 (60)	130 (140)
		Decachlorobiphenyl	1	20	18.4	92		70 (60)	130 (140)
		Tetrachloro-m-xylene	2	20	17.0	85		70 (60)	130 (140)
		Decachlorobiphenyl	2	20	18.1	90		70 (60)	130 (140)
I.BLK-PO106693.D	PIBLK-PO106693.D	Tetrachloro-m-xylene	1	20	24.8	124		70 (60)	130 (140)
		Decachlorobiphenyl	1	20	23.6	118		70 (60)	130 (140)
		Tetrachloro-m-xylene	2	20	22.2	111		70 (60)	130 (140)
		Decachlorobiphenyl	2	20	24.4	122		70 (60)	130 (140)
PB163523BL	PB163523BL	Tetrachloro-m-xylene	1	20	25.8	129		30 (32)	150 (144)
		Decachlorobiphenyl	1	20	24.2	121		30 (32)	150 (175)
		Tetrachloro-m-xylene	2	20	22.7	113		30 (32)	150 (144)
		Decachlorobiphenyl	2	20	25.0	125		30 (32)	150 (175)
PB163523BS	PB163523BS	Tetrachloro-m-xylene	1	20	24.7	124		30 (32)	150 (144)
		Decachlorobiphenyl	1	20	23.4	117		30 (32)	150 (175)
		Tetrachloro-m-xylene	2	20	21.0	105		30 (32)	150 (144)
		Decachlorobiphenyl	2	20	24.2	121		30 (32)	150 (175)
P4103-01	WC-22A	Tetrachloro-m-xylene	1	20	10.3	51		30 (32)	150 (144)
		Decachlorobiphenyl	1	20	14.2	71		30 (32)	150 (175)
		Tetrachloro-m-xylene	2	20	11.1	55		30 (32)	150 (144)
		Decachlorobiphenyl	2	20	12.6	63		30 (32)	150 (175)
P4103-04	WC-14-5-7.5A	Tetrachloro-m-xylene	1	20	11.4	57		30 (32)	150 (144)
		Decachlorobiphenyl	1	20	17.1	85		30 (32)	150 (175)
		Tetrachloro-m-xylene	2	20	13.2	66		30 (32)	150 (144)
		Decachlorobiphenyl	2	20	17.1	85		30 (32)	150 (175)
P4103-07	WC-14-5-7.5B	Tetrachloro-m-xylene	1	20	12.2	61		30 (32)	150 (144)
		Decachlorobiphenyl	1	20	19.2	96		30 (32)	150 (175)
		Tetrachloro-m-xylene	2	20	13.6	68		30 (32)	150 (144)
		Decachlorobiphenyl	2	20	18.6	93		30 (32)	150 (175)
P4103-07MS	WC-14-5-7.5BMS	Tetrachloro-m-xylene	1	20	13.2	66		30 (32)	150 (144)
		Decachlorobiphenyl	1	20	20.5	103		30 (32)	150 (175)
		Tetrachloro-m-xylene	2	20	14.1	70		30 (32)	150 (144)
		Decachlorobiphenyl	2	20	20.3	101		30 (32)	150 (175)
P4103-07MSD	WC-14-5-7.5BMSD	Tetrachloro-m-xylene	1	20	13.3	66		30 (32)	150 (144)
		Decachlorobiphenyl	1	20	20.9	105		30 (32)	150 (175)
		Tetrachloro-m-xylene	2	20	14.0	70		30 (32)	150 (144)
		Decachlorobiphenyl	2	20	20.3	101		30 (32)	150 (175)
P4103-10	WC-14-7.5-10A	Tetrachloro-m-xylene	1	20	10.7	54		30 (32)	150 (144)
		Decachlorobiphenyl	1	20	18.8	94		30 (32)	150 (175)
		Tetrachloro-m-xylene	2	20	11.2	56		30 (32)	150 (144)
		Decachlorobiphenyl	2	20	18.6	93		30 (32)	150 (175)
P4103-13	WC-14-7.5-10B	Tetrachloro-m-xylene	1	20	14.9	74		30 (32)	150 (144)

() = LABORATORY INHOUSE LIMIT

Surrogate Summary

SDG No.: **P4103**

Client: **ENTACT**

Analytical Method: **8082A**

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Rec	Qual	Limits	
								Low	High
P4103-13	WC-14-7.5-10B	Decachlorobiphenyl	1	20	18.2	91		30 (32)	150 (175)
		Tetrachloro-m-xylene	2	20	16.6	83		30 (32)	150 (144)
I.BLK-PO106708.D	PIBLK-PO106708.D	Decachlorobiphenyl	2	20	18.0	90		30 (32)	150 (175)
		Tetrachloro-m-xylene	1	20	24.9	125		70 (60)	130 (140)
		Decachlorobiphenyl	1	20	24.1	120		70 (60)	130 (140)
		Tetrachloro-m-xylene	2	20	22.6	113		70 (60)	130 (140)
		Decachlorobiphenyl	2	20	25.1	126		70 (60)	130 (140)
P4103-16	WC-23	Tetrachloro-m-xylene	1	20	14.4	72		30 (32)	150 (144)
		Decachlorobiphenyl	1	20	18.6	93		30 (32)	150 (175)
		Tetrachloro-m-xylene	2	20	16.9	85		30 (32)	150 (144)
		Decachlorobiphenyl	2	20	19.1	95		30 (32)	150 (175)
I.BLK-PO106715.D	PIBLK-PO106715.D	Tetrachloro-m-xylene	1	20	25.9	130		70 (60)	130 (140)
		Decachlorobiphenyl	1	20	24.7	123		70 (60)	130 (140)
		Tetrachloro-m-xylene	2	20	22.8	114		70 (60)	130 (140)
		Decachlorobiphenyl	2	20	25.5	128		70 (60)	130 (140)