

**Instrument ID:** CV1

**Daily Analysis Runlog For Sequence/QC Batch ID # LB133324**

|                  |                                                 |              |                      |
|------------------|-------------------------------------------------|--------------|----------------------|
| Review By        | Sarabjit Jaswal                                 | Review On    | 11/8/2024 9:50:52 PM |
| Supervise By     | Mohan Bera                                      | Supervise On | 11/8/2024 9:52:01 PM |
| <b>STD. NAME</b> | <b>STD REF.#</b>                                |              |                      |
| ICAL Standard    | MP83107,MP83108,MP83109,MP83110,MP83111,MP83112 |              |                      |
| ICV Standard     | MP83113                                         |              |                      |
| CCV Standard     | MP83115                                         |              |                      |
| ICSA Standard    |                                                 |              |                      |
| CRI Standard     |                                                 |              |                      |
| LCS Standard     |                                                 |              |                      |
| Chk Standard     | MP83114,MP83116,MP83120                         |              |                      |

| Sr# | SampleId | ClientID | QcType | Date           | Comment | Operator | Status |
|-----|----------|----------|--------|----------------|---------|----------|--------|
| 1   | S0       | S0       | CAL1   | 11/07/24 08:03 |         | Mohan    | OK     |
| 2   | S0.2     | S01      | CAL2   | 11/07/24 08:05 |         | Mohan    | OK     |
| 3   | S2.5     | S02      | CAL3   | 11/07/24 08:07 |         | Mohan    | OK     |
| 4   | S5       | S03      | CAL4   | 11/07/24 08:10 |         | Mohan    | OK     |
| 5   | S7.5     | S04      | CAL5   | 11/07/24 08:12 |         | Mohan    | OK     |
| 6   | S10      | S05      | CAL6   | 11/07/24 08:23 |         | Mohan    | OK     |
| 7   | ICV072   | ICV072   | ICV    | 11/07/24 08:26 |         | Mohan    | OK     |
| 8   | ICB072   | ICB072   | ICB    | 11/07/24 08:28 |         | Mohan    | OK     |
| 9   | CCV040   | CCV040   | CCV    | 11/07/24 08:30 |         | Mohan    | OK     |
| 10  | CCB040   | CCB040   | CCB    | 11/07/24 08:33 |         | Mohan    | OK     |
| 11  | P4500-01 | MBHCN1   | SAM    | 11/07/24 08:37 |         | Mohan    | OK     |
| 12  | P4500-02 | MBHCN2   | SAM    | 11/07/24 08:39 |         | Mohan    | OK     |
| 13  | P4500-03 | MBHCN3   | SAM    | 11/07/24 08:42 |         | Mohan    | OK     |
| 14  | P4500-04 | MBHCN4   | SAM    | 11/07/24 08:44 |         | Mohan    | OK     |
| 15  | P4500-05 | MBHCN5   | SAM    | 11/07/24 08:46 |         | Mohan    | OK     |
| 16  | P4500-06 | MBHCN6   | SAM    | 11/07/24 08:49 |         | Mohan    | OK     |
| 17  | P4500-07 | MBHCN7   | SAM    | 11/07/24 08:51 |         | Mohan    | OK     |
| 18  | P4500-08 | MBHCN8   | SAM    | 11/07/24 08:53 |         | Mohan    | OK     |

Instrument ID: CV1

**Daily Analysis Runlog For Sequence/QC Batch ID # LB133324**

|                  |                                                 |              |                      |
|------------------|-------------------------------------------------|--------------|----------------------|
| Review By        | Sarabjit Jaswal                                 | Review On    | 11/8/2024 9:50:52 PM |
| Supervise By     | Mohan Bera                                      | Supervise On | 11/8/2024 9:52:01 PM |
| <b>STD. NAME</b> | <b>STD REF.#</b>                                |              |                      |
| ICAL Standard    | MP83107,MP83108,MP83109,MP83110,MP83111,MP83112 |              |                      |
| ICV Standard     | MP83113                                         |              |                      |
| CCV Standard     | MP83115                                         |              |                      |
| ICSA Standard    |                                                 |              |                      |
| CRI Standard     |                                                 |              |                      |
| LCS Standard     |                                                 |              |                      |
| Chk Standard     | MP83114,MP83116,MP83120                         |              |                      |

|    |            |         |     |                |  |       |    |
|----|------------|---------|-----|----------------|--|-------|----|
| 19 | P4500-09   | MBHCN9  | SAM | 11/07/24 08:55 |  | Mohan | OK |
| 20 | P4500-10   | MBHCP0  | SAM | 11/07/24 08:58 |  | Mohan | OK |
| 21 | P4500-11   | MBHCP0D | DUP | 11/07/24 09:00 |  | Mohan | OK |
| 22 | P4500-12   | MBHCP0S | MS  | 11/07/24 09:02 |  | Mohan | OK |
| 23 | P4500-13   | MBHCP1  | SAM | 11/07/24 09:04 |  | Mohan | OK |
| 24 | P4500-14   | MBHCP2  | SAM | 11/07/24 09:07 |  | Mohan | OK |
| 25 | P4500-15   | MBHCP3  | SAM | 11/07/24 09:09 |  | Mohan | OK |
| 26 | P4500-16   | MBHCP4  | SAM | 11/07/24 09:11 |  | Mohan | OK |
| 27 | P4500-17   | MBHCP5  | SAM | 11/07/24 09:14 |  | Mohan | OK |
| 28 | P4500-18   | MBHCP6  | SAM | 11/07/24 09:16 |  | Mohan | OK |
| 29 | P4500-19   | MBHCP7  | SAM | 11/07/24 09:18 |  | Mohan | OK |
| 30 | CCV041     | CCV041  | CCV | 11/07/24 09:20 |  | Mohan | OK |
| 31 | CCB041     | CCB041  | CCB | 11/07/24 09:23 |  | Mohan | OK |
| 32 | P4500-20   | MBHCX2  | SAM | 11/07/24 09:25 |  | Mohan | OK |
| 33 | P4500-21   | MBHCX3  | SAM | 11/07/24 09:27 |  | Mohan | OK |
| 34 | P4500-22   | MBHCZ6  | SAM | 11/07/24 09:30 |  | Mohan | OK |
| 35 | PB164743BL | PBS743  | MB  | 11/07/24 09:32 |  | Mohan | OK |
| 36 | P4501-01   | MBHCS8  | SAM | 11/07/24 09:34 |  | Mohan | OK |
| 37 | P4501-02   | MBHCS9  | SAM | 11/07/24 09:36 |  | Mohan | OK |
| 38 | P4501-03   | MBHCT0  | SAM | 11/07/24 09:39 |  | Mohan | OK |

Instrument ID: CV1

**Daily Analysis Runlog For Sequence/QC Batch ID # LB133324**

|                  |                                                 |              |                      |
|------------------|-------------------------------------------------|--------------|----------------------|
| Review By        | Sarabjit Jaswal                                 | Review On    | 11/8/2024 9:50:52 PM |
| Supervise By     | Mohan Bera                                      | Supervise On | 11/8/2024 9:52:01 PM |
| <b>STD. NAME</b> | <b>STD REF.#</b>                                |              |                      |
| ICAL Standard    | MP83107,MP83108,MP83109,MP83110,MP83111,MP83112 |              |                      |
| ICV Standard     | MP83113                                         |              |                      |
| CCV Standard     | MP83115                                         |              |                      |
| ICSA Standard    |                                                 |              |                      |
| CRI Standard     |                                                 |              |                      |
| LCS Standard     |                                                 |              |                      |
| Chk Standard     | MP83114,MP83116,MP83120                         |              |                      |

|    |          |         |     |                |         |       |          |
|----|----------|---------|-----|----------------|---------|-------|----------|
| 39 | P4501-04 | MBHCT1  | SAM | 11/07/24 09:41 |         | Mohan | OK       |
| 40 | P4501-05 | MBHCT2  | SAM | 11/07/24 09:43 |         | Mohan | OK       |
| 41 | P4501-06 | MBHCT3  | SAM | 11/07/24 09:45 |         | Mohan | OK       |
| 42 | P4501-07 | MBHCT4  | SAM | 11/07/24 09:48 |         | Mohan | OK       |
| 43 | P4501-08 | MBHCT5  | SAM | 11/07/24 09:50 | HG High | Mohan | Dilution |
| 44 | P4501-09 | MBHCT6  | SAM | 11/07/24 09:52 |         | Mohan | OK       |
| 45 | P4501-10 | MBHCT7  | SAM | 11/07/24 09:54 |         | Mohan | OK       |
| 46 | P4501-11 | MBHCT8  | SAM | 11/07/24 09:57 |         | Mohan | OK       |
| 47 | P4501-12 | MBHCT9  | SAM | 11/07/24 09:59 |         | Mohan | OK       |
| 48 | P4501-13 | MBHCW0  | SAM | 11/07/24 10:01 |         | Mohan | OK       |
| 49 | P4501-14 | MBHCW0D | DUP | 11/07/24 10:04 |         | Mohan | OK       |
| 50 | P4501-15 | MBHCW0S | MS  | 11/07/24 10:06 |         | Mohan | OK       |
| 51 | P4501-16 | MBHCW1  | SAM | 11/07/24 10:08 |         | Mohan | OK       |
| 52 | CCV042   | CCV042  | CCV | 11/07/24 10:10 |         | Mohan | OK       |
| 53 | CCB042   | CCB042  | CCB | 11/07/24 10:13 |         | Mohan | OK       |
| 54 | P4501-17 | MBHCW2  | SAM | 11/07/24 10:15 |         | Mohan | OK       |
| 55 | P4501-18 | MBHCW3  | SAM | 11/07/24 10:17 |         | Mohan | OK       |
| 56 | P4501-19 | MBHCW4  | SAM | 11/07/24 10:19 |         | Mohan | OK       |
| 57 | P4501-20 | MBHCW5  | SAM | 11/07/24 10:22 |         | Mohan | OK       |
| 58 | P4501-21 | MBHCW6  | SAM | 11/07/24 10:24 |         | Mohan | OK       |

Instrument ID: CV1

**Daily Analysis Runlog For Sequence/QC Batch ID # LB133324**

|              |                 |              |                      |
|--------------|-----------------|--------------|----------------------|
| Review By    | Sarabjit Jaswal | Review On    | 11/8/2024 9:50:52 PM |
| Supervise By | Mohan Bera      | Supervise On | 11/8/2024 9:52:01 PM |

| STD. NAME     | STD REF.#                                       |
|---------------|-------------------------------------------------|
| ICAL Standard | MP83107,MP83108,MP83109,MP83110,MP83111,MP83112 |
| ICV Standard  | MP83113                                         |
| CCV Standard  | MP83115                                         |
| ICSA Standard |                                                 |
| CRI Standard  |                                                 |
| LCS Standard  |                                                 |
| Chk Standard  | MP83114,MP83116,MP83120                         |

|    |            |        |     |                |  |       |    |
|----|------------|--------|-----|----------------|--|-------|----|
| 59 | P4501-22   | MBHCW7 | SAM | 11/07/24 10:26 |  | Mohan | OK |
| 60 | PB164742BL | PBS742 | MB  | 11/07/24 10:28 |  | Mohan | OK |
| 61 | CCV043     | CCV043 | CCV | 11/07/24 10:31 |  | Mohan | OK |
| 62 | CCB043     | CCB043 | CCB | 11/07/24 10:33 |  | Mohan | OK |

**Instrument ID:** CV1

**Daily Analysis Runlog For Sequence/QC Batch ID # LB133356**

|              |                 |              |                      |
|--------------|-----------------|--------------|----------------------|
| Review By    | Sarabjit Jaswal | Review On    | 11/8/2024 9:51:01 PM |
| Supervise By | Mohan Bera      | Supervise On | 11/8/2024 9:51:52 PM |

| STD. NAME     | STD REF.#                                       |
|---------------|-------------------------------------------------|
| ICAL Standard | MP83162,MP83163,MP83164,MP83165,MP83166,MP83167 |
| ICV Standard  | MP83168                                         |
| CCV Standard  | MP83170                                         |
| ICSA Standard |                                                 |
| CRI Standard  |                                                 |
| LCS Standard  |                                                 |
| Chk Standard  | MP83169,MP83171,MP83159                         |

| Sr# | SampleId | ClientID | QcType | Date           | Comment | Operator | Status |
|-----|----------|----------|--------|----------------|---------|----------|--------|
| 1   | S0       | S0       | CAL1   | 11/08/24 13:14 |         | Mohan    | OK     |
| 2   | S0.2     | S01      | CAL2   | 11/08/24 13:19 |         | Mohan    | OK     |
| 3   | S2.5     | S02      | CAL3   | 11/08/24 13:22 |         | Mohan    | OK     |
| 4   | S5       | S03      | CAL4   | 11/08/24 13:26 |         | Mohan    | OK     |
| 5   | S7.5     | S04      | CAL5   | 11/08/24 13:28 |         | Mohan    | OK     |
| 6   | S10      | S05      | CAL6   | 11/08/24 13:36 |         | Mohan    | OK     |
| 7   | ICV076   | ICV076   | ICV    | 11/08/24 13:38 |         | Mohan    | OK     |
| 8   | ICB076   | ICB076   | ICB    | 11/08/24 13:43 |         | Mohan    | OK     |
| 9   | CCV054   | CCV054   | CCV    | 11/08/24 13:45 |         | Mohan    | OK     |
| 10  | CCB054   | CCB054   | CCB    | 11/08/24 13:52 |         | Mohan    | OK     |
| 11  | P4496-04 | MBHDO5   | SAM    | 11/08/24 13:57 |         | Mohan    | OK     |
| 12  | P4496-05 | MBHDO6   | SAM    | 11/08/24 13:59 |         | Mohan    | OK     |
| 13  | P4496-06 | MBHZC7   | SAM    | 11/08/24 14:01 |         | Mohan    | OK     |
| 14  | P4496-07 | MBHZC8   | SAM    | 11/08/24 14:04 |         | Mohan    | OK     |
| 15  | P4496-08 | MBHZC8D  | DUP    | 11/08/24 14:06 |         | Mohan    | OK     |
| 16  | P4624-01 | MBH6M7   | SAM    | 11/08/24 14:10 |         | Mohan    | OK     |
| 17  | P4624-02 | MBH6M8   | SAM    | 11/08/24 14:13 |         | Mohan    | OK     |
| 18  | P4624-03 | MBH6M9   | SAM    | 11/08/24 14:15 |         | Mohan    | OK     |

Instrument ID: CV1

**Daily Analysis Runlog For Sequence/QC Batch ID # LB133356**

|                  |                                                 |              |                      |
|------------------|-------------------------------------------------|--------------|----------------------|
| Review By        | Sarabjit Jaswal                                 | Review On    | 11/8/2024 9:51:01 PM |
| Supervise By     | Mohan Bera                                      | Supervise On | 11/8/2024 9:51:52 PM |
| <b>STD. NAME</b> | <b>STD REF.#</b>                                |              |                      |
| ICAL Standard    | MP83162,MP83163,MP83164,MP83165,MP83166,MP83167 |              |                      |
| ICV Standard     | MP83168                                         |              |                      |
| CCV Standard     | MP83170                                         |              |                      |
| ICSA Standard    |                                                 |              |                      |
| CRI Standard     |                                                 |              |                      |
| LCS Standard     |                                                 |              |                      |
| Chk Standard     | MP83169,MP83171,MP83159                         |              |                      |

|    |            |         |     |                |           |       |          |
|----|------------|---------|-----|----------------|-----------|-------|----------|
| 19 | P4624-04   | MBH6N0  | SAM | 11/08/24 14:17 |           | Mohan | OK       |
| 20 | P4624-05   | MBH6N1  | SAM | 11/08/24 14:19 |           | Mohan | OK       |
| 21 | P4624-06   | MBH8M6  | SAM | 11/08/24 14:22 |           | Mohan | OK       |
| 22 | P4624-07   | MBH6N2  | SAM | 11/08/24 14:24 |           | Mohan | OK       |
| 23 | P4624-08   | MBH6N3  | SAM | 11/08/24 14:26 |           | Mohan | OK       |
| 24 | P4624-09   | MBH6N4  | SAM | 11/08/24 14:28 |           | Mohan | OK       |
| 25 | P4624-10   | MBH6N5  | SAM | 11/08/24 14:31 |           | Mohan | OK       |
| 26 | P4624-11   | MBH6N6  | SAM | 11/08/24 14:33 |           | Mohan | OK       |
| 27 | P4624-12   | MBH6N6D | DUP | 11/08/24 14:35 |           | Mohan | OK       |
| 28 | P4624-13   | MBH6N6S | MS  | 11/08/24 14:38 |           | Mohan | OK       |
| 29 | CCV055     | CCV055  | CCV | 11/08/24 14:40 |           | Mohan | OK       |
| 30 | CCB055     | CCB055  | CCB | 11/08/24 14:44 |           | Mohan | OK       |
| 31 | P4624-14   | MBH6N7  | SAM | 11/08/24 14:46 |           | Mohan | OK       |
| 32 | PB164799BL | PBS799  | MB  | 11/08/24 14:49 |           | Mohan | OK       |
| 33 | P4501-08DL | MBHCT5  | SAM | 11/08/24 14:51 | 2X for Hg | Mohan | Confirms |
| 34 | P4496-09   | MBHZC8S | MS  | 11/08/24 14:56 |           | Mohan | OK       |
| 35 | CCV056     | CCV056  | CCV | 11/08/24 14:58 |           | Mohan | OK       |
| 36 | CCB056     | CCB056  | CCB | 11/08/24 15:17 |           | Mohan | OK       |