

Surrogate Summary

SDG No.: **P4602**

Client: **Tetra Tech, EMI**

Analytical Method: **8151A**

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Rec	Qual	Limits	
								Low	High
I.BLK-PS028217.D	PIBLK-PS028217.D	2,4-DCAA	1	500	506	101		39	175
		2,4-DCAA	2	500	449	90		39	175
I.BLK-PS028224.D	PIBLK-PS028224.D	2,4-DCAA	1	500	486	97		39	175
		2,4-DCAA	2	500	559	112		39	175
PB164641BL	PB164641BL	2,4-DCAA	1	500	572	114		10	141
		2,4-DCAA	2	500	568	114		10	141
PB164641BS	PB164641BS	2,4-DCAA	1	500	573	115		10	141
		2,4-DCAA	2	500	572	114		10	141
I.BLK-PS028236.D	PIBLK-PS028236.D	2,4-DCAA	1	500	573	115		39	175
		2,4-DCAA	2	500	538	108		39	175
I.BLK-PS028252.D	PIBLK-PS028252.D	2,4-DCAA	1	500	489	98		39	175
		2,4-DCAA	2	500	491	98		39	175
I.BLK-PS028271.D	PIBLK-PS028271.D	2,4-DCAA	1	500	519	104		39	175
		2,4-DCAA	2	500	516	103		39	175
PB164703BL	PB164703BL	2,4-DCAA	1	500	556	111		10	141
		2,4-DCAA	2	500	520	104		10	141
I.BLK-PS028283.D	PIBLK-PS028283.D	2,4-DCAA	1	500	529	106		39	175
		2,4-DCAA	2	500	516	103		39	175
PB164703BS	PB164703BS	2,4-DCAA	1	500	560	112		10	141
		2,4-DCAA	2	500	535	107		10	141
I.BLK-PS028296.D	PIBLK-PS028296.D	2,4-DCAA	1	500	530	106		39	175
		2,4-DCAA	2	500	508	102		39	175
I.BLK-PS028355.D	PIBLK-PS028355.D	2,4-DCAA	1	500	501	100		39	175
		2,4-DCAA	2	500	503	101		39	175
I.BLK-PS028405.D	PIBLK-PS028405.D	2,4-DCAA	1	500	511	102		39	175
		2,4-DCAA	2	500	472	94		39	175
P4602-43	CC0R8	2,4-DCAA	1	500	301	60		10	141
		2,4-DCAA	2	500	237	47		10	141
P4602-44MS	CC0R8MS	2,4-DCAA	1	500	375	75		10	141
		2,4-DCAA	2	500	253	51		10	141
P4602-45MSD	CC0R8MSD	2,4-DCAA	1	500	379	76		10	141
		2,4-DCAA	2	500	254	51		10	141
P4602-40	CC0R6	2,4-DCAA	1	500	285	57		10	141
		2,4-DCAA	2	500	200	40		10	141
P4602-41MS	CC0R6MS	2,4-DCAA	1	500	390	78		10	141
		2,4-DCAA	2	500	186	37		10	141
P4602-42MSD	CC0R6MSD	2,4-DCAA	1	500	396	79		10	141
		2,4-DCAA	2	500	188	38		10	141
PB164725BL	PB164725BL	2,4-DCAA	1	500	558	112		10	141
		2,4-DCAA	2	500	487	97		10	141
PB164725BS	PB164725BS	2,4-DCAA	1	500	555	111		10	141

Surrogate Summary

SDG No.: **P4602**

Client: **Tetra Tech, EMI**

Analytical Method: **8151A**

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Rec	Qual	Limits	
								Low	High
PB164725BS	PB164725BS	2,4-DCAA	2	500	463	93		10	141
P4602-19	C0P17	2,4-DCAA	1	500	207	41		10	141
		2,4-DCAA	2	500	169	34		10	141
I.BLK-PS028417.D	PIBLK-PS028417.D	2,4-DCAA	1	500	513	103		39	175
		2,4-DCAA	2	500	468	94		39	175
P4602-21	CC0P2	2,4-DCAA	1	500	206	41		10	141
		2,4-DCAA	2	500	168	34		10	141
P4602-22	CC0P4	2,4-DCAA	1	500	216	43		10	141
		2,4-DCAA	2	500	139	28		10	141
P4602-23	CC0P6	2,4-DCAA	1	500	164	33		10	141
		2,4-DCAA	2	500	125	25		10	141
P4602-24	CC0P8	2,4-DCAA	1	500	129	26		10	141
		2,4-DCAA	2	500	105	21		10	141
P4602-25	CC0Q0	2,4-DCAA	1	500	182	36		10	141
		2,4-DCAA	2	500	168	34		10	141
P4602-26	CC0Q2	2,4-DCAA	1	500	233	47		10	141
		2,4-DCAA	2	500	205	41		10	141
P4602-27	CC0Q3	2,4-DCAA	1	500	166	33		10	141
		2,4-DCAA	2	500	140	28		10	141
P4602-29	CC0Q5	2,4-DCAA	1	500	131	26		10	141
		2,4-DCAA	2	500	192	38		10	141
P4602-30	CC0Q7	2,4-DCAA	1	500	161	32		10	141
		2,4-DCAA	2	500	123	25		10	141
I.BLK-PS028429.D	PIBLK-PS028429.D	2,4-DCAA	1	500	508	102		39	175
		2,4-DCAA	2	500	461	92		39	175
P4602-33	CC0R1	2,4-DCAA	1	500	179	36		10	141
		2,4-DCAA	2	500	116	23		10	141
P4602-34	CC0R2	2,4-DCAA	1	500	133	27		10	141
		2,4-DCAA	2	500	110	22		10	141
P4602-31	CC0Q9	2,4-DCAA	1	500	303	61		10	141
		2,4-DCAA	2	500	273	55		10	141
P4602-32	CC0R0	2,4-DCAA	1	500	257	51		10	141
		2,4-DCAA	2	500	255	51		10	141
P4602-37	CC0P0	2,4-DCAA	1	500	156	31		10	141
		2,4-DCAA	2	500	132	26		10	141
P4602-38MS	CC0P0MS	2,4-DCAA	1	500	138	28		10	141
		2,4-DCAA	2	500	108	22		10	141
P4602-39MSD	CC0P0MSD	2,4-DCAA	1	500	137	27		10	141
		2,4-DCAA	2	500	108	22		10	141
I.BLK-PS028441.D	PIBLK-PS028441.D	2,4-DCAA	1	500	509	102		39	175
		2,4-DCAA	2	500	455	91		39	175

Surrogate Summary

SDG No.: P4602

Client: Tetra Tech, EMI

Analytical Method: 8151A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Rec	Qual	Limits	
								Low	High
P4602-28	CC0Q4	2,4-DCAA	1	500	304	61		10	141
		2,4-DCAA	2	500	263	53		10	141
I.BLK-PS028447.D	PIBLK-PS028447.D	2,4-DCAA	1	500	497	99		39	175
		2,4-DCAA	2	500	446	89		39	175