

Ident	Instrument IC-1				Analyst: NF		Method: 300.0 / 9056A		date time	Initial \ Analyst	
	Con F-	Con CL-	Con NO2	Con BR-	Con NO3	Con HPO4	Con SO4	Method name			
STD1	0	0	0	0	0	0	0	0	IC1-101624	10/16/2024 10:55	10 NF/IZ
STD2	0.424	0.632	0.639	2.106	0.535	1.072	3.276	IC1-101624	10/16/2024 11:17	10 NF/IZ	
STD3	0.796	1.206	1.205	4.013	1	2	6.052	IC1-101624	10/16/2024 11:38	10 NF/IZ	
STD4	1.032	1.537	1.531	5.116	1.277	2.557	7.662	IC1-101624	10/16/2024 11:59	10 NF/IZ	
STD5	1.966	2.927	2.921	9.77	2.439	4.871	14.51	IC1-101624	10/16/2024 12:21	10 NF/IZ	
STD6	3.88	5.901	5.91	19.683	4.912	9.831	29.286	IC1-101624	10/16/2024 12:42	10 NF/IZ	
STD7	5.102	7.598	7.593	25.312	6.337	12.67	37.714	IC1-101624	10/16/2024 13:04	10 NF/IZ	
ICV	1.966	2.997	3	9.969	2.497	4.753	14.85	IC1-101624	10/16/2024 13:37	10 NF/IZ	
ICB	0	0	0	0	0	0	0	IC1-101624	10/16/2024 13:59	10 NF/IZ	
CCV	1.988	2.979	2.984	9.93	2.495	4.965	14.701	IC1-101624	11/4/2024 9:38	10 NF/IZ	
CCB	0	0	0	0	0	0	0	IC1-101624	11/4/2024 10:00	10 NF/IZ	
LB133251BLW	0	0	0	0	0	0	0	IC1-101624	11/4/2024 10:21	10 NF/IZ	
LB133251BSW	1.991	2.997	2.992	9.955	2.498	4.962	14.723	IC1-101624	11/4/2024 10:43	10 NF/IZ	
P4673-01	0.411	26.77	0.383	0.335	1.559	0.355	3257.035	IC1-101624	11/4/2024 11:26	10 NF/IZ	
P4673-02	1.198	13.1	0	0.308	1.343	0.726	3367.603	IC1-101624	11/4/2024 11:47	10 NF/IZ	
P4673-03	0.407	11.504	0	0	0.463	0.417	3428.48	IC1-101624	11/4/2024 12:09	10 NF/IZ	
P4673-04	0.488	14.938	0.152	0.297	2.225	0.385	3320.606	IC1-101624	11/4/2024 12:30	10 NF/IZ	
P4673-05	0.474	14.064	0	0.299	2.209	0.361	3234.733	IC1-101624	11/4/2024 13:13	10 NF/IZ	
P4673-06	0.525	15.104	0	0.269	1.651	0.441	3072.351	IC1-101624	11/4/2024 13:35	10 NF/IZ	
P4673-07	0.307	8.987	0	0	0.73	0.324	3255.556	IC1-101624	11/4/2024 13:56	10 NF/IZ	
P4674-01	0.242	0.312	0	0	0.372	0.324	4269.759	IC1-101624	11/4/2024 14:18	10 NF/IZ	
CCV	1.994	3.001	2.991	9.999	2.504	4.931	14.771	IC1-101624	11/4/2024 15:01	10 NF/IZ	
CCB	0	0	0	0	0	0	0	IC1-101624	11/4/2024 15:22	10 NF/IZ	
P4674-01MS	2.02	9.011	1.984	9.966	3.136	2.74	4092.952	IC1-101624	11/4/2024 15:44	10 NF/IZ	
P4674-01MSD	2.108	9.101	1.999	10.195	3.229	2.995	4099.017	IC1-101624	11/4/2024 16:05	10 NF/IZ	
CCV	1.971	2.991	2.997	10.007	2.506	4.882	14.803	IC1-101624	11/4/2024 16:48	10 NF/IZ	
CCB	0	0	0	0	0	0	0	IC1-101624	11/4/2024 17:09	10 NF/IZ	

Clear table


Instrument ID: IC-2 Analyst : IZ Method: 300.0 / 9056A

ident	concentra tion F-	concentratio n CL-	concentrati on NO2	concentrati on BR-	concentrati on NO3	concentratio n HPO4	concentrati on SO4	file name	date time	Initial wt/ Final	Analyst
STD1	0.000	0.000	0.000	0.000	0.000	0.000	0.000	IC1-101624	10/16/2024 10:55	10	NF/IZ
STD2	0.424	0.632	0.639	2.106	0.535	1.072	3.276	IC1-101624	10/16/2024 11:17	10	NF/IZ
STD3	0.796	1.206	1.205	4.013	1.000	2.000	6.052	IC1-101624	10/16/2024 11:38	10	NF/IZ
STD4	1.032	1.537	1.531	5.116	1.277	2.557	7.662	IC1-101624	10/16/2024 11:59	10	NF/IZ
STD5	1.966	2.927	2.921	9.77	2.439	4.871	14.51	IC1-101624	10/16/2024 12:21	10	NF/IZ
STD6	3.88	5.901	5.91	19.683	4.912	9.831	29.286	IC1-101624	10/16/2024 12:42	10	NF/IZ
STD7	5.102	7.598	7.593	25.312	6.337	12.67	37.714	IC1-101624	10/16/2024 13:04	10	NF/IZ

ident	True Value F-	True Value CL-	True Value NO2	True Value BR-	True Value NO3	True Value HPO4	True Value SO4
STD1	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000	0.0000
STD2	0.4000	0.6000	0.6000	2.0000	0.5000	1.0000	3.0000
STD3	0.8000	1.2000	1.2000	4.0000	1.0000	2.0000	6.0000
STD4	1.0000	1.5000	1.5000	5.0000	1.2500	2.5000	7.5000
STD5	2.0000	3.0000	3.0000	10.0000	2.5000	5.0000	15.0000
STD6	4.0000	6.0000	6.0000	20.0000	5.0000	10.0000	30.0000
STD7	5.0000	7.5000	7.5000	25.0000	6.2500	12.5000	37.0000

ident	Relative Error F-	Relative Error CL-	Relative Error NO2	Relative Error BR-	Relative Error NO3	Relative Error HPO4	Relative Error SO4
STD1							
STD2	6.0000	5.3333	6.5000	5.3000	7.0000	7.2000	9.2000
STD3	-0.5000	0.5000	0.4167	0.3250	0.0000	0.0000	0.8667
STD4	3.2000	2.4667	2.0667	2.3200	2.1600	2.2800	2.1600
STD5	-1.7000	-2.4333	-2.6333	-2.3000	-2.4400	-2.5800	-3.2667
STD6	-3.0000	-1.6500	-1.5000	-1.5850	-1.7600	-1.6900	-2.3800
STD7	2.0400	1.3067	1.2400	1.2480	1.3920	1.3600	1.9297

Fluoride (Anions)




Function: $A = -5.06114E-3 + 0.0146570 \times Q$

Relative standard deviation 3.824450 %

Correlation coefficient 0.999229

Sample type	Index	Conc.	Volume	Dilution	Sample amount	Area	Ident	Date	Used
Standard 1	1	0.000	10.0	1.0	1.0	n. d.	STD1	2024-10-16 10:55:46 UTC-4	used
Standard 2	1	0.400	10.0	1.0	1.0	0.057	STD2	2024-10-16 11:17:08 UTC-4	used
Standard 3	1	0.800	10.0	1.0	1.0	0.112	STD3	2024-10-16 11:38:32 UTC-4	used
Standard 4	1	1.000	10.0	1.0	1.0	0.146	STD4	2024-10-16 11:59:57 UTC-4	used
Standard 5	1	2.000	10.0	1.0	1.0	0.283	STD5	2024-10-16 12:21:22 UTC-4	used
Standard 6	1	4.000	10.0	1.0	1.0	0.564	STD6	2024-10-16 12:42:48 UTC-4	used
Standard 7	1	5.000	10.0	1.0	1.0	0.743	STD7	2024-10-16 13:04:15 UTC-4	used

Chloride (Anions)



Function: $A = -4.66155E-3 + 9.86442E-3 \times Q$


Relative standard deviation 2.532075 %

Correlation coefficient 0.999661

Sample type	Index	Conc.	Volume	Dilution	Sample amount	Area	Ident	Date	Used
Standard 1	1	0.000	10.0	1.0	1.0	n. d.	STD1	2024-10-16 10:55:46 UTC-4	used
Standard 2	1	0.600	10.0	1.0	1.0	0.058	STD2	2024-10-16 11:17:08 UTC-4	used
Standard 3	1	1.200	10.0	1.0	1.0	0.114	STD3	2024-10-16 11:38:32 UTC-4	used
Standard 4	1	1.500	10.0	1.0	1.0	0.147	STD4	2024-10-16 11:59:57 UTC-4	used
Standard 5	1	3.000	10.0	1.0	1.0	0.284	STD5	2024-10-16 12:21:22 UTC-4	used
Standard 6	1	6.000	10.0	1.0	1.0	0.577	STD6	2024-10-16 12:42:48 UTC-4	used
Standard 7	1	7.500	10.0	1.0	1.0	0.745	STD7	2024-10-16 13:04:15 UTC-4	used

Nitrite (Anions)

(μS/cm) x min

Function: $A = -0.0158639 + 0.0223869 \times Q$


Relative standard deviation 2.476730 %

Correlation coefficient 0.999680

Sample type	Index	Conc.	Volume	Dilution	Sample amount	Area	Ident	Date	Used
Standard 1	1	0.000	10.0	1.0	1.0	n. d.	STD1	2024-10-16 10:55:46 UTC-4	used
Standard 2	1	0.600	10.0	1.0	1.0	0.127	STD2	2024-10-16 11:17:08 UTC-4	used
Standard 3	1	1.200	10.0	1.0	1.0	0.254	STD3	2024-10-16 11:38:32 UTC-4	used
Standard 4	1	1.500	10.0	1.0	1.0	0.327	STD4	2024-10-16 11:59:57 UTC-4	used
Standard 5	1	3.000	10.0	1.0	1.0	0.638	STD5	2024-10-16 12:21:22 UTC-4	used
Standard 6	1	6.000	10.0	1.0	1.0	1.307	STD6	2024-10-16 12:42:48 UTC-4	used
Standard 7	1	7.500	10.0	1.0	1.0	1.684	STD7	2024-10-16 13:04:15 UTC-4	used

Bromide (Anions)

($\mu\text{S}/\text{cm}$) $\times$ min



Function: $A = -7.83486E-3 + 4.40426E-3 \times Q$


Relative standard deviation 2.425045 %

Correlation coefficient 0.999690

Sample type	Index	Conc.	Volume	Dilution	Sample amount	Area	Ident	Date	Used
Standard 1	1	0.000	10.0	1.0	1.0	n. d.	STD1	2024-10-16 10:55:46 UTC-4	used
Standard 2	1	2.000	10.0	1.0	1.0	0.085	STD2	2024-10-16 11:17:08 UTC-4	used
Standard 3	1	4.000	10.0	1.0	1.0	0.169	STD3	2024-10-16 11:38:32 UTC-4	used
Standard 4	1	5.000	10.0	1.0	1.0	0.217	STD4	2024-10-16 11:59:57 UTC-4	used
Standard 5	1	10.000	10.0	1.0	1.0	0.422	STD5	2024-10-16 12:21:22 UTC-4	used
Standard 6	1	20.000	10.0	1.0	1.0	0.859	STD6	2024-10-16 12:42:48 UTC-4	used
Standard 7	1	25.000	10.0	1.0	1.0	1.107	STD7	2024-10-16 13:04:15 UTC-4	used

Nitrate (Anions)

($\mu\text{S}/\text{cm}$) $\times$ min




Function: $A = -0.0159932 + 0.0247327 \times Q$

Relative standard deviation 2.693966 %

Correlation coefficient 0.999623

Sample type	Index	Conc.	Volume	Dilution	Sample amount	Area	Ident	Date	Used
Standard 1	1	0.000	10.0	1.0	1.0	n. d.	STD1	2024-10-16 10:55:46 UTC-4	used
Standard 2	1	0.500	10.0	1.0	1.0	0.116	STD2	2024-10-16 11:17:08 UTC-4	used
Standard 3	1	1.000	10.0	1.0	1.0	0.231	STD3	2024-10-16 11:38:32 UTC-4	used
Standard 4	1	1.250	10.0	1.0	1.0	0.300	STD4	2024-10-16 11:59:57 UTC-4	used
Standard 5	1	2.500	10.0	1.0	1.0	0.587	STD5	2024-10-16 12:21:22 UTC-4	used
Standard 6	1	5.000	10.0	1.0	1.0	1.199	STD6	2024-10-16 12:42:48 UTC-4	used
Standard 7	1	6.250	10.0	1.0	1.0	1.551	STD7	2024-10-16 13:04:15 UTC-4	used

Phosphate (Anions) $(\mu\text{S}/\text{cm}) \times \text{min}$ Function: $A = -9.96170E-3 + 8.78405E-3 \times Q$


Relative standard deviation 2.667639 %

Correlation coefficient 0.999628

Sample type	Index	Conc.	Volume	Dilution	Sample amount	Area	Ident	Date	Used
Standard 1	1	0.000	10.0	1.0	1.0	n. d.	STD1	2024-10-16 10:55:46 UTC-4	used
Standard 2	1	1.000	10.0	1.0	1.0	0.084	STD2	2024-10-16 11:17:08 UTC-4	used
Standard 3	1	2.000	10.0	1.0	1.0	0.166	STD3	2024-10-16 11:38:32 UTC-4	used
Standard 4	1	2.500	10.0	1.0	1.0	0.215	STD4	2024-10-16 11:59:57 UTC-4	used
Standard 5	1	5.000	10.0	1.0	1.0	0.418	STD5	2024-10-16 12:21:22 UTC-4	used
Standard 6	1	10.000	10.0	1.0	1.0	0.854	STD6	2024-10-16 12:42:48 UTC-4	used
Standard 7	1	12.500	10.0	1.0	1.0	1.103	STD7	2024-10-16 13:04:15 UTC-4	used

Sulfate (Anions)

($\mu\text{S}/\text{cm}$) x min



Function: $A = -0.0322401 + 7.64234E-3 \times Q$

Relative standard deviation 3.651427 %

Correlation coefficient 0.999303

Sample type	Index	Conc.	Volume	Dilution	Sample amount	Area	Ident	Date	Used
Standard 1	1	0.000	10.0	1.0	1.0	n. d.	STD1	2024-10-16 10:55:46 UTC-4	used
Standard 2	1	3.000	10.0	1.0	1.0	0.218	STD2	2024-10-16 11:17:08 UTC-4	used
Standard 3	1	6.000	10.0	1.0	1.0	0.430	STD3	2024-10-16 11:38:32 UTC-4	used
Standard 4	1	7.500	10.0	1.0	1.0	0.553	STD4	2024-10-16 11:59:57 UTC-4	used
Standard 5	1	15.000	10.0	1.0	1.0	1.077	STD5	2024-10-16 12:21:22 UTC-4	used
Standard 6	1	30.000	10.0	1.0	1.0	2.206	STD6	2024-10-16 12:42:48 UTC-4	used
Standard 7	1	37.000	10.0	1.0	1.0	2.850	STD7	2024-10-16 13:04:15 UTC-4	used


Sample data

Ident STD1
Sample type Standard 1
Determination start 2024-10-16 10:55:46 UTC-4
Method IC1-101624
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 12.33 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions




Sample data

Ident STD2
Sample type Standard 2
Determination start 2024-10-16 11:17:08 UTC-4
Method IC1-101624
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 12.56 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions



Peak number	Retention time min	Area (µS/cm) x min	Height µS/cm	Concentration ppm	Component name
1	3.365	0.0571	0.341	0.424	Fluoride
2	6.573	0.0576	0.231	0.632	Chloride
3	8.703	0.1273	0.393	0.639	Nitrite
4	10.758	0.0849	0.219	2.106	Bromide
5	12.458	0.1163	0.261	0.535	Nitrate
6	14.368	0.0842	0.154	1.072	Phosphate
7	15.893	0.2181	0.392	3.276	Sulfate


Sample data

Ident STD3
Sample type Standard 3
Determination start 2024-10-16 11:38:32 UTC-4
Method IC1-101624
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 12.11 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions



Peak number	Retention time min	Area (µS/cm) x min	Height µS/cm	Concentration ppm	Component name
1	3.368	0.1116	0.667	0.796	Fluoride
2	6.560	0.1143	0.459	1.206	Chloride
3	8.673	0.2538	0.788	1.205	Nitrite
4	10.712	0.1689	0.436	4.013	Bromide
5	12.398	0.2314	0.525	1.000	Nitrate
6	14.353	0.1657	0.305	2.000	Phosphate
7	15.922	0.4302	0.775	6.052	Sulfate


Sample data

Ident STD4
Sample type Standard 4
Determination start 2024-10-16 11:59:57 UTC-4
Method IC1-101624
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 11.94 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions



Peak number	Retention time min	Area (µS/cm) x min	Height µS/cm	Concentration ppm	Component name
1	3.367	0.1462	0.882	1.032	Fluoride
2	6.537	0.1469	0.595	1.537	Chloride
3	8.637	0.3270	1.029	1.531	Nitrite
4	10.650	0.2175	0.569	5.116	Bromide
5	12.332	0.2999	0.687	1.277	Nitrate
6	14.337	0.2147	0.401	2.557	Phosphate
7	15.945	0.5533	1.003	7.662	Sulfate


Sample data

Ident STD5
Sample type Standard 5
Determination start 2024-10-16 12:21:22 UTC-4
Method IC1-101624
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 11.99 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions



Peak number	Retention time min	Area (µS/cm) x min	Height µS/cm	Concentration ppm	Component name
1	3.363	0.2831	1.748	1.966	Fluoride
2	6.510	0.2841	1.170	2.927	Chloride
3	8.592	0.6381	2.038	2.921	Nitrite
4	10.580	0.4225	1.121	9.770	Bromide
5	12.252	0.5873	1.364	2.439	Nitrate
6	14.302	0.4179	0.793	4.871	Phosphate
7	15.955	1.0767	1.972	14.510	Sulfate


Sample data

Ident STD6
Sample type Standard 6
Determination start 2024-10-16 12:42:48 UTC-4
Method IC1-101624
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 11.88 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions



Peak number	Retention time min	Area (µS/cm) x min	Height µS/cm	Concentration ppm	Component name
1	3.370	0.5637	3.525	3.880	Fluoride
2	6.500	0.5774	2.404	5.901	Chloride
3	8.568	1.3072	4.215	5.910	Nitrite
4	10.535	0.8590	2.306	19.683	Bromide
5	12.202	1.1988	2.815	4.912	Nitrate
6	14.302	0.8536	1.629	9.831	Phosphate
7	16.002	2.2059	4.054	29.286	Sulfate


Sample data

Ident STD7
Sample type Standard 7
Determination start 2024-10-16 13:04:15 UTC-4
Method IC1-101624
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 12.73 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions



Peak number	Retention time min	Area (µS/cm) x min	Height µS/cm	Concentration ppm	Component name
1	3.373	0.7427	4.687	5.102	Fluoride
2	6.492	0.7448	3.132	7.598	Chloride
3	8.552	1.6841	5.502	7.593	Nitrite
4	10.503	1.1070	3.003	25.312	Bromide
5	12.167	1.5512	3.678	6.337	Nitrate
6	14.297	1.1030	2.137	12.670	Phosphate
7	16.028	2.8500	5.280	37.714	Sulfate


Sample data

Ident ICV
Sample type Check standard 1
Determination start 2024-10-16 13:37:49 UTC-4
Method IC1-101624
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 11.66 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions




Peak number	Retention time min	Area (µS/cm) x min	Height µS/cm	Concentration ppm	Component name
1	3.377	0.2830	1.762	1.966	Fluoride
2	6.498	0.2910	1.209	2.997	Chloride
3	8.558	0.6557	2.109	3.000	Nitrite
4	10.522	0.4312	1.155	9.969	Bromide
5	12.190	0.6017	1.407	2.497	Nitrate
6	14.320	0.4076	0.753	4.753	Phosphate
7	16.035	1.1027	2.020	14.850	Sulfate

Sample data

Ident ICB
Sample type Sample
Determination start 2024-10-16 13:59:18 UTC-4
Method IC1-101624
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 11.60 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions


Sample data

Ident CCV
Sample type Check standard 1
Determination start 2024-11-04 09:38:30 UTC-5
Method IC1-101624
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 12.27 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions



Peak number	Retention time min	Area (µS/cm) x min	Height µS/cm	Concentration ppm	Component name
1	3.377	0.2863	1.892	1.988	Fluoride
2	6.590	0.2892	1.354	2.979	Chloride
3	8.723	0.6521	2.379	2.984	Nitrite
4	10.777	0.4295	1.326	9.930	Bromide
5	12.475	0.6010	1.628	2.495	Nitrate
6	14.120	0.4262	0.918	4.965	Phosphate
7	15.975	1.0912	2.288	14.701	Sulfate


Sample data

Ident CCB
Sample type Sample
Determination start 2024-11-04 10:00:00 UTC-5
Method IC1-101624
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 12.16 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions




Sample data

Ident LB133251BLW
Sample type Sample
Determination start 2024-11-04 10:21:29 UTC-5
Method IC1-101624
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 12.27 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions


Sample data

Ident LB133251BSW
Sample type Check standard 1
Determination start 2024-11-04 10:43:00 UTC-5
Method IC1-101624
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 12.22 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions



Peak number	Retention time min	Area (µS/cm) x min	Height µS/cm	Concentration ppm	Component name
1	3.373	0.2868	1.926	1.991	Fluoride
2	6.577	0.2910	1.389	2.997	Chloride
3	8.702	0.6540	2.424	2.992	Nitrite
4	10.752	0.4306	1.349	9.955	Bromide
5	12.443	0.6019	1.655	2.498	Nitrate
6	14.077	0.4259	0.923	4.962	Phosphate
7	15.932	1.0929	2.313	14.723	Sulfate


Sample data

Ident P4673-01
Sample type Sample
Determination start 2024-11-04 11:26:05 UTC-5
Method IC1-101624
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 12.11 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions



Peak number	Retention time min	Area ($\mu\text{S/cm}$) x min	Height $\mu\text{S/cm}$	Concentration ppm	Component name
1	3.398	0.0552	0.289	0.411	Fluoride
2	3.715	0.0119	0.053	invalid	
3	4.703	0.0212	0.109	invalid	
4	6.802	2.6361	9.738	26.770	Chloride
5	8.753	0.0699	0.264	0.383	Nitrite
6	10.757	0.0069	0.023	0.335	Bromide
7	12.443	0.3695	1.018	1.559	Nitrate
8	14.127	0.0212	0.040	0.355	Phosphate
9	15.563	248.8815	398.987	3257.035	Sulfate


Sample data

Ident P4673-02
Sample type Sample
Determination start 2024-11-04 11:47:38 UTC-5
Method IC1-101624
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 12.16 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions



Peak number	Retention time min	Area (µS/cm) x min	Height µS/cm	Concentration ppm	Component name
1	3.398	0.1706	1.035	1.198	Fluoride
2	4.717	0.0881	0.342	invalid	
3	6.768	1.2876	3.435	13.100	Chloride
4	10.735	0.0057	0.016	0.308	Bromide
5	12.420	0.3162	0.875	1.343	Nitrate
6	14.083	0.0538	0.106	0.726	Phosphate
7	15.533	257.3315	418.307	3367.603	Sulfate


Sample data

Ident P4673-03
Sample type Sample
Determination start 2024-11-04 12:09:11 UTC-5
Method IC1-101624
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 11.99 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions



Peak number	Retention time min	Area (µS/cm) x min	Height µS/cm	Concentration ppm	Component name
1	3.393	0.0546	0.317	0.407	Fluoride
2	3.717	0.0020	0.015	invalid	
3	4.705	0.0337	0.171	invalid	
4	6.803	1.1301	3.264	11.504	Chloride
5	12.425	0.0985	0.272	0.463	Nitrate
6	14.100	0.0266	0.051	0.417	Phosphate
7	15.550	261.9840	424.269	3428.480	Sulfate


Sample data

Ident P4673-04
Sample type Sample
Determination start 2024-11-04 12:30:46 UTC-5
Method IC1-101624
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 11.99 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions



Peak number	Retention time min	Area ($\mu\text{S/cm}$) x min	Height $\mu\text{S/cm}$	Concentration ppm	Component name
1	3.398	0.0665	0.391	0.488	Fluoride
2	4.708	0.0200	0.105	invalid	
3	6.787	1.4689	5.407	14.938	Chloride
4	8.727	0.0181	0.070	0.152	Nitrite
5	10.738	0.0052	0.014	0.297	Bromide
6	12.402	0.5342	1.493	2.225	Nitrate
7	14.085	0.0239	0.046	0.385	Phosphate
8	15.567	253.7398	414.290	3320.606	Sulfate


Sample data

Ident P4673-05
Sample type Sample
Determination start 2024-11-04 13:13:44 UTC-5
Method IC1-101624
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 11.94 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions



Peak number	Retention time min	Area ($\mu\text{S/cm}$) x min	Height $\mu\text{S/cm}$	Concentration ppm	Component name
1	3.388	0.0645	0.388	0.474	Fluoride
2	3.717	0.0023	0.016	invalid	
3	4.697	0.0279	0.132	invalid	
4	6.763	1.3827	5.046	14.064	Chloride
5	10.670	0.0053	0.017	0.299	Bromide
6	12.350	0.5303	1.492	2.209	Nitrate
7	14.073	0.0218	0.042	0.361	Phosphate
8	15.560	247.1772	405.153	3234.733	Sulfate


Sample data

Ident P4673-06
Sample type Sample
Determination start 2024-11-04 13:35:18 UTC-5
Method IC1-101624
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 11.82 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions



Peak number	Retention time min	Area ($\mu\text{S/cm}$) x min	Height $\mu\text{S/cm}$	Concentration ppm	Component name
1	3.400	0.0719	0.374	0.525	Fluoride
2	3.732	0.0209	0.094	invalid	
3	4.003	0.0188	0.059	invalid	
4	4.710	0.0684	0.349	invalid	
5	6.765	1.4852	4.303	15.104	Chloride
6	10.682	0.0040	0.011	0.269	Bromide
7	12.358	0.3923	1.100	1.651	Nitrate
8	14.090	0.0288	0.055	0.441	Phosphate
9	15.613	234.7674	390.770	3072.351	Sulfate


Sample data

Ident P4673-07
Sample type Sample
Determination start 2024-11-04 13:56:45 UTC-5
Method IC1-101624
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 11.82 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions



Peak number	Retention time min	Area (µS/cm) x min	Height µS/cm	Concentration ppm	Component name
1	3.397	0.0400	0.230	0.307	Fluoride
2	4.705	0.0302	0.153	invalid	
3	6.780	0.8819	2.365	8.987	Chloride
4	12.343	0.1646	0.466	0.730	Nitrate
5	14.105	0.0185	0.036	0.324	Phosphate
6	15.602	248.7685	404.794	3255.556	Sulfate


Sample data

Ident P4674-01
Sample type Sample
Determination start 2024-11-04 14:18:18 UTC-5
Method IC1-101624
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 11.82 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions



Peak number	Retention time min	Area (µS/cm) x min	Height µS/cm	Concentration ppm	Component name
1	3.422	0.0304	0.122	0.242	Fluoride
2	4.705	0.0067	0.034	invalid	
3	6.515	0.0261	0.163	0.312	Chloride
4	6.830	0.9756	3.059	invalid	
5	12.352	0.0760	0.213	0.372	Nitrate
6	14.140	0.0185	0.035	0.324	Phosphate
7	15.525	326.2773	483.782	4269.759	Sulfate


Sample data

Ident CCV
Sample type Check standard 1
Determination start 2024-11-04 15:01:12 UTC-5
Method IC1-101624
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 11.77 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions



Peak number	Retention time min	Area (µS/cm) x min	Height µS/cm	Concentration ppm	Component name
1	3.377	0.2871	1.989	1.994	Fluoride
2	6.523	0.2914	1.445	3.001	Chloride
3	8.602	0.6538	2.533	2.991	Nitrite
4	10.590	0.4326	1.416	9.999	Bromide
5	12.265	0.6034	1.733	2.504	Nitrate
6	14.072	0.4232	0.935	4.931	Phosphate
7	16.040	1.0966	2.368	14.771	Sulfate


Sample data

Ident CCB
Sample type Sample
Determination start 2024-11-04 15:22:42 UTC-5
Method IC1-101624
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 11.77 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions




Sample data

Ident P4674-01MS
Sample type Sample
Determination start 2024-11-04 15:44:12 UTC-5
Method IC1-101624
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 11.94 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions



Peak number	Retention time min	Area ($\mu\text{S/cm}$) x min	Height $\mu\text{S/cm}$	Concentration ppm	Component name
1	3.425	0.2909	1.439	2.020	Fluoride
2	4.703	0.0064	0.034	invalid	
3	6.840	0.8843	3.672	9.011	Chloride
4	8.747	0.4283	1.593	1.984	Nitrite
5	10.742	0.4311	1.338	9.966	Bromide
6	12.418	0.7595	2.083	3.136	Nitrate
7	14.063	0.2307	0.433	2.740	Phosphate
8	15.525	312.7652	469.384	4092.952	Sulfate


Sample data

Ident P4674-01MSD
Sample type Sample
Determination start 2024-11-04 16:05:43 UTC-5
Method IC1-101624
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 12.16 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions



Peak number	Retention time min	Area ($\mu\text{S/cm}$) x min	Height $\mu\text{S/cm}$	Concentration ppm	Component name
1	3.422	0.3040	1.495	2.108	Fluoride
2	4.705	0.0065	0.035	invalid	
3	6.855	0.8931	3.664	9.101	Chloride
4	8.778	0.4317	1.594	1.999	Nitrite
5	10.795	0.4412	1.361	10.195	Bromide
6	12.478	0.7826	2.126	3.229	Nitrate
7	14.065	0.2531	0.480	2.995	Phosphate
8	15.508	313.2287	469.462	4099.017	Sulfate


Sample data

Ident CCV
Sample type Check standard 1
Determination start 2024-11-04 16:48:37 UTC-5
Method IC1-101624
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 12.22 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions



Peak number	Retention time min	Area (µS/cm) x min	Height µS/cm	Concentration ppm	Component name
1	3.380	0.2838	1.907	1.971	Fluoride
2	6.580	0.2904	1.382	2.991	Chloride
3	8.702	0.6552	2.419	2.997	Nitrite
4	10.748	0.4329	1.347	10.007	Bromide
5	12.442	0.6039	1.648	2.506	Nitrate
6	14.078	0.4189	0.897	4.882	Phosphate
7	15.953	1.0990	2.301	14.803	Sulfate

Sample data

Ident CCB
Sample type Sample
Determination start 2024-11-04 17:09:52 UTC-5
Method IC1-101624
Operator

Anions

Data source Conductivity detector 1 (Eco IC 1)
Channel Conductivity
Recording time 19.0 min
Integration Automatically
Column type Metrosep A Supp 19 - 150/4.0
Eluent composition not defined
Flow 0.700 mL/min
Maximum flow monitored yes
Pressure 12.27 MPa
Maximum pressure monitored yes
Temperature ---- °C

Anions

