

Instrument ID: FID\_C

**Daily Analysis Runlog For Sequence/QC Batch ID # FC102324AL**

|              |            |                   |                                 |
|--------------|------------|-------------------|---------------------------------|
| Review By    | yogesh     | Review On         | 10/23/2024 4:03:24 PM           |
| Supervise By | Ankita     | Supervise On      | 10/24/2024 9:29:51 AM           |
| SubDirectory | FC102324AL | HP Acquire Method | HP Processing Method FC102324AL |

| STD. NAME                                                          | STD REF.#                               |
|--------------------------------------------------------------------|-----------------------------------------|
| Tune/Reschk<br>Initial Calibration Stds                            | PP23644,PP23646,PP23647,PP23648,PP23649 |
| CCC<br>Internal Standard/PEM                                       | PP23647                                 |
| ICV/I.BLK<br>Surrogate Standard<br>MS/MSD Standard<br>LCS Standard | PP23645,PP23650                         |

| Sr# | SampleID             | ClientID             | Data File Name | Date-Time         | Comment | Operator | Status |
|-----|----------------------|----------------------|----------------|-------------------|---------|----------|--------|
| 1   | MECL2                | MECL2                | FC067504.D     | 23 Oct 2024 13:59 |         | YP/AJ    | Ok     |
| 2   | I.BLK                | I.BLK                | FC067505.D     | 23 Oct 2024 14:35 |         | YP/AJ    | Ok     |
| 3   | 100 PPM ALIPHATIC HC | 100 PPM ALIPHATIC HC | FC067506.D     | 23 Oct 2024 15:12 |         | YP/AJ    | Ok     |
| 4   | 50 PPM ALIPHATIC HC  | 50 PPM ALIPHATIC HC  | FC067507.D     | 23 Oct 2024 15:48 |         | YP/AJ    | Ok     |
| 5   | 20 PPM ALIPHATIC HC  | 20 PPM ALIPHATIC HC  | FC067508.D     | 23 Oct 2024 16:25 |         | YP/AJ    | Ok     |
| 6   | 10 PPM ALIPHATIC HC  | 10 PPM ALIPHATIC HC  | FC067509.D     | 23 Oct 2024 17:02 |         | YP/AJ    | Ok     |
| 7   | 5 PPM ALIPHATIC HC   | 5 PPM ALIPHATIC HC   | FC067510.D     | 23 Oct 2024 17:38 |         | YP/AJ    | Ok     |
| 8   | 20 PPM ALIPHATIC HC  | 20 PPM ALIPHATIC HC  | FC067511.D     | 23 Oct 2024 18:14 |         | YP/AJ    | Ok     |
| 9   | I.BLK                | I.BLK                | FC067512.D     | 23 Oct 2024 19:27 |         | YP/AJ    | Ok     |
| 10  | 20 PPM ALIPHATIC HC  | 20 PPM ALIPHATIC HC  | FC067513.D     | 23 Oct 2024 20:03 |         | YP/AJ    | Ok     |

M : Manual Integration

Instrument ID: FID\_C

**Daily Analysis Runlog For Sequence/QC Batch ID # FC110824AL**

|              |            |                   |                                 |
|--------------|------------|-------------------|---------------------------------|
| Review By    | yogesh     | Review On         | 11/8/2024 1:23:09 PM            |
| Supervise By | Ankita     | Supervise On      | 11/11/2024 9:55:02 AM           |
| SubDirectory | FC110824AL | HP Acquire Method | HP Processing Method FC102324AL |

| STD. NAME                                                          | STD REF.#                               |
|--------------------------------------------------------------------|-----------------------------------------|
| Tune/Reschk<br>Initial Calibration Stds                            | PP23644,PP23646,PP23647,PP23648,PP23649 |
| CCC<br>Internal Standard/PEM                                       | PP23647                                 |
| ICV/I.BLK<br>Surrogate Standard<br>MS/MSD Standard<br>LCS Standard | PP23645,PP23650                         |

| Sr# | SampleID            | ClientID            | Data File Name | Date-Time         | Comment                  | Operator | Status   |
|-----|---------------------|---------------------|----------------|-------------------|--------------------------|----------|----------|
| 1   | MECL2               | MECL2               | FC067667.D     | 08 Nov 2024 10:14 |                          | YP/AJ    | Ok       |
| 2   | I.BLK               | I.BLK               | FC067668.D     | 08 Nov 2024 12:41 |                          | YP/AJ    | Ok       |
| 3   | 20 PPM ALIPHATIC HC | 20 PPM ALIPHATIC HC | FC067669.D     | 08 Nov 2024 13:17 |                          | YP/AJ    | Ok       |
| 4   | PB164790BL          | PB164790BL          | FC067670.D     | 08 Nov 2024 16:19 |                          | YP/AJ    | Ok       |
| 5   | PB164790BS          | PB164790BS          | FC067671.D     | 08 Nov 2024 16:55 |                          | YP/AJ    | Ok       |
| 6   | PB164790BSD         | PB164790BSD         | FC067672.D     | 08 Nov 2024 17:32 |                          | YP/AJ    | Ok       |
| 7   | P4718-01            | WB-307-SB01         | FC067673.D     | 08 Nov 2024 18:08 |                          | YP/AJ    | Ok       |
| 8   | P4718-01D           | P4718-01D           | FC067674.D     | 08 Nov 2024 18:45 |                          | YP/AJ    | Ok       |
| 9   | P4718-01MS          | WB-307-SB01MS       | FC067675.D     | 08 Nov 2024 19:21 | FC067673.D               | YP/AJ    | Ok       |
| 10  | P4718-01MSD         | WB-307-SB01MSD      | FC067676.D     | 08 Nov 2024 19:57 | FC067673.D!FC067675.D    | YP/AJ    | Ok       |
| 11  | P4718-02            | WB-307-SB02         | FC067677.D     | 08 Nov 2024 20:33 |                          | YP/AJ    | Ok       |
| 12  | P4771-05            | P-1                 | FC067678.D     | 08 Nov 2024 21:10 | need 2x dilution         | YP/AJ    | Dilution |
| 13  | P4771-07            | P-2                 | FC067679.D     | 08 Nov 2024 21:46 |                          | YP/AJ    | Ok       |
| 14  | P4721-01            | EX-2-TPH-1          | FC067680.D     | 08 Nov 2024 22:23 | need 50x & 200X dilution | YP/AJ    | Dilution |
| 15  | P4721-02            | EX-2-TPH-2          | FC067681.D     | 08 Nov 2024 23:00 | need 10x & 50X dilution  | YP/AJ    | Dilution |
| 16  | I.BLK               | I.BLK               | FC067682.D     | 08 Nov 2024 23:37 |                          | YP/AJ    | Ok       |
| 17  | 20 PPM ALIPHATIC HC | 20 PPM ALIPHATIC HC | FC067683.D     | 09 Nov 2024 00:50 |                          | YP/AJ    | Ok       |

M : Manual Integration

Instrument ID: FID\_C

**Daily Analysis Runlog For Sequence/QC Batch ID # FC111124AL**

|              |            |                   |                                 |
|--------------|------------|-------------------|---------------------------------|
| Review By    | yogesh     | Review On         | 11/11/2024 11:19:51 AM          |
| Supervise By | Ankita     | Supervise On      | 11/12/2024 10:07:32 AM          |
| SubDirectory | FC111124AL | HP Acquire Method | HP Processing Method FC102324AL |

| STD. NAME                                                          | STD REF.#                               |
|--------------------------------------------------------------------|-----------------------------------------|
| Tune/Reschk<br>Initial Calibration Stds                            | PP23644,PP23646,PP23647,PP23648,PP23649 |
| CCC<br>Internal Standard/PEM                                       | PP23647                                 |
| ICV/I.BLK<br>Surrogate Standard<br>MS/MSD Standard<br>LCS Standard | PP23645,PP23650                         |

| Sr# | SampleID            | ClientID            | Data File Name | Date-Time         | Comment                  | Operator | Status   |
|-----|---------------------|---------------------|----------------|-------------------|--------------------------|----------|----------|
| 1   | MECL2               | MECL2               | FC067684.D     | 11 Nov 2024 07:57 |                          | YP/AJ    | Ok       |
| 2   | I.BLK               | I.BLK               | FC067685.D     | 11 Nov 2024 08:33 |                          | YP/AJ    | Ok       |
| 3   | 20 PPM ALIPHATIC HC | 20 PPM ALIPHATIC HC | FC067686.D     | 11 Nov 2024 09:10 |                          | YP/AJ    | Ok       |
| 4   | P4771-05DL          | P-1DL               | FC067687.D     | 11 Nov 2024 09:48 |                          | YP/AJ    | Ok       |
| 5   | P4721-01DL          | EX-2-TPH-1DL        | FC067688.D     | 11 Nov 2024 10:25 | need further dilution    | YP/AJ    | Dilution |
| 6   | P4721-01DL2         | EX-2-TPH-1DL2       | FC067689.D     | 11 Nov 2024 11:02 |                          | YP/AJ    | Ok       |
| 7   | P4721-02DL          | EX-2-TPH-2DL        | FC067690.D     | 11 Nov 2024 11:39 | need further dilution    | YP/AJ    | Dilution |
| 8   | P4721-02DL2         | EX-2-TPH-2DL2       | FC067691.D     | 11 Nov 2024 12:15 |                          | YP/AJ    | Ok,M     |
| 9   | I.BLK               | I.BLK               | FC067692.D     | 11 Nov 2024 13:29 |                          | YP/AJ    | Ok       |
| 10  | 20 PPM ALIPHATIC HC | 20 PPM ALIPHATIC HC | FC067693.D     | 11 Nov 2024 14:05 |                          | YP/AJ    | Ok       |
| 11  | PB164847BL          | PB164847BL          | FC067694.D     | 11 Nov 2024 15:19 |                          | YP/AJ    | Ok       |
| 12  | PB164847BS          | PB164847BS          | FC067695.D     | 11 Nov 2024 15:56 |                          | YP/AJ    | Ok       |
| 13  | PB164847BSD         | PB164847BSD         | FC067696.D     | 11 Nov 2024 16:33 |                          | YP/AJ    | Ok       |
| 14  | P4718-04            | WB-307-SW           | FC067697.D     | 11 Nov 2024 17:10 |                          | YP/AJ    | Ok       |
| 15  | P4791-02            | HINCHMAN-OILY-WAT   | FC067698.D     | 11 Nov 2024 17:47 | need 20x & 100x dilution | YP/AJ    | Dilution |
| 16  | P4792-04            | MANHOLE-WASTE-DR    | FC067699.D     | 11 Nov 2024 18:23 | need 20x & 100x dilution | YP/AJ    | Dilution |
| 17  | P4792-04            | MANHOLE-WASTE-DR    | FC067699.D     | 11 Nov 2024 18:23 | need 20x & 100x dilution | YP/AJ    | Dilution |
| 18  | I.BLK               | I.BLK               | FC067700.D     | 12 Nov 2024 00:32 |                          | YP/AJ    | Ok       |

Instrument ID: FID\_C

**Daily Analysis Runlog For Sequence/QC Batch ID # FC111124AL**

| Review By                                                          | yogesh                                  | Review On           | 11/11/2024 11:19:51 AM          |
|--------------------------------------------------------------------|-----------------------------------------|---------------------|---------------------------------|
| Supervise By                                                       | Ankita                                  | Supervise On        | 11/12/2024 10:07:32 AM          |
| SubDirectory                                                       | FC111124AL                              | HP Acquire Method   | HP Processing Method FC102324AL |
| STD. NAME                                                          | STD REF.#                               |                     |                                 |
| Tune/Reschk<br>Initial Calibration Stds                            | PP23644,PP23646,PP23647,PP23648,PP23649 |                     |                                 |
| CCC<br>Internal Standard/PEM                                       | PP23647                                 |                     |                                 |
| ICV/I.BLK<br>Surrogate Standard<br>MS/MSD Standard<br>LCS Standard | PP23645,PP23650                         |                     |                                 |
| 19                                                                 | 20 PPM ALIPHATIC HC                     | 20 PPM ALIPHATIC HC | FC067701.D                      |
|                                                                    |                                         |                     | 12 Nov 2024 01:09               |
|                                                                    |                                         |                     | YP/AJ                           |
|                                                                    |                                         |                     | Ok,M                            |

M : Manual Integration

Instrument ID: FID\_D

**Daily Analysis Runlog For Sequence/QC Batch ID # FD102324AR**

|                                         |                                         |                   |                                 |
|-----------------------------------------|-----------------------------------------|-------------------|---------------------------------|
| Review By                               | yogesh                                  | Review On         | 10/23/2024 4:04:27 PM           |
| Supervise By                            | Ankita                                  | Supervise On      | 10/24/2024 9:29:45 AM           |
| SubDirectory                            | FD102324AR                              | HP Acquire Method | HP Processing Method FD102324AR |
| <b>STD. NAME</b>                        | <b>STD REF.#</b>                        |                   |                                 |
| Tune/Reschk<br>Initial Calibration Stds | PP23429,PP23519,PP23520,PP23521,PP23522 |                   |                                 |
| CCC                                     | PP23520                                 |                   |                                 |
| Internal Standard/PEM                   | PP23430,PP23523                         |                   |                                 |
| ICV/I.BLK                               |                                         |                   |                                 |
| Surrogate Standard                      |                                         |                   |                                 |
| MS/MSD Standard                         |                                         |                   |                                 |
| LCS Standard                            |                                         |                   |                                 |

| Sr# | SampleID            | ClientID            | Data File Name | Date-Time         | Comment | Operator | Status |
|-----|---------------------|---------------------|----------------|-------------------|---------|----------|--------|
| 1   | MECL2               | MECL2               | FD048566.D     | 23 Oct 2024 13:59 |         | YP/AJ    | Ok     |
| 2   | I.BLK               | I.BLK               | FD048567.D     | 23 Oct 2024 14:35 |         | YP/AJ    | Ok     |
| 3   | 100 PPM AROMATIC HC | 100 PPM AROMATIC HC | FD048568.D     | 23 Oct 2024 15:12 |         | YP/AJ    | Ok     |
| 4   | 50 PPM AROMATIC HC  | 50 PPM AROMATIC HC  | FD048569.D     | 23 Oct 2024 15:48 |         | YP/AJ    | Ok     |
| 5   | 20 PPM AROMATIC HC  | 20 PPM AROMATIC HC  | FD048570.D     | 23 Oct 2024 16:25 |         | YP/AJ    | Ok     |
| 6   | 10 PPM AROMATIC HC  | 10 PPM AROMATIC HC  | FD048571.D     | 23 Oct 2024 17:02 |         | YP/AJ    | Ok     |
| 7   | 5 PPM AROMATIC HC   | 5 PPM AROMATIC HC   | FD048572.D     | 23 Oct 2024 17:38 |         | YP/AJ    | Ok     |
| 8   | 20 PPM AROMATIC HC  | 20 PPM AROMATIC HC  | FD048573.D     | 23 Oct 2024 18:14 |         | YP/AJ    | Ok     |
| 9   | I.BLK               | I.BLK               | FD048574.D     | 23 Oct 2024 19:27 |         | YP/AJ    | Ok     |
| 10  | 20 PPM AROMATIC HC  | 20 PPM AROMATIC HC  | FD048575.D     | 23 Oct 2024 20:03 |         | YP/AJ    | Ok     |

M : Manual Integration

Instrument ID: FID\_D

**Daily Analysis Runlog For Sequence/QC Batch ID # FD110824AR**

|                                         |            |                                         |                       |  |            |  |  |
|-----------------------------------------|------------|-----------------------------------------|-----------------------|--|------------|--|--|
| Review By                               | yogesh     | Review On                               | 11/8/2024 1:24:18 PM  |  |            |  |  |
| Supervise By                            | Ankita     | Supervise On                            | 11/11/2024 9:55:33 AM |  |            |  |  |
| SubDirectory                            | FD110824AR | HP Acquire Method                       | HP Processing Method  |  | FD102324AR |  |  |
| STD. NAME                               |            | STD REF.#                               |                       |  |            |  |  |
| Tune/Reschk<br>Initial Calibration Stds |            | PP23429,PP23519,PP23520,PP23521,PP23522 |                       |  |            |  |  |
| CCC                                     |            | PP23520                                 |                       |  |            |  |  |
| Internal Standard/PEM                   |            | PP23430,PP23523                         |                       |  |            |  |  |
| ICV/I.BLK                               |            |                                         |                       |  |            |  |  |
| Surrogate Standard                      |            |                                         |                       |  |            |  |  |
| MS/MSD Standard                         |            |                                         |                       |  |            |  |  |
| LCS Standard                            |            |                                         |                       |  |            |  |  |

| Sr# | SampleId           | ClientID           | Data File Name | Date-Time         | Comment                | Operator | Status   |
|-----|--------------------|--------------------|----------------|-------------------|------------------------|----------|----------|
| 1   | MECL2              | MECL2              | FD048667.D     | 08 Nov 2024 10:14 |                        | YP/AJ    | Ok       |
| 2   | I.BLK              | I.BLK              | FD048668.D     | 08 Nov 2024 10:51 |                        | YP/AJ    | Ok       |
| 3   | 20 PPM AROMATIC HC | 20 PPM AROMATIC HC | FD048669.D     | 08 Nov 2024 11:28 |                        | YP/AJ    | Ok,M     |
| 4   | PP23948            | PP23948            | FD048670.D     | 08 Nov 2024 12:05 |                        | YP/AJ    | Ok       |
| 5   | I.BLK              | I.BLK              | FD048671.D     | 08 Nov 2024 12:41 |                        | YP/AJ    | Ok       |
| 6   | 20 PPM AROMATIC HC | 20 PPM AROMATIC HC | FD048672.D     | 08 Nov 2024 13:17 |                        | YP/AJ    | Ok       |
| 7   | PB164790BL         | PB164790BL         | FD048673.D     | 08 Nov 2024 16:19 |                        | YP/AJ    | Ok       |
| 8   | PB164790BS         | PB164790BS         | FD048674.D     | 08 Nov 2024 16:55 |                        | YP/AJ    | Ok,M     |
| 9   | PB164790BSD        | PB164790BSD        | FD048675.D     | 08 Nov 2024 17:32 |                        | YP/AJ    | Ok,M     |
| 10  | P4718-01           | WB-307-SB01        | FD048676.D     | 08 Nov 2024 18:08 |                        | YP/AJ    | Ok       |
| 11  | P4718-01D          | P4718-01D          | FD048677.D     | 08 Nov 2024 18:45 |                        | YP/AJ    | Ok       |
| 12  | P4718-01MS         | WB-307-SB01MS      | FD048678.D     | 08 Nov 2024 19:21 | FD048676.D             | YP/AJ    | Ok,M     |
| 13  | P4718-01MSD        | WB-307-SB01MSD     | FD048679.D     | 08 Nov 2024 19:57 | FD048676.D!FD048678.D  | YP/AJ    | Ok,M     |
| 14  | P4718-02           | WB-307-SB02        | FD048680.D     | 08 Nov 2024 20:33 |                        | YP/AJ    | Ok       |
| 15  | P4721-01           | EX-2-TPH-1         | FD048681.D     | 08 Nov 2024 21:10 | need 20x dilution      | YP/AJ    | Dilution |
| 16  | P4721-02           | EX-2-TPH-2         | FD048682.D     | 08 Nov 2024 21:46 | need 10x dilution      | YP/AJ    | Dilution |
| 17  | P4721-03           | EX-2-TPH-3         | FD048683.D     | 08 Nov 2024 22:23 | need 5x dilution       | YP/AJ    | Dilution |
| 18  | P4721-04           | EX-11-TPH-1        | FD048684.D     | 08 Nov 2024 23:00 | need 5x & 20x dilution | YP/AJ    | Dilution |

Instrument ID: FID\_D

**Daily Analysis Runlog For Sequence/QC Batch ID # FD110824AR**

|              |            |                   |                                 |
|--------------|------------|-------------------|---------------------------------|
| Review By    | yogesh     | Review On         | 11/8/2024 1:24:18 PM            |
| Supervise By | Ankita     | Supervise On      | 11/11/2024 9:55:33 AM           |
| SubDirectory | FD110824AR | HP Acquire Method | HP Processing Method FD102324AR |

| STD. NAME                                                          | STD REF.#                               |
|--------------------------------------------------------------------|-----------------------------------------|
| Tune/Reschk<br>Initial Calibration Stds                            | PP23429,PP23519,PP23520,PP23521,PP23522 |
| CCC<br>Internal Standard/PEM                                       | PP23520                                 |
| ICV/I.BLK<br>Surrogate Standard<br>MS/MSD Standard<br>LCS Standard | PP23430,PP23523                         |

|    |                    |                    |            |                   |  |       |    |
|----|--------------------|--------------------|------------|-------------------|--|-------|----|
| 19 | I.BLK              | I.BLK              | FD048685.D | 08 Nov 2024 23:37 |  | YP/AJ | Ok |
| 20 | 20 PPM AROMATIC HC | 20 PPM AROMATIC HC | FD048686.D | 09 Nov 2024 00:50 |  | YP/AJ | Ok |

M : Manual Integration

Instrument ID: FID\_D

**Daily Analysis Runlog For Sequence/QC Batch ID # FD111124AR**

|              |            |                   |                        |
|--------------|------------|-------------------|------------------------|
| Review By    | yogesh     | Review On         | 11/11/2024 11:20:54 AM |
| Supervise By | Ankita     | Supervise On      | 11/12/2024 10:07:09 AM |
| SubDirectory | FD111124AR | HP Acquire Method | HP Processing Method   |
|              |            |                   | FD102324AR             |

| STD. NAME                                                          | STD REF.#                               |
|--------------------------------------------------------------------|-----------------------------------------|
| Tune/Reschk<br>Initial Calibration Stds                            | PP23429,PP23519,PP23520,PP23521,PP23522 |
| CCC<br>Internal Standard/PEM                                       | PP23520                                 |
| ICV/I.BLK<br>Surrogate Standard<br>MS/MSD Standard<br>LCS Standard | PP23430,PP23523                         |

| Sr# | SampleID           | ClientID           | Data File Name | Date-Time         | Comment                | Operator | Status   |
|-----|--------------------|--------------------|----------------|-------------------|------------------------|----------|----------|
| 1   | MECL2              | MECL2              | FD048687.D     | 11 Nov 2024 07:57 |                        | YP/AJ    | Ok       |
| 2   | I.BLK              | I.BLK              | FD048688.D     | 11 Nov 2024 08:33 |                        | YP/AJ    | Ok       |
| 3   | 20 PPM AROMATIC HC | 20 PPM AROMATIC HC | FD048689.D     | 11 Nov 2024 09:10 |                        | YP/AJ    | Ok       |
| 4   | P4721-06           | EX-11-TPH-3        | FD048690.D     | 11 Nov 2024 09:48 |                        | YP/AJ    | Ok       |
| 5   | P4721-01DL         | EX-2-TPH-1DL       | FD048691.D     | 11 Nov 2024 10:25 |                        | YP/AJ    | Ok,M     |
| 6   | P4721-02DL         | EX-2-TPH-2DL       | FD048692.D     | 11 Nov 2024 11:02 |                        | YP/AJ    | Ok,M     |
| 7   | P4721-03DL         | EX-2-TPH-3DL       | FD048693.D     | 11 Nov 2024 11:39 |                        | YP/AJ    | Ok,M     |
| 8   | P4721-04DL         | EX-11-TPH-1DL      | FD048694.D     | 11 Nov 2024 12:15 | need further dilution  | YP/AJ    | Dilution |
| 9   | P4721-04DL2        | EX-11-TPH-1DL2     | FD048695.D     | 11 Nov 2024 12:52 |                        | YP/AJ    | Ok,M     |
| 10  | I.BLK              | I.BLK              | FD048696.D     | 11 Nov 2024 14:05 |                        | YP/AJ    | Ok       |
| 11  | 20 PPM AROMATIC HC | 20 PPM AROMATIC HC | FD048697.D     | 11 Nov 2024 14:42 |                        | YP/AJ    | Ok       |
| 12  | PB164847BL         | PB164847BL         | FD048698.D     | 11 Nov 2024 15:19 |                        | YP/AJ    | Ok       |
| 13  | PB164847BS         | PB164847BS         | FD048699.D     | 11 Nov 2024 15:56 |                        | YP/AJ    | Ok,M     |
| 14  | PB164847BSD        | PB164847BSD        | FD048700.D     | 11 Nov 2024 16:33 |                        | YP/AJ    | Ok,M     |
| 15  | P4718-04           | WB-307-SW          | FD048701.D     | 11 Nov 2024 17:10 |                        | YP/AJ    | Ok       |
| 16  | P4791-02           | HINCHMAN-OILY-WAT  | FD048702.D     | 11 Nov 2024 17:47 | need 10x dilution      | YP/AJ    | Dilution |
| 17  | P4792-04           | MANHOLE-WASTE-DR   | FD048703.D     | 11 Nov 2024 18:23 | need 2x & 20X dilution | YP/AJ    | Dilution |
| 18  | P4792-04           | MANHOLE-WASTE-DR   | FD048703.D     | 11 Nov 2024 18:23 | need 2x & 20X dilution | YP/AJ    | Dilution |

Instrument ID: FID\_D

**Daily Analysis Runlog For Sequence/QC Batch ID # FD111124AR**

|              |            |                   |                                 |
|--------------|------------|-------------------|---------------------------------|
| Review By    | yogesh     | Review On         | 11/11/2024 11:20:54 AM          |
| Supervise By | Ankita     | Supervise On      | 11/12/2024 10:07:09 AM          |
| SubDirectory | FD111124AR | HP Acquire Method | HP Processing Method FD102324AR |

| STD. NAME                                                          | STD REF.#                               |
|--------------------------------------------------------------------|-----------------------------------------|
| Tune/Reschk<br>Initial Calibration Stds                            | PP23429,PP23519,PP23520,PP23521,PP23522 |
| CCC<br>Internal Standard/PEM                                       | PP23520                                 |
| ICV/I.BLK<br>Surrogate Standard<br>MS/MSD Standard<br>LCS Standard | PP23430,PP23523                         |

|    |                    |                    |            |                   |  |       |      |
|----|--------------------|--------------------|------------|-------------------|--|-------|------|
| 19 | I.BLK              | I.BLK              | FD048704.D | 12 Nov 2024 00:32 |  | YP/AJ | Ok   |
| 20 | 20 PPM AROMATIC HC | 20 PPM AROMATIC HC | FD048705.D | 12 Nov 2024 01:09 |  | YP/AJ | Ok,M |

M : Manual Integration

Instrument ID: FID\_E

**Daily Analysis Runlog For Sequence/QC Batch ID # FE110724AL**

|              |            |                   |                                 |
|--------------|------------|-------------------|---------------------------------|
| Review By    | yogesh     | Review On         | 11/6/2024 3:54:26 PM            |
| Supervise By | Ankita     | Supervise On      | 11/7/2024 8:56:25 AM            |
| SubDirectory | FE110724AL | HP Acquire Method | HP Processing Method FE110724AL |

| STD. NAME                                                          | STD REF.#                               |
|--------------------------------------------------------------------|-----------------------------------------|
| Tune/Reschk<br>Initial Calibration Stds                            | PP23644,PP23646,PP23647,PP23648,PP23649 |
| CCC<br>Internal Standard/PEM                                       | PP23647                                 |
| ICV/I.BLK<br>Surrogate Standard<br>MS/MSD Standard<br>LCS Standard | PP23645,PP23650                         |

| Sr# | SampleID             | ClientID             | Data File Name | Date-Time         | Comment | Operator | Status |
|-----|----------------------|----------------------|----------------|-------------------|---------|----------|--------|
| 1   | MECL2                | MECL2                | FE051073.D     | 06 Nov 2024 15:51 |         | YP\AJ    | Ok     |
| 2   | I.BLK                | I.BLK                | FE051074.D     | 06 Nov 2024 16:21 |         | YP\AJ    | Ok     |
| 3   | 100 PPM ALIPHATIC HC | 100 PPM ALIPHATIC HC | FE051075.D     | 06 Nov 2024 16:52 |         | YP\AJ    | Ok     |
| 4   | 50 PPM ALIPHATIC HC  | 50 PPM ALIPHATIC HC  | FE051076.D     | 06 Nov 2024 17:22 |         | YP\AJ    | Ok     |
| 5   | 20 PPM ALIPHATIC HC  | 20 PPM ALIPHATIC HC  | FE051077.D     | 06 Nov 2024 17:52 |         | YP\AJ    | Ok     |
| 6   | 10 PPM ALIPHATIC HC  | 10 PPM ALIPHATIC HC  | FE051078.D     | 06 Nov 2024 18:22 |         | YP\AJ    | Ok,M   |
| 7   | 5 PPM ALIPHATIC HC   | 5 PPM ALIPHATIC HC   | FE051079.D     | 06 Nov 2024 18:53 |         | YP\AJ    | Ok,M   |
| 8   | 20 PPM ALIPHATIC HC  | 20 PPM ALIPHATIC HC  | FE051080.D     | 06 Nov 2024 19:23 |         | YP\AJ    | Ok     |
| 9   | I.BLK                | I.BLK                | FE051081.D     | 06 Nov 2024 20:23 |         | YP\AJ    | Ok     |
| 10  | 20 PPM ALIPHATIC HC  | 20 PPM ALIPHATIC HC  | FE051082.D     | 06 Nov 2024 20:54 |         | YP\AJ    | Ok     |

M : Manual Integration

Instrument ID: FID\_E

**Daily Analysis Runlog For Sequence/QC Batch ID # FE110924AL**

|              |            |                   |                                 |
|--------------|------------|-------------------|---------------------------------|
| Review By    | yogesh     | Review On         | 11/8/2024 10:15:00 AM           |
| Supervise By | Ankita     | Supervise On      | 11/11/2024 9:52:08 AM           |
| SubDirectory | FE110924AL | HP Acquire Method | HP Processing Method FE110724AL |

| STD. NAME                                                          | STD REF.#                               |
|--------------------------------------------------------------------|-----------------------------------------|
| Tune/Reschk<br>Initial Calibration Stds                            | PP23644,PP23646,PP23647,PP23648,PP23649 |
| CCC<br>Internal Standard/PEM                                       | PP23647                                 |
| ICV/I.BLK<br>Surrogate Standard<br>MS/MSD Standard<br>LCS Standard | PP23645,PP23650                         |

| Sr# | SampleID            | ClientID            | Data File Name | Date-Time         | Comment               | Operator | Status   |
|-----|---------------------|---------------------|----------------|-------------------|-----------------------|----------|----------|
| 1   | MECL2               | MECL2               | FE051104.D     | 08 Nov 2024 04:17 |                       | YPIAJ    | Ok       |
| 2   | I.BLK               | I.BLK               | FE051105.D     | 08 Nov 2024 04:47 |                       | YPIAJ    | Ok       |
| 3   | 20 PPM ALIPHATIC HC | 20 PPM ALIPHATIC HC | FE051106.D     | 08 Nov 2024 05:18 |                       | YPIAJ    | Ok,M     |
| 4   | P4722-03DL          | WC-1(0-6)DL         | FE051107.D     | 08 Nov 2024 08:06 |                       | YPIAJ    | Ok       |
| 5   | P4722-08DL          | WC-2(0-6)DL         | FE051108.D     | 08 Nov 2024 08:36 |                       | YPIAJ    | Ok       |
| 6   | P4722-13DL          | WC-3(0-6)DL         | FE051109.D     | 08 Nov 2024 09:06 |                       | YPIAJ    | Ok       |
| 7   | P4734-01DL          | EO-03-11062024DL    | FE051110.D     | 08 Nov 2024 09:36 |                       | YPIAJ    | Ok       |
| 8   | P4734-02DL          | EO-03-11062024-E2DL | FE051111.D     | 08 Nov 2024 10:06 |                       | YPIAJ    | Ok       |
| 9   | P4739-01DL          | TP-14DL             | FE051112.D     | 08 Nov 2024 10:36 |                       | YPIAJ    | Ok,M     |
| 10  | P4739-02DL          | TP-14-EPHDL         | FE051113.D     | 08 Nov 2024 11:06 | need further dilution | YPIAJ    | Dilution |
| 11  | P4739-02DL2         | TP-14-EPHDL2        | FE051114.D     | 08 Nov 2024 11:36 |                       | YPIAJ    | Ok,M     |
| 12  | I.BLK               | I.BLK               | FE051115.D     | 08 Nov 2024 12:37 |                       | YPIAJ    | Ok       |
| 13  | 20 PPM ALIPHATIC HC | 20 PPM ALIPHATIC HC | FE051116.D     | 08 Nov 2024 13:07 |                       | YPIAJ    | Ok       |
| 14  | P4739-05DL          | BP-G2DL             | FE051117.D     | 08 Nov 2024 13:37 |                       | YPIAJ    | Ok,M     |
| 15  | P4739-09DL          | BP-B2DL             | FE051118.D     | 08 Nov 2024 14:07 |                       | YPIAJ    | Ok       |
| 16  | P4739-10DL          | BP-B2-EPHDL         | FE051119.D     | 08 Nov 2024 14:37 |                       | YPIAJ    | Ok       |
| 17  | P4739-13DL          | TP-11DL             | FE051120.D     | 08 Nov 2024 15:08 |                       | YPIAJ    | Ok       |
| 18  | P4739-14DL          | TP-11-EPHDL         | FE051121.D     | 08 Nov 2024 15:38 |                       | YPIAJ    | Ok       |

Instrument ID: FID\_E

**Daily Analysis Runlog For Sequence/QC Batch ID # FE110924AL**

|                                         |                                         |                   |                                 |
|-----------------------------------------|-----------------------------------------|-------------------|---------------------------------|
| Review By                               | yogesh                                  | Review On         | 11/8/2024 10:15:00 AM           |
| Supervise By                            | Ankita                                  | Supervise On      | 11/11/2024 9:52:08 AM           |
| SubDirectory                            | FE110924AL                              | HP Acquire Method | HP Processing Method FE110724AL |
| <b>STD. NAME</b>                        | <b>STD REF.#</b>                        |                   |                                 |
| Tune/Reschk<br>Initial Calibration Stds | PP23644,PP23646,PP23647,PP23648,PP23649 |                   |                                 |
| CCC                                     | PP23647                                 |                   |                                 |
| Internal Standard/PEM                   |                                         |                   |                                 |
| ICV/I.BLK                               | PP23645,PP23650                         |                   |                                 |
| Surrogate Standard                      |                                         |                   |                                 |
| MS/MSD Standard                         |                                         |                   |                                 |
| LCS Standard                            |                                         |                   |                                 |

|    |                     |                     |            |                   |                           |       |          |
|----|---------------------|---------------------|------------|-------------------|---------------------------|-------|----------|
| 19 | P4739-06DL          | BP-G2-EPHDL         | FE051122.D | 08 Nov 2024 16:08 |                           | YP\AJ | Ok,M     |
| 20 | I.BLK               | I.BLK               | FE051123.D | 08 Nov 2024 17:08 |                           | YP\AJ | Ok       |
| 21 | 20 PPM ALIPHATIC HC | 20 PPM ALIPHATIC HC | FE051124.D | 08 Nov 2024 17:38 |                           | YP\AJ | Ok       |
| 22 | P4721-03            | EX-2-TPH-3          | FE051125.D | 08 Nov 2024 18:38 | need 20x dilution         | YP\AJ | Dilution |
| 23 | P4721-04            | EX-11-TPH-1         | FE051126.D | 08 Nov 2024 19:08 | need 10x & 50x dilution   | YP\AJ | Dilution |
| 24 | P4721-10            | EX-9-TPH-3          | FE051127.D | 08 Nov 2024 19:38 |                           | YP\AJ | Ok       |
| 25 | P4721-06            | EX-11-TPH-3         | FE051128.D | 08 Nov 2024 20:08 | need 2x dilution          | YP\AJ | Dilution |
| 26 | P4721-07            | EX-11-TPH-4         | FE051129.D | 08 Nov 2024 20:38 | need 20x & 100x dilution  | YP\AJ | Dilution |
| 27 | P4721-08            | EX-9-TPH-1          | FE051130.D | 08 Nov 2024 21:08 |                           | YP\AJ | Ok       |
| 28 | P4721-09            | EX-9-TPH-2          | FE051131.D | 08 Nov 2024 21:39 |                           | YP\AJ | Ok       |
| 29 | P4721-05            | EX-11-TPH-2         | FE051132.D | 08 Nov 2024 22:09 | need further 50x dilution | YP\AJ | Dilution |
| 30 | I.BLK               | I.BLK               | FE051133.D | 08 Nov 2024 23:39 |                           | YP\AJ | Ok       |
| 31 | 20 PPM ALIPHATIC HC | 20 PPM ALIPHATIC HC | FE051134.D | 09 Nov 2024 00:09 |                           | YP\AJ | Ok       |

M : Manual Integration

Instrument ID: FID\_E

**Daily Analysis Runlog For Sequence/QC Batch ID # FE111124AL**

|              |            |                   |                                 |
|--------------|------------|-------------------|---------------------------------|
| Review By    | yogesh     | Review On         | 11/11/2024 11:21:52 AM          |
| Supervise By | Ankita     | Supervise On      | 11/12/2024 10:06:35 AM          |
| SubDirectory | FE111124AL | HP Acquire Method | HP Processing Method FE110724AL |

| STD. NAME                                                          | STD REF.#                               |
|--------------------------------------------------------------------|-----------------------------------------|
| Tune/Reschk<br>Initial Calibration Stds                            | PP23644,PP23646,PP23647,PP23648,PP23649 |
| CCC<br>Internal Standard/PEM                                       | PP23647                                 |
| ICV/I.BLK<br>Surrogate Standard<br>MS/MSD Standard<br>LCS Standard | PP23645,PP23650                         |

| Sr# | SampleID            | ClientID            | Data File Name | Date-Time         | Comment               | Operator | Status   |
|-----|---------------------|---------------------|----------------|-------------------|-----------------------|----------|----------|
| 1   | MECL2               | MECL2               | FE051135.D     | 11 Nov 2024 08:10 |                       | YP\AJ    | Ok       |
| 2   | I.BLK               | I.BLK               | FE051136.D     | 11 Nov 2024 08:40 |                       | YP\AJ    | Ok       |
| 3   | 20 PPM ALIPHATIC HC | 20 PPM ALIPHATIC HC | FE051137.D     | 11 Nov 2024 09:11 |                       | YP\AJ    | Ok,M     |
| 4   | P4721-03DL          | EX-2-TPH-3DL        | FE051138.D     | 11 Nov 2024 09:41 |                       | YP\AJ    | Ok,M     |
| 5   | P4721-04DL          | EX-11-TPH-1DL       | FE051139.D     | 11 Nov 2024 10:11 | need further dilution | YP\AJ    | Dilution |
| 6   | P4721-04DL2         | EX-11-TPH-1DL2      | FE051140.D     | 11 Nov 2024 10:41 |                       | YP\AJ    | Ok,M     |
| 7   | P4721-06DL          | EX-11-TPH-3DL       | FE051141.D     | 11 Nov 2024 11:11 |                       | YP\AJ    | Ok,M     |
| 8   | P4721-07DL          | EX-11-TPH-4DL       | FE051142.D     | 11 Nov 2024 11:41 | need further dilution | YP\AJ    | Dilution |
| 9   | P4721-07DL2         | EX-11-TPH-4DL2      | FE051143.D     | 11 Nov 2024 12:11 |                       | YP\AJ    | Ok       |
| 10  | P4721-05DL          | EX-11-TPH-2DL       | FE051144.D     | 11 Nov 2024 12:42 |                       | YP\AJ    | Ok       |
| 11  | I.BLK               | I.BLK               | FE051145.D     | 11 Nov 2024 13:42 |                       | YP\AJ    | Ok       |
| 12  | 20 PPM ALIPHATIC HC | 20 PPM ALIPHATIC HC | FE051146.D     | 11 Nov 2024 14:12 |                       | YP\AJ    | Ok,M     |
| 13  | P4795-01            | LAW-23-00189        | FE051147.D     | 11 Nov 2024 14:42 |                       | YP\AJ    | Ok       |
| 14  | P4796-01            | TP-1                | FE051148.D     | 11 Nov 2024 15:12 |                       | YP\AJ    | Ok       |
| 15  | P4796-02            | TP-1                | FE051149.D     | 11 Nov 2024 15:42 |                       | YP\AJ    | Ok       |
| 16  | P4796-03            | TP-2                | FE051150.D     | 11 Nov 2024 16:12 |                       | YP\AJ    | Ok       |
| 17  | P4796-04            | TP-2                | FE051151.D     | 11 Nov 2024 16:43 |                       | YP\AJ    | Ok       |
| 18  | P4796-05            | TP-3                | FE051152.D     | 11 Nov 2024 17:13 |                       | YP\AJ    | Ok       |

Instrument ID: FID\_E

**Daily Analysis Runlog For Sequence/QC Batch ID # FE111124AL**

| Review By                               | yogesh                                  | Review On         | 11/11/2024 11:21:52 AM          |
|-----------------------------------------|-----------------------------------------|-------------------|---------------------------------|
| Supervise By                            | Ankita                                  | Supervise On      | 11/12/2024 10:06:35 AM          |
| SubDirectory                            | FE111124AL                              | HP Acquire Method | HP Processing Method FE110724AL |
| STD. NAME                               | STD REF.#                               |                   |                                 |
| Tune/Reschk<br>Initial Calibration Stds | PP23644,PP23646,PP23647,PP23648,PP23649 |                   |                                 |
| CCC                                     | PP23647                                 |                   |                                 |
| Internal Standard/PEM                   | PP23645,PP23650                         |                   |                                 |
| ICV/I.BLK                               |                                         |                   |                                 |
| Surrogate Standard                      |                                         |                   |                                 |
| MS/MSD Standard                         |                                         |                   |                                 |
| LCS Standard                            |                                         |                   |                                 |

|    |                     |                     |            |                   |                    |       |          |
|----|---------------------|---------------------|------------|-------------------|--------------------|-------|----------|
| 19 | P4796-06            | TP-3                | FE051153.D | 11 Nov 2024 17:43 |                    | YP\AJ | Ok       |
| 20 | P4796-07            | TP-4                | FE051154.D | 11 Nov 2024 18:13 |                    | YP\AJ | Ok       |
| 21 | P4796-08            | TP-4                | FE051155.D | 11 Nov 2024 18:43 |                    | YP\AJ | Ok       |
| 22 | P4798-01            | MH-6                | FE051156.D | 11 Nov 2024 19:13 | need 5x dilution   | YP\AJ | Dilution |
| 23 | P4798-02            | MH-6-EPH            | FE051157.D | 11 Nov 2024 19:43 | need 5x dilution   | YP\AJ | Dilution |
| 24 | P4792-03            | OILY-STONE-DRUM     | FE051158.D | 11 Nov 2024 20:13 | need 10x dilution  | YP\AJ | Dilution |
| 25 | P4793-01            | M00-24-00345        | FE051159.D | 11 Nov 2024 20:44 | need 500x dilution | YP\AJ | Dilution |
| 26 | I.BLK               | I.BLK               | FE051160.D | 11 Nov 2024 21:44 |                    | YP\AJ | Ok       |
| 27 | 20 PPM ALIPHATIC HC | 20 PPM ALIPHATIC HC | FE051161.D | 11 Nov 2024 22:14 |                    | YP\AJ | Ok       |

M : Manual Integration

Instrument ID: FID\_F

**Daily Analysis Runlog For Sequence/QC Batch ID # FF110724AR**

|              |            |                   |                                 |
|--------------|------------|-------------------|---------------------------------|
| Review By    | yogesh     | Review On         | 11/7/2024 4:37:31 PM            |
| Supervise By | Ankita     | Supervise On      | 11/8/2024 10:01:09 AM           |
| SubDirectory | FF110724AR | HP Acquire Method | HP Processing Method FF110724AR |

| STD. NAME                                                          | STD REF.#                               |
|--------------------------------------------------------------------|-----------------------------------------|
| Tune/Reschk<br>Initial Calibration Stds                            | PP23429,PP23519,PP23520,PP23521,PP23522 |
| CCC<br>Internal Standard/PEM                                       | PP23520                                 |
| ICV/I.BLK<br>Surrogate Standard<br>MS/MSD Standard<br>LCS Standard | PP23430,PP23523                         |

| Sr# | SampleID            | ClientID            | Data File Name | Date-Time         | Comment | Operator | Status |
|-----|---------------------|---------------------|----------------|-------------------|---------|----------|--------|
| 1   | MECL2               | MECL2               | FF014960.D     | 07 Nov 2024 16:53 |         | YP\AJ    | Ok     |
| 2   | I.BLK               | I.BLK               | FF014961.D     | 07 Nov 2024 17:21 |         | YP\AJ    | Ok     |
| 3   | 100 PPM AROMATIC HC | 100 PPM AROMATIC HC | FF014962.D     | 07 Nov 2024 17:49 |         | YP\AJ    | Ok,M   |
| 4   | 50 PPM AROMATIC HC  | 50 PPM AROMATIC HC  | FF014963.D     | 07 Nov 2024 18:18 |         | YP\AJ    | Ok     |
| 5   | 20 PPM AROMATIC HC  | 20 PPM AROMATIC HC  | FF014964.D     | 07 Nov 2024 18:46 |         | YP\AJ    | Ok     |
| 6   | 10 PPM AROMATIC HC  | 10 PPM AROMATIC HC  | FF014965.D     | 07 Nov 2024 19:15 |         | YP\AJ    | Ok     |
| 7   | 5 PPM AROMATIC HC   | 5 PPM AROMATIC HC   | FF014966.D     | 07 Nov 2024 19:43 |         | YP\AJ    | Ok     |
| 8   | 20 PPM AROMATIC HC  | 20 PPM AROMATIC HC  | FF014967.D     | 07 Nov 2024 20:11 |         | YP\AJ    | Ok,M   |
| 9   | I.BLK               | I.BLK               | FF014968.D     | 07 Nov 2024 21:08 |         | YP\AJ    | Ok     |
| 10  | 20 PPM AROMATIC HC  | 20 PPM AROMATIC HC  | FF014969.D     | 07 Nov 2024 21:36 |         | YP\AJ    | Ok,M   |

M : Manual Integration

Instrument ID: FID\_F

**Daily Analysis Runlog For Sequence/QC Batch ID # FF110824AR**

| Review By                               | yogesh                                  | Review On         | 11/8/2024 2:49:13 PM  |            |
|-----------------------------------------|-----------------------------------------|-------------------|-----------------------|------------|
| Supervise By                            | Ankita                                  | Supervise On      | 11/11/2024 9:56:02 AM |            |
| SubDirectory                            | FF110824AR                              | HP Acquire Method | HP Processing Method  | FF110724AR |
| STD. NAME                               | STD REF.#                               |                   |                       |            |
| Tune/Reschk<br>Initial Calibration Stds | PP23429,PP23519,PP23520,PP23521,PP23522 |                   |                       |            |
| CCC                                     | PP23520                                 |                   |                       |            |
| Internal Standard/PEM                   |                                         |                   |                       |            |
| ICV/I.BLK                               | PP23430,PP23523                         |                   |                       |            |
| Surrogate Standard                      |                                         |                   |                       |            |
| MS/MSD Standard                         |                                         |                   |                       |            |
| LCS Standard                            |                                         |                   |                       |            |

| Sr# | SampleId           | ClientID           | Data File Name | Date-Time         | Comment                | Operator | Status   |
|-----|--------------------|--------------------|----------------|-------------------|------------------------|----------|----------|
| 1   | MECL2              | MECL2              | FF014970.D     | 08 Nov 2024 12:20 |                        | YP\AJ    | Ok       |
| 2   | I.BLK              | I.BLK              | FF014971.D     | 08 Nov 2024 12:49 |                        | YP\AJ    | Ok       |
| 3   | 20 PPM AROMATIC HC | 20 PPM AROMATIC HC | FF014972.D     | 08 Nov 2024 14:51 |                        | YP\AJ    | Ok,M     |
| 4   | P4771-05           | P-1                | FF014973.D     | 08 Nov 2024 16:34 | need 2x dilution       | YP\AJ    | Dilution |
| 5   | P4771-07           | P-2                | FF014974.D     | 08 Nov 2024 17:02 |                        | YP\AJ    | Ok       |
| 6   | P4721-15           | EX-9-TPH-8         | FF014975.D     | 08 Nov 2024 17:30 |                        | YP\AJ    | Ok       |
| 7   | P4721-06           | P4721-06           | FF014976.D     | 08 Nov 2024 17:59 | surrogate fail         | YP\AJ    | ReRun    |
| 8   | P4721-07           | EX-11-TPH-4        | FF014977.D     | 08 Nov 2024 18:27 | need 5x & 20x dilution | YP\AJ    | Dilution |
| 9   | P4721-08           | EX-9-TPH-1         | FF014978.D     | 08 Nov 2024 18:55 |                        | YP\AJ    | Ok       |
| 10  | P4721-09           | EX-9-TPH-2         | FF014979.D     | 08 Nov 2024 19:24 |                        | YP\AJ    | Ok       |
| 11  | P4721-10           | EX-9-TPH-3         | FF014980.D     | 08 Nov 2024 19:52 |                        | YP\AJ    | Ok       |
| 12  | P4721-11           | EX-9-TPH-4         | FF014981.D     | 08 Nov 2024 20:20 |                        | YP\AJ    | Ok       |
| 13  | P4721-12           | EX-9-TPH-5         | FF014982.D     | 08 Nov 2024 20:49 |                        | YP\AJ    | Ok,M     |
| 14  | P4721-13           | EX-9-TPH-6         | FF014983.D     | 08 Nov 2024 21:17 |                        | YP\AJ    | Ok       |
| 15  | P4721-14           | EX-9-TPH-7         | FF014984.D     | 08 Nov 2024 21:45 | need 2x & 50x dilution | YP\AJ    | Dilution |
| 16  | P4721-05           | EX-11-TPH-2        | FF014985.D     | 08 Nov 2024 22:14 | further 5x dilution    | YP\AJ    | Dilution |
| 17  | I.BLK              | I.BLK              | FF014986.D     | 08 Nov 2024 22:42 |                        | YP\AJ    | Ok       |
| 18  | 20 PPM AROMATIC HC | 20 PPM AROMATIC HC | FF014987.D     | 08 Nov 2024 23:38 |                        | YP\AJ    | Ok       |

**Instrument ID:** FID\_F

**Daily Analysis Runlog For Sequence/QC Batch ID # FF110824AR**

| Review By                                                          | yogesh                                  | Review On         | 11/8/2024 2:49:13 PM  |            |
|--------------------------------------------------------------------|-----------------------------------------|-------------------|-----------------------|------------|
| Supervise By                                                       | Ankita                                  | Supervise On      | 11/11/2024 9:56:02 AM |            |
| SubDirectory                                                       | FF110824AR                              | HP Acquire Method | HP Processing Method  | FF110724AR |
| STD. NAME                                                          | STD REF.#                               |                   |                       |            |
| Tune/Reschk<br>Initial Calibration Stds                            | PP23429,PP23519,PP23520,PP23521,PP23522 |                   |                       |            |
| CCC<br>Internal Standard/PEM                                       | PP23520                                 |                   |                       |            |
| ICV/I.BLK<br>Surrogate Standard<br>MS/MSD Standard<br>LCS Standard | PP23430,PP23523                         |                   |                       |            |

M : Manual Integration

Instrument ID: FID\_F

**Daily Analysis Runlog For Sequence/QC Batch ID # FF111124AR**

|                                         |                                         |                   |                                 |
|-----------------------------------------|-----------------------------------------|-------------------|---------------------------------|
| Review By                               | yogesh                                  | Review On         | 11/11/2024 11:23:05 AM          |
| Supervise By                            | Ankita                                  | Supervise On      | 11/12/2024 10:05:55 AM          |
| SubDirectory                            | FF111124AR                              | HP Acquire Method | HP Processing Method FF110724AR |
| <b>STD. NAME</b>                        | <b>STD REF.#</b>                        |                   |                                 |
| Tune/Reschk<br>Initial Calibration Stds | PP23429,PP23519,PP23520,PP23521,PP23522 |                   |                                 |
| CCC                                     | PP23520                                 |                   |                                 |
| Internal Standard/PEM                   | PP23430,PP23523                         |                   |                                 |
| ICV/I.BLK                               |                                         |                   |                                 |
| Surrogate Standard                      |                                         |                   |                                 |
| MS/MSD Standard                         |                                         |                   |                                 |
| LCS Standard                            |                                         |                   |                                 |

| Sr# | SampleID           | ClientID           | Data File Name | Date-Time         | Comment               | Operator | Status   |
|-----|--------------------|--------------------|----------------|-------------------|-----------------------|----------|----------|
| 1   | MECL2              | MECL2              | FF014988.D     | 11 Nov 2024 08:39 |                       | YP\AJ    | Ok       |
| 2   | I.BLK              | I.BLK              | FF014989.D     | 11 Nov 2024 09:07 |                       | YP\AJ    | Ok       |
| 3   | 20 PPM AROMATIC HC | 20 PPM AROMATIC HC | FF014990.D     | 11 Nov 2024 09:35 |                       | YP\AJ    | Ok       |
| 4   | P4771-05DL         | P-1DL              | FF014991.D     | 11 Nov 2024 10:22 |                       | YP\AJ    | Ok       |
| 5   | P4721-07DL         | EX-11-TPH-4DL      | FF014992.D     | 11 Nov 2024 10:50 | need further dilution | YP\AJ    | Dilution |
| 6   | P4721-07DL2        | EX-11-TPH-4DL2     | FF014993.D     | 11 Nov 2024 11:18 |                       | YP\AJ    | Ok,M     |
| 7   | P4721-14DL         | EX-9-TPH-7DL       | FF014994.D     | 11 Nov 2024 11:47 | need further dilution | YP\AJ    | Dilution |
| 8   | P4721-14DL2        | EX-9-TPH-7DL2      | FF014995.D     | 11 Nov 2024 12:15 |                       | YP\AJ    | Ok,M     |
| 9   | P4721-05DL         | EX-11-TPH-2DL      | FF014996.D     | 11 Nov 2024 12:43 |                       | YP\AJ    | Ok,M     |
| 10  | I.BLK              | I.BLK              | FF014997.D     | 11 Nov 2024 13:40 |                       | YP\AJ    | Ok       |
| 11  | 20 PPM AROMATIC HC | 20 PPM AROMATIC HC | FF014998.D     | 11 Nov 2024 15:33 |                       | YP\AJ    | Ok       |

M : Manual Integration

Instrument ID: FID\_G

**Daily Analysis Runlog For Sequence/QC Batch ID # FG102124AL**

|              |            |                   |                                 |
|--------------|------------|-------------------|---------------------------------|
| Review By    | yogesh     | Review On         | 10/21/2024 4:07:24 PM           |
| Supervise By | Ankita     | Supervise On      | 10/22/2024 9:03:41 AM           |
| SubDirectory | FG102124AL | HP Acquire Method | HP Processing Method FG102124AL |

| STD. NAME                                                          | STD REF.#                               |
|--------------------------------------------------------------------|-----------------------------------------|
| Tune/Reschk<br>Initial Calibration Stds                            | PP23644,PP23646,PP23647,PP23648,PP23649 |
| CCC<br>Internal Standard/PEM                                       | PP23647                                 |
| ICV/I.BLK<br>Surrogate Standard<br>MS/MSD Standard<br>LCS Standard | PP23645,PP23650                         |

| Sr# | SampleID             | ClientID             | Data File Name | Date-Time         | Comment | Operator | Status |
|-----|----------------------|----------------------|----------------|-------------------|---------|----------|--------|
| 1   | MECL2                | MECL2                | FG014469.D     | 22 Oct 2024 00:00 |         | YP\AJ    | Ok     |
| 2   | I.BLK                | I.BLK                | FG014470.D     | 22 Oct 2024 00:29 |         | YP\AJ    | Ok     |
| 3   | 100 PPM ALIPHATIC HC | 100 PPM ALIPHATIC HC | FG014471.D     | 22 Oct 2024 00:58 |         | YP\AJ    | Ok     |
| 4   | 50 PPM ALIPHATIC HC  | 50 PPM ALIPHATIC HC  | FG014472.D     | 22 Oct 2024 01:27 |         | YP\AJ    | Ok     |
| 5   | 20 PPM ALIPHATIC HC  | 20 PPM ALIPHATIC HC  | FG014473.D     | 22 Oct 2024 01:56 |         | YP\AJ    | Ok     |
| 6   | 10 PPM ALIPHATIC HC  | 10 PPM ALIPHATIC HC  | FG014474.D     | 22 Oct 2024 02:25 |         | YP\AJ    | Ok     |
| 7   | 5 PPM ALIPHATIC HC   | 5 PPM ALIPHATIC HC   | FG014475.D     | 22 Oct 2024 02:54 |         | YP\AJ    | Ok     |
| 8   | 20 PPM ALIPHATIC HC  | 20 PPM ALIPHATIC HC  | FG014476.D     | 22 Oct 2024 03:23 |         | YP\AJ    | Ok     |
| 9   | I.BLK                | I.BLK                | FG014477.D     | 22 Oct 2024 04:21 |         | YP\AJ    | Ok     |
| 10  | 20 PPM ALIPHATIC HC  | 20 PPM ALIPHATIC HC  | FG014478.D     | 22 Oct 2024 04:50 |         | YP\AJ    | Ok     |

M : Manual Integration

Instrument ID: FID\_G

**Daily Analysis Runlog For Sequence/QC Batch ID # FG110824AL**

|                                         |                                         |                   |                                 |
|-----------------------------------------|-----------------------------------------|-------------------|---------------------------------|
| Review By                               | yogesh                                  | Review On         | 11/8/2024 2:48:02 PM            |
| Supervise By                            | Ankita                                  | Supervise On      | 11/11/2024 9:55:44 AM           |
| SubDirectory                            | FG110824AL                              | HP Acquire Method | HP Processing Method FG102124AL |
| <b>STD. NAME</b>                        | <b>STD REF.#</b>                        |                   |                                 |
| Tune/Reschk<br>Initial Calibration Stds | PP23644,PP23646,PP23647,PP23648,PP23649 |                   |                                 |
| CCC                                     | PP23647                                 |                   |                                 |
| Internal Standard/PEM                   |                                         |                   |                                 |
| ICV/I.BLK                               | PP23645,PP23650                         |                   |                                 |
| Surrogate Standard                      |                                         |                   |                                 |
| MS/MSD Standard                         |                                         |                   |                                 |
| LCS Standard                            |                                         |                   |                                 |

| Sr# | SampleID            | ClientID            | Data File Name | Date-Time         | Comment           | Operator | Status   |
|-----|---------------------|---------------------|----------------|-------------------|-------------------|----------|----------|
| 1   | MECL2               | MECL2               | FG014710.D     | 08 Nov 2024 12:20 |                   | YP\AJ    | Ok       |
| 2   | I.BLK               | I.BLK               | FG014711.D     | 08 Nov 2024 12:49 |                   | YP\AJ    | Ok       |
| 3   | 20 PPM ALIPHATIC HC | 20 PPM ALIPHATIC HC | FG014712.D     | 08 Nov 2024 14:51 |                   | YP\AJ    | Ok       |
| 4   | P4721-11            | EX-9-TPH-4          | FG014713.D     | 08 Nov 2024 16:34 | need 10x dilution | YP\AJ    | Dilution |
| 5   | P4721-12            | EX-9-TPH-5          | FG014714.D     | 08 Nov 2024 17:02 |                   | YP\AJ    | Ok       |
| 6   | P4721-13            | EX-9-TPH-6          | FG014715.D     | 08 Nov 2024 17:30 |                   | YP\AJ    | Ok       |
| 7   | P4721-14            | EX-9-TPH-7          | FG014716.D     | 08 Nov 2024 17:59 | need 25x dilution | YP\AJ    | Dilution |
| 8   | P4721-15            | EX-9-TPH-8          | FG014717.D     | 08 Nov 2024 18:27 |                   | YP\AJ    | Ok       |
| 9   | I.BLK               | I.BLK               | FG014718.D     | 08 Nov 2024 19:24 |                   | YP\AJ    | Ok       |
| 10  | 20 PPM ALIPHATIC HC | 20 PPM ALIPHATIC HC | FG014719.D     | 08 Nov 2024 19:52 |                   | YP\AJ    | Ok       |

M : Manual Integration

Instrument ID: FID\_G

**Daily Analysis Runlog For Sequence/QC Batch ID # FG111124AL**

|              |            |                   |                                 |
|--------------|------------|-------------------|---------------------------------|
| Review By    | yogesh     | Review On         | 11/11/2024 11:24:17 AM          |
| Supervise By | Ankita     | Supervise On      | 11/12/2024 10:05:37 AM          |
| SubDirectory | FG111124AL | HP Acquire Method | HP Processing Method FG102124AL |

| STD. NAME                                                          | STD REF.#                               |
|--------------------------------------------------------------------|-----------------------------------------|
| Tune/Reschk<br>Initial Calibration Stds                            | PP23644,PP23646,PP23647,PP23648,PP23649 |
| CCC<br>Internal Standard/PEM                                       | PP23647                                 |
| ICV/I.BLK<br>Surrogate Standard<br>MS/MSD Standard<br>LCS Standard | PP23645,PP23650                         |

| Sr# | SampleID            | ClientID            | Data File Name | Date-Time         | Comment               | Operator | Status |
|-----|---------------------|---------------------|----------------|-------------------|-----------------------|----------|--------|
| 1   | MECL2               | MECL2               | FG014721.D     | 11 Nov 2024 08:39 |                       | YPIAJ    | Ok     |
| 2   | I.BLK               | I.BLK               | FG014722.D     | 11 Nov 2024 09:07 |                       | YPIAJ    | Ok     |
| 3   | 20 PPM ALIPHATIC HC | 20 PPM ALIPHATIC HC | FG014723.D     | 11 Nov 2024 09:35 |                       | YPIAJ    | Ok     |
| 4   | P4721-11DL          | EX-9-TPH-4DL        | FG014724.D     | 11 Nov 2024 10:22 |                       | YPIAJ    | Ok     |
| 5   | P4721-14DL          | EX-9-TPH-7DL        | FG014725.D     | 11 Nov 2024 10:50 |                       | YPIAJ    | Ok     |
| 6   | I.BLK               | I.BLK               | FG014726.D     | 11 Nov 2024 11:47 |                       | YPIAJ    | Ok     |
| 7   | 20 PPM ALIPHATIC HC | 20 PPM ALIPHATIC HC | FG014727.D     | 11 Nov 2024 12:15 |                       | YPIAJ    | Ok     |
| 8   | PB164844BL          | PB164844BL          | FG014728.D     | 11 Nov 2024 12:43 |                       | YPIAJ    | Ok     |
| 9   | PB164844BS          | PB164844BS          | FG014729.D     | 11 Nov 2024 13:12 |                       | YPIAJ    | Ok,M   |
| 10  | PB164844BSD         | PB164844BSD         | FG014730.D     | 11 Nov 2024 13:40 |                       | YPIAJ    | Ok,M   |
| 11  | P4756-01            | BP-B4               | FG014731.D     | 11 Nov 2024 14:08 |                       | YPIAJ    | Ok     |
| 12  | P4756-01D           | P4756-01D           | FG014732.D     | 11 Nov 2024 14:37 |                       | YPIAJ    | Ok     |
| 13  | P4756-01MS          | BP-B4MS             | FG014733.D     | 11 Nov 2024 15:05 | FG014731.D            | YPIAJ    | Ok,M   |
| 14  | P4756-01MSD         | BP-B4MSD            | FG014734.D     | 11 Nov 2024 15:33 | FG014731.D!FG014733.D | YPIAJ    | Ok,M   |
| 15  | P4756-02            | BP-B4-EPH           | FG014735.D     | 11 Nov 2024 16:02 |                       | YPIAJ    | Ok,M   |
| 16  | P4788-01            | BP-G3               | FG014736.D     | 11 Nov 2024 16:30 |                       | YPIAJ    | Ok     |
| 17  | P4788-02            | BP-G3-EPH           | FG014737.D     | 11 Nov 2024 16:59 |                       | YPIAJ    | Ok     |
| 18  | P4789-01            | BP-F21              | FG014738.D     | 11 Nov 2024 17:27 |                       | YPIAJ    | Ok     |

Instrument ID: FID\_G

**Daily Analysis Runlog For Sequence/QC Batch ID # FG111124AL**

| Review By                               | yogesh                                  | Review On         | 11/11/2024 11:24:17 AM          |
|-----------------------------------------|-----------------------------------------|-------------------|---------------------------------|
| Supervise By                            | Ankita                                  | Supervise On      | 11/12/2024 10:05:37 AM          |
| SubDirectory                            | FG111124AL                              | HP Acquire Method | HP Processing Method FG102124AL |
| STD. NAME                               | STD REF.#                               |                   |                                 |
| Tune/Reschk<br>Initial Calibration Stds | PP23644,PP23646,PP23647,PP23648,PP23649 |                   |                                 |
| CCC                                     | PP23647                                 |                   |                                 |
| Internal Standard/PEM                   |                                         |                   |                                 |
| ICV/I.BLK                               | PP23645,PP23650                         |                   |                                 |
| Surrogate Standard                      |                                         |                   |                                 |
| MS/MSD Standard                         |                                         |                   |                                 |
| LCS Standard                            |                                         |                   |                                 |

  

|    |                     |                     |            |                   |  |       |    |
|----|---------------------|---------------------|------------|-------------------|--|-------|----|
| 19 | P4789-02            | BP-F21-EPH          | FG014739.D | 11 Nov 2024 17:55 |  | YP\AJ | Ok |
| 20 | I.BLK               | I.BLK               | FG014740.D | 11 Nov 2024 18:52 |  | YP\AJ | Ok |
| 21 | 20 PPM ALIPHATIC HC | 20 PPM ALIPHATIC HC | FG014741.D | 11 Nov 2024 19:20 |  | YP\AJ | Ok |

M : Manual Integration