

Surrogate Summary

SDG No.: P4799

Client: ENTACT

Analytical Method: 8151A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Rec	Qual	Limits	
								Low	High
I.BLK-PS028487.D	PIBLK-PS028487.D	2,4-DCAA	1	500	548	110		70 (39)	130 (175)
		2,4-DCAA	2	500	547	109		70 (39)	130 (175)
I.BLK-PS028515.D	PIBLK-PS028515.D	2,4-DCAA	1	500	565	113		70 (39)	130 (175)
		2,4-DCAA	2	500	568	114		70 (39)	130 (175)
PB164959BL	PB164959BL	2,4-DCAA	1	500	488	98		70 (39)	130 (175)
		2,4-DCAA	2	500	471	94		70 (39)	130 (175)
PB164959BS	PB164959BS	2,4-DCAA	1	500	537	107		70 (39)	130 (175)
		2,4-DCAA	2	500	506	101		70 (39)	130 (175)
PB164880TB	PB164880TB	2,4-DCAA	1	500	291	58	*	70 (39)	130 (175)
		2,4-DCAA	2	500	286	57	*	70 (39)	130 (175)
P4799-03	WC-TA2-02-C	2,4-DCAA	1	500	398	80		70 (39)	130 (175)
		2,4-DCAA	2	500	361	72		70 (39)	130 (175)
P4799-03MS	WC-TA2-02-CMS	2,4-DCAA	1	500	375	75		70 (39)	130 (175)
		2,4-DCAA	2	500	354	71		70 (39)	130 (175)
P4799-03MSD	WC-TA2-02-CMSD	2,4-DCAA	1	500	376	75		70 (39)	130 (175)
		2,4-DCAA	2	500	355	71		70 (39)	130 (175)
P4799-07	WC-TA2-03-C	2,4-DCAA	1	500	334	67	*	70 (39)	130 (175)
		2,4-DCAA	2	500	313	63	*	70 (39)	130 (175)
P4799-15	WC-TA1-02-C	2,4-DCAA	1	500	335	67	*	70 (39)	130 (175)
		2,4-DCAA	2	500	328	66	*	70 (39)	130 (175)
P4799-19	WC-TA1-03-C	2,4-DCAA	1	500	371	74		70 (39)	130 (175)
		2,4-DCAA	2	500	363	73		70 (39)	130 (175)
I.BLK-PS028527.D	PIBLK-PS028527.D	2,4-DCAA	1	500	571	114		70 (39)	130 (175)
		2,4-DCAA	2	500	577	115		70 (39)	130 (175)
I.BLK-PS028530.D	PIBLK-PS028530.D	2,4-DCAA	1	500	584	117		70 (39)	130 (175)
		2,4-DCAA	2	500	582	116		70 (39)	130 (175)
P4799-11	WC-TA1-01-C	2,4-DCAA	1	500	335	67	*	70 (39)	130 (175)
		2,4-DCAA	2	500	322	64	*	70 (39)	130 (175)
I.BLK-PS028533.D	PIBLK-PS028533.D	2,4-DCAA	1	500	588	118		70 (39)	130 (175)
		2,4-DCAA	2	500	584	117		70 (39)	130 (175)