

Surrogate Summary

SDG No.: **P4799**

Client: **ENTACT**

Analytical Method: **8081B**

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Rec	Qual	Limits	
								Low	High
I.BLK-PL092652.D	PIBLK-PL092652.D	Decachlorobiphenyl	1	20	22.7	114		30 (43)	150 (140)
		Tetrachloro-m-xylene	1	20	21.6	108		30 (77)	150 (126)
		Decachlorobiphenyl	2	20	21.7	109		30 (43)	150 (140)
		Tetrachloro-m-xylene	2	20	20.4	102		30 (77)	150 (126)
I.BLK-PL093038.D	PIBLK-PL093038.D	Decachlorobiphenyl	1	20	22.4	112		30 (43)	150 (140)
		Tetrachloro-m-xylene	1	20	23.4	117		30 (77)	150 (126)
		Decachlorobiphenyl	2	20	22.6	113		30 (43)	150 (140)
		Tetrachloro-m-xylene	2	20	23.5	118		30 (77)	150 (126)
PB164960BL	PB164960BL	Decachlorobiphenyl	1	20	22.2	111		30 (43)	150 (140)
		Tetrachloro-m-xylene	1	20	22.4	112		30 (77)	150 (126)
		Decachlorobiphenyl	2	20	22.8	114		30 (43)	150 (140)
		Tetrachloro-m-xylene	2	20	21.7	109		30 (77)	150 (126)
PB164960BS	PB164960BS	Decachlorobiphenyl	1	20	16.6	83		30 (43)	150 (140)
		Tetrachloro-m-xylene	1	20	18.5	93		30 (77)	150 (126)
		Decachlorobiphenyl	2	20	17.0	85		30 (43)	150 (140)
		Tetrachloro-m-xylene	2	20	17.1	86		30 (77)	150 (126)
PB164880TB	PB164880TB	Decachlorobiphenyl	1	20	22.1	111		30 (43)	150 (140)
		Tetrachloro-m-xylene	1	20	22.1	111		30 (77)	150 (126)
		Decachlorobiphenyl	2	20	23.0	115		30 (43)	150 (140)
		Tetrachloro-m-xylene	2	20	21.8	109		30 (77)	150 (126)
P4799-03	WC-TA2-02-C	Decachlorobiphenyl	1	20	21.8	109		30 (43)	150 (140)
		Tetrachloro-m-xylene	1	20	21.6	108		30 (77)	150 (126)
		Decachlorobiphenyl	2	20	23.1	115		30 (43)	150 (140)
		Tetrachloro-m-xylene	2	20	22.1	110		30 (77)	150 (126)
P4799-03MS	WC-TA2-02-CMS	Decachlorobiphenyl	1	20	22.0	110		30 (43)	150 (140)
		Tetrachloro-m-xylene	1	20	22.6	113		30 (77)	150 (126)
		Decachlorobiphenyl	2	20	23.8	119		30 (43)	150 (140)
		Tetrachloro-m-xylene	2	20	21.9	109		30 (77)	150 (126)
P4799-03MSD	WC-TA2-02-CMSD	Decachlorobiphenyl	1	20	21.5	108		30 (43)	150 (140)
		Tetrachloro-m-xylene	1	20	22.0	110		30 (77)	150 (126)
		Decachlorobiphenyl	2	20	23.7	119		30 (43)	150 (140)
		Tetrachloro-m-xylene	2	20	22.2	111		30 (77)	150 (126)
P4799-07	WC-TA2-03-C	Decachlorobiphenyl	1	20	18.7	93		30 (43)	150 (140)
		Tetrachloro-m-xylene	1	20	21.4	107		30 (77)	150 (126)
		Decachlorobiphenyl	2	20	19.8	99		30 (43)	150 (140)
		Tetrachloro-m-xylene	2	20	20.6	103		30 (77)	150 (126)
P4799-11	WC-TA1-01-C	Decachlorobiphenyl	1	20	21.1	106		30 (43)	150 (140)
		Tetrachloro-m-xylene	1	20	21.5	108		30 (77)	150 (126)
		Decachlorobiphenyl	2	20	22.5	112		30 (43)	150 (140)
		Tetrachloro-m-xylene	2	20	21.2	106		30 (77)	150 (126)
P4799-15	WC-TA1-02-C	Decachlorobiphenyl	1	20	21.6	108		30 (43)	150 (140)

() = LABORATORY INHOUSE LIMIT

Surrogate Summary

SDG No.: P4799

Client: ENTACT

Analytical Method: 8081B

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Rec	Qual	Limits	
								Low	High
P4799-15	WC-TA1-02-C	Tetrachloro-m-xylene	1	20	20.1	101		30 (77)	150 (126)
		Decachlorobiphenyl	2	20	22.6	113		30 (43)	150 (140)
		Tetrachloro-m-xylene	2	20	20.6	103		30 (77)	150 (126)
P4799-19	WC-TA1-03-C	Decachlorobiphenyl	1	20	22.0	110		30 (43)	150 (140)
		Tetrachloro-m-xylene	1	20	19.9	100		30 (77)	150 (126)
		Decachlorobiphenyl	2	20	23.8	119		30 (43)	150 (140)
		Tetrachloro-m-xylene	2	20	19.6	98		30 (77)	150 (126)
		Decachlorobiphenyl	1	20	23.3	117		30 (43)	150 (140)
I.BLK-PL093050.D	PIBLK-PL093050.D	Tetrachloro-m-xylene	1	20	23.4	117		30 (77)	150 (126)
		Decachlorobiphenyl	2	20	24.7	124		30 (43)	150 (140)
		Tetrachloro-m-xylene	2	20	23.7	118		30 (77)	150 (126)
		Decachlorobiphenyl	1	20	23.3	117		30 (43)	150 (140)