

Instrument ID: ECD\_L

**Daily Analysis Runlog For Sequence/QC Batch ID # PL102824**

|                          |                                                                                         |                   |                                    |
|--------------------------|-----------------------------------------------------------------------------------------|-------------------|------------------------------------|
| Review By                | Abdul                                                                                   | Review On         | 10/29/2024 9:09:59 AM              |
| Supervise By             | Ankita                                                                                  | Supervise On      | 10/29/2024 10:02:30 AM             |
| SubDirectory             | PL102824                                                                                | HP Acquire Method | HP Processing Method pl102824 8081 |
| <b>STD. NAME</b>         | <b>STD REF.#</b>                                                                        |                   |                                    |
| Tune/Reschk              | PP23793,PP23517                                                                         |                   |                                    |
| Initial Calibration Stds | PP23673,PP23674,PP23675,PP23676,PP23677,PP23678,PP23679,PP23680,PP23681,PP23682,PP23683 |                   |                                    |
| CCC                      | PP23686,PP23690,PP23695                                                                 |                   |                                    |
| Internal Standard/PEM    |                                                                                         |                   |                                    |
| ICV/I.BLK                | PP23687,PP23693,PP23698                                                                 |                   |                                    |
| Surrogate Standard       |                                                                                         |                   |                                    |
| MS/MSD Standard          |                                                                                         |                   |                                    |
| LCS Standard             |                                                                                         |                   |                                    |

| Sr# | SampleId      | ClientID      | Data File Name | Date-Time         | Comment | Operator | Status |
|-----|---------------|---------------|----------------|-------------------|---------|----------|--------|
| 1   | HEXANE        | HEXANE        | PL092651.D     | 28 Oct 2024 13:41 |         | AR\AJ    | Ok     |
| 2   | I.BLK         | I.BLK         | PL092652.D     | 28 Oct 2024 13:55 |         | AR\AJ    | Ok     |
| 3   | PEM           | PEM           | PL092653.D     | 28 Oct 2024 14:16 |         | AR\AJ    | Ok,M   |
| 4   | RESCHK        | RESCHK        | PL092654.D     | 28 Oct 2024 14:29 |         | AR\AJ    | Ok     |
| 5   | PSTDICC100    | PSTDICC100    | PL092655.D     | 28 Oct 2024 14:43 |         | AR\AJ    | Ok     |
| 6   | PSTDICC075    | PSTDICC075    | PL092656.D     | 28 Oct 2024 14:56 |         | AR\AJ    | Ok     |
| 7   | PSTDICC050    | PSTDICC050    | PL092657.D     | 28 Oct 2024 15:09 |         | AR\AJ    | Ok     |
| 8   | PSTDICC025    | PSTDICC025    | PL092658.D     | 28 Oct 2024 15:23 |         | AR\AJ    | Ok     |
| 9   | PSTDICC005    | PSTDICC005    | PL092659.D     | 28 Oct 2024 15:36 |         | AR\AJ    | Ok,M   |
| 10  | PCHLORICC1000 | PCHLORICC1000 | PL092660.D     | 28 Oct 2024 15:49 |         | AR\AJ    | Ok     |
| 11  | PCHLORICC750  | PCHLORICC750  | PL092661.D     | 28 Oct 2024 16:03 |         | AR\AJ    | Ok     |
| 12  | PCHLORICC500  | PCHLORICC500  | PL092662.D     | 28 Oct 2024 16:16 |         | AR\AJ    | Ok     |
| 13  | PCHLORICC250  | PCHLORICC250  | PL092663.D     | 28 Oct 2024 16:30 |         | AR\AJ    | Ok     |
| 14  | PCHLORICC050  | PCHLORICC050  | PL092664.D     | 28 Oct 2024 16:43 |         | AR\AJ    | Ok     |
| 15  | PTOXICC1000   | PTOXICC1000   | PL092665.D     | 28 Oct 2024 16:56 |         | AR\AJ    | Ok     |
| 16  | PTOXICC750    | PTOXICC750    | PL092666.D     | 28 Oct 2024 17:10 |         | AR\AJ    | Ok     |
| 17  | PTOXICC500    | PTOXICC500    | PL092667.D     | 28 Oct 2024 17:23 |         | AR\AJ    | Ok     |
| 18  | PTOXICC250    | PTOXICC250    | PL092668.D     | 28 Oct 2024 17:37 |         | AR\AJ    | Ok     |

Instrument ID: ECD\_L

**Daily Analysis Runlog For Sequence/QC Batch ID # PL102824**

| Review By                | Abdul                                                                                   | Review On         | 10/29/2024 9:09:59 AM              |
|--------------------------|-----------------------------------------------------------------------------------------|-------------------|------------------------------------|
| Supervise By             | Ankita                                                                                  | Supervise On      | 10/29/2024 10:02:30 AM             |
| SubDirectory             | PL102824                                                                                | HP Acquire Method | HP Processing Method pl102824 8081 |
| STD. NAME                | STD REF.#                                                                               |                   |                                    |
| Tune/Reschk              | PP23793,PP23517                                                                         |                   |                                    |
| Initial Calibration Stds | PP23673,PP23674,PP23675,PP23676,PP23677,PP23678,PP23679,PP23680,PP23681,PP23682,PP23683 |                   |                                    |
| CCC                      | PP23686,PP23690,PP23695                                                                 |                   |                                    |
| Internal Standard/PEM    |                                                                                         |                   |                                    |
| ICV/I.BLK                | PP23687,PP23693,PP23698                                                                 |                   |                                    |
| Surrogate Standard       |                                                                                         |                   |                                    |
| MS/MSD Standard          |                                                                                         |                   |                                    |
| LCS Standard             |                                                                                         |                   |                                    |

|    |              |                  |            |                   |  |       |      |
|----|--------------|------------------|------------|-------------------|--|-------|------|
| 19 | PTOXICC100   | PTOXICC100       | PL092669.D | 28 Oct 2024 17:50 |  | AR\AJ | Ok,M |
| 20 | PSTDICV050   | ICVPL102824      | PL092670.D | 28 Oct 2024 18:03 |  | AR\AJ | Ok   |
| 21 | PCHLORICV500 | ICVPL102824CHLOR | PL092671.D | 28 Oct 2024 18:30 |  | AR\AJ | Ok   |
| 22 | PTOXICV500   | ICVPL102824TOX   | PL092672.D | 28 Oct 2024 18:57 |  | AR\AJ | Ok   |
| 23 | I.BLK        | I.BLK            | PL092673.D | 28 Oct 2024 19:24 |  | AR\AJ | Ok   |
| 24 | PEM          | PEM              | PL092674.D | 28 Oct 2024 19:37 |  | AR\AJ | Ok,M |
| 25 | PSTDCCC050   | PSTDCCC050       | PL092675.D | 28 Oct 2024 19:51 |  | AR\AJ | Ok,M |
| 26 | PB164460BL   | PB164460BL       | PL092676.D | 28 Oct 2024 20:04 |  | AR\AJ | Ok,M |
| 27 | PB164460BS   | PB164460BS       | PL092677.D | 28 Oct 2024 20:17 |  | AR\AJ | Ok,M |
| 28 | P4575-01     | PL-02-102424     | PL092678.D | 28 Oct 2024 20:31 |  | AR\AJ | Ok,M |
| 29 | P4566-01     | HD-01-102524     | PL092679.D | 28 Oct 2024 20:44 |  | AR\AJ | Ok,M |
| 30 | P4567-01     | WC-1             | PL092680.D | 28 Oct 2024 20:57 |  | AR\AJ | Ok,M |
| 31 | P4567-05     | WC-2             | PL092681.D | 28 Oct 2024 21:11 |  | AR\AJ | Ok   |
| 32 | P4567-09     | WC-3             | PL092682.D | 28 Oct 2024 21:24 |  | AR\AJ | Ok,M |
| 33 | P4574-01     | GRAVEL-1         | PL092683.D | 28 Oct 2024 21:38 |  | AR\AJ | Ok,M |
| 34 | P4574-04     | GRAVEL-2         | PL092684.D | 28 Oct 2024 21:51 |  | AR\AJ | Ok,M |
| 35 | P4577-01     | TR-05-102524     | PL092685.D | 28 Oct 2024 22:04 |  | AR\AJ | Ok,M |
| 36 | P4561-01     | BP-F-19          | PL092686.D | 28 Oct 2024 22:18 |  | AR\AJ | Ok,M |
| 37 | P4561-01MS   | BP-F-19MS        | PL092687.D | 28 Oct 2024 22:31 |  | AR\AJ | Ok,M |

Instrument ID: ECD\_L

**Daily Analysis Runlog For Sequence/QC Batch ID # PL102824**

| Review By                | Abdul                                                                                   | Review On         | 10/29/2024 9:09:59 AM              |
|--------------------------|-----------------------------------------------------------------------------------------|-------------------|------------------------------------|
| Supervise By             | Ankita                                                                                  | Supervise On      | 10/29/2024 10:02:30 AM             |
| SubDirectory             | PL102824                                                                                | HP Acquire Method | HP Processing Method pl102824 8081 |
| STD. NAME                | STD REF.#                                                                               |                   |                                    |
| Tune/Reschk              | PP23793,PP23517                                                                         |                   |                                    |
| Initial Calibration Stds | PP23673,PP23674,PP23675,PP23676,PP23677,PP23678,PP23679,PP23680,PP23681,PP23682,PP23683 |                   |                                    |
| CCC                      | PP23686,PP23690,PP23695                                                                 |                   |                                    |
| Internal Standard/PEM    |                                                                                         |                   |                                    |
| ICV/I.BLK                | PP23687,PP23693,PP23698                                                                 |                   |                                    |
| Surrogate Standard       |                                                                                         |                   |                                    |
| MS/MSD Standard          |                                                                                         |                   |                                    |
| LCS Standard             |                                                                                         |                   |                                    |

|    |             |            |            |                   |  |       |      |
|----|-------------|------------|------------|-------------------|--|-------|------|
| 38 | P4561-01MSD | BP-F-19MSD | PL092688.D | 28 Oct 2024 22:44 |  | AR\AJ | Ok,M |
| 39 | P4561-05    | BP-F-18    | PL092689.D | 28 Oct 2024 22:58 |  | AR\AJ | Ok,M |
| 40 | I.BLK       | I.BLK      | PL092690.D | 28 Oct 2024 23:11 |  | AR\AJ | Ok,M |
| 41 | PSTDCCC050  | PSTDCCC050 | PL092691.D | 29 Oct 2024 00:32 |  | AR\AJ | Ok,M |

M : Manual Integration

Instrument ID: ECD\_L

**Daily Analysis Runlog For Sequence/QC Batch ID # PL112124**

|                          |                                                                                         |                   |                               |
|--------------------------|-----------------------------------------------------------------------------------------|-------------------|-------------------------------|
| Review By                | yogesh                                                                                  | Review On         | 11/22/2024 8:48:31 AM         |
| Supervise By             | Ankita                                                                                  | Supervise On      | 11/22/2024 11:29:44 AM        |
| SubDirectory             | PL112124                                                                                | HP Acquire Method | HP Processing Method PL102824 |
| <b>STD. NAME</b>         | <b>STD REF.#</b>                                                                        |                   |                               |
| Tune/Reschk              | PP23793,PP23517                                                                         |                   |                               |
| Initial Calibration Stds | PP23673,PP23674,PP23675,PP23676,PP23677,PP23678,PP23679,PP23680,PP23681,PP23682,PP23683 |                   |                               |
| CCC                      | PP23686,PP23690,PP23695                                                                 |                   |                               |
| Internal Standard/PEM    |                                                                                         |                   |                               |
| ICV/I.BLK                | PP23687,PP23693,PP23698                                                                 |                   |                               |
| Surrogate Standard       |                                                                                         |                   |                               |
| MS/MSD Standard          |                                                                                         |                   |                               |
| LCS Standard             |                                                                                         |                   |                               |

  

| Sr# | SampleId    | ClientID       | Data File Name | Date-Time         | Comment                  | Operator | Status   |
|-----|-------------|----------------|----------------|-------------------|--------------------------|----------|----------|
| 1   | HEXANE      | HEXANE         | PL093190.D     | 21 Nov 2024 09:12 |                          | AR\AJ    | Ok       |
| 2   | I.BLK       | I.BLK          | PL093191.D     | 21 Nov 2024 09:25 |                          | AR\AJ    | Ok       |
| 3   | PEM         | PEM            | PL093192.D     | 21 Nov 2024 09:38 |                          | AR\AJ    | Ok,M     |
| 4   | PSTDCCC050  | PSTDCCC050     | PL093193.D     | 21 Nov 2024 09:51 |                          | AR\AJ    | Ok,M     |
| 5   | P4925-01RE  | MH-741RE       | PL093194.D     | 21 Nov 2024 10:05 | TCMX HIGH in both column | AR\AJ    | Confirms |
| 6   | PB165164BL  | PB165164BL     | PL093195.D     | 21 Nov 2024 10:35 |                          | AR\AJ    | Ok       |
| 7   | PB165164BS  | PB165164BS     | PL093196.D     | 21 Nov 2024 10:48 |                          | AR\AJ    | Ok,M     |
| 8   | PB165060TB  | PB165060TB     | PL093197.D     | 21 Nov 2024 11:02 |                          | AR\AJ    | Ok       |
| 9   | PB165123TB  | PB165123TB     | PL093198.D     | 21 Nov 2024 11:15 |                          | AR\AJ    | Ok       |
| 10  | P4892-03    | WB-310-BOT     | PL093199.D     | 21 Nov 2024 11:28 | contamination            | AR\AJ    | Not Ok   |
| 11  | P4892-03MS  | WB-310-BOTMS   | PL093200.D     | 21 Nov 2024 11:41 |                          | AR\AJ    | Ok,M     |
| 12  | P4892-03MSD | WB-310-BOTMSD  | PL093201.D     | 21 Nov 2024 11:54 |                          | AR\AJ    | Ok,M     |
| 13  | P4921-01    | WC-11-A-202411 | PL093202.D     | 21 Nov 2024 12:08 |                          | AR\AJ    | Ok,M     |
| 14  | PB165154BL  | PB165154BL     | PL093203.D     | 21 Nov 2024 12:21 |                          | AR\AJ    | Ok       |
| 15  | PB165154BS  | PB165154BS     | PL093204.D     | 21 Nov 2024 12:34 |                          | AR\AJ    | Ok,M     |
| 16  | P4938-05    | MH-734         | PL093205.D     | 21 Nov 2024 12:47 |                          | AR\AJ    | Ok,M     |
| 17  | P4936-01    | PL-01-11202024 | PL093206.D     | 21 Nov 2024 13:00 |                          | AR\AJ    | Ok,M     |
| 18  | P4938-01    | MH-732         | PL093207.D     | 21 Nov 2024 14:14 |                          | AR\AJ    | Ok,M     |

Instrument ID: ECD\_L

**Daily Analysis Runlog For Sequence/QC Batch ID # PL112124**

|                          |                                                                                         |                   |                               |
|--------------------------|-----------------------------------------------------------------------------------------|-------------------|-------------------------------|
| Review By                | yogesh                                                                                  | Review On         | 11/22/2024 8:48:31 AM         |
| Supervise By             | Ankita                                                                                  | Supervise On      | 11/22/2024 11:29:44 AM        |
| SubDirectory             | PL112124                                                                                | HP Acquire Method | HP Processing Method PL102824 |
| <b>STD. NAME</b>         | <b>STD REF.#</b>                                                                        |                   |                               |
| Tune/Reschk              | PP23793,PP23517                                                                         |                   |                               |
| Initial Calibration Stds | PP23673,PP23674,PP23675,PP23676,PP23677,PP23678,PP23679,PP23680,PP23681,PP23682,PP23683 |                   |                               |
| CCC                      | PP23686,PP23690,PP23695                                                                 |                   |                               |
| Internal Standard/PEM    |                                                                                         |                   |                               |
| ICV/I.BLK                | PP23687,PP23693,PP23698                                                                 |                   |                               |
| Surrogate Standard       |                                                                                         |                   |                               |
| MS/MSD Standard          |                                                                                         |                   |                               |
| LCS Standard             |                                                                                         |                   |                               |

|    |             |            |            |                   |  |       |      |
|----|-------------|------------|------------|-------------------|--|-------|------|
| 19 | P4938-01MS  | MH-732MS   | PL093208.D | 21 Nov 2024 14:27 |  | AR\AJ | Ok,M |
| 20 | P4938-01MSD | MH-732MSD  | PL093209.D | 21 Nov 2024 14:40 |  | AR\AJ | Ok,M |
| 21 | I.BLK       | I.BLK      | PL093210.D | 21 Nov 2024 14:53 |  | AR\AJ | Ok   |
| 22 | PSTDCCC050  | PSTDCCC050 | PL093211.D | 21 Nov 2024 15:07 |  | AR\AJ | Ok,M |

M : Manual Integration

Instrument ID: ECD\_L

**Daily Analysis Runlog For Sequence/QC Batch ID # PL112224**

|                          |                                                                                         |                   |                               |
|--------------------------|-----------------------------------------------------------------------------------------|-------------------|-------------------------------|
| Review By                | yogesh                                                                                  | Review On         | 11/25/2024 8:55:24 AM         |
| Supervise By             | Ankita                                                                                  | Supervise On      | 11/25/2024 9:35:34 AM         |
| SubDirectory             | PL112224                                                                                | HP Acquire Method | HP Processing Method PL102824 |
| <b>STD. NAME</b>         | <b>STD REF.#</b>                                                                        |                   |                               |
| Tune/Reschk              | PP23793,PP23517                                                                         |                   |                               |
| Initial Calibration Stds | PP23673,PP23674,PP23675,PP23676,PP23677,PP23678,PP23679,PP23680,PP23681,PP23682,PP23683 |                   |                               |
| CCC                      | PP23686,PP23690,PP23695                                                                 |                   |                               |
| Internal Standard/PEM    |                                                                                         |                   |                               |
| ICV/I.BLK                | PP23687,PP23693,PP23698                                                                 |                   |                               |
| Surrogate Standard       |                                                                                         |                   |                               |
| MS/MSD Standard          |                                                                                         |                   |                               |
| LCS Standard             |                                                                                         |                   |                               |

| Sr# | SampleID    | ClientID         | Data File Name | Date-Time         | Comment       | Operator | Status |
|-----|-------------|------------------|----------------|-------------------|---------------|----------|--------|
| 1   | HEXANE      | HEXANE           | PL093212.D     | 22 Nov 2024 09:10 |               | AR\AJ    | Ok     |
| 2   | I.BLK       | I.BLK            | PL093213.D     | 22 Nov 2024 09:24 |               | AR\AJ    | Ok     |
| 3   | PEM         | PEM              | PL093214.D     | 22 Nov 2024 09:37 |               | AR\AJ    | Ok,M   |
| 4   | PSTDCCC050  | PSTDCCC050       | PL093215.D     | 22 Nov 2024 10:03 |               | AR\AJ    | Ok,M   |
| 5   | P4892-03    | WB-310-BOT       | PL093216.D     | 22 Nov 2024 10:27 |               | AR\AJ    | Ok     |
| 6   | PB165184BL  | PB165184BL       | PL093217.D     | 22 Nov 2024 12:10 |               | AR\AJ    | Ok     |
| 7   | PB165184BS  | PB165184BS       | PL093218.D     | 22 Nov 2024 12:23 | recovery high | AR\AJ    | Not Ok |
| 8   | P4948-01    | 337              | PL093219.D     | 22 Nov 2024 12:36 |               | AR\AJ    | Ok,M   |
| 9   | P4948-01MS  | 337MS            | PL093220.D     | 22 Nov 2024 12:49 |               | AR\AJ    | Ok,M   |
| 10  | P4948-01MSD | 337MSD           | PL093221.D     | 22 Nov 2024 13:03 |               | AR\AJ    | Ok,M   |
| 11  | P4954-03    | TR-06-112124     | PL093222.D     | 22 Nov 2024 13:16 |               | AR\AJ    | Ok,M   |
| 12  | P4948-03    | 72-11944         | PL093223.D     | 22 Nov 2024 13:29 |               | AR\AJ    | Ok,M   |
| 13  | P4949-03    | CONTAINMENT-STON | PL093224.D     | 22 Nov 2024 13:42 |               | AR\AJ    | Ok,M   |
| 14  | P4951-01    | AU-05-112124     | PL093225.D     | 22 Nov 2024 13:55 |               | AR\AJ    | Ok,M   |
| 15  | P4954-01    | TR-05-112124     | PL093226.D     | 22 Nov 2024 14:09 |               | AR\AJ    | Ok,M   |
| 16  | I.BLK       | I.BLK            | PL093227.D     | 22 Nov 2024 15:46 |               | AR\AJ    | Ok     |
| 17  | PSTDCCC050  | PSTDCCC050       | PL093228.D     | 22 Nov 2024 16:14 |               | AR\AJ    | Ok,M   |

M : Manual Integration