

SOP ID :	M1311-TCLP-15	Start Prep Date :	11/18/2024	Time :	16:00
SDG No :	N/A	End Prep Date :	11/19/2024	Time :	08:20
Weigh By :	JP	Combination Ratio :	20		
Balance ID :	WC SC-7	ZHE Cleaning Batch :	N/A		
pH Meter ID :	WC PH METER-1	Initial Room Temperature:	24 °C		
Extraction By :	JP	Final Room Temperature:	22 °C		
Filter By :	JP	TCLP Technician Signature :	<i>JP</i>		
Pipette ID :	WC	Supervisor By :	12		
Tumbler ID :	T-1				
TCLP Filter ID :	114771				

Standard Name	MLS USED	STD REF. # FROM LOG
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A

Chemical Used	ML/SAMPLE U	Lot Number
TCLP-FLUID-1	N/A	WP108622
HCL-TCLP,1N	N/A	WP108584
HNO3-TCLP,1N	N/A	WP108585
pH Strips	N/A	W1931,W1934,W2350,W2755
pH Strips	N/A	W1937,W1938,W1939,W1940,W1941,W1942
1 Liter Amber	N/A	23091
120ml Plastic bottle	N/A	21029
1:1 HNO3	N/A	MP83122

Extraction Conformance/Non-Conformance Comments:

Matrix spikes are added after filtration and before preservation. TUMBLER T-1 checked,30 rpm. Particle size reduction is not required. p4910-08 is used for MS-MSD.

Date / Time	Prepped Sample Relinquished By/Location	Received By/Location
11/19/24 10:30	<i>JP</i> / TCLP Room	<i>JP</i> / Met Dig
	Preparation Group	Analysis Group <i>RE</i> / EXT

TCLP EXTRACTION LOGPAGE

PB165060

Sample ID	ClientID	TCLP Vessel ID	Sample Wt (g)	Volume Extraction Fluid #1 (mL)	Multi phasic	Phase Miscible	Phases Combined	Final Leachate PH	Metals Leachate Adj. PH	Prep Pos
P4890-06	D3721	01	100.02	2000	N/A	N/A	N/A	3.0	1.0	T-1
P4892-03	WB-310-BOT	02	100.02	2000	N/A	N/A	N/A	5.8	1.5	T-1
P4893-04	MH-763	03	100.03	2000	N/A	N/A	N/A	6.2	1.0	T-1
P4893-08	MH-762	04	100.04	2000	N/A	N/A	N/A	6.0	1.0	T-1
P4910-04	MH-COTTAGE	05	100.02	2000	N/A	N/A	N/A	7.2	1.5	T-1
P4910-08	MH-759	06	100.03	2000	N/A	N/A	N/A	6.2	1.0	T-1
PB165060TB	LEB060	07	N/A	2000	N/A	N/A	N/A	4.94	1.5	T-1

SampleID	ClientID	Sample Weight (g)	Filter Weight (g)	Filtrate (mL)	Filter + Solid (After 100°C)	% solids	% Dry Solids
P4890-06	D3721	N/A	N/A	N/A	N/A	100	N/A
P4892-03	WB-310-BOT	N/A	N/A	N/A	N/A	100	N/A
P4893-04	MH-763	N/A	N/A	N/A	N/A	100	N/A
P4893-08	MH-762	N/A	N/A	N/A	N/A	100	N/A
P4910-04	MH-COTTAGE	N/A	N/A	N/A	N/A	100	N/A
P4910-08	MH-759	N/A	N/A	N/A	N/A	100	N/A
PB165060TB	LEB060	N/A	N/A	N/A	N/A	N/A	N/A

Hot Block ID : WC S-1 /WC S-2
Thermometer ID : FLASHPOINT

SampleID	ClientID	Sample Weight (g)	Volume DI Water (mL)	PH after 5 min stir	PH after 10 min stir	Extraction Fluid 1 or 2	pH Extraction Fluid
P4890-06	D3721	5.02	96.5	6.2	2.5	#1	4.94
P4892-03	WB-310-BOT	5.03	96.5	8.6	3.5	#1	4.94
P4893-04	MH-763	5.02	96.5	8.6	3.5	#1	4.94
P4893-08	MH-762	5.01	96.5	8.0	3.0	#1	4.94
P4910-04	MH-COTTAGE	5.02	96.5	9.1	4.5	#1	4.94
P4910-08	MH-759	5.03	96.5	8.6	4.0	#1	4.94
PB165060TB	LEB060	N/A	N/A	N/A	N/A	#1	4.94