

SOP ID : M1311-TCLP-15
SDG No : N/A
Weigh By : N/A
Balance ID : N/A
pH Meter ID : WC PH METER-1
Extraction By : N/A
Filter By : JP
Pipette ID : N/A
Tumbler ID : N/A
TCLP Filter ID : 114771

Start Prep Date : N/A Time : N/A
End Prep Date : N/A Time : N/A
Combination Ratio : N/A
ZHE Cleaning Batch : N/A
Initial Room Temperature: N/A
Final Room Temperature: N/A
TCLP Technician Signature : *JP*
Supervisor By : *12*

Standard Name	MLS USED	STD REF. # FROM LOG
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A

Chemical Used	ML/SAMPLE U	Lot Number
N/A	N/A	N/A
N/A	N/A	N/A
HNO3-TCLP,1N	N/A	WP108585
pH Strips	N/A	W1931,W1934,W2350,W2755
pH Strips	N/A	W1937,W1938,W1939,W1940,W1941,W1942
1 Liter Amber	N/A	23091
120ml Plastic bottle	N/A	21029
1:1 HNO3	N/A	MP83122

Extraction Conformance/Non-Conformance Comments:

Matrix spikes are added after filtration and before preservation. p4921-01 is used for MS-MSD.

Date / Time	Prepped Sample Relinquished By/Location	Received By/Location
11/20/22 11:00	<i>JP</i> / TCLP Room	<i>JP</i> / 104014
	Preparation Group	Analysis Group / EXT

TCLP EXTRACTION LOGPAGE

PB165123

Sample ID	ClientID	TCLP Vessel ID	Sample Wt (g)	Volume Extraction Fluid #1 (mL)	Multi phasic	Phase Miscible	Phases Combined	Final Leachate PH	Metals Leachate Adj. PH	Prep Pos
P4921-01	WC-11-A-202411	N/A	N/A	N/A	N/A	N/A	N/A	7.0	1.5	N/A
PB165123TB	LEB123	N/A	N/A	N/A	N/A	N/A	N/A	4.93	1.0	N/A

SampleID	ClientID	Sample Weight (g)	Filter Weight (g)	Filtrate (mL)	Filter + Solid (After 100°C)	% solids	% Dry Solids
P4921-01	WC-11-A-202411	N/A	N/A	N/A	N/A	<0.5	N/A
PB165123TB	LEB123	N/A	N/A	N/A	N/A	N/A	N/A