

Surrogate Summary

SDG No.: P4921

Client: AECOM

Analytical Method: 8082A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Rec	Qual	Limits	
								Low	High
I.BLK-PO108037.D	PIBLK-PO108037.D	Tetrachloro-m-xylene	1	20	21.2	106		60	140
		Decachlorobiphenyl	1	20	23.2	116		60	140
		Tetrachloro-m-xylene	2	20	21.0	105		60	140
		Decachlorobiphenyl	2	20	22.5	112		60	140
I.BLK-PO108153.D	PIBLK-PO108153.D	Tetrachloro-m-xylene	1	20	23.1	116		60	140
		Decachlorobiphenyl	1	20	25.0	125		60	140
		Tetrachloro-m-xylene	2	20	23.2	116		60	140
		Decachlorobiphenyl	2	20	22.0	110		60	140
PB165133BL	PB165133BL	Tetrachloro-m-xylene	1	20	21.7	108		10	157
		Decachlorobiphenyl	1	20	23.8	119		10	173
		Tetrachloro-m-xylene	2	20	21.6	108		10	157
		Decachlorobiphenyl	2	20	20.9	105		10	173
PB165133BS	PB165133BS	Tetrachloro-m-xylene	1	20	22.1	111		10	157
		Decachlorobiphenyl	1	20	23.9	120		10	173
		Tetrachloro-m-xylene	2	20	20.8	104		10	157
		Decachlorobiphenyl	2	20	20.9	104		10	173
PB165133BSD	PB165133BSD	Tetrachloro-m-xylene	1	20	22.2	111		10	157
		Decachlorobiphenyl	1	20	23.7	119		10	173
		Tetrachloro-m-xylene	2	20	20.8	104		10	157
		Decachlorobiphenyl	2	20	20.9	104		10	173
P4921-01	WC-11-A-202411	Tetrachloro-m-xylene	1	20	21.1	106		10	157
		Decachlorobiphenyl	1	20	22.5	113		10	173
		Tetrachloro-m-xylene	2	20	21.4	107		10	157
		Decachlorobiphenyl	2	20	19.8	99		10	173
I.BLK-PO108163.D	PIBLK-PO108163.D	Tetrachloro-m-xylene	1	20	23.4	117		60	140
		Decachlorobiphenyl	1	20	26.2	131		60	140
		Tetrachloro-m-xylene	2	20	23.2	116		60	140
		Decachlorobiphenyl	2	20	23.3	117		60	140