

**Instrument ID:** CV1

**Daily Analysis Runlog For Sequence/QC Batch ID # LB134207**

|                  |                                                 |              |                      |
|------------------|-------------------------------------------------|--------------|----------------------|
| Review By        | Sarabjit Jaswal                                 | Review On    | 1/9/2025 11:42:36 PM |
| Supervise By     | Mohan Bera                                      | Supervise On | 1/9/2025 11:43:26 PM |
| <b>STD. NAME</b> | <b>STD REF.#</b>                                |              |                      |
| ICAL Standard    | MP83951,MP83952,MP83953,MP83954,MP83955,MP83956 |              |                      |
| ICV Standard     | MP83957                                         |              |                      |
| CCV Standard     | MP83959                                         |              |                      |
| ICSA Standard    |                                                 |              |                      |
| CRI Standard     |                                                 |              |                      |
| LCS Standard     |                                                 |              |                      |
| Chk Standard     | MP83958,MP83960,MP83982                         |              |                      |

| Sr# | SampleId   | ClientID | QcType | Date           | Comment | Operator | Status |
|-----|------------|----------|--------|----------------|---------|----------|--------|
| 1   | S0         | S0       | CAL1   | 01/09/25 10:00 |         | Mohan    | OK     |
| 2   | S0.2       | S01      | CAL2   | 01/09/25 10:02 |         | Mohan    | OK     |
| 3   | S2.5       | S02      | CAL3   | 01/09/25 10:04 |         | Mohan    | OK     |
| 4   | S5         | S03      | CAL4   | 01/09/25 10:07 |         | Mohan    | OK     |
| 5   | S7.5       | S04      | CAL5   | 01/09/25 10:12 |         | Mohan    | OK     |
| 6   | S10        | S05      | CAL6   | 01/09/25 10:19 |         | Mohan    | OK     |
| 7   | ICV050     | ICV050   | ICV    | 01/09/25 10:28 |         | Mohan    | OK     |
| 8   | ICB050     | ICB050   | ICB    | 01/09/25 10:30 |         | Mohan    | OK     |
| 9   | CCV020     | CCV020   | CCV    | 01/09/25 10:32 |         | Mohan    | OK     |
| 10  | CCB020     | CCB020   | CCB    | 01/09/25 10:34 |         | Mohan    | OK     |
| 11  | PB165987BL | PBS987   | MB     | 01/09/25 10:37 |         | Mohan    | OK     |
| 12  | P5388-01   | MH2D26   | SAM    | 01/09/25 10:39 |         | Mohan    | OK     |
| 13  | P5388-02   | MH2D29   | SAM    | 01/09/25 10:41 |         | Mohan    | OK     |
| 14  | P5388-03   | MH2D32   | SAM    | 01/09/25 10:43 |         | Mohan    | OK     |
| 15  | P5388-04   | MH2D35   | SAM    | 01/09/25 10:46 |         | Mohan    | OK     |
| 16  | P5388-05   | MH2D38   | SAM    | 01/09/25 10:48 |         | Mohan    | OK     |
| 17  | P5388-06   | MH2D41   | SAM    | 01/09/25 10:50 |         | Mohan    | OK     |
| 18  | P5388-07   | MH2D44   | SAM    | 01/09/25 10:53 |         | Mohan    | OK     |

Instrument ID: CV1

**Daily Analysis Runlog For Sequence/QC Batch ID # LB134207**

|                  |                                                 |              |                      |
|------------------|-------------------------------------------------|--------------|----------------------|
| Review By        | Sarabjit Jaswal                                 | Review On    | 1/9/2025 11:42:36 PM |
| Supervise By     | Mohan Bera                                      | Supervise On | 1/9/2025 11:43:26 PM |
| <b>STD. NAME</b> | <b>STD REF.#</b>                                |              |                      |
| ICAL Standard    | MP83951,MP83952,MP83953,MP83954,MP83955,MP83956 |              |                      |
| ICV Standard     | MP83957                                         |              |                      |
| CCV Standard     | MP83959                                         |              |                      |
| ICSA Standard    |                                                 |              |                      |
| CRI Standard     |                                                 |              |                      |
| LCS Standard     |                                                 |              |                      |
| Chk Standard     | MP83958,MP83960,MP83982                         |              |                      |

|    |            |         |     |                |           |       |          |
|----|------------|---------|-----|----------------|-----------|-------|----------|
| 19 | P5388-08   | MH2D47  | SAM | 01/09/25 10:55 |           | Mohan | OK       |
| 20 | P5388-09   | MH2D50  | SAM | 01/09/25 10:57 |           | Mohan | OK       |
| 21 | P5388-10   | MH2D53  | SAM | 01/09/25 10:59 |           | Mohan | OK       |
| 22 | P5388-11   | MH2D56  | SAM | 01/09/25 11:02 | Hg High   | Mohan | Dilution |
| 23 | P5388-12   | MH2D59  | SAM | 01/09/25 11:04 |           | Mohan | OK       |
| 24 | P5388-13   | MH2D62  | SAM | 01/09/25 11:06 |           | Mohan | OK       |
| 25 | P5388-14   | MH2D65  | SAM | 01/09/25 11:09 |           | Mohan | OK       |
| 26 | P5388-15   | MH2D68  | SAM | 01/09/25 11:11 |           | Mohan | OK       |
| 27 | P5388-16   | MH2D71  | SAM | 01/09/25 11:13 |           | Mohan | OK       |
| 28 | P5388-17   | MH2D74  | SAM | 01/09/25 11:15 |           | Mohan | OK       |
| 29 | P5388-18   | MH2D74D | DUP | 01/09/25 11:18 |           | Mohan | OK       |
| 30 | P5388-19   | MH2D74S | MS  | 01/09/25 11:20 |           | Mohan | OK       |
| 31 | CCV021     | CCV021  | CCV | 01/09/25 11:22 |           | Mohan | OK       |
| 32 | CCB021     | CCB021  | CCB | 01/09/25 11:24 |           | Mohan | OK       |
| 33 | P5388-11DL | MH2D56  | SAM | 01/09/25 11:27 | 5X for Hg | Mohan | Dilution |
| 34 | CCV022     | CCV022  | CCV | 01/09/25 11:29 |           | Mohan | OK       |
| 35 | CCB022     | CCB022  | CCB | 01/09/25 11:31 |           | Mohan | OK       |