

WORKLIST(Hardcopy Internal Chain)

Worklist Name : PB165924

Worklist ID : 186709

Department : Digestion

Date : 12-31-2024 08:03:51

Sample	Customer Sample	Matrix	Test	Preservative	Customer	Raw Sample Storage Location	Collect Date	Method
P5388-01	MH2D26	Solid	Metals CLP MS FULL	Cool 4 deg C	USEP01	C11	10/28/2024	SFAM_MS
P5388-02	MH2D29	Solid	Metals CLP MS FULL	Cool 4 deg C	USEP01	C11	10/28/2024	SFAM_MS
P5388-03	MH2D32	Solid	Metals CLP MS FULL	Cool 4 deg C	USEP01	C11	10/28/2024	SFAM_MS
P5388-04	MH2D35	Solid	Metals CLP MS FULL	Cool 4 deg C	USEP01	C11	10/28/2024	SFAM_MS
P5388-05	MH2D38	Solid	Metals CLP MS FULL	Cool 4 deg C	USEP01	C11	10/28/2024	SFAM_MS
P5388-06	MH2D41	Solid	Metals CLP MS FULL	Cool 4 deg C	USEP01	C11	10/28/2024	SFAM_MS
P5388-07	MH2D44	Solid	Metals CLP MS FULL	Cool 4 deg C	USEP01	C11	10/28/2024	SFAM_MS
P5388-08	MH2D47	Solid	Metals CLP MS FULL	Cool 4 deg C	USEP01	C11	10/28/2024	SFAM_MS
P5388-09	MH2D50	Solid	Metals CLP MS FULL	Cool 4 deg C	USEP01	C11	10/28/2024	SFAM_MS
P5388-10	MH2D53	Solid	Metals CLP MS FULL	Cool 4 deg C	USEP01	C11	10/28/2024	SFAM_MS
P5388-11	MH2D56	Solid	Metals CLP MS FULL	Cool 4 deg C	USEP01	C11	10/28/2024	SFAM_MS
P5388-12	MH2D59	Solid	Metals CLP MS FULL	Cool 4 deg C	USEP01	C11	10/28/2024	SFAM_MS
P5388-13	MH2D62	Solid	Metals CLP MS FULL	Cool 4 deg C	USEP01	C11	10/28/2024	SFAM_MS
P5388-14	MH2D65	Solid	Metals CLP MS FULL	Cool 4 deg C	USEP01	C11	10/28/2024	SFAM_MS
P5388-15	MH2D68	Solid	Metals CLP MS FULL	Cool 4 deg C	USEP01	C11	10/28/2024	SFAM_MS
P5388-16	MH2D71	Solid	Metals CLP MS FULL	Cool 4 deg C	USEP01	C11	10/28/2024	SFAM_MS
P5388-17	MH2D74	Solid	Metals CLP MS FULL	Cool 4 deg C	USEP01	C11	10/28/2024	SFAM_MS
P5388-18	MH2D74D	Solid	Metals CLP MS FULL	Cool 4 deg C	USEP01	C11	10/28/2024	SFAM_MS
P5388-19	MH2D74S	Solid	Metals CLP MS FULL	Cool 4 deg C	USEP01	C11	10/28/2024	SFAM_MS

Date/Time 12/31/24 8:30

Raw Sample Received by: SPB.mch digester

Raw Sample Relinquished by: [Signature]

Date/Time 12/31/24 9:30

Raw Sample Received by: [Signature]

Raw Sample Relinquished by: SPB.mch digester