



TCLP EXTRACTION LOGPAGE

PB166025

SOP ID :	M1311-TCLP-15	Start Prep Date :	N/A	Time :	N/A
SDG No :	N/A	End Prep Date :	N/A	Time :	N/A
Weigh By :	N/A	Combination Ratio :	N/A		
Balance ID :	N/A	ZHE Cleaning Batch :	N/A		
pH Meter ID :	WC PH METER-1	Initial Room Temperature:	N/A		
Extraction By :	JP	Final Room Temperature:	N/A		
Filter By :	JP	TCLP Technician Signature :	JP		
Pipette ID :	N/A	Supervisor By :	IR		
Tumbler ID :	N/A				
TCLP Filter ID :	115525				

Standard Name	MLS USED	STD REF. # FROM LOG
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A

Chemical Used	ML/SAMPLE U	Lot Number
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
pH Strips	N/A	W1931,W1934,W2350,W2755
pH Strips	N/A	W1937,W1938,W1939,W1940,W1941,W1942
N/A	N/A	N/A
120ml Plastic bottle	N/A	405130101
1:1 HNO3	N/A	MP83498

Extraction Conformance/Non-Conformance Comments:

Matrix spikes are added after filtration and before preservation. q1036-02 is used for MS-MSD.

Date / Time	Prepped Sample Relinquished By/Location	Received By/Location
01-15-25 10:00	JP / TCLP Room	SL23. Int. Dig.
	Preparation Group	Analysis Group

Sample ID	ClientID	TCLP Vessel ID	Sample Wt (g)	Volume Extraction Fluid #1 (mL)	Multi phasic	Phase Miscible	Phases Combined	Final Leachate PH	Metals Leachate Adj. PH	Prep Pos
PB166025TB	LEB025	N/A	N/A	N/A	N/A	N/A	N/A	4.94	1.0	N/A
Q1036-02	NP-WS-002	N/A	N/A	N/A	N/A	N/A	N/A	7.0	1.5	N/A

SampleID	ClientID	Sample Weight (g)	Filter Weight (g)	Filtrate (mL)	Filter + Solid (After 100°C)	% solids	% Dry Solids
PB166025TB	LEB025	N/A	N/A	N/A	N/A	N/A	N/A
Q1036-02	NP-WS-002	N/A	N/A	N/A	N/A	<0.5	N/A