

|                     |                                                        |                |                       |
|---------------------|--------------------------------------------------------|----------------|-----------------------|
| Lab Name:           | <u>Alliance Technical Group, LLC</u>                   | Contract:      | <u>68HERH20D0011</u>  |
| Lab Code:           | <u>ACE</u>                                             | Case No.:      | <u>51955</u>          |
|                     |                                                        | MA No.:        | <u>SDG No.: YE8G3</u> |
| Analytical Method:  | <u>ICP-AES</u>                                         | Instrument ID: | <u>P4</u>             |
| Preparation Method: | <u></u>                                                |                |                       |
| Concentration Units | ( $\mu\text{g/L}$ , $\mu\text{g}$ or $\text{mg/kg}$ ): | <u>mg/kg</u>   |                       |

| Analyte   | Wavelength/Mass | MDL  | Date Analyzed |
|-----------|-----------------|------|---------------|
| Aluminum  | 396.15          | 2.5  | 12/17/2024    |
| Calcium   | 373.69          | 14.0 | 12/17/2024    |
| Iron      | 259.837         | 1.6  | 12/17/2024    |
| Magnesium | 279.08          | 18.0 | 12/17/2024    |
| Potassium | 769.896         | 34.0 | 12/17/2024    |
| Sodium    | 818.326         | 67.0 | 12/17/2024    |

FORM 9-IN

METHOD DETECTION LIMIT

|                     |                               |                |               |
|---------------------|-------------------------------|----------------|---------------|
| Lab Name:           | Alliance Technical Group, LLC | Contract:      | 68HERH20D0011 |
| Lab Code:           | ACE                           | Case No.:      | 51955         |
|                     |                               | MA No.:        |               |
|                     |                               | SDG No.:       | YE8G3         |
| Analytical Method:  | ICP-AES                       | Instrument ID: | P4            |
| Preparation Method: |                               |                |               |
| Concentration Units | (µg/L, µg or mg/kg):          | mg/kg          |               |

| Analyte | Wavelength/Mass | MDL | Date Analyzed |
|---------|-----------------|-----|---------------|
|---------|-----------------|-----|---------------|