

Surrogate Summary

SW-846

SDG No.: Q1119

Client: Aramark Uniforms

Analytical Method: 8270E

Lab Sample ID	Client ID	Parameter	Spike (PPM)	Result (PPM)	Recovery (%)	Qual	Limits (%)	
							Low	High
PB166067TB	PB166067TB	2-Fluorophenol	150	129	86		10	139
		Phenol-d6	150	127	85		10	134
		Nitrobenzene-d5	100	88.3	88		49	133
		2-Fluorobiphenyl	100	88.1	88		52	132
		2,4,6-Tribromophenol	150	142	95		44	137
		Terphenyl-d14	100	83.5	83		48	125
PB166115BL	PB166115BL	2-Fluorophenol	150	134	90		10	139
		Phenol-d6	150	132	88		10	134
		Nitrobenzene-d5	100	94.1	94		49	133
		2-Fluorobiphenyl	100	94.3	94		52	132
		2,4,6-Tribromophenol	150	153	102		44	137
		Terphenyl-d14	100	87.9	88		48	125
PB166115BS	PB166115BS	2-Fluorophenol	150	137	92		10	139
		Phenol-d6	150	137	92		10	134
		Nitrobenzene-d5	100	95.8	96		49	133
		2-Fluorobiphenyl	100	95.8	96		52	132
		2,4,6-Tribromophenol	150	158	106		44	137
		Terphenyl-d14	100	97.6	98		48	125
Q1115-02MS	366MS	2-Fluorophenol	150	119	79		10	139
		Phenol-d6	150	114	76		10	134
		Nitrobenzene-d5	100	86.9	87		49	133
		2-Fluorobiphenyl	100	86.0	86		52	132
		2,4,6-Tribromophenol	150	135	90		44	137
		Terphenyl-d14	100	76.2	76		48	125
Q1115-02MSD	366MSD	2-Fluorophenol	150	121	81		10	139
		Phenol-d6	150	114	76		10	134
		Nitrobenzene-d5	100	88.8	89		49	133
		2-Fluorobiphenyl	100	88.4	88		52	132
		2,4,6-Tribromophenol	150	141	94		44	137
		Terphenyl-d14	100	77.8	78		48	125
Q1119-01	FILTER CAKE	2-Fluorophenol	150	112	74		10	139
		Phenol-d6	150	110	74		10	134
		Nitrobenzene-d5	100	77.2	77		49	133
		2-Fluorobiphenyl	100	71.2	71		52	132
		2,4,6-Tribromophenol	150	122	81		44	137
		Terphenyl-d14	100	43.3	43	*	48	125
Q1119-01DL	FILTER CAKEDL	2-Fluorophenol	150	127	85		10	139
		Phenol-d6	150	125	83		10	134
		Nitrobenzene-d5	100	87.4	87		49	133
		2-Fluorobiphenyl	100	83.2	83		52	132
		2,4,6-Tribromophenol	150	146	97		44	137
		Terphenyl-d14	100	54.8	55		48	125