

## Surrogate Summary

SW-846

SDG No.: Q1121

Client: Tetra Tech NUS, Inc.

Analytical Method: 8270-Modified

| Lab Sample ID | Client ID      | Parameter               | Spike (PPM) | Result (PPM) | Recovery (%) | Qual | Limits (%) |      |
|---------------|----------------|-------------------------|-------------|--------------|--------------|------|------------|------|
|               |                |                         |             |              |              |      | Low        | High |
| PB166108BL    | PB166108BL     | 2-Methylnaphthalene-d10 | 0.4         | 0.44         | 109          |      | 30         | 150  |
|               |                | Fluoranthene-d10        | 0.4         | 0.45         | 113          |      | 30         | 150  |
|               |                | Nitrobenzene-d5         | 0.4         | 0.49         | 123          | *    | 55         | 111  |
|               |                | 2-Fluorobiphenyl        | 0.4         | 0.43         | 107          | *    | 53         | 106  |
|               |                | Terphenyl-d14           | 0.4         | 0.51         | 127          |      | 58         | 132  |
| PB166108BS    | PB166108BS     | 2-Methylnaphthalene-d10 | 0.4         | 0.43         | 106          |      | 30         | 150  |
|               |                | Fluoranthene-d10        | 0.4         | 0.44         | 111          |      | 30         | 150  |
|               |                | Nitrobenzene-d5         | 0.4         | 0.53         | 133          | *    | 55         | 111  |
|               |                | 2-Fluorobiphenyl        | 0.4         | 0.43         | 106          |      | 53         | 106  |
|               |                | Terphenyl-d14           | 0.4         | 0.45         | 112          |      | 58         | 132  |
| PB166108BSD   | PB166108BSD    | 2-Methylnaphthalene-d10 | 0.4         | 0.40         | 100          |      | 30         | 150  |
|               |                | Fluoranthene-d10        | 0.4         | 0.38         | 95           |      | 30         | 150  |
|               |                | Nitrobenzene-d5         | 0.4         | 0.51         | 126          | *    | 55         | 111  |
|               |                | 2-Fluorobiphenyl        | 0.4         | 0.43         | 107          | *    | 53         | 106  |
|               |                | Terphenyl-d14           | 0.4         | 0.44         | 111          |      | 58         | 132  |
| Q1121-01      | RW10A-20250116 | 2-Methylnaphthalene-d10 | 0.4         | 0.39         | 97           |      | 30         | 150  |
|               |                | Fluoranthene-d10        | 0.4         | 0.44         | 110          |      | 30         | 150  |
|               |                | Nitrobenzene-d5         | 0.4         | 0.47         | 118          | *    | 55         | 111  |
|               |                | 2-Fluorobiphenyl        | 0.4         | 0.38         | 94           |      | 53         | 106  |
|               |                | Terphenyl-d14           | 0.4         | 0.45         | 112          |      | 58         | 132  |