

Instrument ID: ECD_Q

Daily Analysis Runlog For Sequence/QC Batch ID # PQ013025

Review By	yogesh	Review On	1/31/2025 11:36:54 AM
Supervise By	Ankita	Supervise On	1/31/2025 11:41:16 AM
SubDirectory	PQ013025	HP Acquire Method	HP Processing Method PQ013025
STD. NAME	STD REF.#		
Tune/Reschk Initial Calibration Stds	PP24142,PP24143,PP24144,PP24145,PP24146		
CCC	PP24149		
Internal Standard/PEM			
ICV/I.BLK	PP24147		
Surrogate Standard			
MS/MSD Standard			
LCS Standard			

Sr#	SampleID	ClientID	Data File Name	Date-Time	Comment	Operator	Status
1	HEXANE	HEXANE	PQ069885.D	30 Jan 2025 09:28		YP\AJ	Ok
2	M8011.504.1 0.5 PPB IC	M8011.504.1 0.5 PPB IC	PQ069886.D	30 Jan 2025 09:39		YP\AJ	Ok
3	M8011.504.1 0.25 PPB	M8011.504.1 0.25 PPB	PQ069887.D	30 Jan 2025 09:51		YP\AJ	Ok
4	M8011.504.1 0.1 PPB IC	M8011.504.1 0.1 PPB IC	PQ069888.D	30 Jan 2025 10:03		YP\AJ	Ok
5	M8011.504.1 0.05 PPB	M8011.504.1 0.05 PPB	PQ069889.D	30 Jan 2025 10:15		YP\AJ	Ok
6	M8011.504.1 0.025 PPB	M8011.504.1 0.025 PPB	PQ069890.D	30 Jan 2025 10:27		YP\AJ	Ok
7	M8011.504.1 0.1 PPB IC	CVPQ013025	PQ069891.D	30 Jan 2025 10:38		YP\AJ	Ok
8	RT CHECK	RT CHECK	PQ069892.D	30 Jan 2025 10:50		YP\AJ	Ok
9	M8011.504.1 0.1 PPB CM	M8011.504.1 0.1 PPB CM	PQ069893.D	30 Jan 2025 11:02		YP\AJ	Ok
10	Low-Level LFB-MDL Ch	Low-Level LFB-MDL Ch	PQ069894.D	30 Jan 2025 11:14		YP\AJ	Ok
11	PB166328BL	PB166328BL	PQ069895.D	30 Jan 2025 11:25		YP\AJ	Ok
12	Q1168-07	LOD-MDL-WATER-01-0	PQ069896.D	30 Jan 2025 11:37		YP\AJ	Ok
13	Q1168-08	LOQ-WATER-02-QT1-2	PQ069897.D	30 Jan 2025 11:49		YP\AJ	Ok
14	M8011.504.1 0.1 PPB CM	M8011.504.1 0.1 PPB CM	PQ069898.D	30 Jan 2025 12:01		YP\AJ	Ok
15	PB166329BL	PB166329BL	PQ069899.D	30 Jan 2025 12:12		YP\AJ	Ok
16	PB166329BS	PB166329BS	PQ069900.D	30 Jan 2025 12:24		YP\AJ	Ok
17	PB166329BSD	PB166329BSD	PQ069901.D	30 Jan 2025 12:36		YP\AJ	Ok
18	Q1172-11	PT-EDBCP-WS	PQ069902.D	30 Jan 2025 12:48		YP\AJ	Dilution

Instrument ID: ECD_Q

Daily Analysis Runlog For Sequence/QC Batch ID # PQ013025

Review By	yogesh	Review On	1/31/2025 11:36:54 AM
Supervise By	Ankita	Supervise On	1/31/2025 11:41:16 AM
SubDirectory	PQ013025	HP Acquire Method	HP Processing Method PQ013025

STD. NAME	STD REF.#
Tune/Reschk Initial Calibration Stds	PP24142,PP24143,PP24144,PP24145,PP24146
CCC Internal Standard/PEM	PP24149
ICV/I.BLK Surrogate Standard MS/MSD Standard LCS Standard	PP24147

19	Q1172-11DL	PT-EDBCP-WSDL	PQ069903.D	30 Jan 2025 13:11		YP\AJ	Ok,M
20	M8011.504.1 0.1 PPB	M8011.504.1 0.1 PPB	PQ069904.D	30 Jan 2025 13:23		YP\AJ	Ok

M : Manual Integration