

**Instrument ID:** P7

**Daily Analysis Runlog For Sequence/QC Batch ID # LB134522**

|                  |                                                                         |              |                     |
|------------------|-------------------------------------------------------------------------|--------------|---------------------|
| Review By        | Sarabjit Jaswal                                                         | Review On    | 2/5/2025 7:11:46 AM |
| Supervise By     | Eman Mughal                                                             | Supervise On | 2/5/2025 9:21:10 AM |
| <b>STD. NAME</b> | <b>STD REF.#</b>                                                        |              |                     |
| ICAL Standard    | MP84042,MP84050,MP84049,MP84047,MP84046,MP84045,MP84044,MP84043,MP84073 |              |                     |
| ICV Standard     | MP84051                                                                 |              |                     |
| CCV Standard     | MP84052                                                                 |              |                     |
| ICSA Standard    | MP84053,MP84054                                                         |              |                     |
| CRI Standard     |                                                                         |              |                     |
| LCS Standard     | MP84118                                                                 |              |                     |
| Chk Standard     | MP84057,MP84056                                                         |              |                     |

| Sr# | SampleId   | ClientID | QcType | Date           | Comment                 | Operator | Status |
|-----|------------|----------|--------|----------------|-------------------------|----------|--------|
| 1   | TUNE       | TUNE002  | TUNE   | 02/02/25 13:51 |                         | Jaswal   | OK     |
| 2   | S00        | S0       | CAL1   | 02/02/25 14:25 |                         | Jaswal   | OK     |
| 3   | S02        | S02      | CAL3   | 02/02/25 14:28 |                         | Jaswal   | OK     |
| 4   | S03        | S03      | CAL4   | 02/02/25 14:32 |                         | Jaswal   | OK     |
| 5   | S04        | S04      | CAL5   | 02/02/25 14:35 |                         | Jaswal   | OK     |
| 6   | S05        | S05      | CAL6   | 02/02/25 14:38 |                         | Jaswal   | OK     |
| 7   | S06        | S06      | CAL7   | 02/02/25 14:40 |                         | Jaswal   | OK     |
| 8   | S07        | S07      | CAL8   | 02/02/25 14:43 |                         | Jaswal   | OK     |
| 9   | S08        | S08      | CAL9   | 02/02/25 14:46 |                         | Jaswal   | OK     |
| 10  | ICV002     | ICV002   | ICV    | 02/02/25 14:51 |                         | Jaswal   | OK     |
| 11  | ICB002     | ICB002   | ICB    | 02/02/25 14:55 |                         | Jaswal   | OK     |
| 12  | ICSA002    | ICSA002  | ICSA   | 02/02/25 15:01 | Inst. paused for review | Jaswal   | OK     |
| 13  | ICSAB002   | ICSAB002 | ICSAB  | 02/02/25 15:04 |                         | Jaswal   | OK     |
| 14  | CCV011     | CCV011   | CCV    | 02/02/25 15:07 |                         | Jaswal   | OK     |
| 15  | CCB011     | CCB011   | CCB    | 02/02/25 15:14 | Inst. paused for review | Jaswal   | OK     |
| 16  | PB166384BL | PBS384   | MB     | 02/02/25 15:17 |                         | Jaswal   | OK     |
| 17  | PB166384BS | LCS384   | LCS    | 02/02/25 15:20 |                         | Jaswal   | OK     |
| 18  | PB166471BL | PBS471   | MB     | 02/02/25 15:24 |                         | Jaswal   | OK     |

Instrument ID: P7

**Daily Analysis Runlog For Sequence/QC Batch ID # LB134522**

|                  |                                                                         |              |                     |
|------------------|-------------------------------------------------------------------------|--------------|---------------------|
| Review By        | Sarabjit Jaswal                                                         | Review On    | 2/5/2025 7:11:46 AM |
| Supervise By     | Eman Mughal                                                             | Supervise On | 2/5/2025 9:21:10 AM |
| <b>STD. NAME</b> | <b>STD REF.#</b>                                                        |              |                     |
| ICAL Standard    | MP84042,MP84050,MP84049,MP84047,MP84046,MP84045,MP84044,MP84043,MP84073 |              |                     |
| ICV Standard     | MP84051                                                                 |              |                     |
| CCV Standard     | MP84052                                                                 |              |                     |
| ICSA Standard    | MP84053,MP84054                                                         |              |                     |
| CRI Standard     |                                                                         |              |                     |
| LCS Standard     | MP84118                                                                 |              |                     |
| Chk Standard     | MP84057,MP84056                                                         |              |                     |

|    |            |        |     |                |                                                   |        |          |
|----|------------|--------|-----|----------------|---------------------------------------------------|--------|----------|
| 19 | PB166471BS | LCS471 | LCS | 02/02/25 15:27 |                                                   | Jaswal | OK       |
| 20 | Q1223-05   | A6305  | SAM | 02/02/25 15:30 | Cr High                                           | Jaswal | Dilution |
| 21 | Q1223-06   | A6305D | DUP | 02/02/25 15:33 | Cr High                                           | Jaswal | Dilution |
| 22 | Q1223-05L  | A6305L | SD  | 02/02/25 15:36 | Cr High                                           | Jaswal | Dilution |
| 23 | Q1223-07   | A6305S | MS  | 02/02/25 15:39 | Cr High, MS Fail For Sb(Below RL),Se(Below RL),Zn | Jaswal | Dilution |
| 24 | Q1246-09   | A6342  | SAM | 02/02/25 15:42 |                                                   | Jaswal | OK       |
| 25 | Q1246-10   | A6342D | DUP | 02/02/25 15:46 |                                                   | Jaswal | OK       |
| 26 | Q1246-09L  | A6342L | SD  | 02/02/25 15:49 |                                                   | Jaswal | OK       |
| 27 | Q1246-11   | A6342S | MS  | 02/02/25 15:52 | MS Fail For Se(Below RL)                          | Jaswal | OK       |
| 28 | Q1223-08   | A6306  | SAM | 02/02/25 15:55 | Cr High                                           | Jaswal | Dilution |
| 29 | Q1223-09   | A6307  | SAM | 02/02/25 15:58 |                                                   | Jaswal | OK       |
| 30 | Q1223-10   | A6308  | SAM | 02/02/25 16:02 | Cr High                                           | Jaswal | Dilution |
| 31 | Q1223-11   | A6309  | SAM | 02/02/25 16:05 |                                                   | Jaswal | OK       |
| 32 | Q1223-12   | A6320  | SAM | 02/02/25 16:08 |                                                   | Jaswal | OK       |
| 33 | Q1223-13   | A6321  | SAM | 02/02/25 16:11 | Cr High                                           | Jaswal | Dilution |
| 34 | Q1223-14   | A6322  | SAM | 02/02/25 16:14 | Cr High                                           | Jaswal | Dilution |
| 35 | Q1223-15   | A6323  | SAM | 02/02/25 16:17 | Cr High                                           | Jaswal | Dilution |
| 36 | Q1223-16   | A6324  | SAM | 02/02/25 16:20 |                                                   | Jaswal | OK       |
| 37 | CCV012     | CCV012 | CCV | 02/02/25 16:24 |                                                   | Jaswal | OK       |

Instrument ID: P7

**Daily Analysis Runlog For Sequence/QC Batch ID # LB134522**

|                  |                                                                         |              |                     |
|------------------|-------------------------------------------------------------------------|--------------|---------------------|
| Review By        | Sarabjit Jaswal                                                         | Review On    | 2/5/2025 7:11:46 AM |
| Supervise By     | Eman Mughal                                                             | Supervise On | 2/5/2025 9:21:10 AM |
| <b>STD. NAME</b> | <b>STD REF.#</b>                                                        |              |                     |
| ICAL Standard    | MP84042,MP84050,MP84049,MP84047,MP84046,MP84045,MP84044,MP84043,MP84073 |              |                     |
| ICV Standard     | MP84051                                                                 |              |                     |
| CCV Standard     | MP84052                                                                 |              |                     |
| ICSA Standard    | MP84053,MP84054                                                         |              |                     |
| CRI Standard     |                                                                         |              |                     |
| LCS Standard     | MP84118                                                                 |              |                     |
| Chk Standard     | MP84057,MP84056                                                         |              |                     |

|    |            |        |     |                |            |        |          |
|----|------------|--------|-----|----------------|------------|--------|----------|
| 38 | CCB012     | CCB012 | CCB | 02/02/25 16:29 |            | Jaswal | OK       |
| 39 | Q1246-01   | A6313  | SAM | 02/02/25 16:32 |            | Jaswal | OK       |
| 40 | Q1246-02   | A6314  | SAM | 02/02/25 16:35 |            | Jaswal | OK       |
| 41 | Q1246-03   | A6325  | SAM | 02/02/25 16:38 |            | Jaswal | OK       |
| 42 | Q1246-04   | A6326  | SAM | 02/02/25 16:42 |            | Jaswal | OK       |
| 43 | Q1246-05   | A6331  | SAM | 02/02/25 16:45 |            | Jaswal | OK       |
| 44 | Q1246-06   | A6332  | SAM | 02/02/25 16:48 |            | Jaswal | OK       |
| 45 | Q1246-07   | A6333  | SAM | 02/02/25 16:51 |            | Jaswal | OK       |
| 46 | Q1246-08   | A6334  | SAM | 02/02/25 16:54 |            | Jaswal | OK       |
| 47 | Q1223-01   | A6301  | SAM | 02/02/25 16:57 |            | Jaswal | OK       |
| 48 | Q1223-02   | A6302  | SAM | 02/02/25 17:01 | Cr,Mn High | Jaswal | Dilution |
| 49 | Q1223-03   | A6303  | SAM | 02/02/25 17:04 |            | Jaswal | OK       |
| 50 | Q1223-04   | A6304  | SAM | 02/02/25 17:07 |            | Jaswal | OK       |
| 51 | Q1246-12   | A6343  | SAM | 02/02/25 17:10 |            | Jaswal | OK       |
| 52 | Q1246-13   | A6344  | SAM | 02/02/25 17:13 |            | Jaswal | OK       |
| 53 | Q1246-14   | A6345  | SAM | 02/02/25 17:16 |            | Jaswal | OK       |
| 54 | Q1246-15   | A6346  | SAM | 02/02/25 17:20 |            | Jaswal | OK       |
| 55 | Q1246-16   | A6347  | SAM | 02/02/25 17:23 |            | Jaswal | OK       |
| 56 | Q1223-05DL | A6305  | SAM | 02/02/25 17:26 | 2X For Cr  | Jaswal | Confirms |
| 57 | Q1223-06DL | A6305D | DUP | 02/02/25 17:29 | 2X For Cr  | Jaswal | Confirms |

Instrument ID: P7

**Daily Analysis Runlog For Sequence/QC Batch ID # LB134522**

|                  |                                                                         |              |                     |
|------------------|-------------------------------------------------------------------------|--------------|---------------------|
| Review By        | Sarabjit Jaswal                                                         | Review On    | 2/5/2025 7:11:46 AM |
| Supervise By     | Eman Mughal                                                             | Supervise On | 2/5/2025 9:21:10 AM |
| <b>STD. NAME</b> | <b>STD REF.#</b>                                                        |              |                     |
| ICAL Standard    | MP84042,MP84050,MP84049,MP84047,MP84046,MP84045,MP84044,MP84043,MP84073 |              |                     |
| ICV Standard     | MP84051                                                                 |              |                     |
| CCV Standard     | MP84052                                                                 |              |                     |
| ICSA Standard    | MP84053,MP84054                                                         |              |                     |
| CRI Standard     |                                                                         |              |                     |
| LCS Standard     | MP84118                                                                 |              |                     |
| Chk Standard     | MP84057,MP84056                                                         |              |                     |

|    |             |        |     |                |                                                                                    |        |          |
|----|-------------|--------|-----|----------------|------------------------------------------------------------------------------------|--------|----------|
| 58 | Q1223-05LDL | A6305L | SD  | 02/02/25 17:32 | 10X For Cr                                                                         | Jaswal | Confirms |
| 59 | Q1223-07DL  | A6305S | MS  | 02/02/25 17:35 | 2X For Cr                                                                          | Jaswal | Confirms |
| 60 | Q1223-08DL  | A6306  | SAM | 02/02/25 17:38 | 2X For Cr                                                                          | Jaswal | Confirms |
| 61 | Q1223-10DL  | A6308  | SAM | 02/02/25 17:42 | 25X For Cr                                                                         | Jaswal | Confirms |
| 62 | Q1223-13DL  | A6321  | SAM | 02/02/25 17:45 | 5X For Cr                                                                          | Jaswal | Confirms |
| 63 | Q1223-14DL  | A6322  | SAM | 02/02/25 17:48 | 5X For Cr                                                                          | Jaswal | Confirms |
| 64 | Q1223-15DL  | A6323  | SAM | 02/02/25 17:51 | 5X For Cr                                                                          | Jaswal | Confirms |
| 65 | CCV013      | CCV013 | CCV | 02/02/25 17:54 |                                                                                    | Jaswal | OK       |
| 66 | CCB013      | CCB013 | CCB | 02/02/25 17:57 |                                                                                    | Jaswal | OK       |
| 67 | Q1223-02DL  | A6302  | SAM | 02/02/25 18:00 | 2X For Cr,Mn Still Mn High                                                         | Jaswal | Dilution |
| 68 | Q1223-02DL2 | A6302  | SAM | 02/02/25 18:04 | 5X For Mn                                                                          | Jaswal | Confirms |
| 69 | Q1223-05A   | A6305A | PS  | 02/02/25 18:07 | PS For Sb,Se,Zn(0.01ml<br>mp84317,0.04ml<br>mp84318,0.0044ml<br>m5819-10ml sample) | Jaswal | OK       |
| 70 | Q1246-09A   | A6342A | PS  | 02/02/25 18:10 | PS For Se(0.01ml<br>mp84317-10ml sample)                                           | Jaswal | OK       |
| 71 | CCV014      | CCV014 | CCV | 02/02/25 18:13 |                                                                                    | Jaswal | OK       |
| 72 | CCB014      | CCB014 | CCB | 02/02/25 18:16 |                                                                                    | Jaswal | OK       |

**Instrument ID:** P8

**Daily Analysis Runlog For Sequence/QC Batch ID # LB134612**

|              |                 |              |                      |
|--------------|-----------------|--------------|----------------------|
| Review By    | Sarabjit Jaswal | Review On    | 2/7/2025 9:21:46 AM  |
| Supervise By | Mohan Bera      | Supervise On | 2/7/2025 11:31:39 AM |

| STD. NAME     | STD REF.#                                                               |
|---------------|-------------------------------------------------------------------------|
| ICAL Standard | MP84042,MP84050,MP84049,MP84047,MP84046,MP84045,MP84044,MP84043,MP84073 |
| ICV Standard  | MP84051                                                                 |
| CCV Standard  | MP84052                                                                 |
| ICSA Standard | MP84053,MP84054                                                         |
| CRI Standard  |                                                                         |
| LCS Standard  | MP84118                                                                 |
| Chk Standard  | MP84057,MP84056                                                         |

| Sr# | SampleId  | ClientID | QcType | Date           | Comment                             | Operator | Status |
|-----|-----------|----------|--------|----------------|-------------------------------------|----------|--------|
| 1   | TUNE      | TUNE001  | TUNE   | 02/06/25 11:25 |                                     | Jaswal   | OK     |
| 2   | S00       | S0       | CAL1   | 02/06/25 12:08 |                                     | Jaswal   | OK     |
| 3   | S02       | S02      | CAL3   | 02/06/25 12:12 |                                     | Jaswal   | OK     |
| 4   | S03       | S03      | CAL4   | 02/06/25 12:15 |                                     | Jaswal   | OK     |
| 5   | S04       | S04      | CAL5   | 02/06/25 12:18 |                                     | Jaswal   | OK     |
| 6   | S05       | S05      | CAL6   | 02/06/25 12:21 |                                     | Jaswal   | OK     |
| 7   | S06       | S06      | CAL7   | 02/06/25 12:24 |                                     | Jaswal   | OK     |
| 8   | S07       | S07      | CAL8   | 02/06/25 12:27 |                                     | Jaswal   | OK     |
| 9   | S08       | S08      | CAL9   | 02/06/25 12:29 |                                     | Jaswal   | OK     |
| 10  | ICV001    | ICV001   | ICV    | 02/06/25 12:45 | Inst. paused for Calibration review | Jaswal   | OK     |
| 11  | ICB001    | ICB001   | ICB    | 02/06/25 12:50 |                                     | Jaswal   | OK     |
| 12  | ICSA001   | ICSA001  | ICSA   | 02/06/25 12:54 |                                     | Jaswal   | OK     |
| 13  | ICSAB001  | ICSAB001 | ICSAB  | 02/06/25 12:57 |                                     | Jaswal   | OK     |
| 14  | CCV001    | CCV001   | CCV    | 02/06/25 13:00 |                                     | Jaswal   | OK     |
| 15  | CCB001    | CCB001   | CCB    | 02/06/25 13:07 | Inst. paused for review             | Jaswal   | OK     |
| 16  | Q1177-01A | MH2GW9A  | PS     | 02/06/25 13:10 | PS For Se-0.01ml mp84378            | Jaswal   | OK     |
|     |           |          |        |                |                                     |          |        |

Instrument ID: P8

**Daily Analysis Runlog For Sequence/QC Batch ID # LB134612**

|              |                 |              |                      |
|--------------|-----------------|--------------|----------------------|
| Review By    | Sarabjit Jaswal | Review On    | 2/7/2025 9:21:46 AM  |
| Supervise By | Mohan Bera      | Supervise On | 2/7/2025 11:31:39 AM |

| STD. NAME     | STD REF.#                                                               |
|---------------|-------------------------------------------------------------------------|
| ICAL Standard | MP84042,MP84050,MP84049,MP84047,MP84046,MP84045,MP84044,MP84043,MP84073 |
| ICV Standard  | MP84051                                                                 |
| CCV Standard  | MP84052                                                                 |
| ICSA Standard | MP84053,MP84054                                                         |
| CRI Standard  |                                                                         |
| LCS Standard  | MP84118                                                                 |
| Chk Standard  | MP84057,MP84056                                                         |

|    |            |         |     |                |                                                                                                 |        |          |
|----|------------|---------|-----|----------------|-------------------------------------------------------------------------------------------------|--------|----------|
| 17 | Q1178-01A  | MH2GX0A | PS  | 02/06/25 13:14 | PS For Ba-0.018 ml<br>m6032, Cd-0.056 ml<br>mp84379, Ni-0.0031ml<br>m5748, Se-0.01ml<br>mp84378 | Jaswal | OK       |
| 18 | PB166302BL | PBS302  | MB  | 02/06/25 13:22 |                                                                                                 | Jaswal | OK       |
| 19 | PB166302BS | LCS302  | LCS | 02/06/25 13:26 |                                                                                                 | Jaswal | OK       |
| 20 | Q1159-01   | YE8C9   | SAM | 02/06/25 13:29 | Sr High in MS                                                                                   | Jaswal | Dilution |
| 21 | Q1159-02   | YE8C9D  | DUP | 02/06/25 13:32 | Sr High in MS                                                                                   | Jaswal | Dilution |
| 22 | Q1159-01L  | YE8C9L  | SD  | 02/06/25 13:35 | Sr High in MS                                                                                   | Jaswal | Dilution |
| 23 | Q1159-03   | YE8C9S  | MS  | 02/06/25 13:39 | Sr High , MS Fail in Se<br>Below RI)                                                            | Jaswal | Dilution |
| 24 | PB166317BL | PBW317  | MB  | 02/06/25 13:42 | Not Use                                                                                         | Jaswal | Not Ok   |
| 25 | PB166317BS | LCS317  | LCS | 02/06/25 13:45 | Not Use                                                                                         | Jaswal | Not Ok   |
| 26 | Q1186-01   | ME2948  | SAM | 02/06/25 13:48 | Not Use                                                                                         | Jaswal | Not Ok   |
| 27 | Q1186-02   | ME2949  | SAM | 02/06/25 13:52 |                                                                                                 | Jaswal | OK       |
| 28 | Q1186-03   | ME2955  | SAM | 02/06/25 13:55 |                                                                                                 | Jaswal | OK       |
| 29 | Q1186-04   | ME2956  | SAM | 02/06/25 13:58 |                                                                                                 | Jaswal | OK       |
| 30 | Q1186-05   | ME2957  | SAM | 02/06/25 14:01 |                                                                                                 | Jaswal | OK       |
| 31 | Q1186-06   | ME2960  | SAM | 02/06/25 14:05 |                                                                                                 | Jaswal | OK       |
| 32 | Q1186-07   | ME2961  | SAM | 02/06/25 14:08 |                                                                                                 | Jaswal | OK       |
| 33 | Q1186-08   | ME2962  | SAM | 02/06/25 14:11 |                                                                                                 | Jaswal | OK       |
| 34 | Q1186-09   | ME2959  | SAM | 02/06/25 14:14 |                                                                                                 | Jaswal | OK       |

Instrument ID: P8

**Daily Analysis Runlog For Sequence/QC Batch ID # LB134612**

|                  |                                                                         |              |                      |
|------------------|-------------------------------------------------------------------------|--------------|----------------------|
| Review By        | Sarabjit Jaswal                                                         | Review On    | 2/7/2025 9:21:46 AM  |
| Supervise By     | Mohan Bera                                                              | Supervise On | 2/7/2025 11:31:39 AM |
| <b>STD. NAME</b> | <b>STD REF.#</b>                                                        |              |                      |
| ICAL Standard    | MP84042,MP84050,MP84049,MP84047,MP84046,MP84045,MP84044,MP84043,MP84073 |              |                      |
| ICV Standard     | MP84051                                                                 |              |                      |
| CCV Standard     | MP84052                                                                 |              |                      |
| ICSA Standard    | MP84053,MP84054                                                         |              |                      |
| CRI Standard     |                                                                         |              |                      |
| LCS Standard     | MP84118                                                                 |              |                      |
| Chk Standard     | MP84057,MP84056                                                         |              |                      |

|    |             |         |     |                |         |        |        |
|----|-------------|---------|-----|----------------|---------|--------|--------|
| 35 | Q1186-10    | ME2959D | DUP | 02/06/25 14:18 |         | Jaswal | OK     |
| 36 | Q1186-09L   | ME2959L | SD  | 02/06/25 14:21 |         | Jaswal | OK     |
| 37 | Q1186-11    | ME2959S | MS  | 02/06/25 14:24 |         | Jaswal | OK     |
| 38 | Q1186-12    | ME2963  | SAM | 02/06/25 14:28 |         | Jaswal | OK     |
| 39 | Q1186-13    | ME2967  | SAM | 02/06/25 14:31 |         | Jaswal | OK     |
| 40 | Q1186-14    | ME2965  | SAM | 02/06/25 14:34 |         | Jaswal | OK     |
| 41 | CCV002      | CCV002  | CCV | 02/06/25 14:37 |         | Jaswal | OK     |
| 42 | CCB002      | CCB002  | CCB | 02/06/25 14:40 |         | Jaswal | OK     |
| 43 | Q1186-15    | ME2966  | SAM | 02/06/25 14:44 |         | Jaswal | OK     |
| 44 | Q1186-16    | ME2958  | SAM | 02/06/25 14:47 |         | Jaswal | OK     |
| 45 | Q1186-17    | ME2968  | SAM | 02/06/25 14:50 |         | Jaswal | OK     |
| 46 | Q1186-18    | ME2974  | SAM | 02/06/25 14:53 |         | Jaswal | OK     |
| 47 | Q1186-19    | ME2977  | SAM | 02/06/25 14:57 |         | Jaswal | OK     |
| 48 | Q1186-20    | ME2980  | SAM | 02/06/25 15:00 |         | Jaswal | OK     |
| 49 | Q1186-09RE  | ME2959  | SAM | 02/06/25 15:03 | Not Use | Jaswal | Not Ok |
| 50 | Q1186-10RE  | ME2959D | DUP | 02/06/25 15:06 | Not Use | Jaswal | Not Ok |
| 51 | Q1186-09LRE | ME2959L | SD  | 02/06/25 15:10 | Not Use | Jaswal | Not Ok |
| 52 | Q1186-11RE  | ME2959S | MS  | 02/06/25 15:13 | Not Use | Jaswal | Not Ok |
| 53 | PB166435BL  | PBW435  | MB  | 02/06/25 15:16 |         | Jaswal | OK     |
| 54 | PB166435BS  | LCS435  | LCS | 02/06/25 15:24 |         | Jaswal | OK     |

Instrument ID: P8

**Daily Analysis Runlog For Sequence/QC Batch ID # LB134612**

|                  |                                                                         |              |                      |
|------------------|-------------------------------------------------------------------------|--------------|----------------------|
| Review By        | Sarabjit Jaswal                                                         | Review On    | 2/7/2025 9:21:46 AM  |
| Supervise By     | Mohan Bera                                                              | Supervise On | 2/7/2025 11:31:39 AM |
| <b>STD. NAME</b> | <b>STD REF.#</b>                                                        |              |                      |
| ICAL Standard    | MP84042,MP84050,MP84049,MP84047,MP84046,MP84045,MP84044,MP84043,MP84073 |              |                      |
| ICV Standard     | MP84051                                                                 |              |                      |
| CCV Standard     | MP84052                                                                 |              |                      |
| ICSA Standard    | MP84053,MP84054                                                         |              |                      |
| CRI Standard     |                                                                         |              |                      |
| LCS Standard     | MP84118                                                                 |              |                      |
| Chk Standard     | MP84057,MP84056                                                         |              |                      |

|    |           |         |     |                |  |        |    |
|----|-----------|---------|-----|----------------|--|--------|----|
| 55 | Q1200-01  | ME2964  | SAM | 02/06/25 15:28 |  | Jaswal | OK |
| 56 | Q1200-02  | ME2964D | DUP | 02/06/25 15:32 |  | Jaswal | OK |
| 57 | Q1200-01L | ME2964L | SD  | 02/06/25 15:35 |  | Jaswal | OK |
| 58 | Q1200-03  | ME2964S | MS  | 02/06/25 15:38 |  | Jaswal | OK |
| 59 | Q1200-04  | ME2981  | SAM | 02/06/25 15:41 |  | Jaswal | OK |
| 60 | Q1200-05  | ME2982  | SAM | 02/06/25 15:44 |  | Jaswal | OK |
| 61 | Q1200-06  | ME2983  | SAM | 02/06/25 15:48 |  | Jaswal | OK |
| 62 | Q1200-07  | ME2984  | SAM | 02/06/25 15:51 |  | Jaswal | OK |
| 63 | Q1200-08  | ME2985  | SAM | 02/06/25 15:54 |  | Jaswal | OK |
| 64 | Q1200-09  | ME2986  | SAM | 02/06/25 15:57 |  | Jaswal | OK |
| 65 | Q1200-10  | ME2987  | SAM | 02/06/25 16:01 |  | Jaswal | OK |
| 66 | Q1200-11  | ME2988  | SAM | 02/06/25 16:04 |  | Jaswal | OK |
| 67 | Q1200-12  | ME2992  | SAM | 02/06/25 16:07 |  | Jaswal | OK |
| 68 | Q1200-13  | ME2994  | SAM | 02/06/25 16:10 |  | Jaswal | OK |
| 69 | Q1200-14  | ME2995  | SAM | 02/06/25 16:13 |  | Jaswal | OK |
| 70 | CCV003    | CCV003  | CCV | 02/06/25 16:17 |  | Jaswal | OK |
| 71 | CCB003    | CCB003  | CCB | 02/06/25 16:19 |  | Jaswal | OK |
| 72 | Q1200-15  | ME2999  | SAM | 02/06/25 16:23 |  | Jaswal | OK |
| 73 | Q1200-16  | ME2997  | SAM | 02/06/25 16:26 |  | Jaswal | OK |
| 74 | Q1200-17  | ME2998  | SAM | 02/06/25 16:29 |  | Jaswal | OK |

Instrument ID: P8

**Daily Analysis Runlog For Sequence/QC Batch ID # LB134612**

|                  |                                                                         |              |                      |
|------------------|-------------------------------------------------------------------------|--------------|----------------------|
| Review By        | Sarabjit Jaswal                                                         | Review On    | 2/7/2025 9:21:46 AM  |
| Supervise By     | Mohan Bera                                                              | Supervise On | 2/7/2025 11:31:39 AM |
| <b>STD. NAME</b> | <b>STD REF.#</b>                                                        |              |                      |
| ICAL Standard    | MP84042,MP84050,MP84049,MP84047,MP84046,MP84045,MP84044,MP84043,MP84073 |              |                      |
| ICV Standard     | MP84051                                                                 |              |                      |
| CCV Standard     | MP84052                                                                 |              |                      |
| ICSA Standard    | MP84053,MP84054                                                         |              |                      |
| CRI Standard     |                                                                         |              |                      |
| LCS Standard     | MP84118                                                                 |              |                      |
| Chk Standard     | MP84057,MP84056                                                         |              |                      |

|    |             |            |     |                |                         |        |        |
|----|-------------|------------|-----|----------------|-------------------------|--------|--------|
| 75 | Q1200-18    | ME29A0     | SAM | 02/06/25 16:33 |                         | Jaswal | OK     |
| 76 | Q1200-19    | ME29A1     | SAM | 02/06/25 16:36 |                         | Jaswal | OK     |
| 77 | Q1200-20    | ME29A2     | SAM | 02/06/25 16:39 |                         | Jaswal | OK     |
| 78 | Q1200-21    | ME29A3     | SAM | 02/06/25 16:42 |                         | Jaswal | OK     |
| 79 | WATER-01    | WATER-01   | SAM | 02/06/25 16:46 |                         | Jaswal | OK     |
| 80 | WATER-02RE  | WATER-02RE | SAM | 02/06/25 16:49 | Not Use                 | Jaswal | Not Ok |
| 81 | WATER-02    | WATER-02   | SAM | 02/06/25 16:52 |                         | Jaswal | OK     |
| 82 | CCV004      | CCV004     | CCV | 02/06/25 16:55 |                         | Jaswal | OK     |
| 83 | CCB004      | CCB004     | CCB | 02/06/25 17:02 | Inst. paused for review | Jaswal | OK     |
| 84 | PB166532BL  | PBS532     | MB  | 02/06/25 17:12 | Not Use                 | Jaswal | Not Ok |
| 85 | PB166532BS  | LCS532     | LCS | 02/06/25 17:15 | Not Use                 | Jaswal | Not Ok |
| 86 | Q1271-03    | 3189-3196  | SAM | 02/06/25 17:18 | Not Use                 | Jaswal | Not Ok |
| 87 | Q1271-03DUP | 3189-3196  | DUP | 02/06/25 17:21 | Not Use                 | Jaswal | Not Ok |
| 88 | Q1271-03L   | 3189-3196L | SD  | 02/06/25 17:25 | Not Use                 | Jaswal | Not Ok |
| 89 | Q1271-03MS  | 3189-3196  | MS  | 02/06/25 17:28 | Not Use                 | Jaswal | Not Ok |
| 90 | Q1271-03MSD | 3189-3196  | MSD | 02/06/25 17:31 | Not Use                 | Jaswal | Not Ok |
| 91 | Q1271-03A   | 3189-3196A | PS  | 02/06/25 17:34 | Not Use                 | Jaswal | Not Ok |
| 92 | CCV005      | CCV005     | CCV | 02/06/25 17:39 |                         | Jaswal | OK     |
| 93 | CCB005      | CCB005     | CCB | 02/06/25 17:42 |                         | Jaswal | OK     |
| 94 | PB166282BL  | PBW282     | MB  | 02/06/25 17:45 | Not Use                 | Jaswal | Not Ok |

Instrument ID: P8

**Daily Analysis Runlog For Sequence/QC Batch ID # LB134612**

|                  |                                                                         |              |                      |
|------------------|-------------------------------------------------------------------------|--------------|----------------------|
| Review By        | Sarabjit Jaswal                                                         | Review On    | 2/7/2025 9:21:46 AM  |
| Supervise By     | Mohan Bera                                                              | Supervise On | 2/7/2025 11:31:39 AM |
| <b>STD. NAME</b> | <b>STD REF.#</b>                                                        |              |                      |
| ICAL Standard    | MP84042,MP84050,MP84049,MP84047,MP84046,MP84045,MP84044,MP84043,MP84073 |              |                      |
| ICV Standard     | MP84051                                                                 |              |                      |
| CCV Standard     | MP84052                                                                 |              |                      |
| ICSA Standard    | MP84053,MP84054                                                         |              |                      |
| CRI Standard     |                                                                         |              |                      |
| LCS Standard     | MP84118                                                                 |              |                      |
| Chk Standard     | MP84057,MP84056                                                         |              |                      |

|     |             |                   |     |                |         |        |        |
|-----|-------------|-------------------|-----|----------------|---------|--------|--------|
| 95  | PB166282BS  | LCS282            | LCS | 02/06/25 17:49 | Not Use | Jaswal | Not Ok |
| 96  | Q1193-01    | TAPIAL3-MW04S-012 | SAM | 02/06/25 17:53 | Not Use | Jaswal | Not Ok |
| 97  | Q1193-02    | TAPIAL2-MW01-0124 | SAM | 02/06/25 17:56 | Not Use | Jaswal | Not Ok |
| 98  | Q1193-02DUP | TAPIAL2-MW01-0124 | DUP | 02/06/25 17:59 | Not Use | Jaswal | Not Ok |
| 99  | Q1193-02L   | TAPIAL2-MW01-0124 | SD  | 02/06/25 18:02 | Not Use | Jaswal | Not Ok |
| 100 | Q1193-02MS  | TAPIAL2-MW01-0124 | MS  | 02/06/25 18:06 | Not Use | Jaswal | Not Ok |
| 101 | Q1193-02MSD | TAPIAL2-MW01-0124 | MSD | 02/06/25 18:08 | Not Use | Jaswal | Not Ok |
| 102 | Q1193-02A   | TAPIAL2-MW01-0124 | PS  | 02/06/25 18:11 | Not Use | Jaswal | Not Ok |
| 103 | CCV006      | CCV006            | CCV | 02/06/25 18:14 |         | Jaswal | OK     |
| 104 | CCB006      | CCB006            | CCB | 02/06/25 18:18 |         | Jaswal | OK     |
| 105 | PB166378BL  | PBW378            | MB  | 02/06/25 18:21 | Not Use | Jaswal | Not Ok |
| 106 | PB166378BS  | LCS378            | LCS | 02/06/25 18:25 | Not Use | Jaswal | Not Ok |
| 107 | Q1201-01    | TAPHHA-MW12-0127  | SAM | 02/06/25 18:28 | Not Use | Jaswal | Not Ok |
| 108 | Q1211-01    | TAPHHA-MW01-0128  | SAM | 02/06/25 18:31 | Not Use | Jaswal | Not Ok |
| 109 | Q1211-02    | TAPIAL2-MW03-0128 | SAM | 02/06/25 18:34 | Not Use | Jaswal | Not Ok |
| 110 | Q1211-02DUP | TAPIAL2-MW03-0128 | DUP | 02/06/25 18:37 | Not Use | Jaswal | Not Ok |
| 111 | Q1211-02L   | TAPIAL2-MW03-0128 | SD  | 02/06/25 18:40 | Not Use | Jaswal | Not Ok |
| 112 | Q1211-02MS  | TAPIAL2-MW03-0128 | MS  | 02/06/25 18:44 | Not Use | Jaswal | Not Ok |
| 113 | Q1211-02MSD | TAPIAL2-MW03-0128 | MSD | 02/06/25 18:46 | Not Use | Jaswal | Not Ok |
| 114 | Q1211-02A   | TAPIAL2-MW03-0128 | PS  | 02/06/25 18:49 | Not Use | Jaswal | Not Ok |

Instrument ID: P8

**Daily Analysis Runlog For Sequence/QC Batch ID # LB134612**

|              |                 |              |                      |
|--------------|-----------------|--------------|----------------------|
| Review By    | Sarabjit Jaswal | Review On    | 2/7/2025 9:21:46 AM  |
| Supervise By | Mohan Bera      | Supervise On | 2/7/2025 11:31:39 AM |

| STD. NAME     | STD REF.#                                                               |
|---------------|-------------------------------------------------------------------------|
| ICAL Standard | MP84042,MP84050,MP84049,MP84047,MP84046,MP84045,MP84044,MP84043,MP84073 |
| ICV Standard  | MP84051                                                                 |
| CCV Standard  | MP84052                                                                 |
| ICSA Standard | MP84053,MP84054                                                         |
| CRI Standard  |                                                                         |
| LCS Standard  | MP84118                                                                 |
| Chk Standard  | MP84057,MP84056                                                         |

|     |            |         |     |                |                                       |        |        |
|-----|------------|---------|-----|----------------|---------------------------------------|--------|--------|
| 115 | CCV007     | CCV007  | CCV | 02/06/25 18:52 |                                       | Jaswal | OK     |
| 116 | CCB007     | CCB007  | CCB | 02/06/25 18:58 | Inst. paused for review               | Jaswal | OK     |
| 117 | Q1159-01A  | YE8C9A  | PS  | 02/06/25 19:55 | PS For Se(0.01ML MP84378)             | Jaswal | OK     |
| 118 | Q1186-09A  | ME2959A | PS  | 02/06/25 19:58 | PS For Se(0.01ML MP84378-10ML SAMPLE) | Jaswal | OK     |
| 119 | Q1200-01A  | ME2964A | PS  | 02/06/25 20:01 | PS For Se(0.01ML MP84378-10ML SAMPLE) | Jaswal | OK     |
| 120 | PB166570BL | PBW570  | MB  | 02/06/25 20:04 |                                       | Jaswal | OK     |
| 121 | PB166570BS | LCS570  | LCS | 02/06/25 20:08 |                                       | Jaswal | OK     |
| 122 | Q1135-09   | YE8H7   | SAM | 02/06/25 20:11 | Not Use                               | Jaswal | Not Ok |
| 123 | Q1223-17   | A6310   | SAM | 02/06/25 20:14 |                                       | Jaswal | OK     |
| 124 | CCV008     | CCV008  | CCV | 02/06/25 20:22 |                                       | Jaswal | OK     |
| 125 | CCB008     | CCB008  | CCB | 02/06/25 20:25 |                                       | Jaswal | OK     |