

MDCZN8A

POST-DIGESTION/DISTILLATION SPIKE SAMPLE RECOVERY

Concentration Units (µg/L, mg/L or mg/kg dry weight): ug/L

| Analyte | Control Limit %R | Spiked Sample Result (SSR) Q | | Sample Result (SR) Q | | Spike Added (SA) | %R | Q |
|----------|------------------|------------------------------|--|----------------------|---|------------------|-----|---|
| Selenium | 75 - 125 | 12 | | 5.0 | U | 10.0 | 116 | |