

**SOLID EPH LABORATORY CONTROL SPIKE/LABORATORY CONTROL SPIKE DUPLICATE RECOVERY**

**Lab Name:** Chemtech **Client:** ENTACT  
**Lab Code:** CHEM **Cas No:** Q1308 **SAS No :** Q1308 **SDG No:** Q1308  
**Sample No :** PB166603BS **Datafile:** FC068230.D  
**Client ID :** PB166603BS

| COMPOUND          | SPIKE<br>ADDED<br>mg/kg | LCS/LCSD<br>CONCENTRATION<br>mg/kg | % REC | Qual | QC LIMITS |
|-------------------|-------------------------|------------------------------------|-------|------|-----------|
| Aliphatic C9-C12  | 10                      | 4.57                               | 46    |      | (40-140)  |
| Aliphatic C12-C16 | 6.7                     | 4.09                               | 61    |      | (40-140)  |
| Aliphatic C16-C21 | 10                      | 7.01                               | 70    |      | (40-140)  |
| Aliphatic C21-C28 | 13.3                    | 9.65                               | 72    |      | (40-140)  |
| Aliphatic C28-C40 | 20.0                    | 17.2                               | 86    |      | (40-140)  |