

Surrogate Summary

SDG No.: Q1487

Client: JPCL Engineering

Analytical Method: 8082A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Rec	Qual	Limits	
								Low	High
I.BLK-PP069995.D	PIBLK-PP069995.D	Tetrachloro-m-xylene	1	20	21.9	109		60	140
		Decachlorobiphenyl	1	20	21.8	109		60	140
		Tetrachloro-m-xylene	2	20	22.1	110		60	140
		Decachlorobiphenyl	2	20	21.6	108		60	140
I.BLK-PP070263.D	PIBLK-PP070263.D	Tetrachloro-m-xylene	1	20	20.1	100		60	140
		Decachlorobiphenyl	1	20	18.8	94		60	140
		Tetrachloro-m-xylene	2	20	20.9	105		60	140
		Decachlorobiphenyl	2	20	20.3	102		60	140
PB166985BL	PB166985BL	Tetrachloro-m-xylene	1	20	24.9	125		32	144
		Decachlorobiphenyl	1	20	22.2	111		32	175
		Tetrachloro-m-xylene	2	20	24.1	121		32	144
		Decachlorobiphenyl	2	20	23.8	119		32	175
PB166985BS	PB166985BS	Tetrachloro-m-xylene	1	20	25.4	127		32	144
		Decachlorobiphenyl	1	20	22.5	113		32	175
		Tetrachloro-m-xylene	2	20	24.1	121		32	144
		Decachlorobiphenyl	2	20	23.7	119		32	175
Q1487-01	DN-B-42	Tetrachloro-m-xylene	1	20	23.6	118		32	144
		Decachlorobiphenyl	1	20	18.1	90		32	175
		Tetrachloro-m-xylene	2	20	25.9	130		32	144
		Decachlorobiphenyl	2	20	19.3	97		32	175
I.BLK-PP070278.D	PIBLK-PP070278.D	Tetrachloro-m-xylene	1	20	21.2	106		60	140
		Decachlorobiphenyl	1	20	19.6	98		60	140
		Tetrachloro-m-xylene	2	20	20.7	104		60	140
		Decachlorobiphenyl	2	20	20.6	103		60	140
Q1488-03MS	ENV-101-SB02MS	Tetrachloro-m-xylene	1	20	24.9	125		32	144
		Decachlorobiphenyl	1	20	19.6	98		32	175
		Tetrachloro-m-xylene	2	20	23.8	119		32	144
		Decachlorobiphenyl	2	20	21.1	106		32	175
Q1488-03MSD	ENV-101-SB02MSD	Tetrachloro-m-xylene	1	20	26.5	132		32	144
		Decachlorobiphenyl	1	20	22.3	112		32	175
		Tetrachloro-m-xylene	2	20	24.6	123		32	144
		Decachlorobiphenyl	2	20	24.1	121		32	175
I.BLK-PP070293.D	PIBLK-PP070293.D	Tetrachloro-m-xylene	1	20	21.6	108		60	140
		Decachlorobiphenyl	1	20	20.2	101		60	140
		Tetrachloro-m-xylene	2	20	21.7	109		60	140
		Decachlorobiphenyl	2	20	21.4	107		60	140