

Instrument ID: ECD_Q

Daily Analysis Runlog For Sequence/QC Batch ID # PQ022625

Review By	yogesh	Review On	2/27/2025 9:17:58 AM
Supervise By	Ankita	Supervise On	2/27/2025 9:19:38 AM
SubDirectory	PQ022625	HP Acquire Method	HP Processing Method PQ022625
STD. NAME	STD REF.#		
Tune/Reschk			
Initial Calibration Stds	PP24020,PP24019,PP24018,PP24017,PP24016,PP24025,PP24024,PP24023,PP24022,PP24021,PP24030,PP24029,PP24028,PP24027,PP24026,PP24035,PP24034,PP24033,PP24032,PP24031,PP24040,PP24039,PP24038,PP24037,PP24036,PP24045,PP24044,PP24043,PP24042,PP24041,PP24050,PP24049,PP24048,PP24047,PP24046,PP24055,PP24054,PP24053,PP24052,PP24051		
CCC	PP24018,PP24023,PP24028,PP24033,PP24038,PP24043,PP24048,PP24053		
Internal Standard/PEM			
ICV/I.BLK			
Surrogate Standard			
MS/MSD Standard			
LCS Standard			

Sr#	SampleId	Data File Name	Date-Time	Operator	Status
1	HEXANE	PQ069996.D	26 Feb 2025 12:13	YP\AJ	Ok
2	AIBLK011	PQ069997.D	26 Feb 2025 12:27	YP\AJ	Ok
3	AR1660ICC100	PQ069998.D	26 Feb 2025 12:42	YP\AJ	Ok
4	AR1660ICC200	PQ069999.D	26 Feb 2025 12:56	YP\AJ	Ok
5	AR1660ICC400	PQ070000.D	26 Feb 2025 13:11	YP\AJ	Ok
6	AR1660ICC800	PQ070001.D	26 Feb 2025 13:26	YP\AJ	Ok
7	AR1660ICC1600	PQ070002.D	26 Feb 2025 13:40	YP\AJ	Ok
8	AR1221ICC100	PQ070003.D	26 Feb 2025 13:55	YP\AJ	Ok
9	AR1232ICC100	PQ070004.D	26 Feb 2025 14:09	YP\AJ	Ok
10	AR1242ICC100	PQ070005.D	26 Feb 2025 14:24	YP\AJ	Ok
11	AR1248ICC100	PQ070006.D	26 Feb 2025 14:38	YP\AJ	Ok
12	AR1254ICC100	PQ070007.D	26 Feb 2025 14:53	YP\AJ	Ok
13	AR1262ICC100	PQ070008.D	26 Feb 2025 15:08	YP\AJ	Ok
14	AR1268ICC100	PQ070009.D	26 Feb 2025 15:22	YP\AJ	Ok
15	AIBLK012	PQ070010.D	26 Feb 2025 15:37	YP\AJ	Ok
16	AR1660CCC400	PQ070011.D	26 Feb 2025 15:51	YP\AJ	Ok
17	Q1403-13	PQ070012.D	26 Feb 2025 16:06	YP\AJ	Ok
18	Q1403-14	PQ070013.D	26 Feb 2025 16:20	YP\AJ	Ok
19	AIBLK013	PQ070014.D	26 Feb 2025 16:35	YP\AJ	Ok
20	AR1660CCC400	PQ070015.D	26 Feb 2025 16:50	YP\AJ	Ok

M : Manual Integration

Instrument ID:

Daily Analysis Runlog For Sequence/QCBatch ID #

Review By		Review On	
Supervise By		Supervise On	
SubDirectory		HP Acquire Method	HP Processing Method
STD. NAME	STD REF.#		
Tune/Reschk Initial Calibration Stds			
CCC			
Internal Standard/PEM			
ICV/I.BLK			
Surrogate Standard			
MS/MSD Standard			
LCS Standard			

Sr#	SampleId	Data File Name	Date-Time	Operator	Status

M : Manual Integration