

Instrument ID: MSVOA_L

Daily Analysis Runlog For Sequence/QC Batch ID # VL022525

Review By	Semsettin Yesilyurt	Review On	2/26/2025 8:27:25 AM
Supervise By	Maresh Dadoda	Supervise On	2/26/2025 1:45:21 PM
SubDirectory	VL022525	HP Acquire Method	TO15AIR.M HP Processing Method VL022525AIR.M

STD. NAME	STD REF.#
Tune/Reschk	AP2576
Initial Calibration Stds	AP2582,AP2584,AP2585
CCC	AP2582
Internal Standard/PEM	AP2576
ICV/I.BLK	AP2586
Surrogate Standard	
MS/MSD Standard	
LCS Standard	

Sr#	SampleID	ClientID	Data File Name	Date-Time	Comment	Operator	Status
1	BFB	BFB	VL042040.D	25 Feb 2025 07:56		SY/MD	Ok
2	VSTDICCC010	VSTDICCC010	VL042041.D	25 Feb 2025 08:56		SY/MD	Ok,M
3	VSTDICC002	VSTDICC002	VL042042.D	25 Feb 2025 09:29		SY/MD	Ok,M
4	VSTDICC001	VSTDICC001	VL042043.D	25 Feb 2025 10:00		SY/MD	Ok,M
5	VSTDICC0.5	VSTDICC0.5	VL042044.D	25 Feb 2025 10:31		SY/MD	Ok,M
6	VSTDICC0.1	VSTDICC0.1	VL042045.D	25 Feb 2025 11:03		SY/MD	Ok,M
7	VSTDICC0.03	VSTDICC0.03	VL042046.D	25 Feb 2025 11:34		SY/MD	Ok,M
8	VSTDICC015	VSTDICC015	VL042047.D	25 Feb 2025 12:07		SY/MD	Ok,M
9	VSTDICV010	ICVVL022525	VL042048.D	25 Feb 2025 13:15		SY/MD	Ok,M

M : Manual Integration

Instrument ID: MSVOA_L

Daily Analysis Runlog For Sequence/QC Batch ID # VL031725

Review By	Mahesh Dadoda	Review On	3/18/2025 1:57:03 PM
Supervise By	Semsettin Yesilyurt	Supervise On	3/18/2025 1:57:44 PM
SubDirectory	VL031725	HP Acquire Method	TO15AIR.M
		HP Processing Method	VL022525AIR.M
STD. NAME	STD REF.#		
Tune/Reschk	AP2588		
Initial Calibration Stds	AP2582,AP2584,AP2585		
CCC	AP2582		
Internal Standard/PEM	AP2588		
ICV/I.BLK	AP2586		
Surrogate Standard			
MS/MSD Standard			
LCS Standard			

Sr#	SampleID	ClientID	Data File Name	Date-Time	Comment	Operator	Status
1	BFB	BFB	VL042149.D	17 Mar 2025 09:09		SY/MD	Ok
2	VSTDCCC010	VSTDCCC010	VL042150.D	17 Mar 2025 09:55		SY/MD	Ok,M
3	VL0317ABL01	VL0317ABL01	VL042151.D	17 Mar 2025 10:50		SY/MD	Ok
4	VL0317ABS01	VL0317ABS01	VL042152.D	17 Mar 2025 11:38		SY/MD	Ok,M
5	CAN 10326	CAN 10326	VL042153.D	17 Mar 2025 12:20	carryover	SY/MD	Not Ok
6	Q1586-01	565-SG-1	VL042154.D	17 Mar 2025 12:51	Need 40X	SY/MD	Dilution
7	Q1586-01DL	565-SG-1DL	VL042155.D	17 Mar 2025 13:37	Need 600X	SY/MD	Dilution
8	Q1586-02	557-SG-1	VL042156.D	17 Mar 2025 14:09	Need Straight Run	SY/MD	Not Ok
9	Q1586-03DL	557-SG-2DL	VL042157.D	17 Mar 2025 14:41		SY/MD	Ok
10	Q1586-03	557-SG-2	VL042158.D	17 Mar 2025 15:19	Need 10X	SY/MD	Dilution
11	Q1586-02	557-SG-1	VL042159.D	17 Mar 2025 15:53		SY/MD	Ok,M
12	Q1586-02DUP	557-SG-1DUP	VL042160.D	17 Mar 2025 16:29		SY/MD	Ok,M
13	Q1586-01DL2	565-SG-1DL2	VL042161.D	17 Mar 2025 17:01		SY/MD	Ok
14	Q1586-01DL2	565-SG-1DL2	VL042162.D	17 Mar 2025 17:33	not required	SY/MD	Not Ok

M : Manual Integration