

Surrogate Summary

SDG No.: Q1609

Client: ENTACT

Analytical Method: 8151A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Rec	Qual	Limits	
								Low	High
I.BLK-PS029557.D	PIBLK-PS029557.D	2,4-DCAA	1	500	455	91		70 (39)	130 (175)
		2,4-DCAA	2	500	523	105		70 (39)	130 (175)
I.BLK-PS029564.D	PIBLK-PS029564.D	2,4-DCAA	1	500	438	88		70 (39)	130 (175)
		2,4-DCAA	2	500	523	105		70 (39)	130 (175)
Q1609-03	WC-SCRN-01-C	2,4-DCAA	1	500	470	94		70 (39)	130 (175)
		2,4-DCAA	2	500	483	97		70 (39)	130 (175)
I.BLK-PS029575.D	PIBLK-PS029575.D	2,4-DCAA	1	500	501	100		70 (39)	130 (175)
		2,4-DCAA	2	500	540	108		70 (39)	130 (175)
Q1609-03MS	WC-SCRN-01-CMS	2,4-DCAA	1	500	489	98		70 (39)	130 (175)
		2,4-DCAA	2	500	488	98		70 (39)	130 (175)
Q1609-03MSD	WC-SCRN-01-CMSD	2,4-DCAA	1	500	494	99		70 (39)	130 (175)
		2,4-DCAA	2	500	487	97		70 (39)	130 (175)
I.BLK-PS029587.D	PIBLK-PS029587.D	2,4-DCAA	1	500	518	104		70 (39)	130 (175)
		2,4-DCAA	2	500	541	108		70 (39)	130 (175)
PB167312BL	PB167312BL	2,4-DCAA	1	500	495	99		70 (39)	130 (175)
		2,4-DCAA	2	500	483	97		70 (39)	130 (175)
PB167233TB	PB167233TB	2,4-DCAA	1	500	495	99		70 (39)	130 (175)
		2,4-DCAA	2	500	477	95		70 (39)	130 (175)
I.BLK-PS029599.D	PIBLK-PS029599.D	2,4-DCAA	1	500	528	106		70 (39)	130 (175)
		2,4-DCAA	2	500	545	109		70 (39)	130 (175)
PB167312BS	PB167312BS	2,4-DCAA	1	500	545	109		70 (39)	130 (175)
		2,4-DCAA	2	500	531	106		70 (39)	130 (175)
PB167275TB	PB167275TB	2,4-DCAA	1	500	528	106		70 (39)	130 (175)
		2,4-DCAA	2	500	521	104		70 (39)	130 (175)
I.BLK-PS029611.D	PIBLK-PS029611.D	2,4-DCAA	1	500	529	106		70 (39)	130 (175)
		2,4-DCAA	2	500	554	111		70 (39)	130 (175)