

RETENTION TIMES OF INITIAL CALIBRATION

Contract: WALS01

Lab Code: CHEM Case No.: Q1626 SAS No.: Q1626 SDG NO.: Q1626

Instrument ID: ECD_S Calibration Date(s): 04/02/2025 04/02/2025

Calibration Times: 17:32 20:44

GC Column: RTX-CLP ID: 0.32 (mm)

LAB FILE ID:	RT 200 = <u>PS029657.D</u>	RT 500 = <u>PS029658.D</u>
	RT 750 = <u>PS029659.D</u>	RT 1000 = <u>PS029660.D</u>
		RT 1500 = <u>PS029661.D</u>

COMPOUND	RT 200	RT 500	RT 750	RT 1000	RT 1500	MEAN RT	RT WINDOW	
							FROM	TO
2,4,5-T	9.17	9.17	9.17	9.17	9.17	9.17	9.07	9.27
2,4,5-TP(Silvex)	8.89	8.89	8.89	8.89	8.89	8.89	8.79	8.99
2,4-D	8.04	8.04	8.04	8.04	8.04	8.04	7.94	8.14
2,4-DB	9.73	9.73	9.73	9.73	9.73	9.73	9.63	9.83
2,4-DCAA	6.96	6.96	6.96	6.96	6.96	6.96	6.86	7.06
DICAMBA	7.14	7.14	7.14	7.14	7.14	7.14	7.04	7.24
DICHLORPROP	7.82	7.82	7.82	7.82	7.82	7.82	7.72	7.92
Dinoseb	10.90	10.90	10.90	10.90	10.89	10.90	10.80	11.00

RETENTION TIMES OF INITIAL CALIBRATION

Contract:	<u>WALS01</u>						
Lab Code:	<u>CHEM</u>	Case No.:	<u>Q1626</u>	SAS No.:	<u>Q1626</u>	SDG NO.:	<u>Q1626</u>
Instrument ID:	<u>ECD_S</u>	Calibration Date(s):		<u>04/02/2025</u>	<u>04/02/2025</u>		
		Calibration Times:		<u>17:32</u>	<u>20:44</u>		

GC Column: RTX-CLP2 ID: 0.32 (mm)

LAB FILE ID:	RT 200 =	<u>PS029657.D</u>	RT 500 =	<u>PS029658.D</u>
	RT 750 =	PS029659.D	RT 1000 =	PS029660.D
			RT 1500 =	PS029661.D

[illegible]