

Surrogate Summary

SDG No.: Q1641

Client: ENTACT

Analytical Method: 8082A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Rec	Qual	Limits	
								Low	High
I.BLK-PP070418.D	PIBLK-PP070418.D	Tetrachloro-m-xylene	1	20	22.0	110		70 (60)	130 (140)
		Decachlorobiphenyl	1	20	23.1	116		70 (60)	130 (140)
		Tetrachloro-m-xylene	2	20	23.0	115		70 (60)	130 (140)
		Decachlorobiphenyl	2	20	23.5	117		70 (60)	130 (140)
I.BLK-PP070900.D	PIBLK-PP070900.D	Tetrachloro-m-xylene	1	20	19.0	95		70 (60)	130 (140)
		Decachlorobiphenyl	1	20	20.8	104		70 (60)	130 (140)
		Tetrachloro-m-xylene	2	20	20.4	102		70 (60)	130 (140)
		Decachlorobiphenyl	2	20	17.4	87		70 (60)	130 (140)
PB167324BL	PB167324BL	Tetrachloro-m-xylene	1	20	21.2	106		30 (16)	150 (158)
		Decachlorobiphenyl	1	20	22.7	114		30 (10)	150 (173)
		Tetrachloro-m-xylene	2	20	22.2	111		30 (16)	150 (158)
		Decachlorobiphenyl	2	20	20.2	101		30 (10)	150 (173)
PB167324BS	PB167324BS	Tetrachloro-m-xylene	1	20	21.2	106		30 (16)	150 (158)
		Decachlorobiphenyl	1	20	22.4	112		30 (10)	150 (173)
		Tetrachloro-m-xylene	2	20	21.1	105		30 (16)	150 (158)
		Decachlorobiphenyl	2	20	20.2	101		30 (10)	150 (173)
PB167324BSD	PB167324BSD	Tetrachloro-m-xylene	1	20	20.6	103		30 (16)	150 (158)
		Decachlorobiphenyl	1	20	21.9	110		30 (10)	150 (173)
		Tetrachloro-m-xylene	2	20	21.8	109		30 (16)	150 (158)
		Decachlorobiphenyl	2	20	19.1	95		30 (10)	150 (173)
Q1641-02	NP-WS-003	Tetrachloro-m-xylene	1	20	17.4	87		30 (16)	150 (158)
		Decachlorobiphenyl	1	20	15.8	79		30 (10)	150 (173)
		Tetrachloro-m-xylene	2	20	19.6	98		30 (16)	150 (158)
		Decachlorobiphenyl	2	20	14.4	72		30 (10)	150 (173)
I.BLK-PP070914.D	PIBLK-PP070914.D	Tetrachloro-m-xylene	1	20	18.4	92		70 (60)	130 (140)
		Decachlorobiphenyl	1	20	21.4	107		70 (60)	130 (140)
		Tetrachloro-m-xylene	2	20	19.0	95		70 (60)	130 (140)
		Decachlorobiphenyl	2	20	17.5	88		70 (60)	130 (140)