

Instrument ID: MSVOA_L

Daily Analysis Runlog For Sequence/QC Batch ID # VL032725

Review By	Maresh Dadoda	Review On	4/3/2025 8:14:07 AM
Supervise By	Semsettin Yesilyurt	Supervise On	4/3/2025 8:45:16 AM
SubDirectory	VL032725	HP Acquire Method	TO15AIR.M HP Processing Method VL032725AIR.M

STD. NAME	STD REF.#
Tune/Reschk	AP2589
Initial Calibration Stds	AP2593,AP2595,AP2596
CCC	AP2593
Internal Standard/PEM	AP2589
ICV/I.BLK	AP2594
Surrogate Standard	
MS/MSD Standard	
LCS Standard	

Sr#	SampleID	ClientID	Data File Name	Date-Time	Comment	Operator	Status
1	BFB	BFB	VL042171.D	27 Mar 2025 07:44		SY/MD	Ok
2	VSTDICCC010	VSTDICCC010	VL042172.D	27 Mar 2025 08:26		SY/MD	Ok,M
3	VSTDICC002	VSTDICC002	VL042173.D	27 Mar 2025 08:59		SY/MD	Ok,M
4	VSTDICC001	VSTDICC001	VL042174.D	27 Mar 2025 09:30		SY/MD	Ok,M
5	VSTDICC0.5	VSTDICC0.5	VL042175.D	27 Mar 2025 10:01		SY/MD	Ok,M
6	VSTDICC0.1	VSTDICC0.1	VL042176.D	27 Mar 2025 10:32		SY/MD	Ok,M
7	VSTDICC0.03	VSTDICC0.03	VL042177.D	27 Mar 2025 11:03		SY/MD	Ok,M
8	VSTDICC015	VSTDICC015	VL042178.D	27 Mar 2025 11:36		SY/MD	Ok,M
9	VSTDICV010	ICVVL032725	VL042179.D	27 Mar 2025 12:23		SY/MD	Ok,M
10	VL0327ABL01	VL0327ABL01	VL042180.D	27 Mar 2025 13:50		SY/MD	Ok
11	VL0327ABS01	VL0327ABS01	VL042181.D	27 Mar 2025 14:39		SY/MD	Ok,M
12	QC032625 CAN 10448	QC032625 CAN 10448	VL042182.D	27 Mar 2025 16:41		SY/MD	Ok
13	Q1649-01	SV1	VL042183.D	27 Mar 2025 17:13	rerun 4X,not required	SY/MD	Not Ok
14	Q1649-02	SV2	VL042184.D	27 Mar 2025 17:45	rerun 4X,not required	SY/MD	Not Ok
15	Q1669-01	IA-B9-3-3262025	VL042185.D	27 Mar 2025 18:22	Need 5x	SY/MD	Dilution
16	Q1669-02	IA-B9-2-3262025	VL042186.D	27 Mar 2025 18:59		SY/MD	Ok,M
17	Q1669-03	IA-B9-1-3262025	VL042187.D	27 Mar 2025 19:36		SY/MD	Ok,M
18	Q1669-04	AA-2-3252025	VL042188.D	27 Mar 2025 20:12		SY/MD	Ok,M

Instrument ID: MSVOA_L

Daily Analysis Runlog For Sequence/QC Batch ID # VL032725

Review By	Maresh Dadoda	Review On	4/3/2025 8:14:07 AM
Supervise By	Semsettin Yesilyurt	Supervise On	4/3/2025 8:45:16 AM
SubDirectory	VL032725	HP Acquire Method	TO15AIR.M HP Processing Method VL032725AIR.M
STD. NAME	STD REF.#		
Tune/Reschk	AP2589		
Initial Calibration Stds	AP2593,AP2595,AP2596		
CCC	AP2593		
Internal Standard/PEM	AP2589		
ICV/I.BLK	AP2594		
Surrogate Standard			
MS/MSD Standard			
LCS Standard			

19	Q1669-09	AA-3-3252025	VL042189.D	27 Mar 2025 20:48		SY/MD	Ok,M
20	Q1669-09DUP	AA-3-3252025DUP	VL042190.D	27 Mar 2025 21:24		SY/MD	Ok,M
21	Q1669-05	SSSG-B16-3-3262025	VL042191.D	27 Mar 2025 21:57	Need straight run	SY/MD	Not Ok
22	Q1669-07DL	SSSG-B16-1-3262025	VL042192.D	27 Mar 2025 22:29		SY/MD	Ok,M
23	Q1669-08	SSSG-B16-2-3262025	VL042193.D	27 Mar 2025 23:01	Need straight run	SY/MD	Not Ok
24	VIBLK	VIBLK	VL042194.D	27 Mar 2025 23:34		SY/MD	Ok
25	Q1669-01DL	IA-B9-3-3262025DL	VL042195.D	28 Mar 2025 00:06		SY/MD	Ok,M
26	Q1669-01DUP	IA-B9-3-3262025DUP	VL042196.D	28 Mar 2025 00:44		SY/MD	Ok,M
27	Q1669-05	SSSG-B16-3-3262025	VL042197.D	28 Mar 2025 01:20		SY/MD	Ok,M
28	Q1669-07	SSSG-B16-1-3262025	VL042198.D	28 Mar 2025 01:57	Need 10X	SY/MD	Dilution
29	Q1669-08	SSSG-B16-2-3262025	VL042199.D	28 Mar 2025 02:33		SY/MD	Ok,M

M : Manual Integration

Instrument ID: MSVOA_L

Daily Analysis Runlog For Sequence/QC Batch ID # VL040925

Review By	Semsettin Yesilyurt	Review On	4/10/2025 7:20:08 PM
Supervise By	Maresh Dadoda	Supervise On	4/10/2025 7:21:16 PM
SubDirectory	VL040925	HP Acquire Method	TO15AIR.M
		HP Processing Method	VL032725AIR.M
STD. NAME	STD REF.#		
Tune/Reschk	AP2591		
Initial Calibration Stds	AP2593,AP2595,AP2596		
CCC	AP2593		
Internal Standard/PEM	AP2591		
ICV/I.BLK	AP2594		
Surrogate Standard			
MS/MSD Standard			
LCS Standard			

Sr#	Sampleld	ClientID	Data File Name	Date-Time	Comment	Operator	Status
1	BFB	BFB	VL042285.D	09 Apr 2025 09:38		SY/MD	Ok
2	VSTDCCC010	VSTDCCC010	VL042286.D	09 Apr 2025 10:20		SY/MD	Ok,M
3	VL0409ABL01	VL0409ABL01	VL042287.D	09 Apr 2025 11:03		SY/MD	Ok
4	VL0409ABS01	VL0409ABS01	VL042288.D	09 Apr 2025 11:46	Recovery fail	SY/MD	Not Ok
5	VL0409ABS02	VL0409ABS02	VL042289.D	09 Apr 2025 12:31		SY/MD	Ok,M
6	Q1665-01	AA-DUP-1-3242025	VL042290.D	09 Apr 2025 13:25		SY/MD	Ok,M
7	Q1665-01DUP	AA-DUP-1-3242025DU	VL042291.D	09 Apr 2025 14:02		SY/MD	Ok,M
8	Q1665-04	AA-1-3242025	VL042292.D	09 Apr 2025 14:36		SY/MD	Ok,M
9	Q1665-06	IA-B9-1-3242025	VL042293.D	09 Apr 2025 15:15		SY/MD	Ok,M
10	Q1665-08	IA-DUP-2-3242025	VL042294.D	09 Apr 2025 15:50		SY/MD	Ok,M
11	Q1665-09	IA-B9-2-3242025	VL042295.D	09 Apr 2025 16:23		SY/MD	Ok,M
12	Q1665-20	IA-B9-3-3242025	VL042296.D	09 Apr 2025 16:57		SY/MD	Ok,M
13	Q1757-02	IA1	VL042297.D	09 Apr 2025 17:31		SY/MD	Ok,M
14	Q1757-01	SV1	VL042298.D	10 Apr 2025 07:54		SY/MD	Ok,M

M : Manual Integration

Instrument ID: MSVOA_L

Daily Analysis Runlog For Sequence/QC Batch ID # VL041125

Review By	Semsettin Yesilyurt	Review On	4/16/2025 9:18:47 PM
Supervise By	Maresh Dadoda	Supervise On	4/17/2025 1:40:37 PM
SubDirectory	VL041125	HP Acquire Method	TO15AIR.M
		HP Processing Method	VL032725AIR.M
STD. NAME	STD REF.#		
Tune/Reschk	AP2591		
Initial Calibration Stds	AP2593,AP2595,AP2596		
CCC	AP2593		
Internal Standard/PEM	AP2591		
ICV/I.BLK	AP2594		
Surrogate Standard			
MS/MSD Standard			
LCS Standard			

Sr#	SampleID	ClientID	Data File Name	Date-Time	Comment	Operator	Status
1	BFB	BFB	VL042311.D	11 Apr 2025 08:18		SY/MD	Ok
2	VSTDCCC010	VSTDCCC010	VL042312.D	11 Apr 2025 09:50		SY/MD	Ok,M
3	VL0411ABL01	VL0411ABL01	VL042313.D	11 Apr 2025 10:36		SY/MD	Ok
4	VL0411ABS01	VL0411ABS01	VL042314.D	11 Apr 2025 11:23		SY/MD	Ok,M
5	QC041025CAN 10596	QC041025CAN 10596	VL042315.D	11 Apr 2025 12:08	Contaminated	SY/MD	Not Ok
6	QC041025CAN 10596	QC041025CAN 10596	VL042316.D	11 Apr 2025 12:40		SY/MD	Ok
7	Q1665-03	IA-B10-2-3242025	VL042317.D	11 Apr 2025 13:15	Need 10X	SY/MD	Dilution
8	Q1665-03DUP	IA-B10-2-3242025DUP	VL042318.D	11 Apr 2025 13:52		SY/MD	Ok,M
9	Q1665-10	IA-B10-4-3242025	VL042319.D	11 Apr 2025 14:27	Need 40X	SY/MD	Dilution
10	Q1665-11	IA-B10-1-3242025	VL042320.D	11 Apr 2025 15:01	Need 10X	SY/MD	Dilution
11	Q1665-18	IA-B10-3-3242025	VL042321.D	11 Apr 2025 15:36	Need 10X	SY/MD	Dilution
12	Q1665-03DL	IA-B10-2-3242025DL	VL042322.D	11 Apr 2025 16:08		SY/MD	Ok
13	Q1665-10DL	IA-B10-4-3242025DL	VL042323.D	11 Apr 2025 16:39	Need 100X	SY/MD	Dilution
14	Q1665-11DL	IA-B10-1-3242025DL	VL042324.D	11 Apr 2025 17:10		SY/MD	Ok,M

M : Manual Integration

Instrument ID: MSVOA_L

Daily Analysis Runlog For Sequence/QC Batch ID # VL041525

Review By	Semsettin Yesilyurt	Review On	4/16/2025 2:19:24 PM
Supervise By	Mahesh Dadoda	Supervise On	4/16/2025 3:47:46 PM
SubDirectory	VL041525	HP Acquire Method	TO15AIR.M
		HP Processing Method	VL032725AIR.M
STD. NAME	STD REF.#		
Tune/Reschk	AP2591		
Initial Calibration Stds	AP2593,AP2595,AP2596		
CCC	AP2593		
Internal Standard/PEM	AP2591		
ICV/I.BLK	AP2594		
Surrogate Standard			
MS/MSD Standard			
LCS Standard			

Sr#	SampleID	ClientID	Data File Name	Date-Time	Comment	Operator	Status
1	BFB	BFB	VL042325.D	15 Apr 2025 07:45		SY/MD	Ok
2	VSTDCCC010	VSTDCCC010	VL042326.D	15 Apr 2025 11:30		SY/MD	Ok,M
3	VL0415ABL01	VL0415ABL01	VL042327.D	15 Apr 2025 12:12		SY/MD	Ok
4	VL0415ABS01	VL0415ABS01	VL042328.D	15 Apr 2025 12:54		SY/MD	Ok,M
5	QC041425 CAN 10401	QC041425 CAN 10401	VL042329.D	15 Apr 2025 13:26	carryover	SY/MD	Not Ok
6	QC041425 CAN 10401	QC041425 CAN 10401	VL042330.D	15 Apr 2025 13:57		SY/MD	Ok
7	Q1665-18DL	IA-B10-3-3242025DL	VL042331.D	15 Apr 2025 14:29		SY/MD	Ok
8	Q1665-10DL2	IA-B10-4-3242025DL2	VL042332.D	15 Apr 2025 15:01		SY/MD	Ok,M
9	Q1802-02	IA1	VL042333.D	15 Apr 2025 15:35		SY/MD	Ok,M
10	Q1802-02DUP	IA1DUP	VL042334.D	15 Apr 2025 16:08		SY/MD	Ok,M
11	Q1802-03	OA1	VL042335.D	15 Apr 2025 16:42		SY/MD	Ok,M
12	Q1802-01	SV1	VL042336.D	15 Apr 2025 17:14	Need 40X	SY/MD	Dilution
13	Q1802-01DL	SV1DL	VL042337.D	15 Apr 2025 17:46	Need further dilution	SY/MD	Not Ok
14	Q1802-01DL	SV1DL	VL042338.D	15 Apr 2025 18:21		SY/MD	Ok,M

M : Manual Integration

Instrument ID: MSVOA_L

Daily Analysis Runlog For Sequence/QC Batch ID # VL041725

Review By	Semsettin Yesilyurt	Review On	4/18/2025 4:23:02 PM
Supervise By	Maresh Dadoda	Supervise On	4/18/2025 10:46:12 PM
SubDirectory	VL041725	HP Acquire Method	TO15AIR.M HP Processing Method VL041725AIR.M

STD. NAME	STD REF.#
Tune/Reschk	AP2590
Initial Calibration Stds	AP2598,AP2599,AP2600
CCC	AP2598
Internal Standard/PEM	AP2590
ICV/I.BLK	AP2601
Surrogate Standard	
MS/MSD Standard	
LCS Standard	

Sr#	Sampleld	ClientID	Data File Name	Date-Time	Comment	Operator	Status
1	BFB	BFB	VL042348.D	17 Apr 2025 07:50		SY/MD	Ok
2	VSTDICCC010	VSTDICCC010	VL042349.D	17 Apr 2025 09:44		SY/MD	Ok,M
3	VSTDICC002	VSTDICC002	VL042350.D	17 Apr 2025 10:17		SY/MD	Ok,M
4	VSTDICC001	VSTDICC001	VL042351.D	17 Apr 2025 10:48		SY/MD	Ok,M
5	VSTDICC0.5	VSTDICC0.5	VL042352.D	17 Apr 2025 11:19		SY/MD	Ok,M
6	VSTDICC0.1	VSTDICC0.1	VL042353.D	17 Apr 2025 11:50		SY/MD	Ok,M
7	VSTDICC0.03	VSTDICC0.03	VL042354.D	17 Apr 2025 12:22		SY/MD	Ok,M
8	VSTDICC015	VSTDICC015	VL042355.D	17 Apr 2025 12:54		SY/MD	Ok,M
9	VSTDICV010	ICVVL041725	VL042356.D	17 Apr 2025 13:43		SY/MD	Ok,M
10	VL0417ABL01	VL0417ABL01	VL042357.D	17 Apr 2025 14:23		SY/MD	Ok
11	VL0417ABS01	VL0417ABS01	VL042358.D	17 Apr 2025 15:39		SY/MD	Ok,M
12	Q1665-15	IA-B14-1-3242025	VL042359.D	17 Apr 2025 16:12	Need 10X	SY/MD	Dilution
13	Q1665-15DL	IA-B14-1-3242025DL	VL042360.D	17 Apr 2025 16:44		SY/MD	Ok,M
14	Q1668-01	IA-B1-4-3252025	VL042361.D	17 Apr 2025 17:17	Need 4X	SY/MD	Dilution
15	Q1668-01DUP	IA-B1-4-3252025DUP	VL042362.D	17 Apr 2025 17:51		SY/MD	Ok,M
16	Q1668-03	IA-B1-3-3252025	VL042363.D	17 Apr 2025 18:24	Need 4X	SY/MD	Dilution
17	Q1668-04	IA-B1-1-3252025	VL042364.D	17 Apr 2025 18:58		SY/MD	Ok,M
18	Q1668-02	IA-B1-2-3252025	VL042365.D	17 Apr 2025 19:31		SY/MD	Ok,M

Instrument ID: MSVOA_L

Daily Analysis Runlog For Sequence/QC Batch ID # VL041725

Review By	Semsettin Yesilyurt	Review On	4/18/2025 4:23:02 PM		
Supervise By	Mahesh Dadoda	Supervise On	4/18/2025 10:46:12 PM		
SubDirectory	VL041725	HP Acquire Method	TO15AIR.M	HP Processing Method	VL041725AIR.M
STD. NAME	STD REF.#				
Tune/Reschk	AP2590				
Initial Calibration Stds	AP2598,AP2599,AP2600				
CCC	AP2598				
Internal Standard/PEM	AP2590				
ICV/I.BLK	AP2601				
Surrogate Standard					
MS/MSD Standard					
LCS Standard					

19	Q1726-01	IA-DUP-3-03312025	VL042366.D	17 Apr 2025 20:05	Need 20X	SY/MD	Dilution
20	Q1726-02	IA-B2-6-03312025	VL042367.D	17 Apr 2025 20:38	Need 20X	SY/MD	Dilution
21	Q1726-03	IA-B2-2-03312025	VL042368.D	17 Apr 2025 21:12	Need 20X	SY/MD	Dilution
22	Q1726-04	IA-B2-5-03312025	VL042369.D	17 Apr 2025 21:45	Need 20X	SY/MD	Dilution
23	Q1726-05	IA-B2-1-03312025	VL042370.D	17 Apr 2025 22:19	Need 20X	SY/MD	Dilution
24	Q1726-06	IA-B2-3-03312025	VL042371.D	17 Apr 2025 22:52	Need 20X	SY/MD	Dilution
25	Q1726-07	IA-B2-4-03312025	VL042372.D	17 Apr 2025 23:26	Need 20X	SY/MD	Dilution
26	Q1726-07DL	IA-B2-4-03312025DL	VL042373.D	17 Apr 2025 23:57	Need 40X	SY/MD	Dilution
27	Q1668-01DL	IA-B1-4-3252025DL	VL042374.D	18 Apr 2025 00:28		SY/MD	Ok,M
28	Q1668-03DL	IA-B1-3-3252025DL	VL042375.D	18 Apr 2025 00:59		SY/MD	Ok,M
29	Q1726-01DL	IA-DUP-3-03312025DL	VL042376.D	18 Apr 2025 01:30	Need 40X	SY/MD	Dilution
30	Q1726-02DL	IA-B2-6-03312025DL	VL042377.D	18 Apr 2025 02:01	Need 40X	SY/MD	Dilution
31	Q1726-03DL	IA-B2-2-03312025DL	VL042378.D	18 Apr 2025 02:32	Need 100X	SY/MD	Dilution
32	Q1726-04DL	IA-B2-5-03312025DL	VL042379.D	18 Apr 2025 03:03	Need 40X	SY/MD	Dilution
33	Q1726-05DL	IA-B2-1-03312025DL	VL042380.D	18 Apr 2025 03:34	Need 100X	SY/MD	Dilution
34	Q1726-06DL	IA-B2-3-03312025DL	VL042381.D	18 Apr 2025 04:05	Need 100X	SY/MD	Dilution

M : Manual Integration

Instrument ID: MSVOA_L

Daily Analysis Runlog For Sequence/QC Batch ID # VL042825

Review By	Semsettin Yesilyurt	Review On	4/30/2025 2:24:11 PM
Supervise By	Maresh Dadoda	Supervise On	4/30/2025 2:30:41 PM
SubDirectory	VL042825	HP Acquire Method	TO15AIR.M
		HP Processing Method	VL041725AIR.M

STD. NAME	STD REF.#
Tune/Reschk	AP2602
Initial Calibration Stds	AP2598,AP2599,AP2600
CCC	AP2598
Internal Standard/PEM	AP2602
ICV/I.BLK	AP2601
Surrogate Standard	
MS/MSD Standard	
LCS Standard	

Sr#	SampleID	ClientID	Data File Name	Date-Time	Comment	Operator	Status
1	BFB	BFB	VL042415.D	28 Apr 2025 07:59		SY/MD	Ok
2	VSTDCCC010	VSTDCCC010	VL042416.D	28 Apr 2025 08:44		SY/MD	Ok,M
3	VL0428ABL01	VL0428ABL01	VL042417.D	28 Apr 2025 09:27		SY/MD	Ok
4	VL0428ABS01	VL0428ABS01	VL042418.D	28 Apr 2025 10:10		SY/MD	Ok,M
5	Q1665-05	IA-B2-3-3242025	VL042419.D	28 Apr 2025 10:59		SY/MD	Ok,M
6	Q1665-07	IA-DUP-1-3242025	VL042420.D	28 Apr 2025 11:33	Need 4x	SY/MD	Dilution
7	Q1665-12	IA-B2-6-3242025	VL042421.D	28 Apr 2025 12:17		SY/MD	Ok,M
8	Q1665-16	IA-B2-5-3242025	VL042422.D	28 Apr 2025 12:51		SY/MD	Ok,M
9	Q1665-19	IA-B2-2-3242025	VL042423.D	28 Apr 2025 13:26	Need 4x	SY/MD	Dilution
10	Q1665-21	IA-B2-1-3242025	VL042424.D	28 Apr 2025 14:00	Need 4X	SY/MD	Dilution
11	Q1665-22	IA-B2-4-3242025	VL042425.D	28 Apr 2025 14:34		SY/MD	Ok,M
12	Q1880-02	EBR-P02-SL-IA02-2025	VL042426.D	28 Apr 2025 15:08		SY/MD	Ok,M
13	Q1880-03	EBR-P08-SC-IA01-2025	VL042427.D	28 Apr 2025 15:42		SY/MD	Ok,M
14	Q1880-08	EBR-P02-SC-IA01-2025	VL042428.D	28 Apr 2025 16:16		SY/MD	Ok,M
15	Q1880-10	EBR-P16-SC-AA01-2025	VL042429.D	28 Apr 2025 16:50		SY/MD	Ok,M
16	Q1880-10DUP	EBR-P16-SC-AA01-2025	VL042430.D	28 Apr 2025 17:26		SY/MD	Ok,M
17	Q1880-01	EBR-P11-SC-SS01-2025	VL042431.D	28 Apr 2025 18:00	Need 20X	SY/MD	Dilution
18	Q1880-06	EBR-P02-SC-DUP01-2025	VL042432.D	28 Apr 2025 18:35		SY/MD	Ok,M

Instrument ID: MSVOA_L

Daily Analysis Runlog For Sequence/QC Batch ID # VL042825

Review By	Semsettin Yesilyurt	Review On	4/30/2025 2:24:11 PM
Supervise By	Maresh Dadoda	Supervise On	4/30/2025 2:30:41 PM
SubDirectory	VL042825	HP Acquire Method	TO15AIR.M
		HP Processing Method	VL041725AIR.M

STD. NAME	STD REF.#
Tune/Reschk	AP2602
Initial Calibration Stds	AP2598,AP2599,AP2600
CCC	AP2598
Internal Standard/PEM	AP2602
ICV/I.BLK	AP2601
Surrogate Standard	
MS/MSD Standard	
LCS Standard	

19	Q1880-09	EBR-P16-SC-AA-DUP-	VL042433.D	28 Apr 2025 19:08		SY/MD	Ok,M
20	Q1880-04	EBR-P08-SC-SS01-202	VL042434.D	28 Apr 2025 19:41	rerun straight	SY/MD	Not Ok

M : Manual Integration

Instrument ID: MSVOA_L

Daily Analysis Runlog For Sequence/QC Batch ID # VL042925

Review By	Semsettin Yesilyurt	Review On	4/30/2025 2:23:53 PM
Supervise By	Maresh Dadoda	Supervise On	4/30/2025 2:28:33 PM
SubDirectory	VL042925	HP Acquire Method	TO15AIR.M
		HP Processing Method	VL041725AIR.M

STD. NAME	STD REF.#
Tune/Reschk	AP2602
Initial Calibration Stds	AP2598,AP2599,AP2600
CCC	AP2598
Internal Standard/PEM	AP2602
ICV/I.BLK	AP2601
Surrogate Standard	
MS/MSD Standard	
LCS Standard	

Sr#	SampleID	ClientID	Data File Name	Date-Time	Comment	Operator	Status
1	BFB	BFB	VL042435.D	29 Apr 2025 07:55		SY/MD	Ok
2	VSTDCCC010	VSTDCCC010	VL042436.D	29 Apr 2025 08:39		SY/MD	Ok,M
3	VL0429ABL01	VL0429ABL01	VL042437.D	29 Apr 2025 09:25		SY/MD	Ok
4	VL0429ABS01	VL0429ABS01	VL042438.D	29 Apr 2025 10:08		SY/MD	Ok,M
5	Q1880-04	EBR-P08-SC-SS01-202	VL042439.D	29 Apr 2025 10:59		SY/MD	Ok,M
6	Q1880-05DL	EBR-SC-SS-DUP01-202	VL042440.D	29 Apr 2025 11:30	Need Straight Run	SY/MD	Not Ok
7	Q1880-07DL	EBR-P16-SC-SS01-202	VL042441.D	29 Apr 2025 12:02		SY/MD	Ok,M
8	Q1880-11DL	EBR-P02-SC-SS01-202	VL042442.D	29 Apr 2025 12:33	Need Straight Run	SY/MD	Not Ok
9	Q1908-01	B-1	VL042443.D	29 Apr 2025 13:11		SY/MD	Ok,M
10	Q1908-02	B-2	VL042444.D	29 Apr 2025 13:45	rerun confirmation,4X	SY/MD	Not Ok
11	Q1908-03	B-3	VL042445.D	29 Apr 2025 14:19		SY/MD	Ok,M
12	Q1908-03DL	B-3DL	VL042446.D	29 Apr 2025 14:51	Not Required	SY/MD	Not Ok
13	Q1880-05	EBR-SC-SS-DUP01-202	VL042447.D	29 Apr 2025 15:24		SY/MD	Ok,M
14	Q1880-07	EBR-P16-SC-SS01-202	VL042448.D	29 Apr 2025 15:58	Need 10X	SY/MD	Dilution
15	Q1880-11	EBR-P02-SC-SS01-202	VL042449.D	29 Apr 2025 16:32		SY/MD	Ok,M
16	Q1880-01DL	EBR-P11-SC-SS01-202	VL042450.D	29 Apr 2025 17:04		SY/MD	Ok
17	Q1665-07DL	IA-DUP-1-3242025DL	VL042451.D	29 Apr 2025 17:37		SY/MD	Ok,M
18	Q1665-19DL	IA-B2-2-3242025DL	VL042452.D	29 Apr 2025 18:09		SY/MD	Ok,M

Instrument ID: MSVOA_L

Daily Analysis Runlog For Sequence/QC Batch ID # VL042925

Review By	Semsettin Yesilyurt	Review On	4/30/2025 2:23:53 PM		
Supervise By	Mahesh Dadoda	Supervise On	4/30/2025 2:28:33 PM		
SubDirectory	VL042925	HP Acquire Method	TO15AIR.M	HP Processing Method	VL041725AIR.M
STD. NAME	STD REF.#				
Tune/Reschk	AP2602				
Initial Calibration Stds	AP2598,AP2599,AP2600				
CCC	AP2598				
Internal Standard/PEM	AP2602				
ICV/I.BLK	AP2601				
Surrogate Standard					
MS/MSD Standard					
LCS Standard					

19	Q1665-21DL	IA-B2-1-3242025DL	VL042453.D	29 Apr 2025 18:42		SY/MD	Ok,M
20	Q1908-01DUP	B-1DUP	VL042454.D	29 Apr 2025 19:57		SY/MD	Ok,M
21	Q1908-02DL	B-2DL	VL042455.D	29 Apr 2025 20:30	not required ,not match	SY/MD	Not Ok
22	Q1908-02	B-2	VL042456.D	29 Apr 2025 21:58		SY/MD	Ok,M

M : Manual Integration