

Instrument ID: FID_E

Daily Analysis Runlog For Sequence/QC Batch ID # FE032925AL

Review By	yogesh	Review On	3/28/2025 1:29:33 PM
Supervise By	mohammad	Supervise On	4/1/2025 1:30:29 AM
SubDirectory	FE032925AL	HP Acquire Method	HP Processing Method
			FE032925AL

STD. NAME	STD REF.#
Tune/Reschk Initial Calibration Stds	PP24170,PP24175,PP24176,PP24177,PP24178
CCC Internal Standard/PEM	PP24176
ICV/I.BLK Surrogate Standard MS/MSD Standard LCS Standard	PP24174,PP24179

Sr#	SampleID	ClientID	Data File Name	Date-Time	Comment	Operator	Status
1	MECL2	MECL2	FE053055.D	28 Mar 2025 22:09		YP\AJ	Ok
2	I.BLK	I.BLK	FE053056.D	28 Mar 2025 22:39		YP\AJ	Ok
3	100 PPM ALIPHATIC HC	100 PPM ALIPHATIC HC	FE053057.D	28 Mar 2025 23:09		YP\AJ	Ok
4	50 PPM ALIPHATIC HC	50 PPM ALIPHATIC HC	FE053058.D	29 Mar 2025 00:09		YP\AJ	Ok
5	20 PPM ALIPHATIC HC	20 PPM ALIPHATIC HC	FE053059.D	29 Mar 2025 00:39		YP\AJ	Ok
6	10 PPM ALIPHATIC HC	10 PPM ALIPHATIC HC	FE053060.D	29 Mar 2025 01:09		YP\AJ	Ok
7	5 PPM ALIPHATIC HC	5 PPM ALIPHATIC HC	FE053061.D	29 Mar 2025 01:39		YP\AJ	Ok
8	20 PPM ALIPHATIC HC	20 PPM ALIPHATIC HC	FE053062.D	29 Mar 2025 02:10		YP\AJ	Ok
9	I.BLK	I.BLK	FE053063.D	29 Mar 2025 03:10		YP\AJ	Ok
10	20 PPM ALIPHATIC HC	20 PPM ALIPHATIC HC	FE053064.D	29 Mar 2025 03:40		YP\AJ	Ok

M : Manual Integration

Instrument ID: FID_E

Daily Analysis Runlog For Sequence/QC Batch ID # FE033125AL

Review By	yogesh	Review On	3/31/2025 12:29:00 PM
Supervise By	mohammad	Supervise On	4/2/2025 1:16:46 AM
SubDirectory	FE033125AL	HP Acquire Method	HP Processing Method FE032925AL

STD. NAME	STD REF.#
Tune/Reschk Initial Calibration Stds	PP24170,PP24175,PP24176,PP24177,PP24178
CCC Internal Standard/PEM	PP24176
ICV/I.BLK Surrogate Standard MS/MSD Standard LCS Standard	PP24174,PP24179

Sr#	SampleID	ClientID	Data File Name	Date-Time	Comment	Operator	Status
1	MECL2	MECL2	FE053065.D	31 Mar 2025 08:28		YPIAJ	Ok
2	I.BLK	I.BLK	FE053066.D	31 Mar 2025 08:58		YPIAJ	Ok
3	20 PPM ALIPHATIC HC	20 PPM ALIPHATIC HC	FE053067.D	31 Mar 2025 12:00		YPIAJ	Ok
4	Q1657-01DL	72-11998DL	FE053068.D	31 Mar 2025 12:30		YPIAJ	Ok
5	I.BLK	I.BLK	FE053069.D	31 Mar 2025 13:30		YPIAJ	Ok
6	20 PPM ALIPHATIC HC	20 PPM ALIPHATIC HC	FE053070.D	31 Mar 2025 14:00		YPIAJ	Ok
7	PB167382BL	PB167382BL	FE053071.D	31 Mar 2025 15:01		YPIAJ	Ok
8	PB167382BS	PB167382BS	FE053072.D	31 Mar 2025 15:31		YPIAJ	Ok,M
9	PB167382BSD	PB167382BSD	FE053073.D	31 Mar 2025 16:01		YPIAJ	Ok,M
10	Q1675-01	P1	FE053074.D	31 Mar 2025 16:32		YPIAJ	Ok
11	Q1675-02	P2	FE053075.D	31 Mar 2025 17:02		YPIAJ	Ok
12	Q1675-02MS	P2MS	FE053076.D	31 Mar 2025 17:33	FE053075.D	YPIAJ	Ok
13	Q1675-02MSD	P2MSD	FE053077.D	31 Mar 2025 18:03	FE053075.D!FE053076.D	YPIAJ	Ok
14	Q1675-03	P3	FE053078.D	31 Mar 2025 18:33		YPIAJ	Ok
15	Q1675-05	P5	FE053079.D	31 Mar 2025 19:03	need 20x dilution	YPIAJ	Dilution
16	Q1675-06	P6	FE053080.D	31 Mar 2025 19:34		YPIAJ	Ok
17	Q1675-18	K1P	FE053081.D	31 Mar 2025 20:04	need 2x dilution	YPIAJ	Dilution
18	Q1675-19	K1	FE053082.D	31 Mar 2025 20:34		YPIAJ	Ok

Instrument ID: FID_E

Daily Analysis Runlog For Sequence/QC Batch ID # FE033125AL

Review By	yogesh	Review On	3/31/2025 12:29:00 PM
Supervise By	mohammad	Supervise On	4/2/2025 1:16:46 AM
SubDirectory	FE033125AL	HP Acquire Method	HP Processing Method FE032925AL
STD. NAME	STD REF.#		
Tune/Reschk Initial Calibration Stds	PP24170,PP24175,PP24176,PP24177,PP24178		
CCC	PP24176		
Internal Standard/PEM	PP24174,PP24179		
ICV/I.BLK			
Surrogate Standard			
MS/MSD Standard			
LCS Standard			

19	Q1675-20	K2	FE053083.D	31 Mar 2025 21:04		YP\AJ	Ok
20	Q1675-21	K3	FE053084.D	31 Mar 2025 21:34		YP\AJ	Ok
21	I.BLK	I.BLK	FE053085.D	31 Mar 2025 22:35		YP\AJ	Ok
22	20 PPM ALIPHATIC HC	20 PPM ALIPHATIC HC	FE053086.D	31 Mar 2025 23:05		YP\AJ	Ok

M : Manual Integration

Instrument ID: FID_E

Daily Analysis Runlog For Sequence/QC Batch ID # FE040125AL

Review By	yogesh	Review On	4/1/2025 10:49:14 AM
Supervise By	mohammad	Supervise On	4/3/2025 1:37:54 AM
SubDirectory	FE040125AL	HP Acquire Method	HP Processing Method FE032925AL

STD. NAME	STD REF.#
Tune/Reschk Initial Calibration Stds	PP24170,PP24175,PP24176,PP24177,PP24178
CCC Internal Standard/PEM ICV/I.BLK Surrogate Standard MS/MSD Standard LCS Standard	PP24176 PP24174,PP24179

Sr#	SampleID	ClientID	Data File Name	Date-Time	Comment	Operator	Status
1	MECL2	MECL2	FE053087.D	01 Apr 2025 06:15		YP\AJ	Ok
2	I.BLK	I.BLK	FE053088.D	01 Apr 2025 06:45		YP\AJ	Ok
3	20 PPM ALIPHATIC HC	20 PPM ALIPHATIC HC	FE053089.D	01 Apr 2025 07:15		YP\AJ	Ok
4	Q1675-05DL	P5DL	FE053090.D	01 Apr 2025 08:15		YP\AJ	Ok,M
5	Q1675-18DL	K1PDL	FE053091.D	01 Apr 2025 08:45		YP\AJ	Ok
6	I.BLK	I.BLK	FE053092.D	01 Apr 2025 10:33		YP\AJ	Ok
7	20 PPM ALIPHATIC HC	20 PPM ALIPHATIC HC	FE053093.D	01 Apr 2025 11:03		YP\AJ	Ok
8	PB167403BL	PB167403BL	FE053094.D	01 Apr 2025 15:26		YP\AJ	Ok
9	PB167403BS	PB167403BS	FE053095.D	01 Apr 2025 15:56		YP\AJ	Ok
10	PB167403BSD	PB167403BSD	FE053096.D	01 Apr 2025 16:26		YP\AJ	Ok
11	Q1679-01	TP-10	FE053097.D	01 Apr 2025 16:56		YP\AJ	Ok
12	Q1679-02	TP-10-EPH	FE053098.D	01 Apr 2025 17:26		YP\AJ	Ok
13	Q1679-05	TP-9	FE053099.D	01 Apr 2025 17:57		YP\AJ	Ok
14	Q1679-06	TP-9-EPH	FE053100.D	01 Apr 2025 18:27		YP\AJ	Ok
15	Q1680-01	TP-4	FE053101.D	01 Apr 2025 18:57		YP\AJ	Ok
16	Q1680-02	TP-4-EPH	FE053102.D	01 Apr 2025 19:27		YP\AJ	Ok
17	Q1681-01	TP-12	FE053103.D	01 Apr 2025 19:57		YP\AJ	Ok
18	Q1681-01D	Q1681-01D	FE053104.D	01 Apr 2025 20:27		YP\AJ	Ok

Instrument ID: FID_E

Daily Analysis Runlog For Sequence/QC Batch ID # FE040125AL

Review By	yogesh	Review On	4/1/2025 10:49:14 AM
Supervise By	mohammad	Supervise On	4/3/2025 1:37:54 AM
SubDirectory	FE040125AL	HP Acquire Method	HP Processing Method
			FE032925AL
STD. NAME	STD REF.#		
Tune/Reschk Initial Calibration Stds	PP24170,PP24175,PP24176,PP24177,PP24178		
CCC	PP24176		
Internal Standard/PEM			
ICV/I.BLK	PP24174,PP24179		
Surrogate Standard			
MS/MSD Standard			
LCS Standard			

19	Q1681-01MS	TP-12MS	FE053105.D	01 Apr 2025 20:57	FE053103.D	YPIAJ	Ok
20	Q1681-01MSD	TP-12MSD	FE053106.D	01 Apr 2025 21:28	FE053103.D!FE053105.D	YPIAJ	Ok
21	Q1681-02	TP-12-EPH	FE053107.D	01 Apr 2025 21:58		YPIAJ	Ok
22	I.BLK	I.BLK	FE053108.D	01 Apr 2025 22:58		YPIAJ	Ok
23	20 PPM ALIPHATIC HC	20 PPM ALIPHATIC HC	FE053109.D	01 Apr 2025 23:28		YPIAJ	Ok,M
24	Q1687-01	72-12016	FE053110.D	02 Apr 2025 00:28		YPIAJ	Ok
25	Q1687-02	72-12016-E2	FE053111.D	02 Apr 2025 00:58		YPIAJ	Ok
26	Q1688-01	OR-02-040125	FE053112.D	02 Apr 2025 01:28		YPIAJ	Ok
27	Q1688-02	OR-02-040125-E2	FE053113.D	02 Apr 2025 01:58	Need 2X	YPIAJ	Dilution
28	Q1692-01	NB-08-04012025	FE053114.D	02 Apr 2025 02:29	Need 5X	YPIAJ	Dilution
29	Q1692-02	NB-08-04012025-E2	FE053115.D	02 Apr 2025 02:58		YPIAJ	Ok
30	I.BLK	I.BLK	FE053116.D	02 Apr 2025 03:58		YPIAJ	Ok
31	20 PPM ALIPHATIC HC	20 PPM ALIPHATIC HC	FE053117.D	02 Apr 2025 04:28		YPIAJ	Ok
32	PB167416BL	PB167416BL	FE053118.D	02 Apr 2025 05:29		YPIAJ	Ok
33	PB167416BS	PB167416BS	FE053119.D	02 Apr 2025 05:59		YPIAJ	Ok
34	PB167416BSD	PB167416BSD	FE053120.D	02 Apr 2025 06:29		YPIAJ	Ok
35	Q1687-01D	Q1687-01D	FE053121.D	02 Apr 2025 06:59		YPIAJ	Ok
36	Q1687-01MS	72-12016MS	FE053122.D	02 Apr 2025 07:29	FE053110.D	YPIAJ	Ok,M
37	Q1687-01MSD	72-12016MSD	FE053123.D	02 Apr 2025 07:59	FE053110.D!FE053122.D	YPIAJ	Ok

Instrument ID: FID_E

Daily Analysis Runlog For Sequence/QC Batch ID # FE040125AL

Review By	yogesh	Review On	4/1/2025 10:49:14 AM
Supervise By	mohammad	Supervise On	4/3/2025 1:37:54 AM
SubDirectory	FE040125AL	HP Acquire Method	HP Processing Method FE032925AL

STD. NAME	STD REF.#
Tune/Reschk Initial Calibration Stds	PP24170,PP24175,PP24176,PP24177,PP24178
CCC Internal Standard/PEM	PP24176
ICV/I.BLK Surrogate Standard MS/MSD Standard LCS Standard	PP24174,PP24179

38	I.BLK	I.BLK	FE053124.D	02 Apr 2025 08:29		YP\AJ	Ok
39	20 PPM ALIPHATIC HC	20 PPM ALIPHATIC HC	FE053125.D	02 Apr 2025 08:59		YP\AJ	Ok

M : Manual Integration