

Instrument ID: MSVOA_L

Daily Analysis Runlog For Sequence/QC Batch ID # VL032725

Review By	Maresh Dadoda	Review On	4/3/2025 8:14:07 AM
Supervise By	Semsettin Yesilyurt	Supervise On	4/3/2025 8:45:16 AM
SubDirectory	VL032725	HP Acquire Method	TO15AIR.M HP Processing Method VL032725AIR.M

STD. NAME	STD REF.#
Tune/Reschk	AP2589
Initial Calibration Stds	AP2593,AP2595,AP2596
CCC	AP2593
Internal Standard/PEM	AP2589
ICV/I.BLK	AP2594
Surrogate Standard	
MS/MSD Standard	
LCS Standard	

Sr#	SampleID	ClientID	Data File Name	Date-Time	Comment	Operator	Status
1	BFB	BFB	VL042171.D	27 Mar 2025 07:44		SY/MD	Ok
2	VSTDICCC010	VSTDICCC010	VL042172.D	27 Mar 2025 08:26		SY/MD	Ok,M
3	VSTDICC002	VSTDICC002	VL042173.D	27 Mar 2025 08:59		SY/MD	Ok,M
4	VSTDICC001	VSTDICC001	VL042174.D	27 Mar 2025 09:30		SY/MD	Ok,M
5	VSTDICC0.5	VSTDICC0.5	VL042175.D	27 Mar 2025 10:01		SY/MD	Ok,M
6	VSTDICC0.1	VSTDICC0.1	VL042176.D	27 Mar 2025 10:32		SY/MD	Ok,M
7	VSTDICC0.03	VSTDICC0.03	VL042177.D	27 Mar 2025 11:03		SY/MD	Ok,M
8	VSTDICC015	VSTDICC015	VL042178.D	27 Mar 2025 11:36		SY/MD	Ok,M
9	VSTDICV010	ICVVL032725	VL042179.D	27 Mar 2025 12:23		SY/MD	Ok,M
10	VL0327ABL01	VL0327ABL01	VL042180.D	27 Mar 2025 13:50		SY/MD	Ok
11	VL0327ABS01	VL0327ABS01	VL042181.D	27 Mar 2025 14:39		SY/MD	Ok,M
12	QC032625 CAN 10448	QC032625 CAN 10448	VL042182.D	27 Mar 2025 16:41		SY/MD	Ok
13	Q1649-01	SV1	VL042183.D	27 Mar 2025 17:13	rerun 4X,not required	SY/MD	Not Ok
14	Q1649-02	SV2	VL042184.D	27 Mar 2025 17:45	rerun 4X,not required	SY/MD	Not Ok
15	Q1669-01	IA-B9-3-3262025	VL042185.D	27 Mar 2025 18:22	Need 5x	SY/MD	Dilution
16	Q1669-02	IA-B9-2-3262025	VL042186.D	27 Mar 2025 18:59		SY/MD	Ok,M
17	Q1669-03	IA-B9-1-3262025	VL042187.D	27 Mar 2025 19:36		SY/MD	Ok,M
18	Q1669-04	AA-2-3252025	VL042188.D	27 Mar 2025 20:12		SY/MD	Ok,M

Instrument ID: MSVOA_L

Daily Analysis Runlog For Sequence/QC Batch ID # VL032725

Review By	Maresh Dadoda	Review On	4/3/2025 8:14:07 AM
Supervise By	Semsettin Yesilyurt	Supervise On	4/3/2025 8:45:16 AM
SubDirectory	VL032725	HP Acquire Method	TO15AIR.M HP Processing Method VL032725AIR.M
STD. NAME	STD REF.#		
Tune/Reschk	AP2589		
Initial Calibration Stds	AP2593,AP2595,AP2596		
CCC	AP2593		
Internal Standard/PEM	AP2589		
ICV/I.BLK	AP2594		
Surrogate Standard			
MS/MSD Standard			
LCS Standard			

19	Q1669-09	AA-3-3252025	VL042189.D	27 Mar 2025 20:48		SY/MD	Ok,M
20	Q1669-09DUP	AA-3-3252025DUP	VL042190.D	27 Mar 2025 21:24		SY/MD	Ok,M
21	Q1669-05	SSSG-B16-3-3262025	VL042191.D	27 Mar 2025 21:57	Need straight run	SY/MD	Not Ok
22	Q1669-07DL	SSSG-B16-1-3262025	VL042192.D	27 Mar 2025 22:29		SY/MD	Ok,M
23	Q1669-08	SSSG-B16-2-3262025	VL042193.D	27 Mar 2025 23:01	Need straight run	SY/MD	Not Ok
24	VIBLK	VIBLK	VL042194.D	27 Mar 2025 23:34		SY/MD	Ok
25	Q1669-01DL	IA-B9-3-3262025DL	VL042195.D	28 Mar 2025 00:06		SY/MD	Ok,M
26	Q1669-01DUP	IA-B9-3-3262025DUP	VL042196.D	28 Mar 2025 00:44		SY/MD	Ok,M
27	Q1669-05	SSSG-B16-3-3262025	VL042197.D	28 Mar 2025 01:20		SY/MD	Ok,M
28	Q1669-07	SSSG-B16-1-3262025	VL042198.D	28 Mar 2025 01:57	Need 10X	SY/MD	Dilution
29	Q1669-08	SSSG-B16-2-3262025	VL042199.D	28 Mar 2025 02:33		SY/MD	Ok,M

M : Manual Integration

Instrument ID: MSVOA_L

Daily Analysis Runlog For Sequence/QC Batch ID # VL040725

Review By	Maresh Dadoda	Review On	4/9/2025 11:41:05 AM
Supervise By	Semsettin Yesilyurt	Supervise On	4/9/2025 11:43:39 AM
SubDirectory	VL040725	HP Acquire Method	TO15AIR.M
		HP Processing Method	VL032725AIR.M

STD. NAME	STD REF.#
Tune/Reschk	AP2591
Initial Calibration Stds	AP2593,AP2595,AP2596
CCC	AP2593
Internal Standard/PEM	AP2591
ICV/I.BLK	AP2594
Surrogate Standard	
MS/MSD Standard	
LCS Standard	

Sr#	SampleID	ClientID	Data File Name	Date-Time	Comment	Operator	Status
1	BFB	BFB	VL042227.D	07 Apr 2025 08:46		SY/MD	Ok
2	VSTDCCC010	VSTDCCC010	VL042228.D	07 Apr 2025 09:31		SY/MD	Ok,M
3	VL0407ABL01	VL0407ABL01	VL042229.D	07 Apr 2025 10:23		SY/MD	Ok
4	VL0407ABS01	VL0407ABS01	VL042230.D	07 Apr 2025 11:09		SY/MD	Ok,M
5	Q1725-04	AA-1-03312025	VL042231.D	07 Apr 2025 12:27		SY/MD	Ok,M
6	Q1725-04DUP	AA-1-03312025DUP	VL042232.D	07 Apr 2025 13:00		SY/MD	Ok,M
7	Q1725-01	SSSG-B1-2-04012025	VL042233.D	07 Apr 2025 13:31	Need Straight Run	SY/MD	Not Ok
8	Q1725-01	SSSG-B1-2-04012025	VL042234.D	07 Apr 2025 14:18		SY/MD	Ok,M
9	Q1725-02	SSSG-B10-4-04012025	VL042235.D	07 Apr 2025 14:51	Need 10X	SY/MD	Dilution
10	Q1725-02DL	SSSG-B10-4-04012025	VL042236.D	07 Apr 2025 17:44	Need 40X	SY/MD	Dilution
11	Q1725-02DL2	SSSG-B10-4-04012025	VL042237.D	07 Apr 2025 18:15	Need further dilution	SY/MD	Not Ok
12	Q1725-03DL	SSSG-DUP-2-04012025	VL042238.D	07 Apr 2025 18:46	Not Required	SY/MD	Not Ok
13	Q1725-05DL	SSSG-B14-3-04012025	VL042239.D	07 Apr 2025 19:17		SY/MD	Ok,M
14	Q1725-06DL	SSSG-B10-3-03312025	VL042240.D	07 Apr 2025 19:48	Need further dilution	SY/MD	Not Ok
15	Q1725-07DL	SSSG-B10-1-03312025	VL042241.D	07 Apr 2025 20:19		SY/MD	Ok,M
16	Q1725-08DL	SSSG-B10-2-03312025	VL042242.D	07 Apr 2025 20:50	Need further dilution	SY/MD	Not Ok
17	Q1725-02DL2	SSSG-B10-4-04012025	VL042243.D	07 Apr 2025 21:38		SY/MD	Ok,M
18	Q1725-03DL	SSSG-DUP-2-04012025	VL042244.D	07 Apr 2025 22:09		SY/MD	Ok,M

Instrument ID: MSVOA_L

Daily Analysis Runlog For Sequence/QC Batch ID # VL040725

Review By	Mahesh Dadoda	Review On	4/9/2025 11:41:05 AM		
Supervise By	Semsettin Yesilyurt	Supervise On	4/9/2025 11:43:39 AM		
SubDirectory	VL040725	HP Acquire Method	TO15AIR.M	HP Processing Method	VL032725AIR.M
STD. NAME	STD REF.#				
Tune/Reschk	AP2591				
Initial Calibration Stds	AP2593,AP2595,AP2596				
CCC	AP2593				
Internal Standard/PEM	AP2591				
ICV/I.BLK	AP2594				
Surrogate Standard					
MS/MSD Standard					
LCS Standard					

19	Q1725-05	SSSG-B14-3-04012025	VL042245.D	07 Apr 2025 22:40	Need 10X	SY/MD	Dilution
20	Q1725-05DL	SSSG-B14-3-04012025	VL042246.D	07 Apr 2025 23:11	Not Required	SY/MD	Not Ok
21	Q1725-06	SSSG-B10-3-03312025	VL042247.D	07 Apr 2025 23:41	Need 40X	SY/MD	Dilution
22	Q1725-06DL	SSSG-B10-3-03312025	VL042248.D	08 Apr 2025 00:12		SY/MD	Dilution
23	Q1725-07	SSSG-B10-1-03312025	VL042249.D	08 Apr 2025 00:46	Need 10X	SY/MD	Dilution
24	Q1725-08	SSSG-B10-2-03312025	VL042250.D	08 Apr 2025 01:17	Need 40X	SY/MD	Dilution
25	Q1725-08DL	SSSG-B10-2-03312025	VL042251.D	08 Apr 2025 01:48		SY/MD	Ok,M

M : Manual Integration

Instrument ID: MSVOA_L

Daily Analysis Runlog For Sequence/QC Batch ID # VL040825

Review By	Maresh Dadoda	Review On	4/9/2025 11:42:35 AM
Supervise By	Amit Patel	Supervise On	4/10/2025 11:22:52 AM
SubDirectory	VL040825	HP Acquire Method	TO15AIR.M HP Processing Method VL032725AIR.M

STD. NAME	STD REF.#
Tune/Reschk	AP2591
Initial Calibration Stds	AP2593,AP2595,AP2596
CCC	AP2593
Internal Standard/PEM	AP2591
ICV/I.BLK	AP2594
Surrogate Standard	
MS/MSD Standard	
LCS Standard	

Sr#	SampleID	ClientID	Data File Name	Date-Time	Comment	Operator	Status
1	BFB	BFB	VL042252.D	08 Apr 2025 07:57		SY/MD	Ok
2	VSTDCCC010	VSTDCCC010	VL042253.D	08 Apr 2025 08:39		SY/MD	Ok,M
3	VL0408ABL01	VL0408ABL01	VL042254.D	08 Apr 2025 09:22		SY/MD	Ok
4	VL0408ABS01	VL0408ABS01	VL042255.D	08 Apr 2025 10:05		SY/MD	Ok,M
5	Q1725-03	SSSG-DUP-2-0401202	VL042256.D	08 Apr 2025 10:55	need 4X	SY/MD	Dilution
6	Q1728-03	SSSG-B2-2-04022025	VL042257.D	08 Apr 2025 11:26	Surrogate Fail; Need 40X; Bad chromatograph and RT shifted	SY/MD	Not Ok
7	VIBLK	VIBLK	VL042258.D	08 Apr 2025 12:05	clean up	SY/MD	Ok,M
8	VIBLK	VIBLK	VL042259.D	08 Apr 2025 12:36	clean up	SY/MD	Ok
9	Q1728-06	SSSG-B2-5-04022025	VL042260.D	08 Apr 2025 13:07	Need 40X	SY/MD	Dilution
10	Q1728-06DL	SSSG-B2-5-04022025	VL042261.D	08 Apr 2025 13:52		SY/MD	Ok,M
11	Q1728-07	SSSG-B1-1-04012025	VL042262.D	08 Apr 2025 14:22	Run @ Lower Dilution	SY/MD	Not Ok
12	Q1728-08	SSSG-B2-4-04022025	VL042263.D	08 Apr 2025 14:53	Need 40X	SY/MD	Dilution
13	Q1728-11	SSSG-B1-4-04012025	VL042264.D	08 Apr 2025 15:25	Run @ Lower Dilution	SY/MD	Not Ok
14	Q1728-12	SSSG-B1-3-04012025	VL042265.D	08 Apr 2025 15:56	Run @ Lower Dilution	SY/MD	Not Ok
15	Q1728-09	SSSG-B2-1-04022025	VL042266.D	08 Apr 2025 16:27		SY/MD	Ok,M
16	Q1728-01DL	SSSG-B14-2-04012025	VL042267.D	08 Apr 2025 16:58		SY/MD	Ok,M
17	Q1728-02	SSSG-B1-5-04012025	VL042268.D	08 Apr 2025 17:29	Run @ Lower Dilution	SY/MD	Not Ok

Instrument ID: MSVOA_L

Daily Analysis Runlog For Sequence/QC Batch ID # VL040825

Review By	Maresh Dadoda	Review On	4/9/2025 11:42:35 AM
Supervise By	Amit Patel	Supervise On	4/10/2025 11:22:52 AM
SubDirectory	VL040825	HP Acquire Method	TO15AIR.M
		HP Processing Method	VL032725AIR.M
STD. NAME	STD REF.#		
Tune/Reschk	AP2591		
Initial Calibration Stds	AP2593,AP2595,AP2596		
CCC	AP2593		
Internal Standard/PEM	AP2591		
ICV/I.BLK	AP2594		
Surrogate Standard			
MS/MSD Standard			
LCS Standard			

18	Q1728-04	SSSG-B2-7-04022025	VL042269.D	08 Apr 2025 18:00	Run @ Lower Dilution	SY/MD	Not Ok
19	Q1728-05	SSSG-B2-6-04022025	VL042270.D	08 Apr 2025 18:32	Run @ Lower Dilution	SY/MD	Not Ok
20	Q1728-10	SSSG-B14-4-04022025	VL042271.D	08 Apr 2025 19:03	Run @ Lower Dilution	SY/MD	Not Ok
21	Q1728-07	SSSG-B1-1-04012025	VL042272.D	08 Apr 2025 20:30	no sample to rerun straight	SY/MD	Ok,M
22	Q1728-08DL	SSSG-B2-4-04022025	VL042273.D	08 Apr 2025 21:37		SY/MD	Ok,M
23	Q1728-11	SSSG-B1-4-04012025	VL042274.D	08 Apr 2025 22:11		SY/MD	Ok,M
24	Q1728-12	SSSG-B1-3-04012025	VL042275.D	08 Apr 2025 22:42		SY/MD	Ok,M
25	Q1728-09DUP	SSSG-B2-1-04022025	VL042276.D	08 Apr 2025 23:13		SY/MD	Ok,M
26	Q1728-01	SSSG-B14-2-04012025	VL042277.D	08 Apr 2025 23:44	Need 10X	SY/MD	Dilution
27	Q1728-02	SSSG-B1-5-04012025	VL042278.D	09 Apr 2025 00:15		SY/MD	Ok,M
28	Q1728-04	SSSG-B2-7-04022025	VL042279.D	09 Apr 2025 00:46	Need 40X	SY/MD	Dilution
29	Q1728-04DL	SSSG-B2-7-04022025	VL042280.D	09 Apr 2025 01:18		SY/MD	Ok,M
30	Q1728-05	SSSG-B2-6-04022025	VL042281.D	09 Apr 2025 01:49	Need 40X	SY/MD	Dilution
31	Q1728-05DL	SSSG-B2-6-04022025	VL042282.D	09 Apr 2025 02:20		SY/MD	Ok
32	Q1728-10	SSSG-B14-4-04022025	VL042283.D	09 Apr 2025 02:51		SY/MD	Ok,M
33	Q1728-03	SSSG-B2-2-04022025	VL042284.D	09 Apr 2025 03:22		SY/MD	Ok,M

M : Manual Integration