

SOP ID :	M1312-SPLP-10	
SDG No :	N/A	Start Prep Date : 04/07/2025 Time : 14:00
Weigh By :	JP	End Prep Date : 04/08/2025 Time : 08:10
Balance ID :	WC SC-7	Combination Ratio : 20
pH Meter ID :	WC PH METER-1	ZHE Cleaning Batch : ⁷⁸ N/A
Extraction By :	JP	Initial Room Temperature: 23 °C
Filter By :	JP	Final Room Temperature: 21 °C
Pipette ID :	WC	TCLP Technician Signature : <i>JP</i>
Tumbler ID :	T-1	Supervisor By : <i>12</i>
TCLP Filter ID :	115525	

Standardized Name	MLS USED	STD REF. # FROM LOG
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A

Chemical Used	ML/SAMPLE U	Lot Number
SPLP FLUID	N/A	WP110806
N/A	N/A	N/A
HNO3-TCLP,1N	N/A	WP110804
pH Strips	N/A	W1931,W1934,W3171,W3172
pH Strips	W1941,W1942	W3166,W1938,W1939,W1940,
N/A	N/A	N/A
120ml Plastic bottle	N/A	405130101
1:1 HNO3	N/A	MP84041

Extraction Conformance/Non-Conformance Comments:

Matrix spikes are added after filtration and before preservation. TUMBLER T-1 checked,30 rpm. q1730-20 is used for MS-MSD.

Date / Time	Prepped Sample Relinquished By/Location	Received By/Location
04/08/25 10:30	<i>JP</i> <i>1200A000</i>	<i>SLB.</i> <i>1000000</i>
	Preparation Group	Analysis Group

Sample ID	ClientID	TCLP Vessel ID	Sample Wt (g)	Volume Extraction Fluid #1 (mL)	Multi phasic	Phase Miscible	Phases Combined	Final Leachate PH	Metals Leachate Adj. PH	Prep Pos
PB167493TB	LEB493	11	N/A	2000	N/A	N/A	N/A	4.23	1.5	T-1
Q1730-02	OU4-VSL-15-040325	01	100.02	2000	N/A	N/A	N/A	5.6	1.0	T-1
Q1730-04	OU4-VSL-16-040325	02	100.03	2000	N/A	N/A	N/A	5.5	1.5	T-1
Q1730-06	OU4-VSL-17-040325	03	100.04	2000	N/A	N/A	N/A	5.6	1.0	T-1
Q1730-08	OU4-PCS-TC-21-040325	04	100.02	2000	N/A	N/A	N/A	5.8	1.0	T-1
Q1730-10	OU4-PCS-TC-22-040325	05	100.03	2000	N/A	N/A	N/A	5.6	1.5	T-1
Q1730-12	OU4-PCS-TC-23-040325	06	100.02	2000	N/A	N/A	N/A	5.5	1.0	T-1
Q1730-14	OU4-PCS-TC-24-040325	07	100.01	2000	N/A	N/A	N/A	5.6	1.5	T-1
Q1730-16	OU4-PCS-TC-25-040325	08	100.02	2000	N/A	N/A	N/A	5.6	1.0	T-1
Q1730-18	OU4-PCS-TC-26-040325	09	100.03	2000	N/A	N/A	N/A	5.8	1.5	T-1
Q1730-20	OU4-CF-15-040325	10	100.02	2000	N/A	N/A	N/A	5.6	1.0	T-1

SampleID	ClientID	Sample Weight (g)	Filter Weight (g)	Filtrate (mL)	Filter + Solid (After 100°C)	% solids	% Dry Solids
PB167493TB	LEB493	N/A	N/A	N/A	N/A	N/A	N/A
Q1730-02	OU4-VSL-15-040325	N/A	N/A	N/A	N/A	100	N/A
Q1730-04	OU4-VSL-16-040325	N/A	N/A	N/A	N/A	100	N/A
Q1730-06	OU4-VSL-17-040325	N/A	N/A	N/A	N/A	100	N/A
Q1730-08	OU4-PCS-TC-21-040325	N/A	N/A	N/A	N/A	100	N/A
Q1730-10	OU4-PCS-TC-22-040325	N/A	N/A	N/A	N/A	100	N/A
Q1730-12	OU4-PCS-TC-23-040325	N/A	N/A	N/A	N/A	100	N/A
Q1730-14	OU4-PCS-TC-24-040325	N/A	N/A	N/A	N/A	100	N/A
Q1730-16	OU4-PCS-TC-25-040325	N/A	N/A	N/A	N/A	100	N/A
Q1730-18	OU4-PCS-TC-26-040325	N/A	N/A	N/A	N/A	100	N/A
Q1730-20	OU4-CF-15-040325	N/A	N/A	N/A	N/A	100	N/A

Hot Block ID : WC S-1 /WC S-2
Thermometer ID : FLASHPOINT

SampleID	ClientID	Sample Weight (g)	Volume DI Water (mL)	PH after 5 min stir	PH after 10 min stir	Extraction Fluid 1 or 2	pH Extraction Fluid
PB167493TB	LEB493	N/A	N/A	N/A	N/A	#1	4.23
Q1730-02	OU4-VSL-15-040325	N/A	N/A	N/A	N/A	#1	4.23
Q1730-04	OU4-VSL-16-040325	N/A	N/A	N/A	N/A	#1	4.23
Q1730-06	OU4-VSL-17-040325	N/A	N/A	N/A	N/A	#1	4.23
Q1730-08	OU4-PCS-TC-21-040325	N/A	N/A	N/A	N/A	#1	4.23
Q1730-10	OU4-PCS-TC-22-040325	N/A	N/A	N/A	N/A	#1	4.23
Q1730-12	OU4-PCS-TC-23-040325	N/A	N/A	N/A	N/A	#1	4.23
Q1730-14	OU4-PCS-TC-24-040325	N/A	N/A	N/A	N/A	#1	4.23
Q1730-16	OU4-PCS-TC-25-040325	N/A	N/A	N/A	N/A	#1	4.23
Q1730-18	OU4-PCS-TC-26-040325	N/A	N/A	N/A	N/A	#1	4.23
Q1730-20	OU4-CF-15-040325	N/A	N/A	N/A	N/A	#1	4.23