

SOP ID : M1311-TCLP-15
SDG No : N/A
Weigh By : N/A
Balance ID : N/A
pH Meter ID : WC PH METER-1
Extraction By : N/A
Filter By : JP
Pipette ID : N/A
Tumbler ID : N/A
TCLP Filter ID : 115525

Start Prep Date : N/A Time : N/A
End Prep Date : N/A Time : N/A
Combination Ratio : N/A
ZHE Cleaning Batch : ¹⁰N/A
Initial Room Temperature : N/A
Final Room Temperature : N/A
TCLP Technician Signature : *[Signature]*
Supervisor By : 12

Standard Name	MLS USED	STD REF. # FROM LOG
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A

Chemical Used	ML/SAMPLE U	Lot Number
N/A	N/A	N/A
N/A	N/A	N/A
HNO3-TCLP,1N	N/A	WP110804
pH Strips	N/A	W1931,W1934,W3171,W3172
pH Strips	W1941,W1942	W3166,W1938,W1939,W1940,
1 Liter Amber	N/A	90424-08
120ml Plastic bottle	N/A	405130101
1:1 HNO3	N/A	MP84041

Extraction Conformance/Non-Conformance Comments:

Matrix spikes are added after filtration and before preservation. q1739-02 is used for ms-msd.

Date / Time	Prepped Sample Relinquished By/Location	Received By/Location
04/07/25 12:30	<i>[Signature]</i> 1 CCU Room	5125. RJ / E+H
	Preparation Group	Analysis Group <i>[Signature]</i>