

Surrogate Summary

SDG No.: Q1739

Client: PARSONS Engineering of New York, Inc.

Analytical Method: 8151A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Rec	Qual	Limits	
								Low	High
I.BLK-PS029656.D	PIBLK-PS029656.D	2,4-DCAA	1	500	476	95		39	175
		2,4-DCAA	2	500	485	97		39	175
I.BLK-PS029763.D	PIBLK-PS029763.D	2,4-DCAA	1	500	594	119		39	175
		2,4-DCAA	2	500	582	116		39	175
Q1739-02	WC-LIQUID-20250404	2,4-DCAA	1	500	581	116		39	175
		2,4-DCAA	2	500	589	118		39	175
Q1739-02MS	WC-LIQUID-20250404MS	2,4-DCAA	1	500	589	118		39	175
		2,4-DCAA	2	500	645	129		39	175
Q1739-02MSD	WC-LIQUID-20250404MSD	2,4-DCAA	1	500	593	119		39	175
		2,4-DCAA	2	500	656	131		39	175
PB167536BL	PB167536BL	2,4-DCAA	1	500	559	112		39	175
		2,4-DCAA	2	500	539	108		39	175
PB167536BS	PB167536BS	2,4-DCAA	1	500	556	111		39	175
		2,4-DCAA	2	500	581	116		39	175
PB167488TB	PB167488TB	2,4-DCAA	1	500	626	125		39	175
		2,4-DCAA	2	500	583	117		39	175
I.BLK-PS029772.D	PIBLK-PS029772.D	2,4-DCAA	1	500	618	124		39	175
		2,4-DCAA	2	500	610	122		39	175