

Surrogate Summary

SDG No.: Q1800

Client: ENTACT

Analytical Method: 8151A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Rec	Qual	Limits	
								Low	High
I.BLK-PS029656.D	PIBLK-PS029656.D	2,4-DCAA	1	500	476	95		70 (39)	130 (175)
		2,4-DCAA	2	500	485	97		70 (39)	130 (175)
I.BLK-PS029822.D	PIBLK-PS029822.D	2,4-DCAA	1	500	581	116		70 (39)	130 (175)
		2,4-DCAA	2	500	547	109		70 (39)	130 (175)
PB167608BL	PB167608BL	2,4-DCAA	1	500	681	136	*	70 (39)	130 (175)
		2,4-DCAA	2	500	616	123		70 (39)	130 (175)
PB167587TB	PB167587TB	2,4-DCAA	1	500	638	128		70 (39)	130 (175)
		2,4-DCAA	2	500	603	121		70 (39)	130 (175)
Q1800-03	WC-A4-01-C	2,4-DCAA	1	500	543	109		70 (39)	130 (175)
		2,4-DCAA	2	500	537	107		70 (39)	130 (175)
Q1800-03MS	WC-A4-01-CMS	2,4-DCAA	1	500	519	104		70 (39)	130 (175)
		2,4-DCAA	2	500	500	100		70 (39)	130 (175)
Q1800-03MSD	WC-A4-01-CMSD	2,4-DCAA	1	500	516	103		70 (39)	130 (175)
		2,4-DCAA	2	500	491	98		70 (39)	130 (175)
I.BLK-PS029832.D	PIBLK-PS029832.D	2,4-DCAA	1	500	567	113		70 (39)	130 (175)
		2,4-DCAA	2	500	532	106		70 (39)	130 (175)
I.BLK-PS029835.D	PIBLK-PS029835.D	2,4-DCAA	1	500	558	112		70 (39)	130 (175)
		2,4-DCAA	2	500	522	104		70 (39)	130 (175)
PB167608BS	PB167608BS	2,4-DCAA	1	500	550	110		70 (39)	130 (175)
		2,4-DCAA	2	500	514	103		70 (39)	130 (175)
I.BLK-PS029838.D	PIBLK-PS029838.D	2,4-DCAA	1	500	522	104		70 (39)	130 (175)
		2,4-DCAA	2	500	532	106		70 (39)	130 (175)