

Instrument ID: MSVOA_L

Daily Analysis Runlog For Sequence/QC Batch ID # VL041725

Review By	Semsettin Yesilyurt	Review On	4/18/2025 4:23:02 PM
Supervise By	Maresh Dadoda	Supervise On	4/18/2025 10:46:12 PM
SubDirectory	VL041725	HP Acquire Method	TO15AIR.M HP Processing Method VL041725AIR.M

STD. NAME	STD REF.#
Tune/Reschk	AP2590
Initial Calibration Stds	AP2598,AP2599,AP2600
CCC	AP2598
Internal Standard/PEM	AP2590
ICV/I.BLK	AP2601
Surrogate Standard	
MS/MSD Standard	
LCS Standard	

Sr#	SampleID	ClientID	Data File Name	Date-Time	Comment	Operator	Status
1	BFB	BFB	VL042348.D	17 Apr 2025 07:50		SY/MD	Ok
2	VSTDICCC010	VSTDICCC010	VL042349.D	17 Apr 2025 09:44		SY/MD	Ok,M
3	VSTDICC002	VSTDICC002	VL042350.D	17 Apr 2025 10:17		SY/MD	Ok,M
4	VSTDICC001	VSTDICC001	VL042351.D	17 Apr 2025 10:48		SY/MD	Ok,M
5	VSTDICC0.5	VSTDICC0.5	VL042352.D	17 Apr 2025 11:19		SY/MD	Ok,M
6	VSTDICC0.1	VSTDICC0.1	VL042353.D	17 Apr 2025 11:50		SY/MD	Ok,M
7	VSTDICC0.03	VSTDICC0.03	VL042354.D	17 Apr 2025 12:22		SY/MD	Ok,M
8	VSTDICC015	VSTDICC015	VL042355.D	17 Apr 2025 12:54		SY/MD	Ok,M
9	VSTDICV010	ICVVL041725	VL042356.D	17 Apr 2025 13:43		SY/MD	Ok,M
10	VL0417ABL01	VL0417ABL01	VL042357.D	17 Apr 2025 14:23		SY/MD	Ok
11	VL0417ABS01	VL0417ABS01	VL042358.D	17 Apr 2025 15:39		SY/MD	Ok,M
12	Q1665-15	IA-B14-1-3242025	VL042359.D	17 Apr 2025 16:12	Need 10X	SY/MD	Dilution
13	Q1665-15DL	IA-B14-1-3242025DL	VL042360.D	17 Apr 2025 16:44		SY/MD	Ok,M
14	Q1668-01	IA-B1-4-3252025	VL042361.D	17 Apr 2025 17:17	Need 4X	SY/MD	Dilution
15	Q1668-01DUP	IA-B1-4-3252025DUP	VL042362.D	17 Apr 2025 17:51		SY/MD	Ok,M
16	Q1668-03	IA-B1-3-3252025	VL042363.D	17 Apr 2025 18:24	Need 4X	SY/MD	Dilution
17	Q1668-04	IA-B1-1-3252025	VL042364.D	17 Apr 2025 18:58		SY/MD	Ok,M
18	Q1668-02	IA-B1-2-3252025	VL042365.D	17 Apr 2025 19:31		SY/MD	Ok,M

Instrument ID: MSVOA_L

Daily Analysis Runlog For Sequence/QC Batch ID # VL041725

Review By	Semsettin Yesilyurt	Review On	4/18/2025 4:23:02 PM		
Supervise By	Mahesh Dadoda	Supervise On	4/18/2025 10:46:12 PM		
SubDirectory	VL041725	HP Acquire Method	TO15AIR.M	HP Processing Method	VL041725AIR.M
STD. NAME	STD REF.#				
Tune/Reschk	AP2590				
Initial Calibration Stds	AP2598,AP2599,AP2600				
CCC	AP2598				
Internal Standard/PEM	AP2590				
ICV/I.BLK	AP2601				
Surrogate Standard					
MS/MSD Standard					
LCS Standard					

19	Q1726-01	IA-DUP-3-03312025	VL042366.D	17 Apr 2025 20:05	Need 20X	SY/MD	Dilution
20	Q1726-02	IA-B2-6-03312025	VL042367.D	17 Apr 2025 20:38	Need 20X	SY/MD	Dilution
21	Q1726-03	IA-B2-2-03312025	VL042368.D	17 Apr 2025 21:12	Need 20X	SY/MD	Dilution
22	Q1726-04	IA-B2-5-03312025	VL042369.D	17 Apr 2025 21:45	Need 20X	SY/MD	Dilution
23	Q1726-05	IA-B2-1-03312025	VL042370.D	17 Apr 2025 22:19	Need 20X	SY/MD	Dilution
24	Q1726-06	IA-B2-3-03312025	VL042371.D	17 Apr 2025 22:52	Need 20X	SY/MD	Dilution
25	Q1726-07	IA-B2-4-03312025	VL042372.D	17 Apr 2025 23:26	Need 20X	SY/MD	Dilution
26	Q1726-07DL	IA-B2-4-03312025DL	VL042373.D	17 Apr 2025 23:57	Need 40X	SY/MD	Dilution
27	Q1668-01DL	IA-B1-4-3252025DL	VL042374.D	18 Apr 2025 00:28		SY/MD	Ok,M
28	Q1668-03DL	IA-B1-3-3252025DL	VL042375.D	18 Apr 2025 00:59		SY/MD	Ok,M
29	Q1726-01DL	IA-DUP-3-03312025DL	VL042376.D	18 Apr 2025 01:30	Need 40X	SY/MD	Dilution
30	Q1726-02DL	IA-B2-6-03312025DL	VL042377.D	18 Apr 2025 02:01	Need 40X	SY/MD	Dilution
31	Q1726-03DL	IA-B2-2-03312025DL	VL042378.D	18 Apr 2025 02:32	Need 100X	SY/MD	Dilution
32	Q1726-04DL	IA-B2-5-03312025DL	VL042379.D	18 Apr 2025 03:03	Need 40X	SY/MD	Dilution
33	Q1726-05DL	IA-B2-1-03312025DL	VL042380.D	18 Apr 2025 03:34	Need 100X	SY/MD	Dilution
34	Q1726-06DL	IA-B2-3-03312025DL	VL042381.D	18 Apr 2025 04:05	Need 100X	SY/MD	Dilution

M : Manual Integration

Instrument ID: MSVOA_L

Daily Analysis Runlog For Sequence/QC Batch ID # VL042125

Review By	Semsettin Yesilyurt	Review On	4/22/2025 9:39:52 AM
Supervise By	Maresh Dadoda	Supervise On	4/22/2025 2:00:33 PM
SubDirectory	VL042125	HP Acquire Method	TO15AIR.M
		HP Processing Method	VL041725AIR.M

STD. NAME	STD REF.#
Tune/Reschk	AP2602
Initial Calibration Stds	AP2598,AP2599,AP2600
CCC	AP2598
Internal Standard/PEM	AP2602
ICV/I.BLK	AP2601
Surrogate Standard	
MS/MSD Standard	
LCS Standard	

Sr#	SampleID	ClientID	Data File Name	Date-Time	Comment	Operator	Status
1	BFB	BFB	VL042394.D	21 Apr 2025 07:58		SY/MD	Ok
2	VSTDCCC010	VSTDCCC010	VL042395.D	21 Apr 2025 08:41		SY/MD	Ok,M
3	VL0421ABL01	VL0421ABL01	VL042396.D	21 Apr 2025 09:28		SY/MD	Ok
4	VL0421ABS01	VL0421ABS01	VL042397.D	21 Apr 2025 10:14		SY/MD	Ok,M
5	QC041625 CAN 10329	QC041625 CAN 10329	VL042398.D	21 Apr 2025 11:04	Contaminated	SY/MD	Not Ok
6	QC041625 CAN 10329	QC041625 CAN 10329	VL042399.D	21 Apr 2025 11:36		SY/MD	Ok
7	Q1801-01	SS2R	VL042400.D	21 Apr 2025 12:07		SY/MD	Ok,M
8	Q1801-02	SS3	VL042401.D	21 Apr 2025 12:51		SY/MD	Ok,M
9	Q1844-02	IA1	VL042402.D	21 Apr 2025 13:38	Need 5X	SY/MD	Dilution
10	Q1844-02DUP	IA1DUP	VL042403.D	21 Apr 2025 14:13		SY/MD	Ok,M
11	Q1844-02DL	IA1DL	VL042404.D	21 Apr 2025 14:44		SY/MD	Ok,M
12	Q1843-02	IA1	VL042405.D	21 Apr 2025 15:18		SY/MD	Ok,M
13	Q1843-02DL	IA1DL	VL042406.D	21 Apr 2025 15:49	Not required	SY/MD	Not Ok
14	Q1844-01	SV1	VL042407.D	21 Apr 2025 16:20	Need 40X	SY/MD	Dilution
15	Q1844-01DL	SV1DL	VL042408.D	21 Apr 2025 16:51	Need 600X	SY/MD	Dilution
16	Q1843-01	SV1	VL042409.D	21 Apr 2025 17:23		SY/MD	Ok,M
17	Q1843-01DL	SV1DL	VL042410.D	21 Apr 2025 17:54	Not required	SY/MD	Not Ok
18	Q1843-01DUP	SV1DUP	VL042411.D	21 Apr 2025 19:16		SY/MD	Ok,M

Instrument ID: MSVOA_L

Daily Analysis Runlog For Sequence/QC Batch ID # VL042125

Review By	Semsettin Yesilyurt	Review On	4/22/2025 9:39:52 AM
Supervise By	Maresh Dadoda	Supervise On	4/22/2025 2:00:33 PM
SubDirectory	VL042125	HP Acquire Method	TO15AIR.M
		HP Processing Method	VL041725AIR.M

STD. NAME	STD REF.#
Tune/Reschk	AP2602
Initial Calibration Stds	AP2598,AP2599,AP2600
CCC	AP2598
Internal Standard/PEM	AP2602
ICV/I.BLK	AP2601
Surrogate Standard	
MS/MSD Standard	
LCS Standard	

19	QC041725 CA 27806	QC041725 CA 27806	VL042412.D	21 Apr 2025 19:47	Contaminated	SY/MD	Not Ok
20	QC041725 CA 27806	QC041725 CA 27806	VL042413.D	21 Apr 2025 20:18		SY/MD	Ok
21	Q1844-01DL2	SV1DL2	VL042414.D	22 Apr 2025 07:47		SY/MD	Ok,M

M : Manual Integration