

Instrument ID: MSVOA_L

Daily Analysis Runlog For Sequence/QC Batch ID # VL032725

Review By	Maresh Dadoda	Review On	4/3/2025 8:14:07 AM
Supervise By	Semsettin Yesilyurt	Supervise On	4/3/2025 8:45:16 AM
SubDirectory	VL032725	HP Acquire Method	TO15AIR.M HP Processing Method VL032725AIR.M

STD. NAME	STD REF.#
Tune/Reschk	AP2589
Initial Calibration Stds	AP2593,AP2595,AP2596
CCC	AP2593
Internal Standard/PEM	AP2589
ICV/I.BLK	AP2594
Surrogate Standard	
MS/MSD Standard	
LCS Standard	

Sr#	SampleID	ClientID	Data File Name	Date-Time	Comment	Operator	Status
1	BFB	BFB	VL042171.D	27 Mar 2025 07:44		SY/MD	Ok
2	VSTDICCC010	VSTDICCC010	VL042172.D	27 Mar 2025 08:26		SY/MD	Ok,M
3	VSTDICC002	VSTDICC002	VL042173.D	27 Mar 2025 08:59		SY/MD	Ok,M
4	VSTDICC001	VSTDICC001	VL042174.D	27 Mar 2025 09:30		SY/MD	Ok,M
5	VSTDICC0.5	VSTDICC0.5	VL042175.D	27 Mar 2025 10:01		SY/MD	Ok,M
6	VSTDICC0.1	VSTDICC0.1	VL042176.D	27 Mar 2025 10:32		SY/MD	Ok,M
7	VSTDICC0.03	VSTDICC0.03	VL042177.D	27 Mar 2025 11:03		SY/MD	Ok,M
8	VSTDICC015	VSTDICC015	VL042178.D	27 Mar 2025 11:36		SY/MD	Ok,M
9	VSTDICV010	ICVVL032725	VL042179.D	27 Mar 2025 12:23		SY/MD	Ok,M
10	VL0327ABL01	VL0327ABL01	VL042180.D	27 Mar 2025 13:50		SY/MD	Ok
11	VL0327ABS01	VL0327ABS01	VL042181.D	27 Mar 2025 14:39		SY/MD	Ok,M
12	QC032625 CAN 10448	QC032625 CAN 10448	VL042182.D	27 Mar 2025 16:41		SY/MD	Ok
13	Q1649-01	SV1	VL042183.D	27 Mar 2025 17:13	rerun 4X,not required	SY/MD	Not Ok
14	Q1649-02	SV2	VL042184.D	27 Mar 2025 17:45	rerun 4X,not required	SY/MD	Not Ok
15	Q1669-01	IA-B9-3-3262025	VL042185.D	27 Mar 2025 18:22	Need 5x	SY/MD	Dilution
16	Q1669-02	IA-B9-2-3262025	VL042186.D	27 Mar 2025 18:59		SY/MD	Ok,M
17	Q1669-03	IA-B9-1-3262025	VL042187.D	27 Mar 2025 19:36		SY/MD	Ok,M
18	Q1669-04	AA-2-3252025	VL042188.D	27 Mar 2025 20:12		SY/MD	Ok,M

Instrument ID: MSVOA_L

Daily Analysis Runlog For Sequence/QC Batch ID # VL032725

Review By	Maresh Dadoda	Review On	4/3/2025 8:14:07 AM
Supervise By	Semsettin Yesilyurt	Supervise On	4/3/2025 8:45:16 AM
SubDirectory	VL032725	HP Acquire Method	TO15AIR.M HP Processing Method VL032725AIR.M
STD. NAME	STD REF.#		
Tune/Reschk	AP2589		
Initial Calibration Stds	AP2593,AP2595,AP2596		
CCC	AP2593		
Internal Standard/PEM	AP2589		
ICV/I.BLK	AP2594		
Surrogate Standard			
MS/MSD Standard			
LCS Standard			

19	Q1669-09	AA-3-3252025	VL042189.D	27 Mar 2025 20:48		SY/MD	Ok,M
20	Q1669-09DUP	AA-3-3252025DUP	VL042190.D	27 Mar 2025 21:24		SY/MD	Ok,M
21	Q1669-05	SSSG-B16-3-3262025	VL042191.D	27 Mar 2025 21:57	Need straight run	SY/MD	Not Ok
22	Q1669-07DL	SSSG-B16-1-3262025	VL042192.D	27 Mar 2025 22:29		SY/MD	Ok,M
23	Q1669-08	SSSG-B16-2-3262025	VL042193.D	27 Mar 2025 23:01	Need straight run	SY/MD	Not Ok
24	VIBLK	VIBLK	VL042194.D	27 Mar 2025 23:34		SY/MD	Ok
25	Q1669-01DL	IA-B9-3-3262025DL	VL042195.D	28 Mar 2025 00:06		SY/MD	Ok,M
26	Q1669-01DUP	IA-B9-3-3262025DUP	VL042196.D	28 Mar 2025 00:44		SY/MD	Ok,M
27	Q1669-05	SSSG-B16-3-3262025	VL042197.D	28 Mar 2025 01:20		SY/MD	Ok,M
28	Q1669-07	SSSG-B16-1-3262025	VL042198.D	28 Mar 2025 01:57	Need 10X	SY/MD	Dilution
29	Q1669-08	SSSG-B16-2-3262025	VL042199.D	28 Mar 2025 02:33		SY/MD	Ok,M

M : Manual Integration

Instrument ID: MSVOA_L

Daily Analysis Runlog For Sequence/QC Batch ID # VL041525

Review By	Semsettin Yesilyurt	Review On	4/16/2025 2:19:24 PM
Supervise By	Maresh Dadoda	Supervise On	4/16/2025 3:47:46 PM
SubDirectory	VL041525	HP Acquire Method	TO15AIR.M
		HP Processing Method	VL032725AIR.M
STD. NAME	STD REF.#		
Tune/Reschk	AP2591		
Initial Calibration Stds	AP2593,AP2595,AP2596		
CCC	AP2593		
Internal Standard/PEM	AP2591		
ICV/I.BLK	AP2594		
Surrogate Standard			
MS/MSD Standard			
LCS Standard			

Sr#	SampleID	ClientID	Data File Name	Date-Time	Comment	Operator	Status
1	BFB	BFB	VL042325.D	15 Apr 2025 07:45		SY/MD	Ok
2	VSTDCCC010	VSTDCCC010	VL042326.D	15 Apr 2025 11:30		SY/MD	Ok,M
3	VL0415ABL01	VL0415ABL01	VL042327.D	15 Apr 2025 12:12		SY/MD	Ok
4	VL0415ABS01	VL0415ABS01	VL042328.D	15 Apr 2025 12:54		SY/MD	Ok,M
5	QC041425 CAN 10401	QC041425 CAN 10401	VL042329.D	15 Apr 2025 13:26	carryover	SY/MD	Not Ok
6	QC041425 CAN 10401	QC041425 CAN 10401	VL042330.D	15 Apr 2025 13:57		SY/MD	Ok
7	Q1665-18DL	IA-B10-3-3242025DL	VL042331.D	15 Apr 2025 14:29		SY/MD	Ok
8	Q1665-10DL2	IA-B10-4-3242025DL2	VL042332.D	15 Apr 2025 15:01		SY/MD	Ok,M
9	Q1802-02	IA1	VL042333.D	15 Apr 2025 15:35		SY/MD	Ok,M
10	Q1802-02DUP	IA1DUP	VL042334.D	15 Apr 2025 16:08		SY/MD	Ok,M
11	Q1802-03	OA1	VL042335.D	15 Apr 2025 16:42		SY/MD	Ok,M
12	Q1802-01	SV1	VL042336.D	15 Apr 2025 17:14	Need 40X	SY/MD	Dilution
13	Q1802-01DL	SV1DL	VL042337.D	15 Apr 2025 17:46	Need further dilution	SY/MD	Not Ok
14	Q1802-01DL	SV1DL	VL042338.D	15 Apr 2025 18:21		SY/MD	Ok,M

M : Manual Integration