

Instrument ID: MSVOA\_L

**Daily Analysis Runlog For Sequence/QC Batch ID # VL041725**

|              |                     |                   |                                              |
|--------------|---------------------|-------------------|----------------------------------------------|
| Review By    | Semsettin Yesilyurt | Review On         | 4/18/2025 4:23:02 PM                         |
| Supervise By | Maresh Dadoda       | Supervise On      | 4/18/2025 10:46:12 PM                        |
| SubDirectory | VL041725            | HP Acquire Method | TO15AIR.M HP Processing Method VL041725AIR.M |

| STD. NAME                | STD REF.#            |
|--------------------------|----------------------|
| Tune/Reschk              | AP2590               |
| Initial Calibration Stds | AP2598,AP2599,AP2600 |
| CCC                      | AP2598               |
| Internal Standard/PEM    | AP2590               |
| ICV/I.BLK                | AP2601               |
| Surrogate Standard       |                      |
| MS/MSD Standard          |                      |
| LCS Standard             |                      |

| Sr# | SampleID    | ClientID           | Data File Name | Date-Time         | Comment  | Operator | Status   |
|-----|-------------|--------------------|----------------|-------------------|----------|----------|----------|
| 1   | BFB         | BFB                | VL042348.D     | 17 Apr 2025 07:50 |          | SY/MD    | Ok       |
| 2   | VSTDICCC010 | VSTDICCC010        | VL042349.D     | 17 Apr 2025 09:44 |          | SY/MD    | Ok,M     |
| 3   | VSTDICC002  | VSTDICC002         | VL042350.D     | 17 Apr 2025 10:17 |          | SY/MD    | Ok,M     |
| 4   | VSTDICC001  | VSTDICC001         | VL042351.D     | 17 Apr 2025 10:48 |          | SY/MD    | Ok,M     |
| 5   | VSTDICC0.5  | VSTDICC0.5         | VL042352.D     | 17 Apr 2025 11:19 |          | SY/MD    | Ok,M     |
| 6   | VSTDICC0.1  | VSTDICC0.1         | VL042353.D     | 17 Apr 2025 11:50 |          | SY/MD    | Ok,M     |
| 7   | VSTDICC0.03 | VSTDICC0.03        | VL042354.D     | 17 Apr 2025 12:22 |          | SY/MD    | Ok,M     |
| 8   | VSTDICC015  | VSTDICC015         | VL042355.D     | 17 Apr 2025 12:54 |          | SY/MD    | Ok,M     |
| 9   | VSTDICV010  | ICVVL041725        | VL042356.D     | 17 Apr 2025 13:43 |          | SY/MD    | Ok,M     |
| 10  | VL0417ABL01 | VL0417ABL01        | VL042357.D     | 17 Apr 2025 14:23 |          | SY/MD    | Ok       |
| 11  | VL0417ABS01 | VL0417ABS01        | VL042358.D     | 17 Apr 2025 15:39 |          | SY/MD    | Ok,M     |
| 12  | Q1665-15    | IA-B14-1-3242025   | VL042359.D     | 17 Apr 2025 16:12 | Need 10X | SY/MD    | Dilution |
| 13  | Q1665-15DL  | IA-B14-1-3242025DL | VL042360.D     | 17 Apr 2025 16:44 |          | SY/MD    | Ok,M     |
| 14  | Q1668-01    | IA-B1-4-3252025    | VL042361.D     | 17 Apr 2025 17:17 | Need 4X  | SY/MD    | Dilution |
| 15  | Q1668-01DUP | IA-B1-4-3252025DUP | VL042362.D     | 17 Apr 2025 17:51 |          | SY/MD    | Ok,M     |
| 16  | Q1668-03    | IA-B1-3-3252025    | VL042363.D     | 17 Apr 2025 18:24 | Need 4X  | SY/MD    | Dilution |
| 17  | Q1668-04    | IA-B1-1-3252025    | VL042364.D     | 17 Apr 2025 18:58 |          | SY/MD    | Ok,M     |
| 18  | Q1668-02    | IA-B1-2-3252025    | VL042365.D     | 17 Apr 2025 19:31 |          | SY/MD    | Ok,M     |

Instrument ID: MSVOA\_L

**Daily Analysis Runlog For Sequence/QC Batch ID # VL041725**

| Review By                | Semsettin Yesilyurt  | Review On         | 4/18/2025 4:23:02 PM  |                      |               |
|--------------------------|----------------------|-------------------|-----------------------|----------------------|---------------|
| Supervise By             | Mahesh Dadoda        | Supervise On      | 4/18/2025 10:46:12 PM |                      |               |
| SubDirectory             | VL041725             | HP Acquire Method | TO15AIR.M             | HP Processing Method | VL041725AIR.M |
| STD. NAME                | STD REF.#            |                   |                       |                      |               |
| Tune/Reschk              | AP2590               |                   |                       |                      |               |
| Initial Calibration Stds | AP2598,AP2599,AP2600 |                   |                       |                      |               |
| CCC                      | AP2598               |                   |                       |                      |               |
| Internal Standard/PEM    | AP2590               |                   |                       |                      |               |
| ICV/I.BLK                | AP2601               |                   |                       |                      |               |
| Surrogate Standard       |                      |                   |                       |                      |               |
| MS/MSD Standard          |                      |                   |                       |                      |               |
| LCS Standard             |                      |                   |                       |                      |               |

|    |            |                     |            |                   |           |       |          |
|----|------------|---------------------|------------|-------------------|-----------|-------|----------|
| 19 | Q1726-01   | IA-DUP-3-03312025   | VL042366.D | 17 Apr 2025 20:05 | Need 20X  | SY/MD | Dilution |
| 20 | Q1726-02   | IA-B2-6-03312025    | VL042367.D | 17 Apr 2025 20:38 | Need 20X  | SY/MD | Dilution |
| 21 | Q1726-03   | IA-B2-2-03312025    | VL042368.D | 17 Apr 2025 21:12 | Need 20X  | SY/MD | Dilution |
| 22 | Q1726-04   | IA-B2-5-03312025    | VL042369.D | 17 Apr 2025 21:45 | Need 20X  | SY/MD | Dilution |
| 23 | Q1726-05   | IA-B2-1-03312025    | VL042370.D | 17 Apr 2025 22:19 | Need 20X  | SY/MD | Dilution |
| 24 | Q1726-06   | IA-B2-3-03312025    | VL042371.D | 17 Apr 2025 22:52 | Need 20X  | SY/MD | Dilution |
| 25 | Q1726-07   | IA-B2-4-03312025    | VL042372.D | 17 Apr 2025 23:26 | Need 20X  | SY/MD | Dilution |
| 26 | Q1726-07DL | IA-B2-4-03312025DL  | VL042373.D | 17 Apr 2025 23:57 | Need 40X  | SY/MD | Dilution |
| 27 | Q1668-01DL | IA-B1-4-3252025DL   | VL042374.D | 18 Apr 2025 00:28 |           | SY/MD | Ok,M     |
| 28 | Q1668-03DL | IA-B1-3-3252025DL   | VL042375.D | 18 Apr 2025 00:59 |           | SY/MD | Ok,M     |
| 29 | Q1726-01DL | IA-DUP-3-03312025DL | VL042376.D | 18 Apr 2025 01:30 | Need 40X  | SY/MD | Dilution |
| 30 | Q1726-02DL | IA-B2-6-03312025DL  | VL042377.D | 18 Apr 2025 02:01 | Need 40X  | SY/MD | Dilution |
| 31 | Q1726-03DL | IA-B2-2-03312025DL  | VL042378.D | 18 Apr 2025 02:32 | Need 100X | SY/MD | Dilution |
| 32 | Q1726-04DL | IA-B2-5-03312025DL  | VL042379.D | 18 Apr 2025 03:03 | Need 40X  | SY/MD | Dilution |
| 33 | Q1726-05DL | IA-B2-1-03312025DL  | VL042380.D | 18 Apr 2025 03:34 | Need 100X | SY/MD | Dilution |
| 34 | Q1726-06DL | IA-B2-3-03312025DL  | VL042381.D | 18 Apr 2025 04:05 | Need 100X | SY/MD | Dilution |

M : Manual Integration

Instrument ID: MSVOA\_L

**Daily Analysis Runlog For Sequence/QC Batch ID # VL042825**

|              |                     |                      |                      |
|--------------|---------------------|----------------------|----------------------|
| Review By    | Semsettin Yesilyurt | Review On            | 4/30/2025 2:24:11 PM |
| Supervise By | Mahesh Dadoda       | Supervise On         | 4/30/2025 2:30:41 PM |
| SubDirectory | VL042825            | HP Acquire Method    | TO15AIR.M            |
|              |                     | HP Processing Method | VL041725AIR.M        |

| STD. NAME                | STD REF.#            |
|--------------------------|----------------------|
| Tune/Reschk              | AP2602               |
| Initial Calibration Stds | AP2598,AP2599,AP2600 |
| CCC                      | AP2598               |
| Internal Standard/PEM    | AP2602               |
| ICV/I.BLK                | AP2601               |
| Surrogate Standard       |                      |
| MS/MSD Standard          |                      |
| LCS Standard             |                      |

| Sr# | SampleID    | ClientID              | Data File Name | Date-Time         | Comment  | Operator | Status   |
|-----|-------------|-----------------------|----------------|-------------------|----------|----------|----------|
| 1   | BFB         | BFB                   | VL042415.D     | 28 Apr 2025 07:59 |          | SY/MD    | Ok       |
| 2   | VSTDCCC010  | VSTDCCC010            | VL042416.D     | 28 Apr 2025 08:44 |          | SY/MD    | Ok,M     |
| 3   | VL0428ABL01 | VL0428ABL01           | VL042417.D     | 28 Apr 2025 09:27 |          | SY/MD    | Ok       |
| 4   | VL0428ABS01 | VL0428ABS01           | VL042418.D     | 28 Apr 2025 10:10 |          | SY/MD    | Ok,M     |
| 5   | Q1665-05    | IA-B2-3-3242025       | VL042419.D     | 28 Apr 2025 10:59 |          | SY/MD    | Ok,M     |
| 6   | Q1665-07    | IA-DUP-1-3242025      | VL042420.D     | 28 Apr 2025 11:33 | Need 4x  | SY/MD    | Dilution |
| 7   | Q1665-12    | IA-B2-6-3242025       | VL042421.D     | 28 Apr 2025 12:17 |          | SY/MD    | Ok,M     |
| 8   | Q1665-16    | IA-B2-5-3242025       | VL042422.D     | 28 Apr 2025 12:51 |          | SY/MD    | Ok,M     |
| 9   | Q1665-19    | IA-B2-2-3242025       | VL042423.D     | 28 Apr 2025 13:26 | Need 4x  | SY/MD    | Dilution |
| 10  | Q1665-21    | IA-B2-1-3242025       | VL042424.D     | 28 Apr 2025 14:00 | Need 4X  | SY/MD    | Dilution |
| 11  | Q1665-22    | IA-B2-4-3242025       | VL042425.D     | 28 Apr 2025 14:34 |          | SY/MD    | Ok,M     |
| 12  | Q1880-02    | EBR-P02-SL-IA02-2025  | VL042426.D     | 28 Apr 2025 15:08 |          | SY/MD    | Ok,M     |
| 13  | Q1880-03    | EBR-P08-SC-IA01-2025  | VL042427.D     | 28 Apr 2025 15:42 |          | SY/MD    | Ok,M     |
| 14  | Q1880-08    | EBR-P02-SC-IA01-2025  | VL042428.D     | 28 Apr 2025 16:16 |          | SY/MD    | Ok,M     |
| 15  | Q1880-10    | EBR-P16-SC-AA01-2025  | VL042429.D     | 28 Apr 2025 16:50 |          | SY/MD    | Ok,M     |
| 16  | Q1880-10DUP | EBR-P16-SC-AA01-2025  | VL042430.D     | 28 Apr 2025 17:26 |          | SY/MD    | Ok,M     |
| 17  | Q1880-01    | EBR-P11-SC-SS01-2025  | VL042431.D     | 28 Apr 2025 18:00 | Need 20X | SY/MD    | Dilution |
| 18  | Q1880-06    | EBR-P02-SC-DUP01-2025 | VL042432.D     | 28 Apr 2025 18:35 |          | SY/MD    | Ok,M     |

Instrument ID: MSVOA\_L

**Daily Analysis Runlog For Sequence/QC Batch ID # VL042825**

|              |                     |                      |                      |
|--------------|---------------------|----------------------|----------------------|
| Review By    | Semsettin Yesilyurt | Review On            | 4/30/2025 2:24:11 PM |
| Supervise By | Maresh Dadoda       | Supervise On         | 4/30/2025 2:30:41 PM |
| SubDirectory | VL042825            | HP Acquire Method    | TO15AIR.M            |
|              |                     | HP Processing Method | VL041725AIR.M        |

| STD. NAME                | STD REF.#            |
|--------------------------|----------------------|
| Tune/Reschk              | AP2602               |
| Initial Calibration Stds | AP2598,AP2599,AP2600 |
| CCC                      | AP2598               |
| Internal Standard/PEM    | AP2602               |
| ICV/I.BLK                | AP2601               |
| Surrogate Standard       |                      |
| MS/MSD Standard          |                      |
| LCS Standard             |                      |

|    |          |                     |            |                   |                |       |        |
|----|----------|---------------------|------------|-------------------|----------------|-------|--------|
| 19 | Q1880-09 | EBR-P16-SC-AA-DUP-  | VL042433.D | 28 Apr 2025 19:08 |                | SY/MD | Ok,M   |
| 20 | Q1880-04 | EBR-P08-SC-SS01-202 | VL042434.D | 28 Apr 2025 19:41 | rerun straight | SY/MD | Not Ok |

M : Manual Integration

Instrument ID: MSVOA\_L

**Daily Analysis Runlog For Sequence/QC Batch ID # VL042925**

|              |                     |                      |                      |
|--------------|---------------------|----------------------|----------------------|
| Review By    | Semsettin Yesilyurt | Review On            | 4/30/2025 2:23:53 PM |
| Supervise By | Maresh Dadoda       | Supervise On         | 4/30/2025 2:28:33 PM |
| SubDirectory | VL042925            | HP Acquire Method    | TO15AIR.M            |
|              |                     | HP Processing Method | VL041725AIR.M        |

| STD. NAME                | STD REF.#            |
|--------------------------|----------------------|
| Tune/Reschk              | AP2602               |
| Initial Calibration Stds | AP2598,AP2599,AP2600 |
| CCC                      | AP2598               |
| Internal Standard/PEM    | AP2602               |
| ICV/I.BLK                | AP2601               |
| Surrogate Standard       |                      |
| MS/MSD Standard          |                      |
| LCS Standard             |                      |

| Sr# | SampleID    | ClientID            | Data File Name | Date-Time         | Comment               | Operator | Status   |
|-----|-------------|---------------------|----------------|-------------------|-----------------------|----------|----------|
| 1   | BFB         | BFB                 | VL042435.D     | 29 Apr 2025 07:55 |                       | SY/MD    | Ok       |
| 2   | VSTDCCC010  | VSTDCCC010          | VL042436.D     | 29 Apr 2025 08:39 |                       | SY/MD    | Ok,M     |
| 3   | VL0429ABL01 | VL0429ABL01         | VL042437.D     | 29 Apr 2025 09:25 |                       | SY/MD    | Ok       |
| 4   | VL0429ABS01 | VL0429ABS01         | VL042438.D     | 29 Apr 2025 10:08 |                       | SY/MD    | Ok,M     |
| 5   | Q1880-04    | EBR-P08-SC-SS01-202 | VL042439.D     | 29 Apr 2025 10:59 |                       | SY/MD    | Ok,M     |
| 6   | Q1880-05DL  | EBR-SC-SS-DUP01-202 | VL042440.D     | 29 Apr 2025 11:30 | Need Straight Run     | SY/MD    | Not Ok   |
| 7   | Q1880-07DL  | EBR-P16-SC-SS01-202 | VL042441.D     | 29 Apr 2025 12:02 |                       | SY/MD    | Ok,M     |
| 8   | Q1880-11DL  | EBR-P02-SC-SS01-202 | VL042442.D     | 29 Apr 2025 12:33 | Need Straight Run     | SY/MD    | Not Ok   |
| 9   | Q1908-01    | B-1                 | VL042443.D     | 29 Apr 2025 13:11 |                       | SY/MD    | Ok,M     |
| 10  | Q1908-02    | B-2                 | VL042444.D     | 29 Apr 2025 13:45 | rerun confirmation,4X | SY/MD    | Not Ok   |
| 11  | Q1908-03    | B-3                 | VL042445.D     | 29 Apr 2025 14:19 |                       | SY/MD    | Ok,M     |
| 12  | Q1908-03DL  | B-3DL               | VL042446.D     | 29 Apr 2025 14:51 | Not Required          | SY/MD    | Not Ok   |
| 13  | Q1880-05    | EBR-SC-SS-DUP01-202 | VL042447.D     | 29 Apr 2025 15:24 |                       | SY/MD    | Ok,M     |
| 14  | Q1880-07    | EBR-P16-SC-SS01-202 | VL042448.D     | 29 Apr 2025 15:58 | Need 10X              | SY/MD    | Dilution |
| 15  | Q1880-11    | EBR-P02-SC-SS01-202 | VL042449.D     | 29 Apr 2025 16:32 |                       | SY/MD    | Ok,M     |
| 16  | Q1880-01DL  | EBR-P11-SC-SS01-202 | VL042450.D     | 29 Apr 2025 17:04 |                       | SY/MD    | Ok       |
| 17  | Q1665-07DL  | IA-DUP-1-3242025DL  | VL042451.D     | 29 Apr 2025 17:37 |                       | SY/MD    | Ok,M     |
| 18  | Q1665-19DL  | IA-B2-2-3242025DL   | VL042452.D     | 29 Apr 2025 18:09 |                       | SY/MD    | Ok,M     |

Instrument ID: MSVOA\_L

**Daily Analysis Runlog For Sequence/QC Batch ID # VL042925**

| Review By                | Semsettin Yesilyurt | Review On            | 4/30/2025 2:23:53 PM |                      |               |
|--------------------------|---------------------|----------------------|----------------------|----------------------|---------------|
| Supervise By             | Mahesh Dadoda       | Supervise On         | 4/30/2025 2:28:33 PM |                      |               |
| SubDirectory             | VL042925            | HP Acquire Method    | TO15AIR.M            | HP Processing Method | VL041725AIR.M |
| STD. NAME                |                     | STD REF.#            |                      |                      |               |
| Tune/Reschk              |                     | AP2602               |                      |                      |               |
| Initial Calibration Stds |                     | AP2598,AP2599,AP2600 |                      |                      |               |
| CCC                      |                     | AP2598               |                      |                      |               |
| Internal Standard/PEM    |                     | AP2602               |                      |                      |               |
| ICV/I.BLK                |                     | AP2601               |                      |                      |               |
| Surrogate Standard       |                     |                      |                      |                      |               |
| MS/MSD Standard          |                     |                      |                      |                      |               |
| LCS Standard             |                     |                      |                      |                      |               |

|    |             |                   |            |                   |                         |       |        |
|----|-------------|-------------------|------------|-------------------|-------------------------|-------|--------|
| 19 | Q1665-21DL  | IA-B2-1-3242025DL | VL042453.D | 29 Apr 2025 18:42 |                         | SY/MD | Ok,M   |
| 20 | Q1908-01DUP | B-1DUP            | VL042454.D | 29 Apr 2025 19:57 |                         | SY/MD | Ok,M   |
| 21 | Q1908-02DL  | B-2DL             | VL042455.D | 29 Apr 2025 20:30 | not required ,not match | SY/MD | Not Ok |
| 22 | Q1908-02    | B-2               | VL042456.D | 29 Apr 2025 21:58 |                         | SY/MD | Ok,M   |

M : Manual Integration