

SOP ID : M1311-TCLP-15

SDG No : N/A

Weigh By : N/A

Balance ID : N/A

pH Meter ID : WC PH METER-1

Extraction By : JP

Filter By : JP

Pipette ID : N/A

Tumbler ID : N/A

TCLP Filter ID : 115525

Start Prep Date : N/A Time : N/A

End Prep Date : N/A Time : N/A

Combination Ratio : N/A

ZHE Cleaning Batch : N/A 78

Initial Room Temperature: N/A

Final Room Temperature: N/A

TCLP Technician Signature : 78

Supervisor By : 12

Standardized Name	MLS USED	STD REF. # FROM LOG
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A

Chemical Used	ML/SAMPLE U	Lot Number
N/A	N/A	N/A
N/A	N/A	N/A
HNO3-TCLP, 1N	N/A	WP112799
pH Strips	N/A	W1931, W1934, W3171, W3172
pH Strips	W1940, W1941, W1942	W3166, W1938, W1939,
N/A	N/A	N/A
120ml Plastic bottle	N/A	2738
1:1 HNO3	N/A	MP84041

Extraction Conformance/Non-Conformance Comments:

Matrix spikes are added after filtration and tumbling. Q1913-02 IS USED FOR MS-MSD.

Date / Time	Prepped Sample Relinquished By/Location	Received By/Location
04/20/15 13:00	JA HCL Room	S/S R3/Ext
	Preparation Group	Analysis Group

1/10/14 P. J.