

Instrument ID: BNA_P

Daily Analysis Runlog For Sequence/QC Batch ID # BP050125

Review By	Rahul	Review On	5/2/2025 10:56:34 AM
Supervise By	Jagrut	Supervise On	5/2/2025 5:06:12 PM
SubDirectory	BP050125	HP Acquire Method	BNA_P
		HP Processing Method	bp041425

STD. NAME	STD REF.#
Tune/Reschk	SP6757
Initial Calibration Stds	SP6722,SP6723,SP6724,SP6725,SP6726,SP6727,SP6728,SP6729
CCC	SP6725
Internal Standard/PEM	S12661,10ul/1000ul sample
ICV/I.BLK	SP6770
Surrogate Standard	
MS/MSD Standard	
LCS Standard	

Sr#	SampleID	ClientID	Data File Name	Date-Time	Comment	Operator	Status
1	DFTPP	DFTPP	BP024485.D	01 May 2025 09:57		RC/JU	Ok
2	SSTDCCC040	SSTDCCC040	BP024486.D	01 May 2025 10:37		RC/JU	Ok,M
3	PB167805TB	PB167805TB	BP024487.D	01 May 2025 11:18		RC/JU	Ok
4	Q1905-04	MH-G	BP024488.D	01 May 2025 12:03		RC/JU	Ok
5	Q1913-01	WC-12-S-202504	BP024489.D	01 May 2025 12:44		RC/JU	Ok
6	Q1913-02	WC-12-A-202504	BP024490.D	01 May 2025 13:25		RC/JU	Ok
7	Q1902-02DL	343DL	BP024491.D	01 May 2025 14:05		RC/JU	Ok,M
8	Q1912-08	MH-F	BP024492.D	01 May 2025 14:46		RC/JU	Ok
9	Q1906-12	WC-6	BP024493.D	01 May 2025 15:32		RC/JU	Ok
10	Q1906-16	WC-7	BP024494.D	01 May 2025 16:13		RC/JU	Ok
11	Q1912-04	MH-E	BP024495.D	01 May 2025 16:53		RC/JU	Ok
12	Q1905-08	MH-H	BP024496.D	01 May 2025 17:34		RC/JU	Ok
13	Q1906-04	WC-4	BP024497.D	01 May 2025 18:15		RC/JU	Ok
14	Q1906-08	WC-5	BP024498.D	01 May 2025 18:55		RC/JU	Ok
15	Q1913-03	WC-13-S-202504	BP024499.D	01 May 2025 19:36		RC/JU	Ok
16	Q1910-01	MOO-25-0123-27	BP024500.D	01 May 2025 20:17		RC/JU	Ok
17	Q1912-01	MH-E	BP024501.D	01 May 2025 20:58		RC/JU	Ok,M
18	Q1912-05	MH-F	BP024502.D	01 May 2025 21:38	Need 5X Dilution	RC/JU	Dilution

Instrument ID: BNA_P

Daily Analysis Runlog For Sequence/QC Batch ID # BP050125

Review By	Rahul	Review On	5/2/2025 10:56:34 AM		
Supervise By	Jagrut	Supervise On	5/2/2025 5:06:12 PM		
SubDirectory	BP050125	HP Acquire Method	BNA_P	HP Processing Method	bp041425
STD. NAME	STD REF.#				
Tune/Reschk	SP6757				
Initial Calibration Stds	SP6722,SP6723,SP6724,SP6725,SP6726,SP6727,SP6728,SP6729				
CCC	SP6725				
Internal Standard/PEM	S12661,10ul/1000ul sample				
ICV/I.BLK	SP6770				
Surrogate Standard					
MS/MSD Standard					
LCS Standard					

M : Manual Integration