

Instrument ID: BNA_P

Daily Analysis Runlog For Sequence/QC Batch ID # BP050225

Review By	Rahul	Review On	5/5/2025 10:12:48 AM
Supervise By	Jagrut	Supervise On	5/5/2025 11:50:04 AM
SubDirectory	BP050225	HP Acquire Method	BNA_P
		HP Processing Method	bp041425

STD. NAME	STD REF.#
Tune/Reschk	SP6757
Initial Calibration Stds	SP6722,SP6723,SP6724,SP6725,SP6726,SP6727,SP6728,SP6729
CCC	SP6725
Internal Standard/PEM	S12662,10ul/1000ul sample
ICV/I.BLK	SP6770
Surrogate Standard	
MS/MSD Standard	
LCS Standard	

Sr#	SampleID	ClientID	Data File Name	Date-Time	Comment	Operator	Status
1	DFTPP	DFTPP	BP024503.D	02 May 2025 10:47		RC/JU	Ok
2	SSTDCCC040	SSTDCCC040	BP024504.D	02 May 2025 11:27		RC/JU	Ok
3	PB167823BL	PB167823BL	BP024505.D	02 May 2025 12:08		RC/JU	Ok
4	PB167823BS	PB167823BS	BP024506.D	02 May 2025 12:49		RC/JU	Ok,M
5	Q1916-01	WC-12	BP024507.D	02 May 2025 13:35		RC/JU	Ok
6	Q1916-01MS	WC-12MS	BP024508.D	02 May 2025 14:15		RC/JU	Ok
7	Q1916-01MSD	WC-12MSD	BP024509.D	02 May 2025 14:56		RC/JU	Ok,M
8	Q1917-01	MH-JJ	BP024510.D	02 May 2025 15:38		RC/JU	Ok
9	Q1922-05	MH-S	BP024511.D	02 May 2025 16:19		RC/JU	Ok
10	Q1928-01	MH-Q	BP024512.D	02 May 2025 17:00		RC/JU	Ok
11	Q1874-04MS	RT1491MS	BP024513.D	02 May 2025 17:41	Already Analyzed with OK status	RC/JU	Not Ok
12	Q1874-04MSD	RT1491MSD	BP024514.D	02 May 2025 18:22	Already Analyzed with OK status	RC/JU	Not Ok
13	Q1932-01	COMP-4	BP024515.D	02 May 2025 19:03		RC/JU	Ok
14	Q1932-03	COMP-5	BP024516.D	02 May 2025 19:44		RC/JU	Ok
15	Q1932-05	COMP-6	BP024517.D	02 May 2025 20:25		RC/JU	Ok
16	Q1932-07	COMP-7	BP024518.D	02 May 2025 21:06		RC/JU	Ok
17	Q1912-05DL	MH-FDL	BP024519.D	02 May 2025 21:46		RC/JU	Ok
18	Q1914-07	SS-6	BP024520.D	02 May 2025 22:27		RC/JU	Ok,M

Instrument ID: BNA_P

Daily Analysis Runlog For Sequence/QC Batch ID # BP050225

Review By	Rahul	Review On	5/5/2025 10:12:48 AM		
Supervise By	Jagrut	Supervise On	5/5/2025 11:50:04 AM		
SubDirectory	BP050225	HP Acquire Method	BNA_P	HP Processing Method	bp041425
STD. NAME	STD REF.#				
Tune/Reschk	SP6757				
Initial Calibration Stds	SP6722,SP6723,SP6724,SP6725,SP6726,SP6727,SP6728,SP6729				
CCC	SP6725				
Internal Standard/PEM	S12662,10ul/1000ul sample				
ICV/I.BLK	SP6770				
Surrogate Standard					
MS/MSD Standard					
LCS Standard					

M : Manual Integration