

## Laboratory Control Sample/Laboratory Control Sample Duplicate Summary

SW-846

SDG No.: Q1915

Client: CDM Smith

Analytical Method: 8270E DataFile: BF142255.D

| Lab Sample ID | Parameter             | Spike | Result | Unit | Rec | RPD | Qual | RPD  |     | Limits |     |
|---------------|-----------------------|-------|--------|------|-----|-----|------|------|-----|--------|-----|
|               |                       |       |        |      |     |     |      | Qual | Low | High   | RPD |
| PB167810BS    | Pyridine              | 50    | 37.2   | ug/L | 74  |     |      |      | 29  | 97     |     |
|               | 1,4-Dichlorobenzene   | 50    | 39.6   | ug/L | 79  |     |      |      | 76  | 103    |     |
|               | 2-Methylphenol        | 50    | 41.0   | ug/L | 82  |     |      |      | 69  | 109    |     |
|               | 3+4-Methylphenols     | 50    | 40.8   | ug/L | 82  |     |      |      | 67  | 106    |     |
|               | Hexachloroethane      | 50    | 41.8   | ug/L | 84  |     |      |      | 76  | 118    |     |
|               | Nitrobenzene          | 50    | 45.8   | ug/L | 92  |     |      |      | 58  | 106    |     |
|               | Hexachlorobutadiene   | 50    | 41.5   | ug/L | 83  |     |      |      | 69  | 101    |     |
|               | 2,4,6-Trichlorophenol | 50    | 43.4   | ug/L | 87  |     |      |      | 61  | 110    |     |
|               | 2,4,5-Trichlorophenol | 50    | 44.9   | ug/L | 90  |     |      |      | 70  | 106    |     |
|               | 2,4-Dinitrotoluene    | 50    | 49.8   | ug/L | 100 |     |      |      | 60  | 115    |     |
|               | Hexachlorobenzene     | 50    | 41.6   | ug/L | 83  |     |      |      | 73  | 106    |     |
|               | Pentachlorophenol     | 100   | 93.4   | ug/L | 93  |     |      |      | 47  | 114    |     |