

Instrument ID: MSVOA_L

Daily Analysis Runlog For Sequence/QC Batch ID # VL050125

Review By	Semsettin Yesilyurt	Review On	5/2/2025 7:59:23 AM
Supervise By	Maresh Dadoda	Supervise On	5/5/2025 2:55:39 PM
SubDirectory	VL050125	HP Acquire Method	TO15AIR.M
		HP Processing Method	VL050125AIR.M

STD. NAME	STD REF.#
Tune/Reschk	AP2603
Initial Calibration Stds	AP2606,AP2607,AP2608
CCC	AP2606
Internal Standard/PEM	AP2603
ICV/I.BLK	AP2609
Surrogate Standard	
MS/MSD Standard	
LCS Standard	

Sr#	Sampled	ClientID	Data File Name	Date-Time	Comment	Operator	Status
1	BFB	BFB	VL042457.D	01 May 2025 08:19		SY/MD	Ok
2	VSTDICCC010	VSTDICCC010	VL042458.D	01 May 2025 09:03		SY/MD	Ok,M
3	VSTDICC002	VSTDICC002	VL042459.D	01 May 2025 09:37		SY/MD	Ok,M
4	VSTDICC001	VSTDICC001	VL042460.D	01 May 2025 10:08		SY/MD	Ok,M
5	VSTDICC0.5	VSTDICC0.5	VL042461.D	01 May 2025 10:40		SY/MD	Ok,M
6	VSTDICC0.1	VSTDICC0.1	VL042462.D	01 May 2025 11:12		SY/MD	Ok,M
7	VSTDICC0.03	VSTDICC0.03	VL042463.D	01 May 2025 11:44	Not used,wrong purge	SY/MD	Not Ok
8	VSTDICC0.03	VSTDICC0.03	VL042464.D	01 May 2025 12:15	Not used	SY/MD	Not Ok
9	VSTDICC015	VSTDICC015	VL042465.D	01 May 2025 12:49		SY/MD	Ok,M
10	VSTDICC0.03	VSTDICC0.03	VL042466.D	01 May 2025 14:10		SY/MD	Ok,M
11	VSTDICV010	ICVVL050125	VL042467.D	01 May 2025 15:46		SY/MD	Ok,M

M : Manual Integration

Instrument ID: MSVOA_L

Daily Analysis Runlog For Sequence/QC Batch ID # VL050525

Review By	Semsettin Yesilyurt	Review On	5/6/2025 10:22:30 AM
Supervise By	Mahesh Dadoda	Supervise On	5/6/2025 12:42:58 PM
SubDirectory	VL050525	HP Acquire Method	TO15AIR.M
		HP Processing Method	VL050125AIR.M

STD. NAME	STD REF.#
Tune/Reschk	AP2603
Initial Calibration Stds	AP2606,AP2607,AP2608
CCC	AP2606
Internal Standard/PEM	AP2603
ICV/I.BLK	AP2609
Surrogate Standard	
MS/MSD Standard	
LCS Standard	

Sr#	SampleID	ClientID	Data File Name	Date-Time	Comment	Operator	Status
1	BFB	BFB	VL042468.D	05 May 2025 07:45		SY/MD	Ok
2	VSTDCCC010	VSTDCCC010	VL042469.D	05 May 2025 08:29		SY/MD	Ok,M
3	VL0505ABL01	VL0505ABL01	VL042470.D	05 May 2025 09:19		SY/MD	Ok
4	VL0505ABS01	VL0505ABS01	VL042471.D	05 May 2025 10:07		SY/MD	Ok,M
5	Q1927-03DL	EBR-P04-SC-SS01-202	VL042472.D	05 May 2025 10:50	Not required	SY/MD	Not Ok
6	Q1927-04DL	EBR-P03-SC-SS01-202	VL042473.D	05 May 2025 11:22		SY/MD	Ok,M
7	Q1927-03	EBR-P04-SC-SS01-202	VL042474.D	05 May 2025 13:11		SY/MD	Ok,M
8	Q1927-04	EBR-P03-SC-SS01-202	VL042475.D	05 May 2025 14:44	Need 10X	SY/MD	Dilution
9	Q1927-01	EBR-P06-SC-IA01-202	VL042476.D	05 May 2025 15:18		SY/MD	Ok,M
10	Q1927-01DUP	EBR-P06-SC-IA01-202	VL042477.D	05 May 2025 15:53		SY/MD	Ok,M
11	Q1927-01DL	EBR-P06-SC-IA01-202	VL042478.D	05 May 2025 16:25	Not required	SY/MD	Not Ok
12	Q1927-02	EBR-P03-SC-IA01-202	VL042479.D	05 May 2025 16:59		SY/MD	Ok,M
13	Q1927-02DL	EBR-P03-SC-IA01-202	VL042480.D	05 May 2025 17:31	Not required	SY/MD	Not Ok
14	Q1927-05	EBR-P07-SC-IA01-202	VL042481.D	05 May 2025 18:06		SY/MD	Ok,M
15	Q1927-05DL	EBR-P07-SC-IA01-202	VL042482.D	05 May 2025 18:39	Not required	SY/MD	Not Ok
16	Q1927-06	EBR-P04-SC-IA01-202	VL042483.D	05 May 2025 19:15		SY/MD	Ok,M
17	Q1927-06DL	EBR-P04-SC-IA01-202	VL042484.D	05 May 2025 19:49	Not required	SY/MD	Not Ok
18	QC042825 CAN 10492	QC042825 CAN 10492	VL042485.D	05 May 2025 20:26		SY/MD	Ok

Instrument ID: MSVOA_L

Daily Analysis Runlog For Sequence/QC Batch ID # VL050525

Review By	Semsettin Yesilyurt	Review On	5/6/2025 10:22:30 AM		
Supervise By	Mahesh Dadoda	Supervise On	5/6/2025 12:42:58 PM		
SubDirectory	VL050525	HP Acquire Method	TO15AIR.M	HP Processing Method	VL050125AIR.M
STD. NAME	STD REF.#				
Tune/Reschk	AP2603				
Initial Calibration Stds	AP2606,AP2607,AP2608				
CCC	AP2606				
Internal Standard/PEM	AP2603				
ICV/I.BLK	AP2609				
Surrogate Standard					
MS/MSD Standard					
LCS Standard					

M : Manual Integration