

Surrogate Summary

SDG No.: Q1929

Client: ENTACT

Analytical Method: 8151A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Rec	Qual	Limits	
								Low	High
I.BLK-PS029918.D	PIBLK-PS029918.D	2,4-DCAA	1	500	433	87		70 (39)	130 (175)
		2,4-DCAA	2	500	433	87		70 (39)	130 (175)
I.BLK-PS030051.D	PIBLK-PS030051.D	2,4-DCAA	1	500	541	108		70 (39)	130 (175)
		2,4-DCAA	2	500	561	112		70 (39)	130 (175)
Q1929-03	WC-A4-02-C	2,4-DCAA	1	500	420	84		70 (39)	130 (175)
		2,4-DCAA	2	500	310	62	*	70 (39)	130 (175)
Q1929-07	WC-A1-03-C	2,4-DCAA	1	500	489	98		70 (39)	130 (175)
		2,4-DCAA	2	500	352	70		70 (39)	130 (175)
Q1929-11	WC-A1-04-C	2,4-DCAA	1	500	391	78		70 (39)	130 (175)
		2,4-DCAA	2	500	294	59	*	70 (39)	130 (175)
PB167871BL	PB167871BL	2,4-DCAA	1	500	620	124		70 (39)	130 (175)
		2,4-DCAA	2	500	616	123		70 (39)	130 (175)
PB167871BS	PB167871BS	2,4-DCAA	1	500	555	111		70 (39)	130 (175)
		2,4-DCAA	2	500	562	112		70 (39)	130 (175)
I.BLK-PS030058.D	PIBLK-PS030058.D	2,4-DCAA	1	500	533	107		70 (39)	130 (175)
		2,4-DCAA	2	500	558	112		70 (39)	130 (175)
PB167815TB	PB167815TB	2,4-DCAA	1	500	593	119		70 (39)	130 (175)
		2,4-DCAA	2	500	617	123		70 (39)	130 (175)
Q1901-08MS	B-167-SB01MS	2,4-DCAA	1	500	690	138	*	70 (39)	130 (175)
		2,4-DCAA	2	500	577	115		70 (39)	130 (175)
Q1901-08MSD	B-167-SB01MSD	2,4-DCAA	1	500	691	138	*	70 (39)	130 (175)
		2,4-DCAA	2	500	590	118		70 (39)	130 (175)
I.BLK-PS030069.D	PIBLK-PS030069.D	2,4-DCAA	1	500	537	107		70 (39)	130 (175)
		2,4-DCAA	2	500	563	113		70 (39)	130 (175)