

Instrument ID: ECD_L

Daily Analysis Runlog For Sequence/QC Batch ID # PL042325

Review By	Abdul	Review On	4/24/2025 8:49:02 AM
Supervise By	mohammad	Supervise On	4/26/2025 2:18:30 AM
SubDirectory	PL042325	HP Acquire Method	HP Processing Method pl042325 8081
STD. NAME	STD REF.#		
Tune/Reschk	PP24433,PP24095		
Initial Calibration Stds	PP24260,PP24261,PP24262,PP24269,PP24266,PP24267,PP24268,PP24269,PP24270,PP24271,PP24272,PP24273,PP24274,PP24275,PP24277,PP24278,PP24279,PP24280,PP24281,PP24282,PP24283,PP24284		
CCC	PP24261,PP24273,PP24279,PP24284		
Internal Standard/PEM			
ICV/I.BLK	PP24273,PP24279,PP24284		
Surrogate Standard			
MS/MSD Standard			
LCS Standard			

Sr#	Sampleld	Data File Name	Date-Time	Operator	Status
1	HEXANE	PL095325.D	23 Apr 2025 10:36	AR\AJ	Ok
2	I.BLK	PL095326.D	23 Apr 2025 10:50	AR\AJ	Ok
3	PEM	PL095327.D	23 Apr 2025 11:04	AR\AJ	Ok,M
4	RESCHK	PL095328.D	23 Apr 2025 11:18	AR\AJ	Ok,M
5	PSTDICC100	PL095329.D	23 Apr 2025 11:31	AR\AJ	Ok
6	PSTDICC075	PL095330.D	23 Apr 2025 11:45	AR\AJ	Ok
7	PSTDICC050	PL095331.D	23 Apr 2025 11:59	AR\AJ	Ok
8	PSTDICC025	PL095332.D	23 Apr 2025 12:12	AR\AJ	Ok,M
9	PSTDICC005	PL095333.D	23 Apr 2025 12:26	AR\AJ	Ok,M
10	PCHLORICC1000	PL095334.D	23 Apr 2025 12:40	AR\AJ	Ok
11	PCHLORICC750	PL095335.D	23 Apr 2025 12:53	AR\AJ	Ok
12	PCHLORICC500	PL095336.D	23 Apr 2025 13:07	AR\AJ	Ok
13	PCHLORICC250	PL095337.D	23 Apr 2025 13:21	AR\AJ	Ok
14	PCHLORICC050	PL095338.D	23 Apr 2025 13:51	AR\AJ	Ok,M
15	PTOXICC1000	PL095339.D	23 Apr 2025 14:04	AR\AJ	Ok
16	PTOXICC750	PL095340.D	23 Apr 2025 14:18	AR\AJ	Ok
17	PTOXICC500	PL095341.D	23 Apr 2025 14:32	AR\AJ	Ok
18	PTOXICC250	PL095342.D	23 Apr 2025 14:45	AR\AJ	Ok
19	PTOXICC100	PL095343.D	23 Apr 2025 14:59	AR\AJ	Ok,M
20	PSTDICV050	PL095344.D	23 Apr 2025 15:13	AR\AJ	Ok,M
21	PCHLORICV500	PL095345.D	23 Apr 2025 16:10	AR\AJ	Ok

Instrument ID: ECD_L

Daily Analysis Runlog For Sequence/QC Batch ID # PL042325

Review By	Abdul	Review On	4/24/2025 8:49:02 AM
Supervise By	mohammad	Supervise On	4/26/2025 2:18:30 AM
SubDirectory	PL042325	HP Acquire Method	HP Processing Method pl042325 8081
STD. NAME	STD REF.#		
Tune/Reschk	PP24433,PP24095		
Initial Calibration Stds	PP24260,PP24261,PP24262,PP24269,PP24266,PP24267,PP24268,PP24269,PP24270,PP24271,PP24272,PP24273,PP24274,PP24275,PP24277,PP24278,PP24279,PP24280,PP24281,PP24282,PP24283,PP24284		
CCC	PP24261,PP24273,PP24279,PP24284		
Internal Standard/PEM			
ICV/I.BLK	PP24273,PP24279,PP24284		
Surrogate Standard			
MS/MSD Standard			
LCS Standard			

22	PTOXICV500	PL095346.D	23 Apr 2025 16:24	AR\AJ	Ok
23	I.BLK	PL095347.D	23 Apr 2025 16:38	AR\AJ	Ok
24	PEM	PL095348.D	23 Apr 2025 16:52	AR\AJ	Ok,M
25	PSTDCCC050	PL095349.D	23 Apr 2025 17:05	AR\AJ	Ok,M
26	PB167688BL	PL095350.D	23 Apr 2025 17:19	AR\AJ	Ok
27	PB167688BS	PL095351.D	23 Apr 2025 17:33	AR\AJ	Not Ok
28	Q1848-01	PL095352.D	23 Apr 2025 17:46	AR\AJ	Ok,M
29	Q1848-03	PL095353.D	23 Apr 2025 18:00	AR\AJ	Ok,M
30	Q1848-03MS	PL095354.D	23 Apr 2025 18:14	AR\AJ	Ok,M
31	Q1848-03MSD	PL095355.D	23 Apr 2025 18:28	AR\AJ	Ok,M
32	Q1848-05	PL095356.D	23 Apr 2025 18:41	AR\AJ	Ok,M
33	Q1850-01	PL095357.D	23 Apr 2025 18:55	AR\AJ	Ok,M
34	I.BLK	PL095358.D	23 Apr 2025 19:23	AR\AJ	Ok,M
35	PSTDCCC050	PL095359.D	23 Apr 2025 19:50	AR\AJ	Ok,M

M : Manual Integration

Instrument ID: ECD_L

Daily Analysis Runlog For Sequence/QC Batch ID # PL050625

Review By		Review On	
Supervise By		Supervise On	
SubDirectory	PL050625	HP Acquire Method	HP Processing Method pl042325 8081
STD. NAME	STD REF.#		
Tune/Reschk	PP24433,PP24095		
Initial Calibration Stds	PP24260,PP24261,PP24262,PP24269,PP24266,PP24267,PP24268,PP24269,PP24270,PP24271,PP24272,PP24273,PP24274,PP24275,PP24277,PP24278,PP24279,PP24280,PP24281,PP24282,PP24283,PP24284		
CCC	PP24261,PP24273,PP24279,PP24284		
Internal Standard/PEM			
ICV/I.BLK	PP24273,PP24279,PP24284		
Surrogate Standard			
MS/MSD Standard			
LCS Standard			

Sr#	SampleId	Data File Name	Date-Time	Operator	Status
1	I.BLK	PL095528.D	05 May 2025 08:14	AR\AJ	Ok
2	PEM	PL095529.D	05 May 2025 09:21	AR\AJ	Ok,NS
3	PSTDCCC050	PL095530.D	05 May 2025 09:35	AR\AJ	Ok,NS
4	PB167854BL	PL095531.D	05 May 2025 12:37	AR\AJ	Ok,NS
5	PB167854BS	PL095532.D	05 May 2025 13:25	AR\AJ	Ok,NS
6	Q1946-01	PL095533.D	05 May 2025 13:43	AR\AJ	Ok
7	Q1947-01	PL095534.D	05 May 2025 13:57	AR\AJ	Ok
8	Q1947-05	PL095535.D	05 May 2025 14:10	AR\AJ	Ok
9	Q1951-01	PL095536.D	05 May 2025 14:24	AR\AJ	Ok
10	Q1952-01	PL095537.D	05 May 2025 14:38	AR\AJ	Ok,NS
11	Q1936-01	PL095538.D	05 May 2025 14:52	AR\AJ	Ok
12	Q1936-01MS	PL095539.D	05 May 2025 15:06	AR\AJ	Ok,NS
13	Q1936-01MSD	PL095540.D	05 May 2025 15:19	AR\AJ	Ok,NS
14	Q1936-03	PL095541.D	05 May 2025 16:01	AR\AJ	Ok
15	Q1936-05	PL095542.D	05 May 2025 16:15	AR\AJ	Ok,NS
16	Q1936-07	PL095543.D	05 May 2025 16:28	AR\AJ	Ok
17	Q1937-01	PL095544.D	05 May 2025 16:42	AR\AJ	Ok
18	I.BLK	PL095545.D	05 May 2025 16:56	AR\AJ	Ok,NS
19	PEM	PL095546.D	05 May 2025 18:25	AR\AJ	Ok,NS
20	PSTDCCC050	PL095547.D	05 May 2025 18:39	AR\AJ	Ok,NS
21	Q1937-03	PL095548.D	05 May 2025 19:07	AR\AJ	Ok

Instrument ID: ECD_L

Daily Analysis Runlog For Sequence/QCBatch ID # PL050625

Review By		Review On	
Supervise By		Supervise On	
SubDirectory	PL050625	HP Acquire Method	HP Processing Method pl042325 8081
STD. NAME	STD REF.#		
Tune/Reschk	PP24433,PP24095		
Initial Calibration Stds	PP24260,PP24261,PP24262,PP24269,PP24266,PP24267,PP24268,PP24269,PP24270,PP24271,PP24272,PP24273,PP24274,PP24275,PP24277,PP24278,PP24279,PP24280,PP24281,PP24282,PP24283,PP24284		
CCC	PP24261,PP24273,PP24279,PP24284		
Internal Standard/PEM			
ICV/I.BLK	PP24273,PP24279,PP24284		
Surrogate Standard			
MS/MSD Standard			
LCS Standard			

22	Q1937-05	PL095549.D	05 May 2025 19:21	AR\AJ	Ok
23	Q1937-07	PL095550.D	05 May 2025 19:34	AR\AJ	Ok,NS
24	Q1937-09	PL095551.D	05 May 2025 19:48	AR\AJ	Ok,NS
25	Q1937-11	PL095552.D	05 May 2025 20:02	AR\AJ	Ok,NS
26	Q1938-01	PL095553.D	05 May 2025 20:16	AR\AJ	Ok,NS
27	Q1938-03	PL095554.D	05 May 2025 20:29	AR\AJ	Not Ok
28	Q1938-05	PL095555.D	05 May 2025 20:43	AR\AJ	Ok,NS
29	Q1938-07	PL095556.D	05 May 2025 20:57	AR\AJ	Ok,NS
30	I.BLK	PL095557.D	05 May 2025 21:11	AR\AJ	Ok
31	PSTDCCC050	PL095558.D	05 May 2025 21:25	AR\AJ	Ok,NS

M : Manual Integration