

Instrument ID: FID_C

Daily Analysis Runlog For Sequence/QC Batch ID # FC040925AL

Review By	yogesh	Review On	4/8/2025 3:04:36 PM
Supervise By	mohammad	Supervise On	4/10/2025 5:31:03 AM
SubDirectory	FC040925AL	HP Acquire Method	HP Processing Method FC040925AL
STD. NAME	STD REF.#		
Tune/Reschk Initial Calibration Stds	PP24170,PP24175,PP24176,PP24177,PP24178		
CCC Internal Standard/PEM	PP24176		
ICV/I.BLK Surrogate Standard MS/MSD Standard LCS Standard	PP24174,PP24179		

Sr#	SampleId	Data File Name	Date-Time	Operator	Status
1	MECL2	FC068513.D	08 Apr 2025 13:05	YP/AJ	Ok
2	I.BLK	FC068514.D	08 Apr 2025 13:42	YP/AJ	Ok
3	100 PPM ALIPHATIC HC STD1	FC068515.D	08 Apr 2025 14:19	YP/AJ	Ok
4	50 PPM ALIPHATIC HC STD2	FC068516.D	08 Apr 2025 14:56	YP/AJ	Ok
5	20 PPM ALIPHATIC HC STD3	FC068517.D	08 Apr 2025 15:34	YP/AJ	Ok
6	10 PPM ALIPHATIC HC STD4	FC068518.D	08 Apr 2025 16:12	YP/AJ	Ok
7	5 PPM ALIPHATIC HC STD5	FC068519.D	08 Apr 2025 16:48	YP/AJ	Ok
8	20 PPM ALIPHATIC HC STD ICV	FC068520.D	08 Apr 2025 17:26	YP/AJ	Ok
9	I.BLK	FC068521.D	08 Apr 2025 18:41	YP/AJ	Ok
10	20 PPM ALIPHATIC HC STD	FC068522.D	08 Apr 2025 19:18	YP/AJ	Ok

M : Manual Integration

Instrument ID: FID_C

Daily Analysis Runlog For Sequence/QC Batch ID # FC050225AL

Review By	yogesh	Review On	5/2/2025 10:00:15 AM
Supervise By	mohammad	Supervise On	5/6/2025 4:46:28 AM
SubDirectory	FC050225AL	HP Acquire Method	HP Processing Method FC040925AL
STD. NAME	STD REF.#		
Tune/Reschk Initial Calibration Stds	PP24170,PP24175,PP24176,PP24177,PP24178		
CCC Internal Standard/PEM ICV/I.BLK Surrogate Standard MS/MSD Standard LCS Standard	PP24176 PP24174,PP24179		

Sr#	SampleId	Data File Name	Date-Time	Operator	Status
1	MECL2	FC068740.D	02 May 2025 09:42	YP/AJ	Ok
2	I.BLK	FC068741.D	02 May 2025 10:17	YP/AJ	Ok
3	20 PPM ALIPHATIC HC STD	FC068742.D	02 May 2025 10:53	YP/AJ	Ok
4	PB167840BL	FC068743.D	02 May 2025 11:36	YP/AJ	Ok
5	PB167840BS	FC068744.D	02 May 2025 13:26	YP/AJ	Ok
6	PB167840BSD	FC068745.D	02 May 2025 14:03	YP/AJ	Ok
7	Q1901-06	FC068746.D	02 May 2025 14:40	YP/AJ	Ok
8	I.BLK	FC068747.D	02 May 2025 15:18	YP/AJ	Ok
9	20 PPM ALIPHATIC HC STD	FC068748.D	02 May 2025 15:56	YP/AJ	Ok,M
10	PB167833BL	FC068749.D	02 May 2025 16:33	YP/AJ	Ok
11	PB167833BS	FC068750.D	02 May 2025 17:11	YP/AJ	Ok
12	PB167833BSD	FC068751.D	02 May 2025 17:49	YP/AJ	Ok
13	Q1928-01	FC068752.D	02 May 2025 18:27	YP/AJ	Ok
14	Q1928-01D	FC068753.D	02 May 2025 19:05	YP/AJ	Ok
15	Q1928-01MS	FC068754.D	02 May 2025 19:42	YP/AJ	Ok
16	Q1928-01MSD	FC068755.D	02 May 2025 20:19	YP/AJ	Ok
17	Q1928-02	FC068756.D	02 May 2025 20:57	YP/AJ	ReRun
18	I.BLK	FC068757.D	02 May 2025 22:11	YP/AJ	Ok
19	20 PPM ALIPHATIC HC STD	FC068758.D	02 May 2025 22:47	YP/AJ	Ok

M : Manual Integration

Instrument ID: FID_C

Daily Analysis Runlog For Sequence/QC Batch ID # FC050525AL

Review By	yogesh	Review On	5/5/2025 11:42:23 AM
Supervise By	mohammad	Supervise On	5/7/2025 6:08:25 AM
SubDirectory	FC050525AL	HP Acquire Method	HP Processing Method FC040925AL
STD. NAME	STD REF.#		
Tune/Reschk Initial Calibration Stds	PP24170,PP24175,PP24176,PP24177,PP24178		
CCC Internal Standard/PEM	PP24176		
ICV/I.BLK Surrogate Standard MS/MSD Standard LCS Standard	PP24174,PP24179		

Sr#	SampleId	Data File Name	Date-Time	Operator	Status
1	MECL2	FC068759.D	05 May 2025 08:59	YP/AJ	Ok
2	I.BLK	FC068760.D	05 May 2025 09:35	YP/AJ	Ok
3	20 PPM ALIPHATIC HC STD	FC068761.D	05 May 2025 10:11	YP/AJ	Ok
4	Q1935-02	FC068762.D	05 May 2025 10:48	YP/AJ	Ok
5	Q1928-02	FC068763.D	05 May 2025 11:24	YP/AJ	Ok
6	Q1939-01	FC068764.D	05 May 2025 12:00	YP/AJ	Ok
7	I.BLK	FC068765.D	05 May 2025 12:36	YP/AJ	Ok
8	20 PPM ALIPHATIC HC STD	FC068766.D	05 May 2025 13:13	YP/AJ	Ok
9	PB167855BL	FC068767.D	05 May 2025 13:49	YP/AJ	Ok
10	PB167855BS	FC068768.D	05 May 2025 14:26	YP/AJ	Ok
11	PB167855BSD	FC068769.D	05 May 2025 15:02	YP/AJ	Ok
12	Q1937-02	FC068770.D	05 May 2025 15:38	YP/AJ	Ok,M
13	Q1937-02D	FC068771.D	05 May 2025 16:15	YP/AJ	Ok
14	Q1937-02MS	FC068772.D	05 May 2025 16:51	YP/AJ	Ok
15	Q1937-02MSD	FC068773.D	05 May 2025 17:27	YP/AJ	Ok
16	Q1937-04	FC068774.D	05 May 2025 18:04	YP/AJ	Ok
17	Q1937-06	FC068775.D	05 May 2025 18:40	YP/AJ	Ok
18	Q1937-08	FC068776.D	05 May 2025 19:17	YP/AJ	Ok
19	I.BLK	FC068777.D	05 May 2025 20:30	YP/AJ	Ok
20	20 PPM ALIPHATIC HC STD	FC068778.D	05 May 2025 21:06	YP/AJ	Ok,M

M : Manual Integration

Instrument ID: FID_E

Daily Analysis Runlog For Sequence/QC Batch ID # FE041725AL

Review By	yogesh	Review On	4/16/2025 2:56:14 PM
Supervise By	mohammad	Supervise On	4/18/2025 1:20:00 AM
SubDirectory	FE041725AL	HP Acquire Method	HP Processing Method FE041725AL
STD. NAME	STD REF.#		
Tune/Reschk Initial Calibration Stds	PP24170,PP24175,PP24176,PP24177,PP24178		
CCC Internal Standard/PEM	PP24176		
ICV/I.BLK Surrogate Standard MS/MSD Standard LCS Standard	PP24174,PP24179		

Sr#	SampleId	Data File Name	Date-Time	Operator	Status
1	MECL2	FE053385.D	16 Apr 2025 14:34	YP\AJ	Ok
2	I.BLK	FE053386.D	16 Apr 2025 15:04	YP\AJ	Ok
3	100 PPM ALIPHATIC HC STD1	FE053387.D	16 Apr 2025 15:35	YP\AJ	Ok
4	50 PPM ALIPHATIC HC STD2	FE053388.D	16 Apr 2025 16:35	YP\AJ	Ok
5	20 PPM ALIPHATIC HC STD3	FE053389.D	16 Apr 2025 17:05	YP\AJ	Ok
6	10 PPM ALIPHATIC HC STD4	FE053390.D	16 Apr 2025 17:36	YP\AJ	Ok
7	5 PPM ALIPHATIC HC STD5	FE053391.D	16 Apr 2025 18:06	YP\AJ	Ok
8	20 PPM ALIPHATIC HC STD ICV	FE053392.D	16 Apr 2025 18:36	YP\AJ	Ok
9	I.BLK	FE053393.D	16 Apr 2025 19:36	YP\AJ	Ok
10	20 PPM ALIPHATIC HC STD	FE053394.D	16 Apr 2025 20:06	YP\AJ	Ok

M : Manual Integration

Instrument ID: FID_E

Daily Analysis Runlog For Sequence/QC Batch ID # FE050225AL

Review By	yogesh	Review On	5/2/2025 10:03:08 AM
Supervise By	mohammad	Supervise On	5/6/2025 4:48:01 AM
SubDirectory	FE050225AL	HP Acquire Method	HP Processing Method FE041725AL
STD. NAME	STD REF.#		
Tune/Reschk Initial Calibration Stds	PP24170,PP24175,PP24176,PP24177,PP24178		
CCC Internal Standard/PEM ICV/I.BLK Surrogate Standard MS/MSD Standard LCS Standard	PP24176 PP24174,PP24179		

Sr#	SampleId	Data File Name	Date-Time	Operator	Status
1	MECL2	FE053616.D	02 May 2025 08:40	YP\AJ	Ok
2	I.BLK	FE053617.D	02 May 2025 09:10	YP\AJ	Ok
3	20 PPM ALIPHATIC HC STD	FE053618.D	02 May 2025 10:11	YP\AJ	Ok,M
4	Q1932-01	FE053619.D	02 May 2025 13:52	YP\AJ	Ok
5	Q1932-03	FE053620.D	02 May 2025 14:22	YP\AJ	Dilution
6	Q1932-05	FE053621.D	02 May 2025 14:53	YP\AJ	Ok,M
7	Q1932-07	FE053622.D	02 May 2025 15:23	YP\AJ	Ok
8	Q1933-01	FE053623.D	02 May 2025 15:53	YP\AJ	Ok
9	Q1933-02	FE053624.D	02 May 2025 16:23	YP\AJ	Ok,M
10	Q1935-01	FE053625.D	02 May 2025 16:53	YP\AJ	Ok
11	Q1935-02	FE053626.D	02 May 2025 17:24	YP\AJ	ReRun
12	Q1935-05	FE053627.D	02 May 2025 17:54	YP\AJ	Ok
13	Q1935-06	FE053628.D	02 May 2025 18:24	YP\AJ	Ok,M
14	Q1939-01	FE053629.D	02 May 2025 18:54	YP\AJ	ReRun
15	Q1939-03	FE053630.D	02 May 2025 19:24	YP\AJ	Ok
16	Q1939-05	FE053631.D	02 May 2025 19:55	YP\AJ	ReRun
17	I.BLK	FE053632.D	02 May 2025 20:55	YP\AJ	Ok
18	20 PPM ALIPHATIC HC STD	FE053633.D	02 May 2025 21:25	YP\AJ	Ok

M : Manual Integration

Instrument ID: FID_E

Daily Analysis Runlog For Sequence/QC Batch ID # FE050525AL

Review By	yogesh	Review On	5/5/2025 11:43:54 AM
Supervise By	mohammad	Supervise On	5/7/2025 6:08:32 AM
SubDirectory	FE050525AL	HP Acquire Method	HP Processing Method FE041725AL
STD. NAME	STD REF.#		
Tune/Reschk Initial Calibration Stds	PP24170,PP24175,PP24176,PP24177,PP24178		
CCC Internal Standard/PEM ICV/I.BLK Surrogate Standard MS/MSD Standard LCS Standard	PP24176 PP24174,PP24179		

Sr#	SampleId	Data File Name	Date-Time	Operator	Status
1	MECL2	FE053634.D	05 May 2025 09:10	YP\AJ	Ok
2	I.BLK	FE053635.D	05 May 2025 09:40	YP\AJ	Ok
3	20 PPM ALIPHATIC HC STD	FE053636.D	05 May 2025 11:14	YP\AJ	Ok,M
4	Q1939-05	FE053637.D	05 May 2025 11:44	YP\AJ	Ok
5	Q1932-03DL	FE053638.D	05 May 2025 12:14	YP\AJ	Not Ok
6	I.BLK	FE053639.D	05 May 2025 13:14	YP\AJ	Ok
7	20 PPM ALIPHATIC HC STD	FE053640.D	05 May 2025 13:44	YP\AJ	Ok
8	Q1936-02	FE053641.D	05 May 2025 14:14	YP\AJ	Ok
9	Q1936-04	FE053642.D	05 May 2025 14:45	YP\AJ	Ok
10	Q1936-06	FE053643.D	05 May 2025 15:14	YP\AJ	Ok
11	Q1936-08	FE053644.D	05 May 2025 15:45	YP\AJ	Ok
12	Q1937-10	FE053645.D	05 May 2025 16:15	YP\AJ	Ok
13	Q1937-12	FE053646.D	05 May 2025 16:45	YP\AJ	Ok
14	Q1938-02	FE053647.D	05 May 2025 17:16	YP\AJ	Ok
15	Q1938-04	FE053648.D	05 May 2025 17:46	YP\AJ	Ok
16	Q1938-06	FE053649.D	05 May 2025 18:16	YP\AJ	Ok
17	Q1938-08	FE053650.D	05 May 2025 18:46	YP\AJ	Ok
18	Q1946-01	FE053651.D	05 May 2025 19:16	YP\AJ	Ok
19	Q1946-02	FE053652.D	05 May 2025 19:46	YP\AJ	Ok
20	Q1947-01	FE053653.D	05 May 2025 20:16	YP\AJ	Ok
21	Q1947-02	FE053654.D	05 May 2025 20:47	YP\AJ	Ok

Instrument ID: FID_E

Daily Analysis Runlog For Sequence/QC Batch ID # FE050525AL

Review By	yogesh	Review On	5/5/2025 11:43:54 AM
Supervise By	mohammad	Supervise On	5/7/2025 6:08:32 AM
SubDirectory	FE050525AL	HP Acquire Method	HP Processing Method FE041725AL
STD. NAME	STD REF.#		
Tune/Reschk Initial Calibration Stds	PP24170,PP24175,PP24176,PP24177,PP24178		
CCC Internal Standard/PEM	PP24176		
ICV/I.BLK Surrogate Standard MS/MSD Standard LCS Standard	PP24174,PP24179		

22	Q1947-05	FE053655.D	05 May 2025 21:17	YP\AJ	Ok
23	Q1947-06	FE053656.D	05 May 2025 21:47	YP\AJ	Ok
24	I.BLK	FE053657.D	05 May 2025 22:47	YP\AJ	Ok
25	20 PPM ALIPHATIC HC STD	FE053658.D	05 May 2025 23:18	YP\AJ	Ok,M

M : Manual Integration