

Surrogate Summary

SDG No.: Q1956

Client: CDM Smith

Analytical Method: 8151A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Rec	Qual	Limits	
								Low	High
I.BLK-PS029918.D	PIBLK-PS029918.D	2,4-DCAA	1	500	433	87		39	175
		2,4-DCAA	2	500	433	87		39	175
I.BLK-PS030072.D	PIBLK-PS030072.D	2,4-DCAA	1	500	513	103		39	175
		2,4-DCAA	2	500	522	104		39	175
PB167891BL	PB167891BL	2,4-DCAA	1	500	619	124		10	141
		2,4-DCAA	2	500	603	121		10	141
PB167891BS	PB167891BS	2,4-DCAA	1	500	565	113		10	141
		2,4-DCAA	2	500	555	111		10	141
I.BLK-PS030080.D	PIBLK-PS030080.D	2,4-DCAA	1	500	508	102		39	175
		2,4-DCAA	2	500	510	102		39	175
Q1956-01	SB1-3-4	2,4-DCAA	1	500	399	80		10	141
		2,4-DCAA	2	500	390	78		10	141
Q1956-02	SB2-4-5	2,4-DCAA	1	500	407	81		10	141
		2,4-DCAA	2	500	377	75		10	141
Q1956-03MS	SB2-4-5MS	2,4-DCAA	1	500	509	102		10	141
		2,4-DCAA	2	500	503	101		10	141
Q1956-04MSD	SB2-4-5MSD	2,4-DCAA	1	500	478	96		10	141
		2,4-DCAA	2	500	471	94		10	141
Q1956-06	SB91-3-4	2,4-DCAA	1	500	446	89		10	141
		2,4-DCAA	2	500	450	90		10	141
PB167912BL	PB167912BL	2,4-DCAA	1	500	572	114		39	175
		2,4-DCAA	2	500	554	111		39	175
PB167912BS	PB167912BS	2,4-DCAA	1	500	572	114		39	175
		2,4-DCAA	2	500	576	115		39	175
PB167912BSD	PB167912BSD	2,4-DCAA	1	500	580	116		39	175
		2,4-DCAA	2	500	580	116		39	175
Q1956-07	FB-05022025	2,4-DCAA	1	500	587	117		39	175
		2,4-DCAA	2	500	574	115		39	175
I.BLK-PS030092.D	PIBLK-PS030092.D	2,4-DCAA	1	500	511	102		39	175
		2,4-DCAA	2	500	522	104		39	175