

Surrogate Summary

SDG No.: Q1956

Client: CDM Smith

Analytical Method: 8151A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Rec	Qual	Limits	
								Low	High
I.BLK-PS030124.D	PIBLK-PS030124.D	2,4-DCAA	1	500	419	84		39	175
		2,4-DCAA	2	500	418	84		39	175
I.BLK-PS030137.D	PIBLK-PS030137.D	2,4-DCAA	1	500	422	84		39	175
		2,4-DCAA	2	500	432	86		39	175
PB167968BL	PB167968BL	2,4-DCAA	1	500	540	108		39	175
		2,4-DCAA	2	500	511	102		39	175
PB167968BS	PB167968BS	2,4-DCAA	1	500	517	103		39	175
		2,4-DCAA	2	500	523	105		39	175
PB167851TB	PB167851TB	2,4-DCAA	1	500	825	165		39	175
		2,4-DCAA	2	500	462	92		39	175
Q1956-05	COMP1	2,4-DCAA	1	500	805	161		39	175
		2,4-DCAA	2	500	432	86		39	175
Q1956-05MS	COMP1MS	2,4-DCAA	1	500	779	156		39	175
		2,4-DCAA	2	500	439	88		39	175
Q1956-05MSD	COMP1MSD	2,4-DCAA	1	500	759	152		39	175
		2,4-DCAA	2	500	434	87		39	175
I.BLK-PS030145.D	PIBLK-PS030145.D	2,4-DCAA	1	500	432	86		39	175
		2,4-DCAA	2	500	440	88		39	175