

Surrogate Summary

SDG No.: Q1982

Client: CDM Smith

Analytical Method: 8151A

Lab Sample ID	Client ID	Parameter	Column	Spike	Result	Rec	Qual	Limits	
								Low	High
I.BLK-PS030124.D	PIBLK-PS030124.D	2,4-DCAA	1	500	419	84		39	175
		2,4-DCAA	2	500	418	84		39	175
I.BLK-PS030223.D	PIBLK-PS030223.D	2,4-DCAA	1	500	491	98		39	175
		2,4-DCAA	2	500	503	101		39	175
Q1982-01	TP-1	2,4-DCAA	1	500	372	74		10	141
		2,4-DCAA	2	500	362	72		10	141
Q1982-02	TP-2B	2,4-DCAA	1	500	377	75		10	141
		2,4-DCAA	2	500	354	71		10	141
Q1982-03	TP-3	2,4-DCAA	1	500	321	64		10	141
		2,4-DCAA	2	500	305	61		10	141
Q1982-04	TP-4	2,4-DCAA	1	500	300	60		10	141
		2,4-DCAA	2	500	288	58		10	141
Q1982-05	TP-5	2,4-DCAA	1	500	320	64		10	141
		2,4-DCAA	2	500	312	62		10	141
Q1982-06	TP-6	2,4-DCAA	1	500	356	71		10	141
		2,4-DCAA	2	500	324	65		10	141
Q1982-07	TP-8	2,4-DCAA	1	500	300	60		10	141
		2,4-DCAA	2	500	302	60		10	141
Q1982-08	TP-9	2,4-DCAA	1	500	386	77		10	141
		2,4-DCAA	2	500	369	74		10	141
I.BLK-PS030235.D	PIBLK-PS030235.D	2,4-DCAA	1	500	451	90		39	175
		2,4-DCAA	2	500	462	92		39	175
Q1982-08MS	TP-9MS	2,4-DCAA	1	500	493	99		10	141
		2,4-DCAA	2	500	458	92		10	141
Q1982-08MSD	TP-9MSD	2,4-DCAA	1	500	529	106		10	141
		2,4-DCAA	2	500	499	100		10	141
I.BLK-PS030247.D	PIBLK-PS030247.D	2,4-DCAA	1	500	519	104		39	175
		2,4-DCAA	2	500	531	106		39	175
I.BLK-PS030275.D	PIBLK-PS030275.D	2,4-DCAA	1	500	471	94		39	175
		2,4-DCAA	2	500	469	94		39	175
PB167996BL	PB167996BL	2,4-DCAA	1	500	473	95		10	141
		2,4-DCAA	2	500	450	90		10	141
PB167996BS	PB167996BS	2,4-DCAA	1	500	554	111		10	141
		2,4-DCAA	2	500	539	108		10	141
I.BLK-PS030287.D	PIBLK-PS030287.D	2,4-DCAA	1	500	463	93		39	175
		2,4-DCAA	2	500	462	92		39	175