

SOP ID : M1312-SPLP-10

SDG No : N/A

Weigh By : JP

Balance ID : WC SC-7

pH Meter ID : WC PH METER-1

Extraction By : JP

Filter By : JP

Pipette ID : WC

Tumbler ID : T-2

TCLP Filter ID : 115525

Start Prep Date : 05/09/2025 **Time :** 17:20

End Prep Date : 05/10/2025 **Time :** 11:35

Combination Ratio : 20

ZHE Cleaning Batch : ⁷⁸ N/A

Initial Room Temperature: 23 °C

Final Room Temperature: 22 °C

TCLP Technician Signature : ⁷⁸

Supervisor By : ¹²

Standardized Name	MLS USED	STD REF. # FROM LOG
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A

Chemical Used	ML/SAMPLE U	Lot Number
SPLP FLUID	WP112802	N/A
N/A	N/A	N/A
HNO3-TCLP, 1N	N/A	WP112799
pH Strips	N/A	W1931, W1934, W3171, W3172
pH Strips	N/A	W3166, W1938, W1939,
N/A	N/A	N/A
120ml Plastic bottle	N/A	2738
1:1 HNO3	N/A	MP84041

Extraction Conformance/Non-Conformance Comments:

Matrix spikes are added after filtration and before preservation. TUMBLER T-2 checked, 30 rpm. q1984-16 is used for MS-MSD.

Date / Time	Prepped Sample Relinquished By/Location	Received By/Location
05/12/25 11:30	⁸⁰ / Lab Room	S123. / INC
	Preparation Group	Analysis Group ¹²