

SOP ID :	M1312-SPLP-10	
SDG No :	N/A	Start Prep Date : 05/09/2025 Time : 17:20
Weigh By :	JP	End Prep Date : 05/10/2025 Time : 11:35
Balance ID :	WC SC-7	Combination Ratio : 20
pH Meter ID :	WC PH METER-1	ZHE Cleaning Batch : ⁷⁸ N/A
Extraction By :	JP	Initial Room Temperature: 23 °C
Filter By :	JP	Final Room Temperature: 22 °C
Pipette ID :	WC	TCLP Technician Signature : ⁷⁸
Tumbler ID :	T-2	Supervisor By : ¹²
TCLP Filter ID :	115525	

Standardized Name	MLS USED	STD REF. # FROM LOG
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A
N/A	N/A	N/A

Chemical Used	ML/SAMPLE U	Lot Number
SPLP FLUID	WP112802	N/A
N/A	N/A	N/A
HNO3-TCLP,1N	N/A	WP112799
pH Strips	N/A	W1931,W1934,W3171,W3172
pH Strips	N/A	W3166,W1938,W1939,
N/A	N/A	N/A
120ml Plastic bottle	N/A	2738
1:1 HNO3	N/A	MP84041

Extraction Conformance/Non-Conformance Comments:

Matrix spikes are added after filtration and before preservation. TUMBLER T-2 checked,30 rpm. q1984-16 is used for MS-MSD.

Date / Time	Prepped Sample Relinquished By/Location	Received By/Location
05/12/25 11:30	⁸⁰ / Lab Room	S125. / 124
	Preparation Group	Analysis Group ¹²⁴

Sample ID	ClientID	TCLP Vessel ID	Sample Wt (g)	Volume Extraction Fluid #1 (mL)	Multi phasic	Phase Miscible	Phases Combined	Final Leachate PH	Metals Leachate Adj. PH	Prep Pos
PB167933TB	LEB933	19	N/A	2000	N/A	N/A	N/A	4.24	1.0	T-2
Q1984-02	OU4-PCS-TC-33-050725	11	100.02	2000	N/A	N/A	N/A	5.6	1.5	T-2
Q1984-04	OU4-PCS-TC-34-050725	12	100.03	2000	N/A	N/A	N/A	5.5	1.0	T-2
Q1984-06	OU4-PCS-TC-35-050725	13	100.02	2000	N/A	N/A	N/A	5.6	1.0	T-2
Q1984-08	OU4-TS-24-050725	14	100.02	2000	N/A	N/A	N/A	4.5	1.5	T-2
Q1984-10	OU4-TS-25-050725	15	100.03	2000	N/A	N/A	N/A	4.0	1.5	T-2
Q1984-12	OU4-TS-26-050725	16	100.02	2000	N/A	N/A	N/A	4.5	1.0	T-2
Q1984-14	OU4-TS-27-050725	17	100.02	2000	N/A	N/A	N/A	4.0	1.5	T-2
Q1984-16	OU4-TS-28-050725	18	100.03	2000	N/A	N/A	N/A	4.5	1.0	T-2

SampleID	ClientID	Sample Weight (g)	Filter Weight (g)	Filtrate (mL)	Filter + Solid (After 100°C)	% solids	% Dry Solids
PB167933TB	LEB933	N/A	N/A	N/A	N/A	N/A	N/A
Q1984-02	OU4-PCS-TC-33-050725	N/A	N/A	N/A	N/A	100	N/A
Q1984-04	OU4-PCS-TC-34-050725	N/A	N/A	N/A	N/A	100	N/A
Q1984-06	OU4-PCS-TC-35-050725	N/A	N/A	N/A	N/A	100	N/A
Q1984-08	OU4-TS-24-050725	N/A	N/A	N/A	N/A	100	N/A
Q1984-10	OU4-TS-25-050725	N/A	N/A	N/A	N/A	100	N/A
Q1984-12	OU4-TS-26-050725	N/A	N/A	N/A	N/A	100	N/A
Q1984-14	OU4-TS-27-050725	N/A	N/A	N/A	N/A	100	N/A
Q1984-16	OU4-TS-28-050725	N/A	N/A	N/A	N/A	100	N/A

Hot Block ID : WC S-1 / WC S-2

Thermometer ID : FLASHPOINT

SampleID	ClientID	Sample Weight (g)	Volume DI Water (mL)	PH after 5 min stir	PH after 10 min stir	Extraction Fluid 1 or 2	pH Extraction Fluid
PB167933TB	LEB933	N/A	N/A	N/A	N/A	#1	4.24
Q1984-02	OU4-PCS-TC-33-050725	N/A	N/A	N/A	N/A	#1	4.24
Q1984-04	OU4-PCS-TC-34-050725	N/A	N/A	N/A	N/A	#1	4.24
Q1984-06	OU4-PCS-TC-35-050725	N/A	N/A	N/A	N/A	#1	4.24
Q1984-08	OU4-TS-24-050725	N/A	N/A	N/A	N/A	#1	4.24
Q1984-10	OU4-TS-25-050725	N/A	N/A	N/A	N/A	#1	4.24
Q1984-12	OU4-TS-26-050725	N/A	N/A	N/A	N/A	#1	4.24
Q1984-14	OU4-TS-27-050725	N/A	N/A	N/A	N/A	#1	4.24
Q1984-16	OU4-TS-28-050725	N/A	N/A	N/A	N/A	#1	4.24